EL CASO PARA PRO-FLEX® PIVOT

Less load, more dynamics™

INFORME DE ESTUDIOS CLÍNICOS - PRO-FLEX® PIVOT

Less load, more dynamics™ - Menos Carga, Mayor Dinámica

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

1. Struyf, Pieter A., et al. "The prevalence of osteoarthritis of the intact hip and knee among traumatic leg amputees." *Archives of physical medicine and rehabilitation* 90.3 (2009): 440-446.
2. Engsborg, J. R., et al. "Normative ground reaction force data for able-bodied and below-knee-amputee children during walking." *Journal of pediatric orthopedics* 13.2 (1992): 169-173.
3. Chen, A., et al. "The global economic cost of osteoarthritis: how the UK compares." *Arthritis* 2012 (2012).
4. Heitzmann DWW, et al; A novel prosthetic foot leads to increased ankle power and reduced sound side loads in trans-tibial amputees; Abstract, Oral Presentation at the AOPA National Assembly San Antonio, TX, USA, October 7-10, 2015;
5. HCUP: Healthcare Cost and Utilization Project, June 2015. Agency for Healthcare Research and Quality, Rockville, MD, U.S. Department of Health and Human Services. <http://hcupnet.ahrq.gov/HCU/Net.jsp> Agency for Healthcare Research and Quality
6. Ziegler-Graham, Kathryn, et al. "Estimating the prevalence of limb loss in the United States: 2005 to 2050." *Archives of physical medicine and rehabilitation* 89.3 (2008): 422-429. Miller, William C., et al. "The influence of falling, fear of falling, and balance confidence on prosthetic mobility and social activity among individuals with a lower extremity amputation." *Archives of physical medicine and rehabilitation* 82.9 (2001): 1238-1244.
7. Burger, Helena, C. R. T. Marincek, and Eli Isakov. "Mobility of persons after traumatic lower limb amputation." *Disability & Rehabilitation* 19.7 (1997): 272-277.
8. Geertzen JH, Bosmans JC, Van der Schans CP. Claimed walking distance of lower limb amputees. *Disabil Rehabil* 2005;27:101-4.
9. Nolan L, Wit A, Dudzinski K, Lees A, Lake M, Wychowanski M. Adjustments in gait symmetry with walking speed in trans-femoral and trans-tibial amputees. *Gait Posture*.2003;17(2):142-51 prosthetic knee. *Arch Phys Med Rehabil* 2007;88:207-17.
10. Burke MJ, Roman V, Wright V. Bone and joint changes in lower limb amputees. *Ann Rheum Dis*. 1978;37(3): 252-54.
11. Gailey R, Allen K, Castles J, Kucharik J, Roeder M. Review of secondary physical conditions associated with lower-limb amputation and long-term prosthesis use. *J Rehabil Res Dev* 2008;45(1):15-29.
12. Kulkarni J, Adams J, Thomas E, Silman A. Association between amputation, arthritis and osteopenia in British male war veterans with major lower limb amputations. *Clin. Rehabil.*, 12 (4) (1998), pp. 348-353
13. Melzer I, Yekutieli M, Sukenik S. Comparative study of osteoarthritis of the contralateral knee joint of male amputees who do and do not play volleyball. *J. Rheumatol.*, 28 (1) (2001), pp. 169-172
14. Nolan L, Wit A, Dudzinski K, Lees A, Lake M, Wychowanski M. Adjustments in gait symmetry with walking speed in trans-femoral and trans-tibial amputees. *Gait Posture*.2003;17(2):142-51 prosthetic knee. *Arch Phys Med Rehabil* 2007;88:207-17.
15. Hungerford D, Cockin J. Fate of the retained lower limb joints in World War II amputees. *J. Bone Jt. Surg.*, 57 (1975), p. 111
16. Norvell DC, Czerniecki JM, Reiber GE, Maynard C, Pecoraro JA, Weiss NS. The prevalence of knee pain and symptomatic knee osteoarthritis among veteran traumatic amputees and nonamputees. *Arch Phys Med Rehabil* 2005;86(3):487-93.
17. Lemaire ED, Fisher FR. Osteoarthritis elderly amputee gait. *Arch Phys Med Rehabil* 1994;75(10):1094-9.
18. Breakley J. Gait of unilateral trans-tibial amputees. *Orthot Prosthet.* 1976;30:17-24.
19. Murray MP, Mollinger LA, Sepic SB, Gardner GM, Linder MT. Gait patterns in above-knee amputee patients: Hydraulic swing control vs constant-friction knee components. *Arch Phys Med Rehabil.* 1983;64(8):339-45.
20. Engsborg JR, Lee AG, Tedford KG, Harder JA. Normative ground reaction force data for able-bodied and below knee amputee children during walking. *J Pediatr Orthop.* 1993;13(2):169-73.
21. Zernicke RF, Hoy MG, Whiting WC. Ground reaction forces and center of pressure patterns in the gait of children with amputation: Preliminary report. *Arch Phys Med Rehabil.* 1985;66(11):736-41.
22. Menard MR, McBride ME, Sanderson DJ, Murray D. Comparative biomechanical analysis of energy-storing prosthetic feet. *Arch Phys Med Rehabil.* 1992;73(5):451-58.
23. Schneider K, Hart T, Zernicke RF, Setoguchi Y, Oppenheim W. Dynamics of below-knee amputee child gait: SACH foot versus Flex foot. *J Biomech.* 1993;26(10): 1191-1204.
24. Suzuki K. Force plate study on the artificial limb gait. *J Jpn Orthop Assoc.* 1972;46:503-16.
25. Engsborg JR, Lee AG, Patterson JL, Harder JA. External loading comparisons between able-bodied and below knee amputee children during walking. *Arch Phys Med Rehabil.* 1991;72(9): 657-61
26. Centers for Disease Control and Prevention (CDC). Prevalence and impact of chronic joint symptoms—seven states, 1996. *MMWR Morb. Mortal. Wkly Rep.* 47, 345-351 (1998).
27. Dunlop, D. D., Manheim, L. M., Song, J. & Chang, R. W. Arthritis prevalence and activity limitations in older adults. *Arthritis Rheum.* 44, 212-221 (2001).
28. Hunter, D. J. Lower extremity osteoarthritis management needs a paradigm shift. *Br. J. Sports Med.* 45, 283-288 (2011).
29. Hootman, J. M. & Helmick, C. G. Projections of US prevalence of arthritis and associated activity limitations. *Arthritis Rheum.* 54, 226-229 (2006).
30. Perruccio, A. V., Power, J. D. & Badley, E. M. Revisiting arthritis prevalence projections—it's more than just the aging of the population. *J. Rheumatol.* 33, 1856-1862 (2006).
31. Hunter, David J., Deborah Schofield, and Emily Callander. "The individual and socioeconomic impact of osteoarthritis." *Nature Reviews Rheumatology* 10.7 (2014): 437-441.
32. Nho, Shane J., Steven M. Kymes, John J. Callaghan, and David T. Felson 2013, The Burden of Hip Osteoarthritis in the United States: Epidemiologic and Economic Considerations. The Journal of the American Academy of Orthopaedic Surgeons 21 Suppl 1: S1-6.
33. E. Levy, A. Ferme, D. Perocheau, and I. Bono, "Socioeconomic costs of osteoarthritis in France," *Revue du Rhumatisme*, vol. 60, no. 6, pp. 635-675, 1993.
34. Murray, C. J. et al. Disability-adjusted life years (DALYs) for 291 diseases and injuries in 21 regions, 1990-2010: a systematic analysis for the Global Burden of Disease Study 2010. *Lancet* 380, 2197-2223 (2013).
35. Centers for Disease Control and Prevention (CDC). Racial disparities in total knee replacement among Medicare enrollees— United States, 2000-2006. *MMWR Morb Mortal Wkly Rep* 2009;58(6):1338.
36. Osteoarthritis kneebracing – A health economic evaluation – USA, 2012. On file at Össur
37. Snyder, R. D., et al., The effect of five prosthetic feet on the gait and loading of the sound limb in dysvascular below-knee amputees. *J Rehabil Res Dev*, 1995. 32(4): p. 309-15.
38. Lehmann JF, Price R, Boswell-Bessette S, Dralle A, Questad K. Comprehensive analysis of dynamic elastic response feet: Seattle Ankle/Lite Foot versus SACH foot. *Archives of Physical Medicine & Rehabilitation* 1993;74(8):853-61.
39. Heitzmann DWW, et al; Evaluation of a novel prosthetic foot while walking on level ground, ascending and descending a ramp; *Gait & Posture* 42 (2015): S94-S95. Abstract, Oral Presentation at the ESMAC 24th annual Meeting Heidelberg, Germany, September 10-12, 2015;
40. Powers, Christopher M., et al. "Influence of prosthetic foot design on sound limb loading in adults with unilateral below-knee amputations." *Archives of physical medicine and rehabilitation* 75.7 (1994): 825-829.
41. Segal, Ava D., et al. "The effects of a controlled energy storage and return prototype prosthetic foot on transtibial amputee ambulation." *Human movement science* 31.4 (2012): 918-931.
42. Kuo, Arthur D. "The six determinants of gait and the inverted pendulum analogy: A dynamic walking perspective." *Human movement science* 26.4 (2007): 617-656.
43. Kuo, Arthur D., J. Maxwell Donelan, and Andy Ruina. "Energetic consequences of walking like an inverted pendulum: step-to-step transitions." *Exercise and sport sciences reviews* 33.2 (2005): 88-97.
44. Morgenroth, David C., et al. "The effect of prosthetic foot push-off on mechanical loading associated with knee osteoarthritis in lower extremity amputees." *Gait & posture* 34.4 (2011): 502-507.

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