ENGLISH

USES

Post surgery, shoulder, wrist, forearm, and elbow injuries, humeral fractures, and to support arm casts.

CAUTION

Patients should be instructed to inspect their shoulder, elbow, hand and this product before and after each use. Modifications to this product are the sole responsibility of the prescribing physician.

PHOTO CAPTIONS

- A. Shoulder strap
- B. Sling Pocket
- C. Contact closures
- D. D-ring

FITTING INSTRUCTIONS

1. Open the sling.
2. Place injured arm into the sling. Place elbow deep into the pocket of the sling at a 90° angle when possible. Ensure that the open end of sling supports the hand and fingers.
3. Place strap behind the neck and over the opposite shoulder of the injured arm. Position strap and adjust padding on shoulder for comfort.
4. Undo contact closure on the end of strap. Loop through d-rings and secure.
5. Adjust length of strap using buckle (circled) to give proper support. Adjustments may be necessary as swelling patterns change.

WASHING INSTRUCTIONS

1. The sling can be hand washed using mild detergent and cold water. Do not wash in washing machine. Do not use bleach.
2. Lay sling out flat on a towel and air dry. Do not place in dryer. Avoid direct heat and sun.

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www.ossur.com ©Copyright Össur 2019 IFU 1450_001 Rev. 2
Össur Deutschland GmbH