## DISPOSAL

The device and packaging must be disposed of in accordance with respective local or national environmental regulations.

## LIABILITY TEXT

Össur does not assume liability for the following:

- Device not maintained as instructed by the instructions for use.
- Device assembled with components from other manufacturers.
- Device used outside of recommended use condition, application, or environment.

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**MD** Medical Device

### GENERAL SAFETY INFORMATION

**WARNING**: In case of functional change or functional loss, the patient should stop using the device and contact a healthcare professional.

**WARNING**: If the device shows signs of damage or wear hindering its normal functions, the patient should stop using the device and contact a healthcare professional.

### APPLICATION INSTRUCTIONS

1. Position sling against rib cage. Bring straps straight up and over the patient's shoulders. Place affected arm through the sling (**Figure 1**).
2. Criss cross both straps against the patient's back (**Figure 2**).
3. Attach hook straps to the front of the sling in a criss cross motion supporting the affected arm (**Figure 3**).
4. Pull the end of the straps to adjust the elbow to 90° position.
5. Secure arm against the body by attaching the swathe strap to the sling. Place one end of the strap to the sling. Encircle the affected arm and bring around the patient's back. Secure loose end to the front of the sling (**Figure 4**).

**WASHING INSTRUCTIONS**: Join hook and loop. Machine wash in cold water (30°C/86°F) with mild detergent. Air dry.

### REPORT A SERIOUS INCIDENT

**Important notice to users and/or patients established in Europe**:

The user and/or patient must report any serious incident that has occurred in relation to the device to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

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©Copyright Össur 2020
**Össur Deutschland GmbH**