Instructions for Use

TECHFORM[®] PREMIUM
CASTING TAPE

EN		Instructions for Use	3
DE		Gebrauchsanweisung	6
FR		Notice d'utilisation	9
ES		Instrucciones para el uso	12
IT		Istruzioni per l'uso	15
NO		Bruksanvisning	18
DA		Brugsanvisning	21
SV		Bruksanvisning	24
FI		Käyttöohjeet	27
PT		Instruções de Utilização	30

ENGLISH

DESCRIPTION OF SYMBOLS

MD

Medical Device

Please read these instructions before use. The healthcare professional applying this casting material is responsible for the proper care of the patient. Only responsible, qualified, and experienced healthcare professionals should apply casting tape to patients.

INDICATIONS FOR USE:

TECHFORM® casting materials are to be used for casting, splinting, and external immobilization.

MATERIAL DESCRIPTION:

TECHFORM® casting materials consist of a knitted synthetic substrate impregnated with a water activated polyurethane resin.

PRODUCT STORAGE:

This product should be stored at room temperature. Storage at high temperatures will result in a reduced shelf life and tacky handling characteristics. In addition, care should be taken to avoid puncturing or tearing the pouch in storage as this will cause premature hardening in the package.

FLAMMABILITY:

Uncured polyurethane tape will burn over an open flame. Finished casts using TECHFORM® cast tape will not support combustion under normal conditions as defined by the NFPA Standard #702-1975 and the Flammable Fabrics Act (CS-191-53).

PRE-APPLICATION PRECAUTIONS:

- The healthcare professional applying this casting material is responsible for the
- proper care of the patient. Only responsible, qualified, and experienced healthcare
- practitioners should apply casting material to patients.
- The polyurethane resin will adhere to unprotected skin or clothing. Apply gloves
- before handling the casting material.
- Exposure of the resin to ambient air or water will begin the activation process and the
- tape will become rigid. Therefore, it is important that the cast tape pouch remain free
- from punctures until it is time to submerge the cast tape in water and apply the cast.
- As the material begins to set, it releases heat (exotherm).

The following are conditions that can increase heat to uncomfortable/injurious levels:

1. The more layers applied to the cast, the higher the exotherm, therefore, care should be taken when applying multiple layers of cast tape.

2. Activating the cast tape by submerging the roll in water over 80° Fahrenheit (27° Celsius) may result in a higher exotherm temperature.
 3. Applying any coverings such as wraps, dressings, and/or padding over a setting cast, or splint, or lying the cast on any surface such as a pillow during the setting process should be avoided. In addition, proper care should be exercised with children and patients with sensitive skin.
- Do not insert any materials in between the cast and the patient or padding.
 - Inform the patient that they should not get the cast wet. Skin maceration and irritation may result from prolonged and/or frequent wetting of the cast.

CAST APPLICATION:

- Immerse the cast tape roll in room temperature water. (68° F (20° C) to 75° F (24°C)) and firmly squeeze the roll three times while it is submerged in the water. Only open one pouch at a time.
- Traditional casting techniques are to be exercised when applying the casting material. 3-6 layers of cast tape are needed for proper immobilization. Weight bearing casts should consist of a minimum of four layers. (See heat precautions)
- The product should be rubbed well during application for ideal lamination of the cast tape layers. Set times vary between 3 to 5 minutes when the material is activated with room temperature water.
- Allow 30 minutes before allowing the patient to bear weight or loads on the cast. Removal, alteration, uni/bivalving, and repairing should be instituted with care. The applicator is responsible for any alteration and/or removal of the cast or splint.

CAST REMOVAL:

Use a standard cast saw for cast removal. Care should be taken to avoid bony prominences, as well as, excessive heat build up from the cast saw blades. Avoid exposed edges of the cast to prevent scratches.

PATIENT INSTRUCTIONS:

Instruct the patient of the following information to help insure patient comfort and safety, as well as, care for the cast integrity.

- Patients should be notified by their healthcare professional as to proper care for their injury, cast instructions, mobility instructions, and how to notify their healthcare professional should problems with the cast occur.
- Patients should be encouraged to not get the cast wet. The cast may be cleaned with a damp cloth if necessary. Should the cast become wet, use a towel to remove excess moisture and then use a blower dryer (such as a hair dryer on a cool setting) to sufficiently dry the undercast materials.
- Skin maceration and irritation may result from frequent wetting of the cast without proper drying.
- Do not insert any materials into the cast. Avoid activities that may result in dirt, sand, or other foreign articles getting under the cast.
- If necessary, rough edges of the cast may be reduced with a nail file.

- The patient should immediately contact their healthcare professional if they experience; painful pressure or rubbing from the cast, their toes or fingers become numb or bluish in color, or if the cast feels too tight. In addition, the patient should notify their healthcare professional if the cast softens or breaks at any time.

GENERAL SAFETY INFORMATION

WARNING: In cases of functional change or functional loss, the patient should stop using the device and contact a healthcare professional.

WARNING: If the device shows signs of damage or wear hindering its normal functions, the patient should stop using the device and contact a healthcare professional.

REPORT A SERIOUS INCIDENT

Important notice to users and/or patients established in Europe:

The user and/or patient must report any serious incident that has occurred in relation to the device to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

DISPOSAL

The device and packaging must be disposed of in accordance with respective local or national environmental regulations.

LIABILITY TEXT

Össur does not assume liability for the following:

- Device not maintained as instructed by the instructions for use.
- Device assembled with components from other manufacturers.
- Device used outside of recommended use condition, application, or environment.

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