A BRACE, ONE OF THOSE THINGS WILL DO.

Timothy Burke - cPO



Acknowledgement







Basics Classifications Clinical Considerations Putting into practice

The Basics

- Main aims for spinal orthoses
- Prevent uncontrolled movement
 - Control gross flexion
 - Maintain or prompt extension
 - Limit rotation
 - Prevent further fracture loss of height (LOS)



The Basics

- Main aims for spinal orthoses
- Prevent uncontrolled movement
 - Control gross flexion
 - Maintain or prompt extension
 - Limit rotation
 - Prevent further fracture loss of height (LOS)

- Influence Pain
- Prevent further loss of condition / function
- Enable function
 - Global mobility
 - Bed mobility
 - Compensate for weakness



The Basics

There multiple answers to the same question and no single correct answer

Evidence based practice

- Paucity of evidence
- Conflicted findings

The Evidence Gap for Assistive Devices

- Multiple known challenges to knowledge development and mobilisation (knowledge translation)
- Lemaire, E. D. (2016). "Mobilizing knowledge: The evidence gap for assistive devices." <u>Technology</u> <u>Innovation Management Review</u> 6(9).



Spinal Columns

- Three main columns of the spine
- Anterior, Middle and Posterior



Compression Fractures

- Mechanism of injury determines the presentation
 - Eg: Anterior or lateral flexion
 - Failure of anterior column with middle column acting as hinge.

Stability it dependent upon

- Number of columns involved
 - 1 column = minor injury
 - 3 column = significant injury
- Involvement of ligaments
- Loss of height and bony fragments



Spinal Columns

- Three main columns of the spine
- Anterior, Middle and Posterior

Stability it dependent upon

- Number of columns involved
 - 1 column = minor injury
 - 3 column = significant injury
- Involvement of ligaments
- Loss of height and bony fragments





Spinal Columns

- Three main columns of the spine
- Anterior, Middle and Posterior



Stability it dependent upon

- Number of columns involved
 - 1 column = minor injury
 - 3 column = significant injury
- Involvement of ligaments
- Loss of height and bony fragments

Compression Fractures

- Mechanism of injury determines the presentation
 - Eg: Anterior or lateral flexion
 - Failure of anterior column with middle column acting as hinge.



Stable Fractures without neurology

- Conservative management
- Pain management
- No management



Stable Fractures without neurology

- Conservative management
- Pain management
- No management

Unstable Fractures with neurology

Surgical management and internal fixation



Stable Fractures without neurology

- Conservative management
- Pain management
- No management

Stable Fractures / Unstable Fractures without neurology

- Surgical management
 - Involvement of 2/3 vertebral columns
 - >50% LOS
 - Angulation at thoraco-lumbar junction >20°
 - Multi-adjacent compression fractures

Conservative management

Unstable Fractures with neurology

Surgical management and internal fixation

Controversies in Spinal Trauma

Controversies in Spinal Trauma and Evolution of Care 🕮

James S. Harrop, MD, George N. Rymarczuk, MD, Alexander R. Vaccaro, MD, PhD, Michael P. Steinmetz, MD, Lindsay A. Tetreault, HBSc, Michael G. Fehlings, MD, PhD

Neurosurgery, Volume 80, Issue 3S, March 2017, Pages S23-S32,

https://doi.org/10.1093/neuros/nyw076

Published: 21 February 2017 Article history •

🔺 PDF 🛛 💵 Split View 🛛 🖌 Cite 🎾 Permissions 🛛 🤜 Share 🔻



Controversies in Spinal Trauma

- 1. corticosteroid therapy in acute spinal trauma,
- 2. odontoid fractures,
- 3. timing of surgery in central cord syndrome,
- 4. thoracolumbar burst fracture treatment,
- 5. thoracolumbar vertebral compression fractures (VCFs),
- 6. timing of surgery in spinal trauma,
- 7. blood pressure management goal,
- 8. timing of magnetic resonance imaging (MRI) in spinal trauma,
- 9. hypothermia in spinal trauma,
- 10. effect of global sagittal balance on clinical outcomes,
- 11. lumbar subarachnoid drainage,
- 12. minimally invasive surgery (MIS) strategies for thoracolumbar trauma,
- 13. diffusion MRI,
- 14. classification schemes in spinal trauma.

(Harrop, Rymarczuk et al. 2017)



Controversies in Spinal Trauma

1. corticosteroid therapy in acute spinal trauma,

2. odontoid fractures,

3. timing of surgery in central cord syndrome,

4. thoracolumbar burst fracture treatment,

5. thoracolumbar vertebral compression fractures (VCFs),

6. timing of surgery in spinal trauma,

7. blood pressure management goal,

8. timing of magnetic resonance imaging (MRI) in spinal trauma,

9. hypothermia in spinal trauma,

10. effect of global sagittal balance on clinical outcomes,

11. lumbar subarachnoid drainage,

12. minimally invasive surgery (MIS) strategies for thoracolumbar trauma,

13. diffusion MRI,

14. classification schemes in spinal trauma.

(Harrop, Rymarczuk et al. 2017)





www.aospine.aofoundation.org

Upper Cervical Injury Classification System	~
Subaxial Injury Classification System	~
Thoracolumbar Injury Classification System	~
Sacral Injury Classification System	~



- AO classification of spinal injuries
 - Three separate components to every fracture are considered, with only the first fully assessable on imaging alone.
 - **1**. morphology of the fracture
 - 2. presence of neurological signs
 - 3. presence of ligamentous injuries or co-morbid conditions (referred to as modifiers)



Morphology (A, B or C)

- Injuries are broadly categorized into three groups:
 - A. compression injuries
 - B. distraction injuries
 - C. displacement or dislocation



- Neurology (N0 Nx / +)
 - Neurology is divided into 7 categories
 - NO. Neurological intact
 - N1. Transcient neurologic deficit
 - N2. Radicular symptoms
 - N3. Incomplete spinal cord injury or any degree
 - of cauda equina injury
 - N4. Complete spinal cord injury

Nx. Cannot be examined

+. Continued spinal cord compression



Modifers (M1 – M4)

- Modifiers are spinal region specific
 - Cervical
 - M4. Vascular Injury
 - Thoracic
 - M2. Patient specific comorbidity (ankylosing spondylitis



AO Spine Thoracolumbar Injury Classification System



SPINE

© 2020 AO Foundation, AO Spine his work is licensed under the Creative Commons Attribution-NonCommercial-NoDerivatives 4.0 International license o view a copy of this license, visit http://creativecommons.org/licenses/by-nc-nd/4.0/

Further information: www.aospine.org/classification



Type B Distraction Injuries















B3 Hyperextension injury through the disc or vertebral body leading to a hyperextended position of the spinal column.







Commonly seen in ankylotic disorders. Anterior

A3 Incomplete burst

Fracture with any involvement of the posterior wall; only a single endplate fractured. Vertical fracture of the lamina is usually present and does not constitute a tension band failure.



A4 Complete burst

Fracture with any involvement of the posterior wall and both endplates. Vertical fracture of the lamina is usually present and does not constitute a tension band failure.



Type C Translation Injuries

Displacement or dislocation There are no subtypes because various configurations are possible due to dissociation/ dislocation. Can be combined with subtypes of A or B.



AO Spine Thoracolumbar Injury SPINE Classification System

Algorithm for morphologic classification



Neurology



Modifiers

Туре	Description
M1	This modifier is used to designate fractures with an indeterminate injury to the tension band based on spinal imaging with or without MRI. This modifier is important for designating those injuries with stable injuries from a bony standpoint for which ligamentous insufficiency may help determine whether operative stabilization is a consideration.
M2	Is used to designate a patient-specific comorbidity, which might argue either for or against surgery for patients with relative surgical indications. Examples of an M2 modifier include ankylosing spondylitis or burns affecting the skin overlying the injured spine.

Classification Nomenclature



Disclaimen

1. Vaccaro, A. R., C. Oner, C. K. Kepler, M. Dvorak, K. Schnake, C. Bellabarba, M. Reinhold, B. Aarabi, F. Kandziora, J. Chapman, R. Shanmuganathan, M. Fehlings, L. Vialle, A. O. S. C. Injury and F. Trauma Knowledge (2013). "AOSpine thoracolumbar spine injury classification system: fracture description, neurological status, and key modifiers." <u>Scine (Phila Pa 1976)</u> 38(23): 2028-2037.

 Kepler, C. K., A. R. Vaccaro, J. D. Koemer, M. F. Dvorak, F. Kandziora, S. Rajasekaran, B. Aarabi, L. R. Vialle, M. G. Fehlings, G. D. Schroeder, M. Reinhold, K. J. Schnake, C. Bellabarba and F. Cumhur Oner (2016), "Reliability analysis of the AOSpine thoracolumbar spine injury classification system by a worldwide group of naive spinal surgeons," <u>Eur</u> <u>Science 1</u>(e-pub)



© 2020 AO Foundation, AO Spine This work is licensed under the Creative Commons Attribution-NonCommercial-NoDerivatives 4.0 International License

To view a copy of this license, visit http://creativecommons.org/licenses/by-nc-nd/4.0/

Further information: www.aospine.org/classification

Know Your Goal







Prescription matrix

"Clinicians need to know the best type of [orthoses] to prescribe, for whom they should be prescribed, the optimal time to prescribe one, how long they should be used, the adverse effects, and the factors influencing acceptability and adherence to their use." (Tyson and Kent 2013)

	Collar	СТО	HALO	Antiflexion	TLSO
Flex / Extension	Mild - Mod	Mod - Sig	Sig	Mild - Mod	Mod - Sig
Rotation	Mild	Mild - Mod	Sig	N/A - Sig	Mod - Sig
Skin Care	Mild	Mod	Sig	N/A - Mild	Sig
Trache Care	Mod	Mod	Mild	N/A	N/A
Unique bits	Occian backs	Ability to add / subtract sections	Captain Complex	Captain Dependable	Custom vs Prefab



N/A - Not applicable, Mild - Mild, Mod - moderate, Sig - Significant.

Location, Location

- Who decides on prescription and how much say do they have?
 - O&P Department, Physiotherapy Department, Nursing etc...
 - Surgeons (Neurosurgeon vs Orthopaedic Surgeon)
- Historical prescription trends
 - Confidence and knowledge in products OR Cost?
- Facility priorities
 - Length of Stay (LOS), Pressure Injury prevention and Complex Care



Location, Location

- How much support does the facility have and for which aspect of care?
 - Acute vs Subacute / Rehab
 - Post Discharge vs Outpatient
- Is the facility linked to another organisation?
 - Officially / Non-officially
- Governance
 - Internal vs External oversight
 - Clinical guidelines



Clinical Guidelines

- Organisation Specific Guidelines
 - Organisations specific guidelines around spinal management, depending on the size and specialty of the service
 - These can cover device indications, timelines for management and clinical expectations
- Product Guidelines
 - Internal Ossur 'guidance' documents that indicated what fracture types / levels devices are appropriate for



Cervical Spine Immobilisation

Orthoses	Use	Timeframe
Patriot Extrication	Emergency Immobilisation	< 6/24
Philadelphia Collar	Short term Immobilisation	< 3/7
Miami J / Miami J Select	Medium / Long term Immobilisation	+/- 12/52

Additional considerations:

- Ongoing investigations
- Progression of immobilisation
- Education and carer support
- Skin care



Cervico-Thoraic Spine Immobilisation

Orthoses	Use	Timeframe
Miami JTO / Minerva	Long term Immobilisation	+/- 12/52
Halo Thoracic Orthosis	Long term Immobilisation	+/- 12/52
CTLSO	Long term Immobilisation	+/- 12/52

- Additional considerations:
 - Long term care options and follow up
 - Education and carer support
 - Skin care
 - Pressure injury risk level



Thoraco-Lumbar Spine Immobilisation

Orthoses	Use	Timeframe
Anti-flexion: Jewett/Cash	Single Anterior Column Fractures	+/- 12/52
Boston Overlap Brace*	Multiple Column Fractures	+/- 12/52
Miami LSO / TLSO*	Multiple Column Fractures	+/- 12/52

Additional considerations:

- Ability to don/doff independently
- Education and carer support
- Skin care and pressure injury risk level
- Burden of care and discharge destination





Personal Care

Clarify showering restriction at the outset, on for showering or sit to fit / off for showering*

Cervical Collars:

- Shower daily to 2nd/3rd daily activity dependent
- Change liners following shower
- Shave following shower with patient supine and head held
- Clean/wash liners daily but they must be dry before being re-used

Halo's and CTO:

- Shower and liner change during Orthotic review only (fortnightly)
- Very tiring for patients, so consider dry liner change if patients do not have necessary physical endurance.



Personal Care

Clarify showering restriction at the outset, on for showering or sit to fit / off for showering*

TLSO/LSO:

- Shower daily to 2nd/3rd daily
- Typically will not have liners and can be donned directly on skin of showers
- Aim to shower in the evening to enable the straps to dry overnight
- Consider provision of second device for showering



To Who, When, How

- Education is a core part of spinal management
- Things that seem straight forward, might not also be
 - The arrow points up

Discharge

- Where is the patient going? Who we be there?
- Who else needs to know?
 - Secondary discharge



To Who, When, How

- Education is a core part of spinal management
- Things that seem straight forward, might not also be



- To Who, When, How
 - Education is a core part of spinal management
 - Things that seem straight forward, might not also be
 - The arrow points UP!



Image source: Netflix (2018). New Amsterdam



- To Who, When, How
 - Education is a core part of spinal management
 - Things that seem straight forward, might not also be
 - The arrow points up

Discharge

- Where is the patient going? Who we be there?
- Who else needs to know?
 - Secondary discharge



Image source: https://www.hartmann.info/en-cn/-/media/



Complications





Complications

"Was it that colour before?"





Complications







"The room smelt funny and we couldn't figure out why..."

Complications — when things aren't right

Pain

- Un-resolving pain
- Pain 'spiking' to an abnormal level during mobility

Neurology

- Previously documented?
- Previously investigated?
- Changing neurology

What to do?

- Escalate to treating team
- Liaise with treating spinal team



Reason for admission:

 Patient painting 2nd storey of house, ladder slipped (propped against house from carport), patient fell in gap. Estimated fall 3-5 meters.













XR - Left Knee AP



Reason for admission:

 Patient painting 2nd storey of house, ladder slipped (propped against house from carport), patient fell in gap. Estimated fall 3-5 meters.

Injury:

- Fracture L1, Superior End Plate (SEP) AO: L1 A1, NO
- Open Left fracture proximal tibia / fibular

Subjective:

Prolonged NWB Left L/L.



Reason for admission:

 Patient painting 2nd storey of house, ladder slipped (propped against house from carport), patient fell in gap. Estimated fall 3-5 metres.

Injury:

- Fracture L1, Superior End Plate (SEP) AO: L1 A1, (NO)
- Open Left fracture proximal tibia / fibular

Management:

TLSO - Antiflexion, Jewett





Reason for admission:

 Pt sitting on stool (two step) and fell onto floor. Direct axial load and radicular pain for last 3/7. Presented to ED due to pain and reduced mobility.







CT - AP View

CT - Sagittal View



Reason for admission:

 Pt sitting on stool (two step) and fell onto floor. Direct axial load and radicular pain for last 3/7. Presented to ED due to pain and reduced mobility.

Injury:

- Acute/subacute Fracture T10 vertebral body with widening and gas within fracture defect AO: T10 A3 (NO; M2)
- Multiple additional subacute rib fractures

Subjective / Objective:

- Third hospital admission in 6 months, previous presentations for decreasing mobility, SOB and poor glycemic control (BSL 2.7 on past admission)
- Previously home alone





CT - AP Scout







CASE STUDY 2 - 74YO MALE, 101KG

Reason for admission:

 Pt sitting on stool (two step) and fell onto floor. Direct axial load and radicular pain for last 3/7. Presented to ED due to pain and reduced mobility.

Injury:

- Acute/subacute Fracture T10 vertebral body with widening and gas within fracture defect AO: T10 A3 (NO; M2)
- Multiple additional subacute rib fractures
- Large body habitus

Management:

- Modified Bi-valved BOB with custom anterior section.
- 3-4 x assist to don = Subacute / Slow Stream Rehab given de-conditioning and poor base line.







Anterior Section

Midline Overlap



Reason for referral:

 Previous admission for acute fracture L4 / L5 intervertebral disc following fall from step ladder. Fell while vacuuming her walls at home. Re-referred 6 months later for long term pain management.

Injury:

- Fracture L4/L5 3 column AO: L4-L5 B3 (L4: A0; N0; M1)
- Intervertebral disc fracture with associated pars fracture.





CT - Sagittal View



CT - AP View



Reason for referral:

 Previous admission for acute fracture L4 / L5 intervertebral disc following fall from step ladder. Fell while vacuuming her walls at home. Re-referred 6 months later for long term pain management.

Injury:

- Fracture L4/L5-3 column
- Intervertebral disc fracture with associated pars fracture.
- Large body habitus

Management:

- Previous managed in BOB, moderate compliance but requesting something less rigid.
- Primary request for pain relief



Reason for referral:

 Previous admission for acute fracture L4 / L5 intervertebral disc following fall from step ladder. Fell while vacuuming her walls at home. Re-referred 6 months later for long term pain management.

Injury:

- Fracture L4/L5 3 column
- Intervertebral disc fracture with associated pars fracture.
- Large body habitus

Management:

- Miami Lumbar posterior panel +/- anterior panel
- ? No management



https://www.ossur.com.au



References

Ackland, H. M., et al. (2007). "Factors predicting cervical collar-related decubitus ulceration in major trauma patients." Spine **32**(4): 423-428. Ackland, H. M., et al. (2016). "Danger at every rung: epidemiology and outcomes of ICU-admitted ladder-related trauma." Injury 47(5): 1109-1117. Al-Khalifa, F. K., et al. (2005). "Patterns of collapse in thoracolumbar burst fractures." Clinical Spine Surgery 18(5): 410-412. Bakhsheshian, J., et al. (2014). "Evidence-based management of traumatic thoracolumbar burst fractures: a systematic review of nonoperative management." Neurosurgical focus 37(1): E1. Cassell, E. and A. Clapperton (2006). "Consumer product-related injury (3): Injury related to the use of ladders." Hazard, Edition 63: 4. Chang, V. and L. T. Holly (2014). "Bracing for thoracolumbar fractures." Neurosurgical focus 37(1): E3. Faergemann, C. and L. B. Larsen (2000). "Non-occupational ladder and scaffold fall injuries." Accident Analysis & Prevention 32(6): 745-750. Harrop, J. S., et al. (2017). "Controversies in spinal trauma and evolution of care." Neurosurgery 80(3S): S23-S32. Hayward, G. (1996). "Risk of injury per hour of exposure to consumer products." Accident Analysis & Prevention 28(1): 115-121. Lemaire, E. D. (2016). "Mobilizing knowledge: The evidence gap for assistive devices." Technology Innovation Management Review 6(9) Kortebein, P. (2009). "Rehabilitation for hospital-associated deconditioning." American Journal of Physical Medicine & Rehabilitation 88(1): 66-77. Newman, M., et al. (2016). "Spinal orthoses for vertebral osteoporosis and osteoporotic vertebral fracture: a systematic review." Archives of physical medicine and rehabilitation 97(6): 1013-1025. Rubery, P. T., et al. (2013). "Stabilization of 2-column thoracolumbar fractures with orthoses: a cadaver model." Spine 38(5): E270-E275. Shen, J., et al. (2015). "Risk factors for the failure of spinal burst fractures treated conservatively according to the Thoracolumbar Injury Classification and Severity Score (TLICS): a retrospective cohort trial." PloS one 10(8): e0135735. Skoch, J., et al. (2016). "Bracing after surgical stabilization of thoracolumbar fractures: a systematic review of evidence, indications, and practices." World neurosurgery 93: 221-228. Takasaki, H. and T. Miki (2017). "The impact of continuous use of lumbosacral orthoses on trunk motor performance: a systematic review with meta-analysis." The Spine Journal 17(6): 889-900. Tyson, S. F. and R. M. Kent (2013). "Effects of an ankle-foot orthosis on balance and walking after stroke: a systematic review and pooled meta-analysis." Archives of physical medicine and rehabilitation 94(7): 1377-1385. Urquhart, J. C., et al. (2017). "Treatment of thoracolumbar burst fractures: extended follow-up of a randomized clinical trial comparing orthosis versus no orthosis." Journal of Neurosurgery: Spine 27(1): 42-47. Wood, K., et al. (2003). "Operative compared with nonoperative treatment of a thoracolumbar burst fracture without neurological deficit: a prospective, randomized study." JBJS 85(5): 773-781.



OUESTIONS?