

**50**  
1971-2021

 **ÖSSUR**®  
LIFE WITHOUT LIMITATIONS

## ICEFORM® Liners & Sleeve

Rob Hodgson – Clinical Specialist (Prosthetics)

Webinar - July 2021



# Objectives

- Foundation of liner technology
- Introduction
  - What is TPE
- User Profile
- Product Characteristics
  - Liners
  - Sleeve
- Sizing
- Best Practices



Iceform® Locking




Iceform® Cushion



Iceform® Sleeve



# Foundations of Liner Technology




## Liner Materials

Thermoplastic elastomer (TPE)

Polyurethane (TPU / PUR)

Silicon elastomer (Silicon)




## Liner Materials

Thermoplastic elastomer (TPE)

Polyurethane (TPU / PUR)

Silicone elastomer (Silicone)



TPE - Thermoplastic Elastomer, AKA: GEL

(Styrene) plastics

Mineral oil



+



=



Stiffness modified by *level of mineral oil* incorporated in the polymer matrix

## Benefits

- Gel Flow = High to Low areas of pressure
- Cushioning
- Even Pressure Distribution
- Reduced Friction
- Provide Suspension





## Liner Materials

Thermoplastic elastomer (TPE)

Polyurethane (TPU / PUR)

Silicone elastomer (Silicone)



Silicone elastomer, Iceross

Dimethylsiloxane-chains

SkinCare ingredients



Stiffness modified by chain length, level of crosslinking & skincare additives

## Benefits

- Many suspension possibilities, all activity levels
- Effective comfortable suspension
- Improved function
- Active skin care
- Ideal for total weight bearing socket solution



## Durometer Shore Hardness Scale




- Durometer is a measure of the hardness (flexibility) used typically with polymers, elastomers, and rubbers
- **The purpose of the scale is to provide individuals with a common point of reference for the observed material or similar product**
- Developed in the 1920's by Albert Ferdinand Shore
- The hardness of a material can be defined as the resistance to permanent indentation of that material (Softer materials = lower numbers)



# ICEROSS - DUROMETER OPTIONS



	Original • 60
	Transfemoral • 50
	Sport • 50/25
	Activa • 50/25
	Synergy • 50/25
	Iceform TPE • 37
	Comfort • 30
	Dermo • 25



# What is TPE

TPE = thermoplastic elastomer or “gel”

- Polystyrene – microscopic crosslinks that can be rearranged with heat, but reform during cooling



(Styrene) plastics

+



Mineral oil

=



- Stiffness of material modified by level of mineral oil incorporated into polymer mix





# What is TPE


- Gel flows from areas of high pressure to areas of lower pressure
- Helps to reduce friction against the skin that may occur while walking
- Conform to shape of residual limb
- Wear-and-tear of daily activities
- Liners vary in thickness and stiffness and usually have some type of fabric on one or both sides



# User Profile

- TT users
- K1-K2
- Low to moderate impact levels
- Balance Solution

User Information	
Amputation Level:	Transtibial
Impact Level:	Low to Moderate



# Iceform Liners and Sleeve



Iceform Locking



Iceform Cushion



Iceform Sleeve



# Iceform liners


TPE material

Locking / Cushion

Seamless textile

Indicators for anterior  
midline

10cm matrix



To be used with  
- Icelock series  
- Iceform Sleeve

Shape capture

Posterior thickness  
3 mm

Anterior thickness  
3, 6 or 9 mm

Uniform liner profile



# Iceform sleeve

TPE material

To be used with  
Iceform cushion  
liner

Expulsion valve  
recommended




Seamless textile

Shape capture

Profile


Wave feature





# Sizing


## Available sizes



# Sizing

## Sizing liners


- Measure circumference 4 cm and 30 cm from distal end residual limb



## Item numbers

## Liners

Part#	Profile	Suspension Method
T-0103XX	Iceform Uniform 3 mm	Locking
T-0106XX	Iceform Uniform 6 mm	Locking
T-0109XX	Iceform Uniform 9 mm	Locking
T-CL03XX	Iceform Uniform 3 mm	Cushion
T-CL06XX	Iceform Uniform 6 mm	Cushion
T-CL09XX	Iceform Uniform 9 mm	Cushion



# Sizing

## Sizing sleeve

- Measure circumference 15 cm above mid-patella to choose the correct size.

Part#	Size	Measurements
S-201135	35	35-41 cm
S-201141	41	41-49 cm
S-201149	49	49-55 cm
S-201155	55	55-62 cm
S-201162	62	62-76 cm



# Best Practices

## Custom Shape Capture:

After < 2 weeks of use liner conformed to shape assisting with liner orientation

- Automatic customisation
- Cold Capture
  - Add a protective layer to the clean positive model
  - Apply Iceform liner to positive model
  - Leave for a minimum time of one hour
- Hot Capture
  - Add a protective layer to the clean dry positive model
  - Place positive model into preheated oven at 65°C for 60 minutes
  - Apply Iceform liner to positive model and place into the oven
  - After maximum of 10 minutes remove from oven







# Best Practices

## Maintenance, Cleaning, and Storage:

- Turn liner inside out
- Wash the liner with pH-balanced dye-free soap
- Rinse the liner with approx 40°C warm water
- Pat dry the liner with lint-free cloth



# Best Practices

## Casting:

- Cast and modification for Total Surface Bearing (TSB) concept
  - Pressure distribution
  - Liner life expectancy
- Recommended to create check socket first



## Conclusion

- Össur TPE solution
- Locking, cushion and sleeve
- Trans-tibial amputee with low to moderate impact level

Providing a safe, secure and comfortable suspension for those patients that need/require a TPE liner



Rob Hodgson: [rhodgson@ossur.com](mailto:rhodgson@ossur.com) (NSW, ACT, SA)

Sarah Mulroy-Lang: [smulroy-lang@ossur.com](mailto:smulroy-lang@ossur.com) (VIC, TAS)

Lizzie Carey: [ecarey@ossur.com](mailto:ecarey@ossur.com) (NZ)

Laura Roberts: [lroberts@ossur.com](mailto:lroberts@ossur.com) (QLD, WA, NT)

