

**50**  
1971-2021

 **ÖSSUR®**  
LIFE WITHOUT LIMITATIONS

## Össur Legs: Optimising Mobility Outcomes with Complete Solutions

Sarah Mulroy-Lang – Senior Clinical Specialist

November 2021

# Össur Legs - Complete Solutions



## ÖSSUR LEGS

COMPLETE SOLUTIONS | CLINICALLY FLEXIBLE | SIMPLE & EASY

- Commencing in January 2022 we will be launching the Össur legs program
- Össur legs program allows you to quick and easily configure a complete and clinically appropriate solution to meet your clients needs
- Össur Legs program incorporates all of Össur's Mechanical and Bionic products, including Direct Socket and Connect TF socket solutions
- It features digital order forms for a fast and efficient ordering service



# Why Össur Legs?: Complete solutions addressing the needs of clients and clinicians



## TRANSTIBIAL LEG SOLUTIONS




- Complete range of complementary solutions
- Uncompromised quality
- Expedited fitting capability
- Flexible and fast ordering & service processes


## TRANSFEMORAL LEG SOLUTIONS



## For Clients



## Superior Outcomes



## For Clinicians & Workshops



# Complete Solutions: Impact Levels

## IMPACT LEVELS

### LOW

Daily activities involving gentle, steady walking with the use of a walking aid.

Example: Moving around at home, modest walking in the community.

### MODERATE

Daily activities involving average walking with the ability to vary speed or walking pattern.

Example: Going to the shops, confident outdoor walking.

### HIGH


Daily activities involving fast walking, jogging and climbing stairs.

Example: Heavy-lifting, manual labor, recreational sports.

### EXTREME


Activities involving running, track and field, sprinting, and long-distance running.

Example: Track and field sports.



# Complete Solutions – The Journey

**Aim:** Improve clinical outcomes, increase efficiency, and patient satisfaction



## Low Active User Profile

Low active prosthetic users have different needs and wants from their prosthesis compared to the active amputee

**Comfort**


**Safety**

**Mobility**

- They also have a different set of mobility goals
- Getting that balance right can mean
  - greater independence
  - increased confidence
  - increase in both the frequency & safety of transfers and ambulation




## Low Active User – Challenges



# Medical Necessity – Internal Research

## Amputee Compensatory Actions for ADL's





## Low Active User - Challenges



# Low Active User - Challenges



Low activity user

Moderate activity user

## Low Activity Solutions - Transtibial

Dave has type 2 diabetes and has spent the last 2 years hoping that the ulcer on his foot would heal, but the wound persisted.

He's been using a wheelchair for the past 8 months since a partial foot amputation which also refused to heal, and his fitness has declined

Dave recently underwent a transtibial amputation



**Dave**

"I'm keen to get back to bowling and seeing my friends"

## Low Activity Solutions - Transtibial

- 54 year old male
- Diabetic Ulcer
- Fishing, gardening, outdoor bowling, active fundraiser
- Former military, very eager, ulcer and pain in sound foot,
- Goals: Reduce his wheelchair use, get back to bowling, K2-K3 potential





## Low Activity Solutions – Direct Socket



Direct Socket is a **standardised, efficient** and **repeatable** socket solution for TT amputees

- Reduces patient visits
- Reduces labour hours
- Ability to use **all suspension options**
- Suitable for all **activity levels**
- Rated to **166kg**

## Low Activity Solutions – Direct Socket



### Basalt, Glass and Carbon options

- Fabricate in under 2 hours
- **Accurate anatomical fitting** utilising Icecast
- Consistent and repeatable results between clinicians
- User is seated when manufactured
  - **Safe & efficient process** done directly on client.
  - Fit same day possible – ideal for patients with volume fluctuations
- Compatible with, Locking, Cushion, Seal –In and Unity Elevated Vacuum

## Low Activity Solutions – Iceross Dermo Locking



### Iceross Dermo Locking Liner

- Softest durometer of silicone
- Ideal for vascular users or those with sensitive skin
- Features:
  - Active Skin care
  - Silken Inner surface
  - Conforms well to complex shapes
  - Ultra-Elastic textile for improved fit

## Low Activity Solutions – Icelock 214



### Icelock 214

- Attach to the Direct Socket
- Clutch mechanism provide easy and secure wind-down donning and easy release
- Can be donned in sitting
- Durable and strong, tested to tolerate at least 300kg of pull





## Low Activity Solutions – Balance Foot S Torsion



### Balance Foot S Torsion

C-Shaped Glass Fibre Foot



Cushioning heel foam



Smooth Roll Over



Gentle energy return



Balance and Stability



WATERPROOF

Waterproof design

## Low Active Solutions – Balance Foot S Torsion



Torsion version provides shock absorption and rotation



Reduces shear forces and stress on residual limb

### *PERMITTED ROTATION AND SHOCK MOVEMENTS REDUCE SHEAR STRESS*

- Twiste M. Transverse rotation and longitudinal translation during prosthetic gait – A literature review. Journal of Rehabilitation, 2003; 40:9-18

### *LOWER SHEAR STRESSES → REDUCED DISCOMFORT*

- Segal AD, et al. Transtibial amputee joint rotation moments during straight-line walking and a common turning task with and without a torsion adapter. Journal of Rehabilitation R&D. 2009 46:375-384

### *SELF PERCEIVED PAIN REDUCTION*

- Segal AD, et al. Does a Torsion Adapter Improve Functional Mobility, Pain and Fatigue in Patients with Transtibial Amputation? Clin Orthop Relat Res 2014 472:3085–3092

### *REDUCED SHEAR STRESS*

- Heitzmann D, et al. Functional effects of a prosthetic torsion adapter in transtibial amputees during unplanned spin and step turns. Prosthet Orthot Int 2015; 1-8.

## Low Activity Solutions - Summary



- Direct Socket provides TSB socket and manufactured quickly
- Iceross Dermo locking liner provides cushioning and protection
- Icelock 214 provides secure mechanical lock
- Balance Foot S Torsion combines stability, shock and rotation and a smooth roll over, providing balance and safety for the low active user



## Low Activity Solutions: Alternative Prescription



- Iceross Dermo Cushion Liner
    - Provides cushioning and protection for sensitive residual limbs
  - Iceform sleeve
    - Easy Donning for users with compromised hand dexterity
  - Can be used with Direct Socket, with 511 or 544 expulsion valves
- 
- Suitable for interim use, simple to manage volume changes with socks



## Low Activity Solutions: Alternative Prescription



C-Shaped Glass Fibre Foot



Cushioning heel foam



Smooth Roll Over



Gentle energy return



Balance and Stability



WATERPROOF

Waterproof design

## Alternative Foot Considerations: Low Activity Solutions



### Balance Foot

- Balance Foot is lightweight, multi-axial foot for low activity users
- Provides balance comfort and stability
- Single, slow speed walking, with step to gait



### DP Flexion

- Hydraulic ankle with adjustable resistance to plantar and dorsiflexion.
- Range of 18°, for adaptability on uneven terrain
- Provides 4° swing dorsiflexion for clearance in swing, increasing toe clearance to reduce trips and falls



### Balance Foot J

- J shaped Carbon keel provides good roll over characteristics
- Polyurethane cushioned heel strike benefits users struggling with heel strike and getting into midstance
- New amputees and slow speed walkers



### Flex-Foot Assure


- J Shaped Carbon keel provides full length toe lever
- Active Carbon heel for shock absorption
- Heel wedge enables a slower energy return to accommodate slower walking speeds

## Low Activity Solutions - Transfemoral

Bill has been unwell for a while and has recently undergone a transfemoral amputation due to PVD

His wife, Kathryn is his main carer and helps Bill with most daily tasks since his recent amputation.

Bill wants to be able to make Kathryn cups of tea and slow dance again




**Bill**

"I want to get out of the house and see my grandchildren"

## Low Activity Solutions - Transfemoral

- 78 year old male
- Vascular Disease
- Meeting friends at the RSL, Community outings
- Lacks confidence, low energy
- Goals: Limited ambulation, play with grandchildren, K2



**Bill**

"I want to get out of the house and see my grandchildren"



## Low Activity Solutions – Connect TF



### Connect TF

- Adjustable prosthetic socket solution
- Can be custom fit in a single session
- Height and volume adjustable
- Can be donned and doffed in sitting
- Designed to be comfortable whilst sitting and standing

## Low Activity Solutions – Connect TF



### Connect TF

#### Safety

- Don in sitting
- Mechanical clutch lock suspension

#### Comfortable

- Flexible socket brim
- Very comfortable when sitting
- Tension handle can be released in sitting

#### Ease of Use

- Easy for user or carer to don prosthesis
- Easy to clean - waterproof

## Low Activity Solutions – Connect TF



### Connect TF – Benefits:

- Beneficial for use in early rehab
  - Can be adjusted to accommodate changes in volume
- K1 users who may not be considered suitable for a prosthesis
  - Can be fit with a Connect TF socket whilst in sitting
  - Can try a comfortable socket solution and assess prosthetic suitability

## Low Activity Solutions – Iceross TF Locking Liner



- **Iceross TF Locking Liner**
- Provides firm suspension and soft tissue stabilisation
- Durable fabric cover allows for radial stretch and comfortable elasticity
- Recommended for use with Connect TF or Icelock 200 series

## Low Activity Solutions – Balance Knee OFM1



### • Balance Knee OFM1

- K-Level 1 and 2
- Maximum user weight 136Kg
- 4 bar geometry for stance phase stability
- Mid-swing shortening for toe clearance
- Optional locking mechanism

### Indications:

- Household ambulators
- Limited community ambulators
- Users in transition between these levels
- Long residual limbs
- Benefit from shock absorption from optional Stance Flexion adapter



## Low Activity Solutions – Balance Knee OFM1



- **Balance Knee OFM1**

- Optional locking mechanism-
  - users with low voluntary control or use in early rehab
- Adjustable integrated extension spring
  - Optimise extension speed for the user
  - Adjustable axis friction setting

### Contraindications:

- Users capable of varying their walking speeds
- Unable to fully extend the knee

# Low Activity Solutions – Balance Foot S



LOW  
ACTIVITY

SECURE | COMFORTABLE | EASY DONNING



CONNECT® TF



ICEROSS® TF LOCKING LINER



BALANCE™ KNEE OFM1



BALANCE™ FOOT S



C-Shaped Glass Fibre Foot



Cushioning heel foam



Smooth Roll Over



Gentle energy return



Balance and Stability




WATERPROOF

Waterproof design

## Low Activity Solutions - Transfemoral



- Connect TF is adjustable, can be donned in sitting and provides a comfortable socket solution
- Iceross TF Locking liner provides soft tissue stabilisation and secure suspension
- Balance Knee OFM1 provides geometric stability in stance, mid swing shortening for toe clearance and an optional lock for use in early rehab
- Balance Foot S combines stability, and a smooth roll over, providing balance and safety for the low active user




## Low Activity Solutions – Alternative Prescription



### Balance Knee OFM2

- Features:
  - Adjustable weight activated brake
  - Individual adjustable stance safety
  - Knee will remain stiff even when loading a flexed knee
  - Once fully unloaded the brake mechanism is released and the knee will swing freely
  - Optional locking mechanism for use in early rehab

- Brake sensitivity adjustment (B)
  - Clockwise (-): Decrease sensitivity
  - Anti-clockwise (+): Increase sensitivity
- (P) Brake play screw



## Low Activity Solutions – Alternative Prescription



### Balance Knee OFM2

Indications:

- Users with low voluntary control
- Elderly/Weak/Co-morbidities
- Low activity user utilising gait aids
- New amputees (interim use)
- Single, slow speed ambulators



## Alternative Knee Considerations: Low Activity Solutions



### Locking Knee

- K1 level users requiring highest level of stance control
- Single speed slow walkers
- Transfers
- Elderly/weak users unable to utilise a polycentric or friction brake knee



### Balance Knee OM8

- K2 Level
- 4 bar geometry for stance phase stability
- Adjustable centroid for to fine tune stance release
- Mid-swing shortening for toe clearance
- Household and limited community ambulators
- Manual lock not required



### OP4

- K2-K3 Level
- Adjustable weight activated friction brake knee
- Pneumatic swing phase control for users capable of variable cadence
- Don't require a manual lock, but require stance phase stability from weight activated friction brake




### Total Knee 1900

- K2-K3 Level
- Enhanced stability from Geometric Lock
- Mid-swing shortening for toe clearance
- Incorporates adjustable stance flexion for shock absorption
- Polymer friction for variable cadence
- Low build height

## Moderate Activity Solutions - Transtibial

Stewart is a family man, who underwent a traumatic transtibial amputation 1 year ago following an accident at work

He loves to bushwalk and go fishing with his family however he is still on his interim prosthesis and has issues with socket discomfort, which limits how long he can wear his prosthesis



**Stewart**

“I love being outdoors and fishing with my family”

## Moderate Activity Solutions - Transtibial

- 41 year old male
- Accident at work, lead to traumatic amputation
- Young family, bush walking, fishing
- 1 year post amputation, motivated, active
- Goals: Walk without discomfort, fishing K3



**Stewart**

"I love being outdoors and fishing with my family"

## Moderate Activity Solutions – Direct Socket



- Moderate active users can benefit from a total surface bearing socket
- Accurate anatomical fit utilising Icecast
- Manufactured quickly
- Can easily change components to try different suspension methods
- Rated to **166kg**



## User experience– Direct Socket



- ‘... from fabrication through to completion of my new legs, the process was **simple and exceptionally quick**’
- ‘The major benefit I found was **knowing how my new sockets would feel on using the pressure test**. This allowed me to identify any pressure areas or discomfort prior to the resin setting stage’
- ‘The fit is perfect and they’re so **light and comfortable** that I don’t feel like I’m wearing anything compared to my previous prosthetics’
- ‘The Direct Socket System provides the **custom comfort, support, and sense of freedom** I never had in my previous sockets’
- ‘**I can walk longer distances and even jog without any pain or discomfort**. They also look so awesome that I’m more confident and happier to show them off. I love my legs!’





## Moderate Activity Solutions – Seal-In X Locking



### Seal-In X Locking Liner

- Dual suspension silicone liner
- Provides mechanical suspension in combination with Seal-In or elevated vacuum
- Minimises pistoning and rotation within the socket
- Choice of Seal Rings, Classic, Grip or Volume
- Provides confidence and security for the user



### Icelock 562 Hybrid Lock

Provides mechanical suspension in combination with Seal-In or elevated vacuum

- Unity elevated vacuum provides:
  - Very firm suspension
  - Manages volume fluctuations
  - Reduces pistoning and movement providing increased comfort
- 562 Hybrid Lock Features:
  - Low build height
  - Waterproof
  - Rated to 166Kg
  - Unity and Direct Socket compatible

\*Must be used in combination with Seal-In X Locking



## Advantages of Vacuum Suspension

- Very firm suspension providing excellent security and improved **proprioception** and **balance**<sup>1</sup>
- Elevated vacuum manage limb volume fluctuation<sup>2</sup>. A **constant limb volume**, decreases the need for sock management
- Elevated vacuum **assists with wound healing** by improving circulation through the residual limb<sup>3</sup>
- Provides good distal comfort for bony and sensitive distal ends; as long as the socket fits properly with good **volume and length matching**
- The combination of **reduced pistoning** and maintenance of volume is thought to account for the **more symmetrical gait** observed with vacuum<sup>1,4</sup>

## Moderate Activity Solutions – Pro-Flex Pivot



Unique pivot technology  
+  
Interaction of 3 carbon fibre blades  
↓  
ROM of 27 degrees  
↓  
Powerful push off  
↓  
Physiological progression of COP



Weatherproof

## Moderate Activity Solutions – Pro-Flex Pivot



- Pro-Flex Pivot adapts very well on uneven ground, suitable for:
  - Bush walking
  - Farms environment
  - Slopes/Ramps
  - Stairs
  - Active users that do a lot of walking



## Moderate Activity Solutions – Pro-Flex Pivot



- Pro-Flex Pivot delivers functional improvements and benefits to users
- This combination of functional features delivers overall medical benefit:
- Reduced impact of 13% on the sound side
- Reduced varus moment of 19% on the sound side
- Reimbursement Document available

\* Data compared to Vari-Flex foot

# Moderate Activity Solutions – Pro-Flex Pivot



## User Experience – Pro-Flex Pivot



... 'the Pro-Flex Pivot allows my gait to be as natural as possible'

... 'Everyone comments on how natural my gait is'

'Gives me practical movement that allows me to perform movements that amputees usually struggle with, balance particularly!'

... 'Prior to Pro-Flex Pivot, I would be fatigued by the afternoon. Now nothing stops me!'


<https://www.ossur.com/en-au/prosthetics/client-experiences>



## Moderate Activity Solutions



- Direct Socket provides TSB socket
- Iceross Seal-In X Locking provides dual suspension for security and confidence
- Icelock 562 hybrid lock provides secure mechanical lock with Unity and Seal-In options
- Pro-Flex Pivot combines a large ROM, anatomical roll over and powerful push off, helping to reduce the impact on the sound limb



## Moderate Activity Solution: Alternative Prescription



Weatherproof



Össur Logic connectivity

Terrain adaptation and  
4° mid-swing dorsiflexion



Integrated battery:  
18 - 36 hours battery life  
depending on usage

33° ROM

Excellent stance dynamics from  
Pro-Flex LP foot blade




# Moderate Activity Solution: PROPRIO FOOT



Core functions:

Swing dorsiflexion




Ankle alignment

Stair adaptation



Ramp adaptation

Relax/Chair exit



# PROPRIO FOOT - Medical Necessity

- Amputees are more at risk of falling than their peers
- Reduction to Minimum Toe Clearance (MTC) from amputation could increase the incidence of trips and falls



PROPRIO FOOT provides toe clearance

- Reduces risk of trips and falls
- Increases safety



# PROPRIO FOOT: Medical Necessity



- 16 TTA + 16 non-amputees
- Neutral ankle angle vs. 4 degree adaptation

Knee flexion is restricted  
because of limited (ankle)  
dorsal flexion

Both stair ascent and descent  
improvements of knee  
kinematics and kinetics

Increased knee flexion and  
increased knee moment

More physiological knee flexion  
during stair ascent and descent



PROPRIO FOOT: Medical Necessity



**MODERATE  
ACTIVITY**

ADAPTABLE | STABLE | COMFORTABLE







- 'The PROPRIO FOOT is very reliable, and I find I'm less tired throughout the day'
- 'I feel more confident now I'm using the PROPRIO FOOT due to the ankle being used in combination with a great energy storing carbon fibre foot'
- 'The PROPRIO FOOT has improved my comfort by enabling me to walk more naturally'
- 'I find it easier to navigate uneven surfaces such as building sites and mine sites. This makes me feel more comfortable in my socket and I'm able to walk further than before'
- 'I find I'm less fatigued on slopes and ramps as this foot adjusts to accommodate the surface angles and makes it more comfortable for me to walk on'
- 'The PROPRIO FOOT has positively impacted my life, giving me more flexibility at work and at home to utilise the prosthesis for longer periods'



# Moderate Activity Solution: PROPRIO FOOT



### MODERATE ACTIVITY

ADAPTABLE | STABLE | COMFORTABLE

PROPRIO FOOT

## PROPRIO FOOT®

Reimbursement Support

### FOLLOW UP CONSIDERATIONS/MAINTENANCE REQUIREMENTS

PROPRIO FOOT limited warranty period: 24 months. Extended warranty is available for purchase, contact Customer Service for options and prices.  
Total warranty period is 5 years (must be purchased within one year of original purchase date).  
Includes support to the customer with loaner units if any repairs or maintenance are required.

### REFERENCES

Active dorsiflexing prostheses may reduce trip-related fall risk in people with transtibial amputation. Tenenblatt, Noah J., et al. Journal of Rehabilitation Research and Development 51.8 (2014): 1229-1242.  
Risk factors and costs associated with accidental falls among adults with above-knee amputations: a population based study. Kaufman, K. American Orthotic and Prosthetic Association 2016. (May 2016). <http://www.aopa.net.org/resources/research>  
Asymmetry in external work (SEW): A novel method of quantifying gait differences between prosthetic and natural feet. Agrawal, Vibhor, et al. Prosthetics and orthotics international 33.2 (2009): 148-156.  
Assessment of the effects of carbon fibre and bionic foot during overground and treadmill walking in transtibial amputees. Delussu, AS, et al. Gait and Posture (2013).  
Validation of self-report and performance-based outcome measures to determine functional differences between four categories of prosthetic feet. Gailey, R. et al. Journal of Rehabilitation Research and Development (2012): pp.597-612.  
Mechanical analysis of ramp ambulation of transtibial amputees with an adaptive ankle foot prosthesis. Fradet L, Alimusaj M, Braatz F, Wolf SI. Gait & Posture. 2010; 32(2): 191 - 198.  
Preliminary results of trans-femoral amputees walking with a microprocessor controlled prosthetic foot. Gait & Posture. June 2012; 36 Supplement 1:59. Heitzmann DWW, Alimusaj M, Braatz F, Wolf SI.  
Comparison of compensation mechanisms in transfemoral amputees fitted to a conventional energy storing foot versus microprocessor controlled energy storing foot. Lechler K. 2008.  
Gait and kinetics with an adaptive ankle foot system during stair ambulation of trans-tibial amputees. Gait & Posture. Alimusaj M, Fradet L, Braatz F, Gerner HJ, Wolf SI. 2009; 30:3356-363.  
Effects of an increased prosthetic ankle range of motion for individuals with a trans-tibial amputation walking with a new prosthetic foot. Gait & posture. 2018 Jul 1;64:174-80. Heitzmann DW, Gerner HJ, De Asha AR, Block J, Putz C, Wolf SI, Alimusaj M.  
Increasing prosthetic foot energy return affects whole-body mechanics during walking on level ground and slopes. Scientific reports. 2018 Mar 29;8(1):3354. Childers WL, Takahashi KZ.

ÖSSUR ON

ÖSSUR.COM/AU

Össur Australia - Prosthetics

TEL: +61 2 9636 2600

FAX: +61 2 9475 1114

[info@ossur.com.au](mailto:info@ossur.com.au)

[Download Here](#)

## Alternative Foot Considerations: Moderate Activity Solutions



Vari-Flex

- Offers very high energy return in a lightweight design
- Forward tibial progression to assist in propelling the body forward
- Shock absorption at heel strike
- Full-length toe lever
- Split heel and keel to accommodate uneven terrain



Talux

- Ideal for community ambulator, for effortless walking and moderate impact recreational sports
- Multi-axial core providing inversion/eversion compliance on uneven ground
- Features Sandal toe




Pro-Flex LP

- K2-K4 Level
- Greater ankle range of motion than other low profile feet
- Anatomical sole blade and multi-axial compliance
- Smooth Roll over for everyday use
- Great solution when space is compromised

## Moderate Activity Solutions - Transfemoral

Sue works in marketing, and loves socialising with friends and trying new restaurants

Sue underwent a transfemoral amputation 7 years ago, due to Osteosarcoma, and has been getting by as best she can, however she sometimes utilises a walking aid when outdoors, for security.




**Sue**

“I’d love to have the confidence to go out without my walking stick”



## Moderate Activity Solutions - Transfemoral

- 28 year old female
- Osteosarcoma, resulting in TF Amputation
- Foodie, socialising, marketing role
- 7 years post amputation, motivated, active
- Goals: Walk without discomfort, K3



**Sue**

"I'd love to have the confidence to go out without my walking stick"





## Moderate Activity Solutions – Direct Socket TF



### Direct Socket TF

- Easy to Donn & Doff as compared to traditional socket
- Superior Functional Outcomes
  - Improved ROM/Suspension/Stability in all planes & Superior femur control
- Superior Comfort
  - Sitting, standing & walking
- Ultra-Lightweight & Low-Profile
- True Total Contact

Direct Socket TF: Clinical Study



UNIVERSITY OF TORONTO  
LIBRARIES

CPOJ  
ISSN: 2561-987X

RESEARCH ARTICLE

# Canadian Prosthetics & Orthotics Journal

All articles are permanently available online to the public without restrictions or subscription fees. All articles are free to be used, cited, and distributed, on condition that appropriate acknowledgment is included. Authors are the copyright holders of their original contributions and grant the Canadian Online Publication Group (COPG) a license to publish the article and identify itself as the original publisher. CPOJ articles are licensed under the Creative Commons Attribution 4.0 International License.

CPOJ Website: <https://jpls.library.utoronto.ca/index.php/cpoj/index>  
Editorial Office: [cpoj@online-publication.com](mailto:cpoj@online-publication.com)  
ISSN: 2561-987X  
CPOJ is a member of, and subscribes to the principles of, the Committee on Publication Ethics (COPE).

Volume 3, Issue 2  
2020

Marable W.R., Smith C., Sigurjónsson B.D., Atlason I.F., Jóhannesson G.A. Transfemoral socket fabrication method using direct casting: outcomes regarding patient satisfaction with device and services. Canadian Prosthetics & Orthotics Journal. 2020;Volume 3, Issue 2, No 6. <https://doi.org/10.33137/cpoj.v3i2.34672>

Journal Homepage: <https://jpls.library.utoronto.ca/index.php/cpoj/index>  
Volume 3, Issue 2, Article No 6, 2020

## FABRICATION METHOD USING DIRECT CASTING: OUTCOMES REGARDING PATIENT SATISFACTION WITH DEVICE AND SERVICES

Marable W.R.<sup>1</sup>, Atlason I.F.<sup>1</sup>, Jóhannesson G.A.<sup>2,4\*</sup>

**ABSTRACT**  
The Direct Socket (DS-TF) prosthetic user is a novel method of fabricating a prosthetic socket. It requires different training, production method and service delivery. This method and model may improve patient satisfaction and delivery in one visit.

**KEYWORDS**  
Transfemoral amputation, Prosthetic, Socket, Interface, Outcome measure, Satisfaction, Direct casting

**ARTICLE INFO**  
Received: August 3, 2020  
Accepted: November 16, 2020  
Published: November 23, 2020

**CITATION**  
Marable W.R., Smith C., Sigurjónsson B.D., Atlason I.F., Jóhannesson G.A. Transfemoral socket fabrication method using direct casting: outcomes regarding patient satisfaction with device and services. Canadian Prosthetics & Orthotics Journal. 2020;Volume 3, Issue 2, No 6. <https://doi.org/10.33137/cpoj.v3i2.34672>

**KEYWORDS**  
Transfemoral amputation, Prosthetic, Socket, Interface, Outcome measure, Satisfaction, Direct casting

**INTRODUCTION**  
The Direct Socket (DS-TF) prosthetic user is a novel method of fabricating a prosthetic socket. It requires different training, production method and service delivery. This method and model may improve patient satisfaction and delivery in one visit.

**CONCLUSION**  
The Direct Socket (DS-TF) prosthetic user is a novel method of fabricating a prosthetic socket. It requires different training, production method and service delivery. This method and model may improve patient satisfaction and delivery in one visit.

Marable et al. 2020  
TRANSFEMORAL SOCKET FABRICATION METHOD USING DIRECT CASTING

## TRANSFEMORAL SOCKET FABRICATION METHOD USING DIRECT CASTING

The problem with this 'definition' is that a transtibial or toe prosthesis could also be referred to as sub-ischial socket. This lack of detail and consistency in definitions prevents objective comparison.<sup>27</sup> In an effort to provide detail and resolve the inconsistency, ISO standard 13405-2:2015 contains recommendations on how to systematically describe an interface and provides a base for comparison between different interface design and function.<sup>28</sup>

A method of fabricating a finished laminated TT interface directly on the patient's residual limb has been on the market since 1996. First introduced as ICEX by Össur HF of Iceland, ICEX was based on the inventor's philosophy of pressure casting<sup>29</sup> and was included in a variety of studies related to user satisfaction, cost and function.<sup>30-32</sup> Since then ICEX was improved to Modular Socket System (MSS) in 2005,<sup>31</sup> and finally to Direct Socket in 2018.<sup>33</sup> This system enables a prosthetist to fabricate a custom-made interface directly on the TT residual limb in a single visit. The lesson learned from the use of DS-TT in Scandinavia for over two decades and the finding related to the process, e.g. shorter rehabilitation time, demanded a solution for TF level.<sup>34</sup>

A TF version of Direct Socket (DS-TF) began testing in 2016 in select Scandinavian clinics.<sup>35,31</sup> The proximal portion of a DS-TF interface design differs significantly from the proximal portion of sockets typically called Ischial Containment or Inter Ramus Containment sockets, as the proximal part of the DS-TF includes a size-specific silicone brim. The method of fabricating directly on the residual limb requires a different approach to prosthetist training and fabrication compared to the socket fabrication process most prosthetists apply today.<sup>35</sup>

The primary aim of this study was to collect data on prosthetic users satisfaction regarding DS-TF interface in terms of both interface function and the clinic service model (i.e. one patient visit to the clinic for fabrication of custom TF interface, prosthesis assembly, alignment, gait analysis, and delivery).

Marable et al. 2020  
TRANSFEMORAL SOCKET FABRICATION METHOD USING DIRECT CASTING

## TRANSFEMORAL SOCKET FABRICATION METHOD USING DIRECT CASTING

creating axial and transverse stabilization. During the rest of the gait (pre-, initial-, mid-, and late swing) the brim is only following the hip movement.

The selection criteria for study principle investigators included a Certified Prosthetist (CP) with more than 5-years clinical experience in serving TF amputees with interest in improving patient outcomes and satisfaction, willingness to follow a defined novel fabrication protocol, and commitment to document outcome measures at defined time intervals. The inclusion criteria (rationale) is listed in Table 1.

Figure 1: Direct Socket

Table 1: Amputee inclusion criteria (rationale)

Amputee inclusion criteria (rationale)
• 50Kg+ body weight < 160Kg (the ISO validated weight limit of the DS-TF)
• Cognitive ability to understand all instructions and questionnaires in the study
• Patients who have undergone a TF amputation > 1-year post amputation (this was to avoid postoperative problems and/or adjustments related the initial prosthetic fitting of a new amputee)
• Older than 18 years
• Willing and able to participate in the study and follow the protocol
• Circular dimension of 40-65 cm at the crotch (limited to available silicone brim sizes)
• Residual limb length at least 20 cm from ischium to distal end (fabrication limitation of the DS-TF)
• Currently using a prosthetic liner (this was to avoid potential confounding influence from transitioning an amputee from a skin fitting interface (i.e. without a liner), to an interface with a liner)
• Willing to use a silicone prosthetic liner as called for in Direct Socket Instructions For Use. <sup>34</sup>

Document **patient satisfaction** of new Direct Socket TF (DS-TF) interface (socket) compared to amputees' current interface (socket) in terms of:

- **Comfort & function**
- **Prosthetic Service**  
Including Prosthetist training using a standardized training protocol

n=38	Socket satisfaction (DS TF vs previous)	Service satisfaction (DS TF vs previous)	DS TF Preference (DS TF vs previous)
6 months follow-up	<b>+29,8%</b>	<b>+14,8%</b>	<b>93.6%</b>

Investigate patient outcomes of new Direct Socket TF (DS-TF) interface (socket) compared to amputees' current interface (socket) in terms of:

- Quality of Life (QoL)
- Comfort
- Mobility

- **Subjective:** EQ-5D-5L®, PLUS-M™, CLASS, ABC, (\*OPUS-CSD and CSS),
- **Objective:** TUG and AMPPRO

K Level	Subjects
1	4
2	11
3	21
4	11

Outcome Measure	EQ-5D-5L®	PLUS-M™	CLASS	ABC	TUG	AMPPRO
Initial	0.75	46	74%	63%	14.7s	38
6 weeks	0.82	54	86%	75%	13.0	40
6 months	0.84	61	86%	78%	12.8s	40



### Iceross Seal-In HSM

- Features full length matrix for soft tissue stabilisation
- Good compliance with complex limb shapes and bony residual limbs
- Active skin care
- Can be used with Seal-In Socks
  - Maintenance of suspension for volume fluctuations
- Can be used in conjunction with Icelock 544 for seal in and Unity elevated vacuum to provide a firm suspension



## Moderate Activity Solutions – PASO Knee



### PASO Knee

- Closing Geometry for stability in stance
- **Auto-adaptive swing phase control**
  - No valves to adjust
- Adapts to variable walking speed and smooth and powerful running
- Mid swing shortening (15mm) – helps to prevent trips and falls

Moderate Activity Solutions – PASO Knee



MODERATE  
ACTIVITY

ADAPTABLE | STABLE | COMFORTABLE

DIRECT SOCKET TF

ICEROSS SEAL-IN<sup>®</sup> HSM  
WITH ICELOCK<sup>™</sup> 544

PASO KNEE<sup>™</sup>

PRO-FLEX<sup>®</sup> LP ALIGN



Moderate Activity Solutions – PASO Knee



**MODERATE  
ACTIVITY**

ADAPTABLE | STABLE | COMFORTABLE







Moderate Activity Solutions – PASO Knee



MODERATE  
ACTIVITY

ADAPTABLE | STABLE | COMFORTABLE





## Moderate Activity Solutions – PASO Knee





## User Experience: PASO KNEE



- 'The **PASO Knee maintains the same speed as my natural leg** no matter how fast or slow I am walking. This removes the stress of having a prosthetic knee that either won't keep up with me or has too much terminal impact which is uncomfortable for my hip'
- '...I feel more confident and stable'
- '...with the **polycentric geometry** of the knee and the short build height I find the PASO Knee folds away really well, making sitting on the bus and fitting into an aircraft seat a lot easier than before'
- 'When using the Pro-Flex XC Torsion and PASO Knee I can also **run with absolute ease**'
- '...I now look forward to walking, no matter how far or how technical the walk is'

## Moderate Activity Solutions – Pro-Flex LP Align



Greater ROM than other feet with a low profile design



Functional joint centre is closer to anatomical



Improved physiological gait



Increased function and user satisfaction

- Hydraulic housing, allows alignment adjustment of the foot to match different heel heights.

**0 – 7cm**


## Moderate Activity Solutions – Pro-Flex LP Align



### PRO-FLEX® LP ALIGN

Stay aligned

Up to 2 ¾" Heel Height  
Adjustability



Malalignment leads to:

Excessive pressures on residual limb:

- Residual limb problems
- Skin breakdown

Disturbed roll-over:

- Asymmetrical gait
- Low back pain

Inefficient gait:

- Higher energy expenditure

## Moderate Activity Solutions



- Direct Socket TF provides a total contact, comfortable socket, manufactured in a single appointment
- Iceross Seal-In HSM can be used in conjunction with Icelock 544 for seal in and Unity elevated vacuum to provide a firm suspension
- PASO Knee provides closing geometry for safety in stance, and automatically adapts to variable walking speeds
- Pro-Flex LP Align combines a wide range of motion, smooth roll over and hydraulic housing, for adjustment of heel height from 0 – 7cm







### RHEO KNEE


- Stance and Swing controlled MPK, providing adjustable resistance utilising a Magnetorheologic actuator
- Sensors measure angle, load and movement to detect its environment and interaction between the device and the user

#### Features:

- Kinematic Sensor
  - Creates safer transitions from stance to swing
  - Identifies quiet standing
  - Prevents inadvertent release when backwards walking



## Alternative Prescription - Moderate Active Solutions



- Extension Assist
  - Faster extension rate
  - Lowers energy consumption/exertion
- Manual Extension Lock
  - Traversing extreme terrain
  - Exercise related activities
- Enhanced Stability Control Settings
  - **Loading Response**
  - **Swing Initiation**
- Stable options are aimed at insecure, lower active users or primary users



### PROPRIO FOOT

- Provides 4° of mid swing dorsiflexion during swing, for toe clearance
  - Beneficial for users who like to bush walk
- Adapts to uneven terrain, and stairs, providing stability, comfort and safety
- Relax and Chair exit mode
- Pro-Flex LP foot blade provides stance dynamics and a smooth roll over

## Alternative Knee Considerations: Moderate Activity Solutions



OHP3

- K2 – K3 Level
- Enhanced stability from closing geometry – Provides safety in stance
- High pressure pneumatics for faster speed walkers
- Mid swing shortening for toe clearance in swing, helps to prevent trips and falls
- Optional stance flexion adapter (IKF) for shock absorption



Total Knee 2000

- K3 – K4 Level
- Enhanced stability from Geometric Lock
- Mid-swing shortening for toe clearance
- Incorporates adjustable stance flexion for shock absorption
- 3 phase hydraulic swing phase control for faster speed walking
- Low build height




MAUCH Knee

- K3 – 4 Level
- Hydraulic cylinder for control of swing and stance phase
- Hydraulic swing phase control for a wide range of cadence variation, with adjustable flexion and extension resistance
- Hydraulic stance control provides yielding resistance for stair and ramp descent
- Stumble control for safety

## High Activity Solutions - Transtibial

Karu works in finance, however outside of work he very into fitness, participating in cycling and cross fit, and going to the gym daily.

Karu became a transtibial amputee following a motorbike accident, and recently has been experiencing socket discomfort due to volume changes in his residual limb




**Karu**

"I am experiencing discomfort in my socket when exercising"



## High Activity Solutions - Transtibial

- 31 year old male
- Motorbike accident
- Cycling, Gym work, Cross Fit
- 3 years post amputation, highly active
- Goals: Cycle and Gym work without limitations



**Karu**

"I am experiencing discomfort in my socket when exercising"



## High Activity Solutions – Direct Socket



- High active users can benefit from a total surface bearing socket
- Accurate anatomical fit utilising Icecast
- Can easily change components to try different suspension methods
- Lightweight and strong - rated to **166kg**

## High Activity Solutions – Iceross Synergy Cushion



### Iceross Synergy Cushion

- Stability, cushioning and impact protection
  - Soft silicone inner layer provides cushioning
  - Firm silicone outer layer provides tissue stabilisation
- Wave feature improves comfort during knee flexion
- Features:
  - Double layer of silicon
  - Active Skin care
  - Silken surface
  - Wave feature available

Iceross sleeve features wave technology for ease of knee flexion, whilst generating firm suction

## High Activity Solutions – Icelock 544



### Icelock 544 expulSION valve

- For use with cushion liner and sleeve set up, where expulsion only is desired
- Valve sits distal to the socket, protecting it from external damage

# High Activity Solutions – Pro-Flex XC



Anatomical sole blade



Smooth roll over



Balance and Stability



Multiaxiality on uneven ground



## High Activity Solutions – Pro-Flex XC



**Waterproof**

C-shaped design provides 10mm vertical compression



Shock absorption for high impact activities



C-shaped design provides progressive stiffening for a more natural roll over




Great ankle ROM and reduced sound side loading



## High Activity Solutions: Summary



- Direct Socket is a lightweight and strong socket (rated to 166kg)
- Iceross Sleeve and Synergy cushion liner provide firm suspension, cushioning and protection
- IceLock 544 Expulsion valve sits distal to the socket, protecting it from external damage
- Pro-Flex XC provides a smooth roll over, shock absorption for high impact activities, and is waterproof (Salt and Chlorine)



**Karu**

"I can now exercise without discomfort – nothing limits me now"

## High Activity Solutions – Alternative Prescription



- High active user can benefit from the dual suspension Seal-In X Locking liner
- Use combination with the Icelock 562 Hybrid lock with Unity elevated vacuum to minimise movement within the socket during reactional activities
- Pro-Flex XC Torsion provides all the benefits of the Pro-Flex XC, with integrated Torsion Shock
- Compensates for physiological rotation and vertical shock absorption of lost joints
- Benefits users who participate in high impact activities such as basketball or running

## Alternative Foot Considerations: High Activity Solutions



### Pro-Flex LP Torsion

- Combines the dynamics of the Pro-Flex LP with torsion and shock
- Provides comfort and reduces shear stresses on the residual limb, for co-morbidities, scarring and Osseointegration
- Everyday use and recreational activities
- Rated to 147Kg



### Re-Flex Rotate

- K3-K4 Level
- J-Shaped carbon fiber foot
- Integrated torsion shock adapter
- Reduces shear stress on residual limb for high performance individuals
- Available with male or female pyramid options



### Re-Flex Shock


- K3-K4 Level
- J-Shaped Carbon fiber foot
- Integrated shock- absorbing adapter
- Designed for active individuals taking part in sport or high impact manual labour
- Require a foot for everyday use and recreation

## High Activity Solutions - Transfemoral

Tim is adventurous, loves to travel, and participates in extreme sports such as snow boarding.

Tim was born with a congenital deficiency which led to a transfemoral amputation at a very early age. He has always lived life to the fullest and doesn't want his amputation to stop him from doing anything.

He has been on a mechanical knee for years and wants to know what other options are available, to keep up with his active lifestyle.




**Tim**

"I do my best be as active as possible, but I get frustrated when my prosthesis limits me "



# High Activity Solutions - Transfemoral

- 28 year old male
- Congenital – fibular hemimelia
- Travelling, Cycling, Snowboarding, Adventurous
- Highly active, freelance web designer
- Goals: Mobilise without discomfort and perform recreational activities



**Tim**

“I do my best be as active as possible, but I get frustrated when my prosthesis limits me “



## High Activity Solutions – Direct Socket TF



### Direct Socket TF

- Easy to Donn & Doff as compared to traditional socket
- Superior Functional Outcomes
  - Improved ROM/Suspension/Stability in all planes & Superior femur control
- Superior Comfort
  - Sitting, standing & walking
- Ultra-Lightweight & Low-Profile
- True Total Contact

## High Activity Solutions – Iceross Seal-In X5



### Iceross Seal-In X5

- Features five integrated seals
- Provides security, minimises pistoning and provides enhanced rotational control
- Ideal for active amputees looking for security during recreational activities

## High Activity Solutions – RHEO KNEE XC



### RHEO KNEE XC

- High active users looking for the increased stability in stance
- Adapts to variable cadence, for fast paced walking and running
- Negotiate ramps/stairs – capable of ascending stairs step over step
- Automatic Cycling recognition – for efficient power transfer to the pedals, allows free movement without resistance

## User Experiences: RHEO KNEE XC



'The microprocessor in the RHEO KNEE XC provides **safety and prevents falls**, and has made the biggest difference to my mobility. '

'...The RHEO KNEE XC **features artificial intelligence** and is thinking for you the whole time, measuring your speed as you go and **constantly adapting for you**'

'...My gait and my walking have also improved, and people have commented on how **normal my walking looks**'

'The RHEO KNEE XC provides resistance which allows me to **feel very comfortable and confident going down stairs and down slopes too**. I feel you can really get your resistance into the slope and it's quite nice ambling down the hill'

'...I have also enjoyed the cycling function because when you get on the bicycle, the **RHEO KNEE XC automatically enters the cycling mode** which provides the free resistance that you need for cycling'



## High Activity Solutions – Pro-Flex XC Torsion



C-shaped design and Torsion unit provide 16mm vertical compression



Shock absorption for high impact activities



C-shaped design provides progressive stiffening for a more natural roll over



Great ankle ROM and reduced sound side loading



Torsion version provides shock absorption and rotation



Reduces shear forces and stress on residual limb

### *PERMITTED ROTATION AND SHOCK MOVEMENTS REDUCE SHEAR STRESS*

- Twiste M. Transverse rotation and longitudinal translation during prosthetic gait – A literature review. Journal of Rehabilitation, 2003; 40: 9-18

### *LOWER SHEAR STRESSES → REDUCED DISCOMFORT*

- Segal AD, et al. Transtibial amputee joint rotation moments during straight-line walking and a common turning task with and without a torsion adapter. Journal of Rehabilitation R&D. 2009 46:375-384

### *SELF PERCEIVED PAIN REDUCTION*

- Segal AD, et al. Does a Torsion Adapter Improve Functional Mobility, Pain and Fatigue in Patients with Transtibial Amputation? Clin Orthop Relat Res 2014 472:3085–3092


### *REDUCED SHEAR STRESS*

- Heitzmann D, et al. Functional effects of a prosthetic torsion adapter in transtibial amputees during unplanned spin and step turns. Prosthet Orthot Int 2015; 1-8.

## High Activity Solutions



- Direct Socket TF provides superior comfort in sitting, standing and walking
- Iceross Seal-In X5 provides security during recreational activities
- RHEO KNEE XC provides high level of stability in stance, whilst keeping up with active users, looking to ascend stairs step-over-step, run and cycle
- Pro-Flex XC Torsion provides shock absorption and rotation for users participating in high impact activities and recreational sports





### Mauch Knee Plus

#### Features:

- Hydraulic swing control:
  - wide range of cadence variation and independent control of flexion and extension resistance
- Hydraulic stance control:
  - Yielding resistance allows foot over foot stair and ramp descent.
  - Stumble recovery feature improves security
- Hydraulics create a 'yielding' brake for control of knee flexion during stance
  - Stair and ramp descent
  - Sitting
  - Allows for a more natural approach to stairs and sitting





### Mauch Knee Plus

- Swing Only Mode:
  - Allows for cycling, rowing, or use of various exercise machines.
- Manual Lock Mode:
  - Allows for ultimate security when descending steep or slippery slopes, climbing ladders, or standing for extended periods
- Durable, suitable for high impact activities
  - Weight limit: 166 kg
  - Knee flexion range: 125°
  - 4 hole adapter connections for low build height
  - Heavy duty use
  - Recreational sports

## Alternative Knee Considerations: High Activity Solutions



### OH7

- K3 – K4 Level – rated to 136Kg
- Enhanced stability from closing geometry – Provides safety in stance
- 3 phase hydraulic swing phase accommodates changes in walking speed
- Capable of fast paced walking and running – up to 10km/hr
- Mid swing shortening for toe clearance in swing, helps to prevent trips and falls



### Total Knee 2100

- K3 – K4 Level
- Enhanced stability from Geometric Lock
- Mid-swing shortening for toe clearance
- Incorporates adjustable stance flexion for shock absorption
- 3 phase hydraulic swing phase control for faster speed walking
- Low build height



### Cheetah Knee

- K3 – 4 Level
- For jogging and distance running
- Enhanced Hydraulics for a faster swing speed
- 4-bar geometry with stable stance control for controlled deceleration
- Extension bumper for softer terminal impact
- Can be combined with the Flex Run for running and jogging

# Mechanical Knee Prescription Guide



## Mechanical Knee Prescription Guide




### PASO KNEE™

The PASO Knee features an innovative autoadaptive pneumatic swing phase control which requires no valve adjustments. The knee will automatically respond to any change in cadence up to 7m/hr with no adjustments required from the prosthetist. The PASO Knee is suitable for a wide range of users due to its stability and responsiveness.

Unique 'closing geometry' creates an inherently more stable system during stance than a standard 4-bar design. It ensures optimal stance phase stability from loading response through terminal stance until weight is applied on the toe in conjunction with a hip flexion moment to release the knee.

#### Unique Features

- Low proximal build height ideal for long residual limbs and knee disarticulation users
- Auto-adaptive pneumatic swing phase control adapts automatically
- Closing geometry (refer to page 22)
- Choice of 4 socket connection options, providing a variety of shift, rotation, tilt and stance flexion adjustments



#### MOBILITY GOALS AND REQUIREMENTS

- Requires enhanced knee stability beyond basic 4-bar during early stance
- Ability or potential for moderate - high cadence variation
- Regularly traversing different terrain: level and uneven ground, shallow slopes
- Higher loading or impact requirements

#### RECOMMENDED FEET



#### TECHNICAL SPECIFICATION

<b>WEIGHT LIMIT</b>	<b>KNEE FLEXION</b>
• 136kg	• 150°
<b>WEIGHT</b>	<b>SWING CONTROL</b>
• 1050g	• Auto-adaptive Pneumatic
<b>EFFECTIVE BUILD</b>	<b>STANCE CONTROL</b>
• 158mm	• 4-Bar Closing Geometry (↑ Stability) - adjustable

Part #	Description
1728180	PASO Knee with pyramid adapter
1728181	PASO Knee with IKF adapter
1728182	PASO Knee with loop adapter
1728183	PASO Knee with 3-arm adapter

K2 - K4 ACTIVITY LEVEL

11



# Össur Legs - Complete Solutions




## ÖSSUR LEGS

COMPLETE SOLUTIONS | CLINICALLY FLEXIBLE | SIMPLE & EASY



### For Clients



### For Clinicians & Clinics



# Össur Legs - Complete Solutions



## ÖSSUR LEGS

COMPLETE SOLUTIONS | CLINICALLY FLEXIBLE | SIMPLE & EASY



- Össur has a complete range of solutions to suit low, moderate, high users
- Commencing in January 2022 we will be launching the Össur legs program
- Please contact your local Össur Clinical Specialist for further information on Össur Legs

# Össur Legs – Clinical Specialist Contacts



## ÖSSUR LEGS

COMPLETE SOLUTIONS | CLINICALLY FLEXIBLE | SIMPLE & EASY



- Sarah Mulroy-Lang [smulroy-lang@ossur.com](mailto:smulroy-lang@ossur.com)
- (Senior Clinical Specialist –Prosthetics)
  - Area Manager for VIC, TAS and APC Clinics



- Lizzie Carey [ecarey@ossur.com](mailto:ecarey@ossur.com)
- (Clinical Specialist- Prosthetics)
  - Area Manager for NZ, NSW and ACT



- Laura Roberts [lroberts@ossur.com](mailto:lroberts@ossur.com)
- (Clinical Specialist –Prosthetics)
  - Area Manager for QLD, WA, NT and SA

# References

- Avita HersHKoviz et al, Rehabilitation outcome of post-acute lower limb geriatric amputee. Disability and rehabilitation. February 2013. Vol. 35, No. 3, Pages 221-227
- Spruit v E, Harmen van der Linde, Bianca Buijck, Alexander Geurts, Sytse Zuideman „Predicting prosthetic use in elderly patients after major lower limb amputation, Prosthetic Orthot Int March 2012 vol 36 no. 1 45-52
- Johannesson, A., G. U. Larsson, et al (2010) „Outcomes of a standardized surgical and rehabilitation program in transtibial amputation for peripheral vascular disease: a prospective cohort study“ Am J Phys Med Rehabil 89 (4): 293-303
- Rajiv Kumar Singh (2016) Long-term mortality after lower limb amputation. Prosthet Orthot Int 2016; Vol. 40(5): 545-551
- Rosenblatt, Noah J., et al. "Active dorsiflexing prostheses may reduce trip-related fall risk in people with transtibial amputation." J Rehabil Res Dev 51.8 (2014): 1229-1242.
- Alimusaj M, Fradet L, Braatz F, Gerner HJ, Wolf SI. Kinematics and kinetics with an adaptive ankle foot system during stair ambulation of trans-tibial amputees. Gait & Posture. 2009; 30:3:356-363.
- Twiste M. Transverse rotation and longitudinal translation during prosthetic gait – A literature review. Journal of Rehabilitation, 2003; 40: 9-18
- Segal AD, et al. Transtibial amputee joint rotation moments during straight-line walking and a common turning task with and without a torsion adapter. Journal of Rehabilitation R&D. 2009 46:375-384
- Segal AD, et al. Does a Torsion Adapter Improve Functional Mobility, Pain and Fatigue in Patients with Transtibial Amputation? Clin Orthop Relat Res 2014 472:3085–3092
- Heitzmann D, et al. Functional effects of a prosthetic torsion adapter in transtibial amputees during unplanned spin and step turns. Prosthet Orthot Int 2015; 1-8.
- Rosenblatt, Noah J., et al. "Active dorsiflexing prostheses may reduce trip-related fall risk in people with transtibial amputation." J Rehabil Res Dev 51.8 (2014): 1229-1242.
- Alimusaj M, Fradet L, Braatz F, Gerner HJ, Wolf SI. Kinematics and kinetics with an adaptive ankle foot system during stair ambulation of trans-tibial amputees. Gait & Posture. 2009; 30:3:356-363.

# References

- Ludviksdottir A, Gruben K, Gunnsteinsson K, Ingvarsson Th, Nicholls M. Effects on user mobility and safety when changing from a carbon fiber prosthetic foot to a bionic prosthetic foot. Presented at Orthopadie&Reha-Technik Congress, Leipzig, May 2012.
- Miller, William C., Mark Speechley, and Barry Deathe. "The prevalence and risk factors of falling and fear of falling among lower extremity amputees." Archives of physical medicine and rehabilitation 82.8 (2001): 1031-1037.
- K. Kaufman, B. Mundell, S. Visscher, H. M. Kremers, D. Larson, and J. Ransom, "Risk factors and costs associated with accidental falls among adults with above-knee amputations: a population-based study," Mayo Clinic, Rochester, MN, Apr. 2015.
- Jonkergouw, N., Prins, M. R., Buis, A. W., & van der Wurff, P. (2016). The effect of alignment changes on unilateral transtibial amputee's gait: a systematic review. PloS one, 11(12), e0167466.
- Wolf, S. I., Alimusaj, M., Fradet, L., Siegel, J., & Braatz, F. (2009). Pressure characteristics at the stump/socket interface in transtibial amputees using an adaptive prosthetic foot. Clinical Biomechanics, 24(10), 860-865.
- Ebbeling, C. J., Hamill, J., & Crussemeyer, J. A. (1994). Lower extremity mechanics and energy cost of walking in high-heeled shoes. Journal of Orthopaedic & Sports Physical Therapy, 19(4), 190-196.
- Anton, Johannesson. (2020). TRANSFEMORAL SOCKET FABRICATION METHOD USING DIRECT CASTING OUTCOMES REGARDING PATIENT SATISFACTION WITH DEVICE AND SERVICES. Canadian Prosthetics & Orthotics Journal. 3. 6. 10.33137/cpoj.v3i.
- Walker, Joel & Marable, W. & Smith, Christian & Sigurjónsson, Benedikt & Atlason, Ingi & Anton, Johannesson. (2021). CLINICAL OUTCOME OF TRANSFEMORAL DIRECT SOCKET INTERFACE (PART 2). CANADIAN PROSTHETICS & ORTHOTICS JOURNAL. 4. 10.33137/cpoj.v4i1.36065.



**50**  
**1971-2021**



**ÖSSUR®**  
LIFE WITHOUT LIMITATIONS