




Instructions for Use

AFO DYNAMIC

		3
EN	Instructions for Use	4
DE	Gebrauchsanweisung	6
FR	Notice d'utilisation	9
ES	Instrucciones para el uso	12
IT	Istruzioni per l'uso	15
NO	Bruksanvisning	18
DA	Brugsanvisning	20
SV	Bruksanvisning	22
EL	Οδηγίες Χρήσης	24
FI	Käyttöohjeet	27
NL	Gebruiksaanwijzing	30
PT	Instruções de Utilização	33
PL	Instrukcja użytkownika	36
CS	Návod k použití	39
HR	Upute za uporabu	42
JA	取扱説明書	44
ZH	中文说明书	46
KO	사용 설명서	48
RO	Indicații de utilizare	51
AR	تعليمات الاستخدام	54

ENGLISH

MD Medical Device

INTENDED USE

The device is intended to support drop-foot.

The device must be fitted and adjusted by a healthcare professional.

Indications for use

Drop foot.

Contraindications

Loss of sensation, spasticity of the calf muscles, foot neuropathy

Warnings and Cautions:

- Customization beyond that specified in this IFU is not recommended.

GENERAL SAFETY INSTRUCTIONS

The healthcare professional should inform the patient about everything in this document that is required for safe use of this device.

Any serious incident in relation to the device must be reported to the manufacturer and relevant authorities.

The patient should stop using the device and contact a healthcare professional:

- If there is a change or loss in device functionality, or if the device shows signs of damage or wear hindering its normal functions.
- If any pain, skin irritation, or unusual reaction occurs with the use of the device.

The device is for single patient – multiple use.

FITTING INSTRUCTIONS

To obtain the best effect of the device, proper footwear is critically important to the overall success. Shoes should be well constructed to include:

- Firm heel counter, for proper control of the rear foot.
- The foot and the device should be secured by the shoe.
- Rocker-type sole at toe end for smoother transition from 2nd to 3rd rocker.
- Use shoes that allow for adjustable compressive support at mid-foot.
- Rubber sole, to minimize the chance of slipping on wet surfaces.
- Removable insole, to allow space for the footplate with modifications.
- Adequate height and adjustability to avoid pressure from the upper leather on the upper side of the foot. The heel height should be approximately 1,7 cm, wedges can be used as needed for adjustment.

Device Adjustments

Foot plate:

- If present remove the insole from the shoe and use as a template for adjusting the footplate shape (A).
- The foot plate can be ground using a grinding machine to fit the shoe shape. All edges should be smoothed so no sharp edges are left.

- Avoid grinding excessive material from the lateral side of the footplate, ensuring sufficient space between the medial strut (B) and the medial arch (**Fig. 1**).
- To prevent movement of the device in the shoe the footplate should cover the whole length of the sole (**Fig. 2**).
- If possible, replace the insert on top of the footplate. Foot orthoses can be used on top of the footplate if required to support the foot function.

Tibial shell:

- Apply padding and straps (C, D) as needed. Additional padding can be added to the medial strut if required.
- Shift the anterior/posterior position of the device to obtain a close fit between the anterior tibial shell (E) and the lower leg (**Fig. 3**).
- Attention should be paid to the position of the medial strut relative to the medial malleolus, especially for a patient with a pronated foot position. Up to 1 cm of the strut may be ground off providing it is at least 7 cm from the footplate (**Fig. 4**).
- The lower lateral support (F) can be removed if less control is required (**Fig. 5**).

Note: The device cannot be reshaped by heating. Heating will cause delamination and negatively alter the device's dynamics.

Device Application

- Insert the device into the shoe.
- The patient should be seated with the knee flexed to 90° and foot in a neutral position.
- Align the anterior tibial shell on the tibia and fasten the lower strap (**Fig. 6**).
- Fasten the upper calf strap, the distance between the two anterior support members can be adjusted with the strap to fit the tibial contouring (**Fig. 7, 8**).

Accessories and Replacement Parts

Please refer to the Össur catalog for a list of available replacement parts or accessories.

USAGE

Cleaning and care

- Hand-wash using mild detergent and rinse thoroughly.
- Air dry.

Note: Do not machine-wash, tumble dry, iron, bleach, or wash with fabric softener.

Note: Avoid contact with salt water or chlorinated water. In case of contact, rinse with fresh water and air dry.

DISPOSAL

The device and packaging must be disposed of in accordance with respective local or national environmental regulations.

LIABILITY

Össur does not assume liability for the following:

- Device not maintained as instructed by the instructions for use.
- Device assembled with components from other manufacturers.
- Device used outside of recommended use condition, application, or environment.

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