


Instructions for Use



---

AFO Light

		3
EN	Instructions for Use	5
DE	Gebrauchsanweisung	8
FR	Notice d'utilisation	11
ES	Instrucciones para el uso	14
IT	Istruzioni per l'uso	17
NO	Bruksanvisning	20
DA	Brugsanvisning	23
SV	Bruksanvisning	26
EL	Οδηγίες Χρήσης	29
FI	Käyttöohjeet	33
NL	Gebruiksaanwijzing	36
PT	Instruções de Utilização	39
PL	Instrukcja użytkowania	42
HR	Upute za uporabu	45
ZH	中文说明书	48
KO	사용 설명서	51
BG	Указания за употреба	54
AR	تعليمات الاستخدام	57

## **INTENDED USE**

The device is intended to support drop-foot.

The device must be fitted and adjusted by a healthcare professional.

## ***Indications for use***

Drop foot.

## ***Contraindications***

Loss of sensation, spasticity of the calf muscles, and foot neuropathy.

## **Warnings and Cautions:**

- Make sure the strut doesn't reach the back of the knee in a seated or slightly kneeled position and does not compress the perineal artery or nerve.
- Patients with foot neuropathy should not use this device.

## **GENERAL SAFETY INSTRUCTIONS**

The healthcare professional should inform the patient about everything in this document that is required for safe use of this device.

Any serious incident in relation to the device must be reported to the manufacturer and relevant authorities.

The patient should stop using the device and contact a healthcare professional:

- If there is a change or loss in device functionality, or if the device shows signs of damage or wear hindering its normal functions.
- If any pain, skin irritation, or unusual reaction occurs with the use of the device.

The device is for single patient – multiple use.

This device is a supportive brace only and is not intended or guaranteed to prevent injury.

## **ASSEMBLY INSTRUCTIONS**

### **FITTING INSTRUCTIONS**

While carrying out the following instructions, please refer to the overview figure for locating components mentioned in the text (**Fig. 1**).

To obtain the best effect of the device, proper footwear is critically important to the overall success. Shoes should be well constructed to include:

- Firm heel counter, for proper control of the rear foot.
- The foot and the device should be secured by the shoe.
- Rocker-type sole at toe end for smoother transition from 2nd to 3rd rocker.
- Use shoes that allow for adjustable compressive support at mid-foot.
- Rubber sole, to minimize the chance of slipping on wet surfaces.
- Removable insole, to allow space for the footplate with modifications.
- Adequate height and adjustability to avoid pressure from the upper leather on the upper side of the foot. The heel height should be approximately 1,7 cm, wedges can be used as needed for adjustment.

## **Device Adjustments**

**Note:** The device cannot be reshaped by heating. Heating will cause delamination and negatively alter the dynamics of the device.

### **Foot Plate Customization:**

1. If present, remove the insole from the shoe and use it as a template for adjusting the Foot Plate (A) shape.

**Note:** The Foot Plate can be ground using a grinding machine to fit the shoe shape. All edges should be smoothed such that no sharp edges are left.

**Note:** Do not grind the device at the heel area as this may jeopardize product performance and strength.

2. As applicable, re-insert the insole on top of the Foot Plate. Foot orthoses can be used on top of the Foot Plate if required to support the foot function.

### **AFO Height Customization and Assembly:**

3. Slide the Calf Piece (B) over the Posterior Strut (C) (**Fig. 2**).
4. Verify the height of the Calf Piece (**Fig. 3**), the top of the Calf Piece should be approximately 3 cm below the tibial head. To reduce the length of the posterior strut, remove the Calf Piece and grind Posterior Strut to the appropriate level.

**Warning:** Make sure the Posterior Strut doesn't reach the back of the knee in a seated or slightly kneeled position and does not compress the perineal artery or nerve.

5. Slide the Calf Piece over the Posterior Strut (**Fig. 2**).
6. Apply slight tension on the Elastic Connection Strap (D) and attach the PSA Hook (E).
7. Connect the Elastic Connection Strap using the loop (**Fig. 4**).
8. Apply the tapered part (rounded edge) of the Posterior Strut Pad (F) to the distal part of the Posterior Strut and have it cover the entire length of the Posterior Strut as well as the Elastic Connection Strap.
9. Attach the Calf Pad (G) and the Calf Strap (H) to the Calf Piece.  
**Note:** The padding is slightly larger than the Calf Piece to avoid contact between the piece and the leg.
10. Insert the device into the shoe (**Fig. 5**) and allow the patient to walk around for a while to verify the fit of the device. Decrease or increase padding as needed.  
**Note:** In case of point pressure at the heel bone, it is advised to apply padding under the existing padding just proximal to the pressure point (Distal calf). Remove the padding and glue additional padding to the Posterior Strut before re-attaching the padding.

When final fit has been confirmed:

11. Remove the upper part of the padding of the Posterior Strut (**Fig. 6**).
12. Remove the Elastic Connection Strap from the PSA Hook on the Posterior Strut.
13. Remove the Calf Piece from the Posterior Strut.
14. Slide the Shrink Tube (I) over the Posterior Strut and Posterior Strut Pad until it covers the heel area.
15. Use a heat gun to shrink the tube to the proper size. Evenly distribute the heat to avoid burning and wrinkles (**Fig. 7**). First apply the heat to the inside of the AFO, then to the outside of the AFO.
16. Re-apply the Calf Piece, Elastic Connection Strap, and padding.

17. If necessary, remove excessive strap material with a pair of scissors.  
**Note:** For patients with lower activity levels, where the Calf Piece doesn't move, it is advised to permanently fixate the Calf Piece by application of strong glue (Loctite 495.) into the Calf Piece before attaching it to the Posterior Strut (**Fig. 8**). The Calf Piece can also be riveted to the Posterior Strut. The Elastic Connection Strap and the PSA Hook may be removed in this case.

### ***Device Application***

1. Insert the device into the shoe.
2. The patient should be seated with the knee flexed to 90° and the foot in a neutral position.  
**Note:** When donning the device, use the posterior strut of the AFO as a shoe-horn (**Fig. 9**).
3. Secure the Calf Strap around the calf (**Fig. 10**).

### ***Accessories and Replacement Parts***

Please refer to the Össur catalog for a list of available replacement parts or accessories.

## **USAGE**

### ***Cleaning and care***

- Hand-wash using mild detergent and rinse thoroughly.
- Air dry.

**Note:** Do not machine-wash, tumble dry, iron, bleach, or wash with fabric softener.

**Note:** Avoid contact with salt water or chlorinated water. In case of contact, rinse with fresh water and air dry.

## **DISPOSAL**

The device and packaging must be disposed of in accordance with respective local or national environmental regulations.

## **LIABILITY**

Össur does not assume liability for the following:

- Device not maintained as instructed by the instructions for use.
- Device assembled with components from other manufacturers.
- Device used outside of recommended use condition, application, or environment.

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