





Instructions for Use

CTI[®] CUSTOM
CTI[®] OA


		3
EN	Instructions for Use	6
DE	Gebrauchsanweisung	10
FR	Notice d'utilisation	14
ES	Instrucciones para el uso	19
IT	Istruzioni per l'uso	23
NO	Bruksanvisning	27
DA	Brugsanvisning	31
SV	Bruksanvisning	35
EL	Οδηγίες Χρήσης	39
FI	Käyttöohjeet	44
NL	Gebruiksaanwijzing	48
PT	Instruções de Utilização	52
PL	Instrukcja użytkowania	57
CS	Návod k použití	61
JA	取扱説明書	65
SK	Indikácie použitia	69
HU	Használati javallat	73
RO	Indicații de utilizare	77

20

INTENDED USE

CTi Custom

The device is intended for external support, stabilization and protection of the knee.

CTi OA

The device is intended for unicompartmental unloading of the knee

The device must be fitted and adjusted by a healthcare professional.

Indications for use

- For knee conditions that may benefit from increased AP and ML stability around the knee, such as: ACL, MCL, LCL, PCL, rotary, and combined instabilities.
- The OA model also treats mild to severe unicompartmental osteoarthritis.

No known contraindications.

Warnings and Cautions:

- The use of the device may increase the risk of deep vein thrombosis and pulmonary embolism.

GENERAL SAFETY INSTRUCTIONS

The healthcare professional should inform the patient about everything in this document that is required for safe use of this device.

Any serious incident in relation to the device must be reported to the manufacturer and relevant authorities.

The patient should stop using the device and contact a healthcare professional:

- If there is a change or loss in device functionality, or if the device shows signs of damage or wear hindering its normal functions.
- If any pain, skin irritation, excessive pressure, or unusual reaction occurs with the use of the device.

The device is for single patient – multiple use.

This device is not guaranteed to prevent injury.

FITTING INSTRUCTIONS

While carrying out the following instructions, please refer to the overview figure for locating components mentioned in the text (**Fig. 1**).

Fitting the CTi Custom and The CTi OA is the same procedure, except the CTi OA frame alignment is adjusted at manufacture to provide unloading of the affected compartment of the knee.

Device Application

1. Unfasten all straps from the outer side of the device.
2. Ask the patient to bend the knee at a 90° angle, with the foot flat on the floor (**Fig. 2**).
3. Ensure that the CTi logo goes on the outside of the leg. Place the device on the patient's leg such that the patella is centered between the Condyle Pads (A) (**Fig. 3**).

Ensure proper alignment of the device on the leg:

4. Height positioning: Align the center of the Hinge slightly above the middle of the patella (**Fig. 3**).
5. The device is pre-assembled with thick Condyle Pads. There are also thin and medium Condyle Pads, as well as a hook-and-loop, shim included with the device. Use any combination of the Condyle Pads and shim to widen or narrow the width of the device at the knee.
6. The supplied tibial wedge pads can be placed between liner to adjust the fit around the tibia (**Fig. 11**). They can also be used to position the hinge center forward slightly or removed to position the hinge further backwards (**Fig. 12**).
Fasten the straps in the following order. All straps can be adjusted and/or trimmed to length.
7. Fit the Lower Cruciate Strap (B) first, located directly below the knee (**Fig. 4**). This strap is located above the calf muscle and aids in the correct positioning of the device as it helps create the device-to-bone contact required to help keep the device in place.
8. Next, snugly fasten the Bottom Strap (C) by inserting it through the D-Ring.
9. Continue with the Upper Cruciate Strap (D), located just above the knee.
10. Fasten the Top Strap (E). Make sure that the strap hook tab does not overlap the elastic insert on the strap.

Note: Overtightening the top two straps can cause the device to migrate. The following modifications come pre-installed, if selected during the ordering process of the device.

Device Application with ACL Cable System

1. Apply the device according to the device application procedure. The four straps should be fastened, but the secondary cable strap should still be loose.
2. Once the device is firmly in place, pass the cable strap through the smaller D-Ring.
3. Tighten the cable strap to place the desired amount of pressure on the tibia.

Device Application with PCL Strap Kit

1. Remove the Lower Cruciate Strap pad, attach the PCL pad onto the straps below the knee.
2. Apply the device according to the device application procedure.
3. Tighten the PCL strap with the leg in full extension so that the strap and pad are comfortably snug. The pad should be centered within the frame. If the PCL strap is too long, trim the strap on the alligator tab side. Replace the alligator tab and fasten. Make sure that the alligator tab is under the pull tab.

Device Application with Ski Boot Attachment

1. Apply the device according to device application procedure.
2. Apply ski boots, ask the patient to bend his knees and lean forward as though skiing. Position the bottom bearing of the aluminum rod on the hard plastic shell of the boot and mark the bent leg position of the bearing (**Fig. 15**).
Note: this mark should be positioned to allow for at least ½" (13mm) of plastic shell surrounding it on all sides.
3. Remove the boot and device and drill a 3/16th" (5mm) hole at the spot marked on the boot. Use care so as not to drill the tongue and the shell (**Fig. 16**).

4. Remove clevis pin from the bottom bearing and press it from the inside out through the newly drilled hole. Slip the bearing on the rod portion of the attachment over the clevis pin and secure it with the stay ring (**Fig. 17**).
5. Ensure the rod portion is attached to the device. Slip the fabric loop over the tibial extension of the device with the pad on the inside of the device.
6. Re-apply ski boot, slide the tube and rod portions of the attachment together by lowering the device down on the boot, be sure the tube portion of the attachment passes through the fabric loop between the buckles, (**Fig. 18**) and firmly place the device over the leg and apply it according to the standard fitting instructions (**Fig. 19**).

Note: For non-skiing usage, the rod portion of the attachment can be unbolted from the device using the 1/8" (3mm) hex key provided. The buckle system is attached via a hook-and-loop closure for easy removal. A separate pad is provided with hook-and-loop fastener for attachment to the brace in place of the buckle system. To adjust the amount of internal or external rotation allowed by the ski boot attachment, loosen or tighten respective buckles.

Device Application with Flexion Stop Kit

1. Set flexion angle using the included goniometer. The holes are positioned in approximately 10–14° increments with "A" at approximately 0° and "I" approximately equal to 90°.
2. Place the slotted screw on outside of plate in the desired flexion angle and attach to the roller on the inside of the hinge plate and tighten. If needed use one or more washers to ensure that the roller sits in position on the frame.
3. Repeat the same process for the other side.
4. Check the angle and that the rollers are engaging nicely with the hinge arms. If required the flexion stop plate can be contoured to optimize positioning.
5. Apply the device according to the device application procedure.

Device Adjustments

Range of Motion (ROM) Adjustment

The device comes with a 5-10° extension stop installed and includes the following additional steps:

0°, 5–10° (installed), 10–20°, 20–30°, and 30–40° (**Fig. 20**)

1. Remove the Condyle Pad from the hinge.
2. Use a short size 2 Philips head screwdriver to remove the screw securing the extension stop (**Fig. 13**).
3. Remove the extension stop out of the hinge.
4. Select the desired extension stop and place the stop into the hinge from the back with the small lip facing down (**Fig. 14**).
5. Position the hole in the extension stop with the screw hole in the hinge plate. Re-insert the screw and fasten. Making sure to not overtighten.
6. Replace the condyle pad.
7. Repeat this procedure for the opposite hinge.

Note: Ensure the same extension stop is used for both hinges.

How to Replace or Move the Field-Serviceable D-Rings

The field-serviceable D-Rings can be changed from internal to external configuration in case of discomfort or replaced if a D-Ring is broken. If the D-Rings are riveted onto the device, please reach out to Customer Service or CTi distributor for assistance.

1. Use a size 1 Philips head screwdriver to undo the central screw (**Fig. 5**). Remove the D-Ring and both posts from the device (**Fig. 6**).
2. Select a same sized D-Ring. Place the female post inside the D-Ring (**Fig. 8**).
3. The default D-Ring placement is on the inside of the device. However, the D-Rings can be moved to the outside of the frame, if preferred (**Fig. 9**).
4. Place the D-Ring with the female post into the hole on the frame.
5. Place the male post through the outside of the frame and slide the two together. We suggest using the longer male post on all Pro Sport models or when the frame is too thick for the short male post.
6. Re-insert the central screw and fasten. Making sure not to over tighten.

Device Removal

1. Ask the patient to sit down with the leg bent 80° with the foot flat on the floor.
2. Unfasten all straps from the outer side of the device and attach the straps back to themselves. This will make it easier to put the device on next time and will also prolong the life of the straps.

Accessories and Replacement Parts

Please refer to the Össur catalog for a list of available replacement parts or accessories.

USAGE

Refurbishing

The device should be regularly sent to Össur for refurbishing. The frequency of which depends on amount of activity by patient. Your healthcare professional should help determine frequency of service; the recommended minimum frequency is yearly. Refurbishment is provided at the cost of the customer. There is a usual turnaround time of three to five working days, plus shipping time.

Cleaning and care

- Hand-wash using mild detergent and rinse thoroughly.
- Air dry.

Note: Do not machine-wash, tumble dry, iron, bleach, or wash with fabric softener.

Note: If the device is used in salt water or chlorinated water, rinse with fresh water and air dry.

Hinge

- Remove foreign materials (e.g., dirt or grass) and clean using fresh water.
- Lubricate with silicone spray or dry graphite if needed.

DISPOSAL

The device and packaging must be disposed of in accordance with respective local or national environmental regulations.

LIABILITY

Össur does not assume liability for the following:

- Device not maintained as instructed by the instructions for use.
- Device assembled with components from other manufacturers.
- Device used outside of recommended use condition, application, or environment.

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