








Instructions for Use










---

CHEETAH<sup>®</sup> XPLORE

		3
EN	Instructions for Use	5
DE	Gebrauchsanweisung	11
FR	Notice d'utilisation	17
ES	Instrucciones para el uso	23
IT	Istruzioni per l'uso	29
NO	Bruksanvisning	35
DA	Brugsanvisning	41
SV	Bruksanvisning	47
EL	Οδηγίες Χρήσης	53
FI	Käyttöohjeet	59
NL	Gebruiksaanwijzing	65
PT	Instruções de Utilização	71
PL	Instrukcja użytkownika	77
TR	Kullanım Talimatları	83
RU	Инструкция по использованию	89
JA	取扱説明書	95
ZH	中文说明书	101
KO	사용 설명서	106

### DESCRIPTION

The device is a prosthetic foot attached to the posterior wall of the socket. It is composed of a cheetah inspired toe blade and a heel blade secured using heel bolts.

The device consists of the following components (**Fig. 1**):

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This device must be used with a Foot Cover and a Spectra Sock.

### INTENDED USE

The device is intended as a part of a prosthetic system that replaces the foot and ankle function of a missing lower limb.

Suitability of the device for the prosthesis and the patient must be evaluated by a healthcare professional.

The device must be fitted and adjusted by a healthcare professional.

### *Indications For Use and Target Patient Population*

- Lower limb loss, amputation, or deficiency
- No known contraindications

The device is for moderate to extreme impact use, e.g., walking, running, and sports.

### GENERAL SAFETY INSTRUCTIONS

**Warning:** Use of a lower limb prosthetic device carries an inherent risk of falling which may lead to injury.

The healthcare professional should inform the patient about everything in this document that is required for safe use of this device.

**Warning:** If there is a change or loss in device functionality, or if the device shows signs of damage or wear hindering its normal functions, the patient should stop using the device and contact a healthcare professional.

The device is for single patient use.

**Caution:** Epoxy can cause allergic reactions when in direct contact with body.

### REQUIRED COMPONENTS

The following components are required to use the device:

- Lamination plate
- Attachment hardware

### DEVICE SELECTION

Verify that selected variant of the device is suitable for the impact level and weight limit according to the following table.

**Warning:** Do not exceed weight limit. Risk of device failure.

Incorrect category selection may also result in poor device function.

Weight kg	45-52	53-59	60-68	69-77	78-88	89-100	101-116	117-130	131-147
Weight lbs	99-115	116-130	131-150	151-170	171-194	195-220	221-256	257-287	288-324
Moderate Impact Level	1	2	3	4	5	6	7	8	9
High Impact Level	2	3	4	5	6	7	8	9	N/A
Extreme Impact Level	3	4	5	6	7	8	9	N/A	N/A

## ASSEMBLY INSTRUCTIONS

### **Foot Cover and Spectra Sock (Fig. 2)**

**Caution:** To avoid pinching fingers, always use a shoehorn.

1. Put the foot into the Spectra Sock.
2. Remove the foot cover attachment.
3. Use the applicable (straight) end of a shoehorn to put the foot with the Spectra Sock into the Foot Cover.
4. Move the shoehorn up to fully push the foot into the Foot Cover.
5. Put the foot cover attachment back onto the Foot Cover.
6. After alignment is complete, fix the Spectra Sock to the prosthesis to seal against dust and dirt.

**Note:** The Spectra Sock must be pulled up to prevent it from interfering with moving parts of the foot.

If required doff the Foot Cover as follows:

1. Remove the foot cover attachment.
2. Insert the applicable (edged) end of a shoehorn behind the foot.
3. Push the shoehorn down and pull the foot out of the Foot Cover.
4. Fully remove the Spectra Sock.

### **Prosthesis**

Assemble prosthesis with applicable devices.

**Warning:** Risk of structural failure. Components from other manufacturers have not been tested and may cause excessive load on the device.

**Warning:** Ensure proper attachment by following the applicable device assembly instructions.

### **Adjusting Foot Size**

Grind foot to length, with help from the grinding sticker placed on the top of the foot plate (Fig. 4).

**Note:** When using a Foot Cover the size of the foot needs to be 2 cm shorter when grinding the toe to length. This requirement has to be taken into account when grinding the sole plate.

### **Lamination**

The foot may be connected with either Lamination Connector or direct lamination.

#### **Lamination Connector**

The Lamination Connector must be laminated directly into the posterior wall of the socket.

**Warning:** Before starting make sure to orient the Lamination Connector the correct way. The shoulders of the inserts (marked with arrow) should be oriented towards the socket (Fig. 5).

1. Put the Lamination Connector in the middle third of the socket to make height adjustment easier.
2. Make sure the socket surface is flat. A flat surface below the Lamination Connector distributes the load to the connected foot.
3. Use a friction pad between the Lamination Connector and the foot.
4. Apply Loctite 410 and torque bolts to 45 Nm for bolt size 9,5 mm.

### **Cutting Foot Module**

- Check for clearance in swing phase.
- Trim foot module to the top end of the laminating connector.
- Cut foot module using a hack saw or band saw.

### **Direct Lamination**

Before you start you need a well fitting conventional prosthesis. Either an existing socket or alternatively a prosthesis with a check socket on conventional prosthesis can be used.

#### **1. Laminate the Socket (Fig. 6)**

- a. Laminate the socket using carbon braid. Add extra layers of carbon in the upper posterior surface where the foot will be attached.
- b. Roughen the upper posterior surface to improve the adhesion to the socket when the foot module is attached (shown by grey area).

#### **2. Trace around the Original Walking Prosthesis (Fig. 7)**

- a. Ensure that the prosthetic foot is resting firmly on the shoes insole and fasten the laces tightly.
- b. Lay the original walking prosthesis over the grid on the A1 alignment sheet (included)
- c. Mark the lowest points of the anterior and posterior brims as a reference.
- d. With the shoe on, trace a line around the entire prosthesis. Then remove the shoe and trace the outline of the foot cover, ensuring the prosthesis remains aligned with the original tracing.

#### **3. Fill in the Patient Information Table on the Alignment Sheet (Fig. 8)**

- a. The socket adduction can be observed by looking at the original walking prosthesis in the frontal plane.
- b. Use the angle guide on the A1 alignment sheet to help measure the foot rotation in relation to the socket.

#### **4. Attach Sole and Heel Counter (Fig. 9)**

- a. Glue the heel counter onto the upper surface of the heel, aligning the center line with the split in the carbon heel.
- b. Glue the sole onto the underside of the foot, again aligning the center line with the splits in the carbon fiber.
- c. Grind the sole and heel into a shape that gives a good fit in the patient's shoe. Use the sole template on the back of the product box.

#### **5. Cut the Pylon to Length (Fig. 10)**

- a. Align the new socket and the foot with the tracing of the original prosthesis.
- b. Mark the pylon at the height of the posterior brim of the socket using the drawing of the original prosthesis.
- c. Cut the pylon at this mark and round off the end, removing any sharp edges.

#### **6. Mark the Length of the Attachment Surface on the Foot Pylon (Fig. 11)**

- d. Mark the pylon at 90 mm from the top (cut end).
- e. Roughen all the surfaces of the pylon above this mark (shown by grey area).

## **7. Align the Foot and the Socket with the Tracing of the Original Walking Prosthesis (Fig. 12)**

- a. The sole should be lined up with the outline of the foot cover sole.
- b. The new socket should be positioned over the tracing of the original socket with additional 5-7° of knee extension in the socket.
- c. Use a foam spacer/s (or other spacer) placed between the pylon and the socket to maintain the correct extension. Spacers should be placed at the end of the roughened attachment surface. Ensure the minimum clearance of 5 mm is maintained between the distal posterior area of the socket and the foot pylon.
- d. Any spacer used at the top of the pylon can have maximum thickness of 5 mm.
- e. Fasten the foam to the foot using tape.
- f. Mark a clear height on the socket at height of a foam spacer on the pylon.

## **8. Align the Socket according to the Patient Information (Fig. 13)**

- a. Put the shoe on the foot and place it upright in the work surface.
- b. Hold the socket in position on the pylon according to:
  - The height
  - The recorded socket adduction
  - The recorded foot rotation in relation to socket (using angle guides on A1 sheet)
- c. The foam spacers will hold the socket in the correct extension.

## **9. Mark the Attachment Surface onto the Socket (Fig. 14)**

- a. Holding the socket in the desired position, trace around the end of the pylon onto the socket above the height mark.

## **10. Fix the Foot to the Socket (Fig. 15)**

- a. Use Composite 1 adhesive that is provided with the foot kit, read all information on the adhesive packaging.  
**Note:** Adhesive working time is 1 minute.
- b. Cut the tip of the adhesive mixer at 16 mm from the end.
- c. Apply the adhesive generously onto the marked attachment surface on the socket.
- d. Attach the foot onto the marked attachment area and make any small adjustment to ensure the alignment is correct.
- e. Smooth the surface of any adhesive that has been pushed out around the outside of the pylon.
- f. Hold the socket in position for 1 minute, then leave for a further 5 minutes for the adhesive to set.

## **11. Apply Adhesive around the outside of the Pylon (Fig. 16)**

- a. Apply adhesive generously onto the outside roughed surface of the pylon and spread outwards onto the roughened surface of the socket.
- b. Smooth the finish of the adhesive and allow to set for 5 minutes.

## **12. Test Alignment before Final Lamination (Fig. 17)**

- a. Remove any foam spacers and tape before testing the alignment on the patient.
- b. If needed, minor adjustments can be made by shaping the underside of the sole.
- c. If the alignment needs further adjustment, grind/cut the foot off the socket and make the desired changes before repeating steps 10 and 11.
- d. If the alignment is good, proceed to laminate over the attachment area onto the socket using carbon braid and a hand lay up.

## **ALIGNMENT INSTRUCTIONS**

### ***Bench Alignment (Fig. 3)***

#### **Sagittal Plane**

Determine vertical socket position in alignment jig.

Adjust to appropriate heel height.

Introduce appropriate socket angles flexion.

The load line should fall at the junction of the posterior and the middle third portions as shown in **Fig. 3**.

#### **Coronal Plane**

Introduce appropriate socket abduction/adduction using alignment jig.

Determine appropriate internal or external socket rotation.

## **USAGE**

### ***Cleaning and care***

Clean with a damp cloth and a mild soap. Dry with a cloth after cleaning.

### ***Environmental Conditions***

The device is Waterproof.

A Waterproof device can be used in a wet or humid environment and submerged in up to 3-meter-deep water for a maximum of 1 hour.

It can tolerate contact with: Salt water, chlorinated water, perspiration, urine, and mild soaps.

It can also tolerate occasional exposure to sand, dust, and dirt.

Continuous exposure is not allowed.

Dry with a cloth after contact with fresh water or humidity.

Clean with fresh water after exposure to other liquids, chemicals, sand, dust, or dirt and dry with a cloth.

## **MAINTENANCE**

The device and the overall prosthesis should be examined by a healthcare professional. Interval should be determined based on patient activity.

### ***Noise from Foot***

Noise may occur if sand or debris is present in device. In that case, the healthcare professional should doff the foot, clean it with the help of compressed air and replace the Spectra sock if it is damaged.

## **REPORT OF SERIOUS INCIDENT**

Any serious incident in relation to the device must be reported to the manufacturer and relevant authorities.

## **DISPOSAL**

The device and packaging must be disposed of in accordance with respective local or national environmental regulations.

## **LIABILITY**

Össur prosthetic devices are designed and verified to be safe and compatible in combination with each other and custom-made prosthetic sockets with Össur adapters, and when used in accordance with their intended use.

Össur does not assume liability for the following:

- Device not maintained as instructed by the instructions for use.
- Device assembled with components from other manufacturers.
- Device used outside of recommended use condition, application, or environment.

### Compliance

This device has been tested according to ISO 10328 standard to two million load cycles.

Depending on patient activity this may correspond to 2-3 years of use.

ISO 10328 - Label		
Category	Weight (Kg) Extreme Impact	Label Text
1	52	ISO 10328-P3-52 kg
2	59	ISO 10328-P3-59 kg
3	68	ISO 10328-P4-68 kg
4	77	ISO 10328-P5-77 kg
5	88	ISO 10328-P6-88 kg
6	100	ISO 10328-P7-100 kg
7	116	ISO 10328-P8-116 kg
8	130	ISO 10328-P8-130 kg
9	147	ISO 10328-P8-147 kg

\*Actual test loads reflect maximum body mass

ISO 10328 - "P" - "m"kg \*)

\*) Body mass limit not to be exceeded!

For specific conditions and limitations of use see manufacturer's written instructions on intended use!

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