




Instructions for Use

CHEETAH[®] KNEE

		3
EN	Instructions for Use	4
DE	Gebrauchsanweisung	8
FR	Notice d'utilisation	13
ES	Instrucciones para el uso	18
IT	Istruzioni per l'uso	23
NO	Bruksanvisning	28
DA	Brugsanvisning	32
SV	Bruksanvisning	36
EL	Οδηγίες Χρήσης	40
FI	Käyttöohjeet	45
NL	Gebruiksaanwijzing	49
PT	Instruções de Utilização	54
PL	Instrukcja użytkowania	59
TR	Kullanım Talimatları	64
RU	Инструкция по использованию	68
JA	取扱説明書	73
ZH	中文说明书	77
KO	사용 설명서	81

DESCRIPTION

The device is a polycentric (four-bar) prosthetic knee with a 3-phase hydraulic swing control. It is intended to extend the patient's activity to sport-specific applications such as long-distance running and sprinting. Proximal adapter options are available, and the device has a Male Pyramid distal connection.

INTENDED USE

The device is intended as a part of a prosthetic system that replaces knee function of a missing lower limb.

Suitability of the device for the prosthesis and the patient must be evaluated by a healthcare professional.

The device must be fitted and adjusted by a healthcare professional.

Indications For Use and Target Patient Population

- Lower limb loss, amputation, or deficiency
- No known contraindications

The device is for moderate to extreme impact use, e.g., walking, running, and sports.

The weight limit for the device is 100 kg.

GENERAL SAFETY INSTRUCTIONS

Warning: Use of a lower limb prosthetic device carries an inherent risk of falling which may lead to injury.

The healthcare professional should inform the patient about everything in this document that is required for safe use of this device.

Warning: If there is a change or loss in device functionality, or if the device shows signs of damage or wear hindering its normal functions, the patient should stop using the device and contact a healthcare professional.

Warning: Avoid placing hands or fingers near moving joints.

Caution: Do not adjust other screws than described in these instructions. The device is for single patient use.

ASSEMBLY INSTRUCTIONS

Warning: Risk of structural failure. Components from other manufacturers have not been tested and may cause excessive load on the device.

Proximal Adapter Options

The following proximal adapter options are available for the device (**Fig. 1**):

- 4-Prong Socket Adapter (1A),
- Female Pyramid Adapter (1B),
- Euro 4-Hole Adapter (1C), or
- Male Pyramid Adapter (1D).

Caution: Only use the specified Total Knee proximal adapters, please refer to Össur catalogue for more information.

1. Fully engage the threads of the selected proximal adapter.

Caution: Adapters from other manufacturers may not fully engage all the threads and compromise the strength of the prosthesis.

2. Align the Pinch Bolt (1E) of the Male Pyramid Adapter (1D), 4-Prong Socket Adapter (1A), and the Euro 4-hole Adapter (1C) perpendicular

to the anterior/posterior axis of the device, except for the Female Pyramid Adapter (1B), which should be angled 45°.

Caution: Maximum connection strength will not be established unless the adapter is correctly aligned.

3. For secure fastening to the device, tighten the adapter's Pinch Bolt according to its instructions for use.

Caution: Do not tighten the proximal adapter against the housing of the device.

Caution: When prong adapters are used, resin may get under the adapter and hinder full thread engagement. Carefully remove resin without damaging the adapter so that adapter threads can be fully engaged when mounted on the device.

ALIGNMENT INSTRUCTIONS

Bench Alignment (Fig. 2)

Alignment Goal

Alignment reference line (B) should:

- pass through midpoint of socket at the ischial tuberosity level (D)
- pass through pivot axis (A)
- fall at the 1/3 mark on the inside of the foot cover.

Note: Prioritize knee alignment over foot alignment if there is a mismatch.

Alignment Instructions

1. Position the foot so that the alignment reference line (B) falls at the 1/3 mark on the inside of the foot cover (with the foot cover and shoe on). Consider the external rotation of the foot.
2. Use the applicable adapters to connect the knee to the foot and establish the correct knee center height.
3. Position the knee so that the alignment reference line passes through pivot axis (A)
4. On the lateral side of the socket, make a first mark at the midpoint of the socket at the ischial tuberosity level (D). Make a second mark at the midpoint of the socket distally (E). Draw a line through both marks.
5. Position the socket so the alignment reference line (B) goes through first mark at the midpoint of the socket at the ischial tuberosity level (D).
6. Adjust the socket flexion to 5° in addition to the existing position (i. e., hip flexion contracture) and set the height of the full prosthesis.
7. Use the applicable adapters to connect the knee to the socket.

Warning: After adjustments, all screws must be secured with a medium strength threadlocker and tightened with the correct torque.

Tighten the screws with the following torque:

- Proximal Adapter Pinch Bolt (**Fig. 1:E**): Please refer to the respective proximal adapter's instructions for use.

Static Alignment

- Make sure the patient stands with equal weight on both legs.
- Check for correct prosthesis length.
- Check internal / external rotation.
- Check for correct load on toe and heel.

Dynamic Alignment

Ensure the patient is familiar with the functioning of the device.

Ensure the patient can activate and deactivate the geometric lock.

Carrying out the initial dynamic assessment is important for getting the patient to feel more confident in relying on the prosthesis. The best way

to acquire a dynamic alignment is by having a symmetric gait pattern. Any deviations will be more pronounced in running gait.

1. Evaluate the patient's level of voluntary control, the patient should be able to maintain knee stability.
2. Use hip extension to maintain stability into stance.
3. Train the patient to maintain normal step length on the sound side.

Caution: For safety, please make initial assessments and adjustments with the patient standing between parallel bars.

Swing Control (Fig. 3)

The hydraulic system is supplied in its basic setting.

In the basic setting, Valves "F" and "H" are both opened by $\frac{3}{4}$ turn. Valve "E" is opened by $1\frac{1}{2}$ turns.

Adjusting the Flexion and Extension valves

The following valves can be adjusted:

- Valve "F": Affects swing flexion resistance above 60°.
- Valve "H": Affects swing flexion resistance from 0° - 60°.
- Valve "E": Affects swing extension resistance.

Start with adjusting the flexion with the patient walking slowly and at faster walking speeds.

For each adjustment, turn the valve in small increments of about a $\frac{1}{4}$ of a turn. Always check the result right after each adjustment step.

Possible observations and actions to take are outlined in the following:

- If the patient is walking fast and you observe excessive heel rise, increase the swing flexion resistance by turning Valve "F" to the right until heel rise is normalized.
- To smooth flexion resistance changes between 0° and 60°, turn Valve "H" to the right or left to increase or decrease the resistance, respectively.

Afterwards, adjust the extension to harmonize the gait.

- Increase the swing extension resistance by turning Valve "E" to the right until terminal impact at full extension is reduced. Patient should feel a slight bump at full extension.
- Fine tune Valves "F" and "E" until a smooth and secure gait pattern is achieved for slow and fast walking speeds.

Caution: The valves should never be fully closed. Over tightening valves or flexing knee with all valves completely closed may damage valves.

Warning: Valves adjustment must always allow for complete extension of the knee. Excessive swing extension resistance will prevent knee from moving into full extension and may compromise stance security.

USAGE

Cleaning and care

Wipe the device with a soft cloth. Do not use solvents.

Environmental Conditions

The device should not come in contact with fresh water, salt water or chlorinated water.

Caution: Device should not be used in dusty environment. Exposure to sand, talcum, or similar should be avoided.

The device can be used in temperatures between -15°C to 50°C.

MAINTENANCE

The device and the overall prosthesis should be examined by a healthcare professional. Interval should be determined based on patient activity.

Recommended interval is every 6 months.

Check for damage, excessive wear and dirt.

If the device is subjected to excessive moisture or a corrosive environment, it is recommended to clean and lubricate the device frequently.

Two sets of roller bearings are visible on the device's left or right side (**Fig. 4: A**). Place few drops of general purpose oil or sewing machine oil on the roller bearings. Move the knee several times and wipe off excess oil with a soft cloth. No other parts need external lubrication.

Wipe the device with a soft cloth moistened with small amount of general purpose oil or sewing machine oil.

Caution: Do not use compressed air to clean device. Air forces pollutants into bearings and may cause malfunctions and wear.

Caution: Do not use talcum powder for lubrication.

REPORT OF SERIOUS INCIDENT

Any serious incident in relation to the device must be reported to the manufacturer and relevant authorities.

DISPOSAL

The device and packaging must be disposed of in accordance with respective local or national environmental regulations.

LIABILITY

Össur prosthetic devices are designed and verified to be safe and compatible in combination with each other and custom-made prosthetic sockets with Össur adapters, and when used in accordance with their intended use.

Össur does not assume liability for the following:

- Device not maintained as instructed by the instructions for use.
- Device assembled with components from other manufacturers.
- Device used outside of recommended use condition, application, or environment.

Compliance

This device has been tested according to ISO 10328 standard to three million load cycles.

Depending on patient activity, this may correspond to 3-5 years of use.

ISO 10328 - P5 - 100 kg *)

*) Body mass limit not to be exceeded!

For specific conditions and limitations of use see manufacturer's written instructions on intended use!

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