



Instructions for Use

CHOPART

		3
EN	Instructions for Use	4
DE	Gebrauchsanweisung	8
FR	Notice d'utilisation	13
ES	Instrucciones para el uso	18
IT	Istruzioni per l'uso	23
NO	Bruksanvisning	28
DA	Brugsanvisning	32
SV	Bruksanvisning	36
EL	Οδηγίες Χρήσης	40
FI	Käyttöohjeet	45
NL	Gebruiksaanwijzing	49
PT	Instruções de Utilização	53
PL	Instrukcja użytkowania	57
TR	Kullanım Talimatları	62
RU	Инструкция по использованию	66
JA	取扱説明書	71
ZH	中文说明书	75
KO	사용 설명서	79

DESCRIPTION

The device consists of a single carbon blade for direct lamination (**Fig. 1**). This device must be used with a Foot Cover and a Spectra Sock.

INTENDED USE

The device is intended as a part of a prosthetic system that replaces the foot and ankle function of a missing lower limb.

Suitability of the device for the prosthesis and the patient must be evaluated by a healthcare professional.

The device must be fitted and adjusted by a healthcare professional.

Indications For Use and Target Patient Population

- Lower limb loss, amputation, or deficiency
- No known contraindications

The device is for low to high impact use, e.g., walking and occasional running.

The device is intended for transtibial and partial foot prosthesis.

GENERAL SAFETY INSTRUCTIONS

Warning: Use of a lower limb prosthetic device carries an inherent risk of falling which may lead to injury.

The healthcare professional should inform the patient about everything in this document that is required for safe use of this device.

Warning: If there is a change or loss in device functionality, or if the device shows signs of damage or wear hindering its normal functions, the patient should stop using the device and contact a healthcare professional.

The device is for single patient use.

DEVICE SELECTION

Verify that selected variant of the device is suitable for the impact level and weight limit according to the following table.

Warning: Do not exceed weight limit. Risk of device failure.

Incorrect category selection may also result in poor device function.

Weight kg	45-52	53-59	60-68	69-77	78-88	89-100	101-116	117-130	131-147
Weight lbs	99-115	116-130	131-150	151-170	171-194	195-220	221-256	257-287	288-324
Low Impact Level	1	1	2	3	4	5	6	7	8
Moderate Impact Level	1	2	3	4	5	6	7	8	N/A
High Impact Level	2	3	4	5	6	7	8	N/A	N/A

ASSEMBLY INSTRUCTIONS

Foot Cover and Spectra Sock

Caution: To avoid pinching fingers, always use a shoehorn.

1. Put the foot into the Spectra Sock.
2. Remove the foot cover attachment.

3. Use the applicable (straight) end of a shoehorn to put the foot with the Spectra Sock into the Foot Cover.
4. Move the shoehorn up to fully push the foot into the Foot Cover.
5. Put the foot cover attachment back onto the Foot Cover.
6. After alignment is complete, fix the Spectra Sock to the prosthesis to seal against dust and dirt.

Note: The Spectra Sock must be pulled up to prevent it from interfering with moving parts of the foot.

If required doff the Foot Cover as follows:

1. Remove the foot cover attachment.
2. Insert the applicable (edged) end of a shoehorn behind the foot.
3. Push the shoehorn down and pull the foot out of the Foot Cover.
4. Fully remove the Spectra Sock.

Prosthesis

Assemble prosthesis with applicable devices.

Warning: Risk of structural failure. Components from other manufacturers have not been tested and may cause excessive load on the device.

Warning: Ensure proper attachment by following the applicable device assembly instructions.

ALIGNMENT INSTRUCTIONS

Bench Alignment (Fig. 2)

- Adjust to appropriate heel height (using the shoe).
- Introduce appropriate socket angles flexion/extension and abduction/adduction.
- The load line taken from the center of the socket should fall through the midline of the foot module in a neutral M-L position. Depending on proximal or distal bearing socket design, use appropriate load line.
- Set appropriate toe rotation.
- Divide the foot into 3 equal portions.
- The load line should fall at the junction of the posterior and the middle third portions.

Initial Fitting Procedure

Preparation

- From Bench Alignment, Chopart foot plate should be bonded to the distal socket with potting epoxy or Siegelharz.
- Roughen carbon fiber surfaces (upper plate and socket) and clean with acetone to ensure best adhesion of epoxy or lamination.
- While in alignment jig, raise socket on the plate and lower unit back in place.
- Fill in around the outside edges with potting epoxy to create a smooth transition from the socket to the plate.
- Wrap with synthetic casting tape.
- Carry out dynamic alignment and make changes as necessary.
- During dynamic alignment you may cut off the tip of the heel and replace with soft foam material (**Fig. 3**).

Dynamic Alignment

- The roll over action of the foot can be influenced by the heel and the shape of the sole.
- Secure achieved dynamic alignment using a jig.
- Remove sole materials and synthetic casting tape and proceed to direct lamination.

Lamination Procedure

- The lay-up sequence involves posterior and M-L positioned carbon fiber cloth. Refer to table below for suggested number of layers.
- Apply on posterior layer, then follow with a M-L layer.
- Continue to alternate using the required number of layers. The long posterior layer is folded in half over the M-L layer and reflected backup to the socket to form a T layup, locking the layers together and producing a clean lamination edge (**Fig. 4**).
- Seal the not laminated part of the plates with tape.
- Pull a PVA bag over the prostheses and laminate with acrylic or epoxy resin.
- After lamination has set, sand off any protruding edges.

Note: Do not cut into the lamination.

Caution: Epoxy can cause allergic reactions when in direct contact with body.

Suggested Reinforcement				
		Min # of Layers*		
Layup Chart Sequence	Cloth Size	Cat 1-3	Cat 4-6	Cat 7-8
M-L Layers	76 mm x 355 mm	3	4	6
Posterior Layers	76 mm x 355 mm	2	2	3
Circumferencial Wrap	114 mm	2	2	3

*The appropriate number of layers depends on user activity and weight and must be evaluated by the CPO.

USAGE

Cleaning and care

Clean with a damp cloth and a mild soap. Dry with a cloth after cleaning.

Environmental Conditions

The device is Waterproof.

A Waterproof device can be used in a wet or humid environment and submerged in up to 3-meter-deep water for a maximum of 1 hour.

It can tolerate contact with: Salt water, chlorinated water, perspiration, urine, and mild soaps.

It can also tolerate occasional exposure to sand, dust, and dirt.

Continuous exposure is not allowed.

Clean with fresh water after exposure to other liquids, chemicals, sand, dust, or dirt and dry with a cloth.

Dry with a cloth after contact with fresh water or humidity.

MAINTENANCE

The device and the overall prosthesis should be examined by a healthcare professional. Interval should be determined based on patient activity.

Noise from Foot

Noise may occur if sand or debris is present in device. In that case, the healthcare professional should doff the foot, clean it with the help of compressed air and replace the Spectra sock if it is damaged.

REPORT OF SERIOUS INCIDENT

Any serious incident in relation to the device must be reported to the manufacturer and relevant authorities.

DISPOSAL

The device and packaging must be disposed of in accordance with respective local or national environmental regulations.

LIABILITY

Össur prosthetic devices are designed and verified to be safe and compatible in combination with each other and custom-made prosthetic sockets with Össur adapters, and when used in accordance with their intended use.

Össur does not assume liability for the following:

- Device not maintained as instructed by the instructions for use.
- Device assembled with components from other manufacturers.
- Device used outside of recommended use condition, application, or environment.

Compliance

This device has been tested according to ISO 10328 standard to two million load cycles.

Depending on patient activity this may correspond to 2-3 years of use.

ISO 10328 - Label		
Category	Weight (Kg) Moderate Impact	Label Text
1	52	ISO 10328 - P3*-52 kg
2	59	ISO 10328 - P3*-59 kg
3	68	ISO 10328 - P3-68 kg
4	77	ISO 10328 - P3-77 kg
5	88	ISO 10328 - P4-88 kg
6	100	ISO 10328 - P5-100 kg
7	116	ISO 10328 - P6-116 kg
8	130	ISO 10328 - P7-130 kg
*Actual test loads reflect maximum body mass		

ISO 10328 - "P" - "m"kg *)

*) Body mass limit not to be exceeded!

For specific conditions and limitations of use
see manufacturer's written instructions on
intended use!

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