

Instructions for Use

CONNECT[®] TF








		3
EN	Instructions for Use	8
DE	Gebrauchsanweisung	14
FR	Notice d'utilisation	20
ES	Instrucciones para el uso	26
IT	Istruzioni per l'uso	32
NO	Bruksanvisning	38
DA	Brugsanvisning	44
SV	Bruksanvisning	50
EL	Οδηγίες Χρήσης	56
FI	Käyttöohjeet	62
NL	Gebruiksaanwijzing	68
PT	Instruções de Utilização	74
PL	Instrukcja użytkowania	80
CZ	Návod k použití	86
TR	Kullanım Talimatları	92
RU	Инструкция по использованию	98
JA	取扱説明書	104
ZH	中文说明书	110
KO	사용 설명서	116







1







2

MD Medical Device

DESCRIPTION

The device is an above-knee socket made for lower active and impaired patients. Four socket shells can be adjusted in height, angle, and tension. This lets the device be adapted to the shape and volume of a patient residual limb. A handle opens and closes the device.

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INTENDED USE

The device is a pre-assembled, size adjustable socket intended as a part of a prosthetic system that replaces a missing lower limb.

Suitability of the device for the prosthesis and the patient must be evaluated by a healthcare professional.

The device must be fitted and adjusted by a healthcare professional.

Indications For Use and Target Patient Population

- Lower limb amputation and/or congenital deficiency
- No known contraindications

The device is for low impact use, e.g., gentle walking.

The weight limit for the device is 125 kg.

GENERAL SAFETY INSTRUCTIONS

The healthcare professional should inform the patient about everything in this document that is required for safe use of this device.

Warning: If there is a change or loss in device functionality, or if the device shows signs of damage or wear hindering its normal functions, the patient should stop using the device and contact a healthcare professional.

Warning: Risk of structural failure. Components from other manufacturers have not been tested and may cause excessive load on the device.

The device is for single patient use.

Warning: Do not use the device if the socket is too tight or too loose.

Warning: Any changes made by the patient can cause damage.

Warning: Do not use heat guns, ovens, and heaters on the device.

DEVICE SELECTION

Refer to Össur catalog for more information, e.g., sizing.

The device must be used together with the Iceross® Transfemoral Locking liner with a locking pin. A release lock secures the pin of the liner. The patient opens the device with the handle and then pulls the knob to unlock the liner to pull out the residual limb.

If the device is used with a locking knee, then it must be used together with the Locking Knee Lanyard Kit.

Size Selection

The patient should be in a seated position.

1. Measure the circumference of the residual limb at 4 cm from the distal end (**Fig. 2A**). It is important to take this measurement with the patient sitting.
2. Measure the length of the residual limb from the distal end to the perineum (**Fig. 2B**).
3. Select the correct device size from the chart (**Fig. 2C**). If the patient is between sizes, refer to the advanced selection chart (**Fig. 20**).

FITTING INSTRUCTIONS

Tools

An adjustable torque screwdriver with:

- Range from 0.6 Nm to 3 Nm
- 0.1 Nm increments
- a 4 mm hex bit

Donning Instructions

The patient should be in a seated position on a flat hard surface, such as a chair.

1. Fully invert the liner and roll it on the residual limb (**Fig. 3**).
Note: Make sure to align the pin correctly.
2. Put the device onto the liner (**Fig. 4**).
3. Push the pin on the liner into the hole. The Knob will turn until the pin is fully in the lock (**Fig. 5**).
4. If needed, turn the Knob to pull the pin fully into the lock.
5. Fully push the handle down. This will close the device to the residual limb (**Fig. 6**).

Device Adjustments

Adjusting the Height

The height of the shells can be adjusted to fit the length of the residual limb (**Fig. 7**).

Prior to changing the height ensure:

1. distal screws are loose
2. bobbins are not tensioned, and lace has some slack
3. handle being in the 'up' position
4. Loosen the Medial telescoping screw.
5. Adjust the height of the Medial proximal shell first (**Fig. 7A**).

Note:

SHORT version of the device uses the 0–40 mm marking on the medial strut. 0 mm corresponds to a residual limb length of 200 mm when the strut is vertical.

STANDARD version of the device uses the 40–80 mm marking on the medial strut. 40 mm corresponds to a residual limb length of 240 mm when the strut is vertical.

6. Loosen the Lateral telescoping screw.
7. Adjust the height of the Lateral proximal shell (**Fig. 7B**). The recommendation is that the Lateral proximal shell is a very small distance above the Medial proximal shell.
8. After dynamic alignment of the device, torque the telescoping screws to 5 Nm (**Fig. 7C**).
9. Note down the final heights of medial and lateral struts in your clinical notes.

Adjusting the Tension

The shells can be adjusted to fit the shape of the residual limb (**Fig. 8**).

Laces are used to tighten or loosen a shell. Four bobbins above the handle adjust the tension of the laces as follows:

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1. Ensure that all laces are very loose and that the shells are significantly larger in volume than the patient's residual limb.
2. Ensure that distal end screws are loose.
3. Tension should best be adjusted in a seated position.
4. Close the handle before tensioning (**Fig. 8A**).
5. Push and turn each bobbin to remove excess lace length first (remove slack) (**Fig. 8B**). When you push the bobbin, this releases it from its locking teeth, allowing it to rotate.
Note: When tensioning, it is important not to get lace tangled up inside of the bobbin. To prevent this, ensure that when turning a bobbin there is always some tension on the lace.
6. Use a torque screwdriver set to 1 Nm and tighten the bobbins in increments, doing multiple rounds across all 4 bobbins (do not tighten 1 bobbin straight to 1 Nm when the other 3 are still slack).
Note: When you tension the bobbins, the angle of the struts will rotate in and align to the shape of the residual limb.
7. When 1 Nm is achieved on all 4 bobbins, open and close the handle. This evenly distributes tension throughout the whole tensioning system. Reapply 1 Nm torque again.
8. 1 Nm is an initial guideline only. Clinical experience will determine final torque.
9. Note down the final torque value of all 4 bobbins in your clinical notes.

Adjusting the Distal End

The distal end can be adjusted to let the strut angles align with the residual limb.

1. Move the distal end to the medial or lateral side until it is perpendicular to the axis of the residual limb (**Fig. 9A**). This will make it easier for the pin to engage the lock.
2. Torque the two distal end screws to 20 Nm (**Fig. 9B**).
3. Note the final angle of medial and lateral struts in your clinical notes.

Dynamic Alignment

Attach the device to the rest of the prosthetic system.

You can adjust the device for better fit during dynamic alignment. For example:

1. Loosen the two telescoping screws to tune the height of the proximal shells.
2. Tune the tension of the laces. 1 Nm is an initial guideline for low active patients. Use your clinical judgment to decide if more or less tension is required.
3. Strut angle/distal end orientation can also be fine-tuned.

Secure fastening

When you get good dynamic alignment, do as follows:

1. Torque telescoping screw to 5 Nm.
2. Torque the distal end screws to 20 Nm.

No threadlocker is required.

Attaching Cosmetic Covers (Fig. 10)

1. Carefully tighten the lateral proximal and distal cover screws until the lateral cover is securely in place.
2. Carefully tighten the medial proximal and distal cover screws until the medial cover is securely in place.

Doffing Instructions

The patient should be in a seated position on a flat hard surface, such as a chair.

1. Fully pull the handle up. This will open the socket (Fig. 11).
2. Pull up and hold up the knob (Fig. 11). Pull the residual limb out of the socket.

USAGE

Cleaning and care

Clean with a damp cloth and a mild soap.

- Or use disinfectant wipes to clean the device.
- Air dry, do not expose to direct heat.

Daily cleaning of the residual limb is essential. Use of a mild, pH balance, 100% fragrance- and dye-free liquid soap is recommended. If the skin is dry, apply a pH balance, 100% fragrance- and dye-free lotion to nourish and soften the skin.

Caution: Skin health should be monitored daily. In case of abnormal symptoms, discontinue use of the device and contact the healthcare professional immediately.

Environmental Conditions

The device is Weatherproof.

A Weatherproof device can be used in a wet or humid environment and can tolerate being splashed by fresh water (e.g., rain), no submersion is allowed.

No contact with salt water or chlorinated water is allowed.

Clean with fresh water in case of accidental exposure to other liquids, chemicals, sand, dust, or dirt and dry with a cloth.

Dry with a cloth after contact with fresh water or humidity.

The device can be used in temperatures between 10°C to 50°C.

OPTIONAL FITTING

A left device can be made into a right device (and vice versa). To do this the handle needs to be changed and the distal end needs to be rotated 180°.

Changing the handle as follows (Fig. 13):

1. Ensure the laces are loose with lots of slack.
2. Remove the handle cover from the handle.
3. Remove the handle screw (Fig. 13A).
4. Replace with the handle for left- or right-side use (Fig. 13B).
5. Replace the handle screw.
6. Torque the handle screw to 3 Nm (Fig. 13C).
7. Replace the handle cover.

Rotate the distal end as follows (Fig. 14):

8. Remove the distal end screws (Fig. 14A).
9. Rotate the distal end. The knob must face towards the anterior (Fig. 14B).
10. Turn the distal end screws back into the distal end (Fig. 14C).

Attaching the Locking Knee Lanyard Kit

Attaching the Lever and Guides:

1. Remove the lateral cover screws (leave the lateral cover in place).
2. Put the locking knee lanyard washer over the proximal hole on the lateral cover.
3. Attach the lanyard lever, lanyard washer, and lanyard screw to the strut (**Fig. 15A**).
4. Attach the proximal guide to the strut, below the lanyard lever (**Fig. 15B**).
5. Attach the distal guide to the posterior side of the strut (**Fig. 15C**).

Attaching the Lanyard Cord to the Proximal End:

6. Attach the spring, nylon lanyard, lanyard collar, and heat-shrink tubing to the lever.
7. Attach the assembled end of the lanyard cord to the spring (**Fig. 16**).
8. Heat the heat-shrink tube. **Note:** When heating the heat-shrink tube be careful to heat it only for a short duration away from the device.
9. Put the tube with the nylon lanyard through the two guides (**Fig. 17**).

Attaching the Lanyard Cord at the Distal End:

10. At the distal end, loop the end of the nylon lanyard through the lanyard collar (**Fig. 18**).
11. Examine the function of the locking function to find the correct length of the nylon lanyard at the distal end.
12. Fasten the lanyard collar around the nylon lanyard (**Fig. 18**).
13. Attach the loop to the knee lock (**Fig. 19**).
14. Cut off excess nylon lanyard.

Distal Elevation

To allow for fitting of a shorter limb, change the funnel from 0cm to 2cm

1. Unscrew the 2 funnel screws.
2. Remove funnel 0cm.
3. Put in place funnel 2 cm, the guide pin will help locate the correct position.
4. Screw in the 2 funnel screws to 5 Nm.

Note: A longer locking pin might be required.

Best Practice

1. The key to a good device fitting is to spend time on iterative adjustments to fine-tune the socket fit after the initial fit. As device adjustments are fast and easy, it is possible to make real-time rapid changes and observe and discuss with the patient their effectiveness. In this way improvements are '*dialed in*'.
2. Do not be afraid to adjust, it is easy to '*undo*' the socket and start over.
3. When fine-tuning the tension, work in 0.2 Nm increments.
4. To ease donning:
 - The patient should be in a seated position on a flat hard surface, such as a chair.
 - Lift up the residual limb or move to the edge of the seat and lean backward.

MAINTENANCE

The device and the overall prosthesis should be examined by a healthcare professional. Interval should be determined based on patient activity.

REPORT OF SERIOUS INCIDENT

Any serious incident in relation to the device must be reported to the manufacturer and relevant authorities.

DISPOSAL

The device and packaging must be disposed of in accordance with respective local or national environmental regulations.

LIABILITY

Össur does not assume liability for the following:

- Device not maintained as instructed by the instructions for use.
- Device assembled with components from other manufacturers.
- Device used outside of recommended use condition, application, or environment.

Compliance

This device has been tested according to ISO 10328 standard to three million load cycles. Depending on patient activity, this may correspond to up to 2 years of use.

ISO 10328 - P7 - 125 kg *)

*) Body mass limit not to be exceeded!

For specific conditions and limitations of use
see manufacturer's written instructions on
intended use!

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