








Instructions for Use

---

## CUSHION & LOCKING

		3
EN	Instructions for Use	4
CS	Návod k použití	11
SK	Indikácie použitia	18
RO	Indicații de utilizare	25
HR	Upute za uporabu	32
HU	Használati javallat	39
BG	Указания за употреба	46
SL	Indikácie pre použitie	53

## DESCRIPTION

The device is a silicone liner that covers the residual limb and provides an interface between the residual limb and the socket of the prosthesis. The device provides a secure suspension to the prosthesis, protection of the residual limb, and comfort.

A locking liner either has a fixed or an interchangeable distal attachment. A cushion liner has no distal attachment and needs to be used together with a sleeve that holds vacuum.

## INTENDED USE

The device is a prosthetic interface with suspension properties intended to be used as part of a system that replaces a missing lower limb.

Suitability of the device for the prosthesis and the patient must be evaluated by a healthcare professional.

The device must be fitted and adjusted by a healthcare professional.

### **Indications For Use and Target Patient Population**

- Lower limb loss, amputation, or deficiency
- No known contraindications

The impact level for the device is listed in the Table below.

Transtibial Liners	Impact Level
Iceross Dermo® Locking	The device is for low to moderate impact use, e.g., walking.
Iceross Dermo® Uniform Locking	
Iceross Dermo® Cushion	
Iceross Dermo® Uniform Cushion	
Iceross Comfort® Locking	
Iceross Comfort® Cushion	
Protect Locking	
Protect Cushion	
First	
Sensitive Locking	
Iceross® Original Locking	The device is for low to high impact use, e.g., walking and occasional running.
Iceross Stabilo® Junior Locking	
Iceross Dermo® Junior Locking	
Iceross® Synergy Locking	The device is for moderate to extreme impact use, e.g., walking, running, and sports.
Iceross® Synergy Cushion	
Iceross® Activa Cushion	
Iceross® Sport Locking	The device is for high to extreme impact use, e.g., fast walking, running, and sports.

Transtibial Liners	Impact Level
Soft C Locking	The device is for low to extreme impact use, e.g., walking, running, and sports.
Relax Locking	
Relax Cushion	
Transfemoral Liners	Impact Level
Iceross® Transfemoral	The device is for low to extreme impact use, e.g., walking, running, and sports.
AKOS	
Relax Transfemoral	

## GENERAL SAFETY INSTRUCTIONS

**Warning:** Use of a lower limb prosthetic device carries an inherent risk of falling which may lead to injury.

The healthcare professional should inform the patient about everything in this document that is required for safe use of this device.

**Warning:** If there is a change or loss in device functionality, or if the device shows signs of damage or wear hindering its normal functions, the patient should stop using the device and contact a healthcare professional.

**Warning:** Misuse of the device may result in loss of suspension.

**Caution:** Potential side-effects are skin irritation/reaction problems.

**Caution:** The device must not come in contact with glass, carbon fibers, or other foreign particles. Washing alone may not be sufficient to eliminate the problem. If the device is exposed to foreign substances or chemicals, it should be returned to the healthcare professional.

The device is for single patient use.

## DEVICE SELECTION

For a transtibial liner, the correct size is determined by bending the knee slightly, allowing the soft tissue to hang down.

To determine the correct size:

- Measure the circumference of the residual limb 4 cm from the distal end (**Fig. 1**)
- Choose the corresponding device size. If the measurement taken does not have a corresponding device size, round the measurement down.

**Note:** Do not choose a device size which exceeds the measurement, i.e., do not round up.

**Caution:** A too small/tight device may result in pistoning, numbness, or blistering.

**Caution:** A too large/loose device may result in increased perspiration and movement of the limb inside the device, which may lead to blisters and rashes.

If any of the above symptoms, or any other indications of improper fit are experienced, the patient should contact the healthcare professional immediately.

**Note:** For a transtibial locking liner, make sure that the liner matrix covers at least 4 cm of bony structure. Do not allow the matrix to extend higher than the fibula head.

## FITTING INSTRUCTIONS

**Caution:** The device must be thoroughly inspected for defects before fitting.

### **Trimming Instructions**

The device may be trimmed to increase range of motion, or for patient preference (**Fig. 2**). The Iceross Curvemaster is a recommended tool, specifically designed to produce a smooth, rounded edge when trimming. **Caution:** Do not cut the device below the level of the matrix or below the proximal trim-lines of the socket. Excessive trimming may reduce the suction between the device and limb, and compromise suspension.

### **Socket Fitting**

Achieving the intended function of the device is dependent on the fit of the socket. Appropriate socket design is essential to achieve suspension, control, and comfort.

The following checklist should be used when checking the socket fit:

- Patient can don the socket easily without the use of lubricant aids.
- Full contact between the liner and the socket inner wall.
- No air pockets develop when the patient loads the prosthesis.
- Minimized pistoning in the socket thus, no visible movement when pulling on the prosthesis.
- Full movement possible without compromising the suspension.

**Caution:** If feedback indicates any of the above conditions are not met, the socket will need to be adjusted or remade accordingly.

**Caution:** For a cushion liner, the socket must be airtight in order to maintain vacuum suspension when a sleeve is applied.

### **Locking Liner**

**Warning:** Failure to adhere to the following instructions may cause the pin to detach from the liner resulting in loss of suspension, which may lead to injury.

It is recommended that the device is used in conjunction with a compatible Össur lock.

Iceross Dermo Junior and Iceross Stabilo Junior are only compatible with the Icelock 700 series lock.

**Caution:** Do not use a flangeless pin with the device (**Fig. 3**). A flangeless pin can push through the distal attachment.

### **Fixed distal attachment**

If a lanyard suspension system is used: Turn the connecting screw into the distal attachment of the liner. Refer to the instructions for use for the lanyard suspension system for correct instructions.

If a locking pin is used: Apply Loctite 410/411 or an equivalent compound to the threads of the pin and tighten the attachment pin to the fixed distal attachment, torque to 4 Nm. Refer to the instructions for use for the lock used for correct instructions.

### **Interchangeable distal attachment (only applicable for Iceross Sport)**

Tighten the distal attachment to the liner to 4 Nm. To fasten the attachment pin to the distal attachment, apply Loctite 410/411 or an equivalent compound to the threads of the pin and torque to 4 Nm.

## USAGE

### **Donning**

**Caution:** Open wounds or damaged skin should be covered with a bandage or other appropriate covering to prevent direct contact between the wound and device.

**Caution:** Ensure the inside of the device is clean, dry, and free from any foreign objects that may cause skin irritation.

**Caution:** Do not apply lotion to the residual limb immediately before donning the device. The residual limb should be clean and dry to prevent device damage.

1. Grip the top of the device from the inside and slide it over the hand until the inner surface is fully exposed (**Fig. 4**). Take care not to damage the device with fingernails.
2. After exposing as much of the distal end of the device as possible, position it against the distal end of the residual limb (**Fig. 5**).  
**Note:** Be careful that the pin on the locking liner is positioned centrally on the distal end of the residual limb (**Fig. 7**).
3. For a transtibial liner: Bend the knee slightly.
4. With light compression, roll the device all the way upwards onto the limb (**Fig. 6**).  
**Caution:** Do not tug or pull the upper end of the liner when rolling the device upwards onto the limb. This can result in tension on the skin, causing blisters or rashes.
5. Check that no air pockets are present. If air pockets are found, reapply the device.

The Iceross Dermo Uniform Locking and Cushion Liners must be orientated with the external labeling at the front of the residual limb. For Iceross Activa the TibiaGuard with the wave structure must be aligned along the front side of the tibia.

### Socket Donning

**Note:** Do not use alcohol spray as a donning aid for a device without a textile cover (with Easy Glide) as this can destroy the low-friction coating and reduce device functionality.

**Cushion Liner:** Use together with a sleeve. Refer to the instructions for use for the sleeve for correct use.

Step into the socket using body weight to achieve full engagement before using the sleeve to maintain vacuum suspension.

**Locking Liner:** Step into the socket using bodyweight to fully engage the distal attachment into the lock until a secure connection is achieved.

**Warning:** Ensure that the connection between the distal attachment and lock is secure and without play.

### **Cleaning and care**

Daily cleaning of the residual limb is essential. Use of a mild, pH balance, 100% fragrance- and dye-free liquid soap is recommended. If the skin is dry, apply a pH balance, 100% fragrance- and dye-free lotion to nourish and soften the skin.

Clean the device daily after use and before first use.

- Turn the device inside out (**Fig. 4**).

- Wash with pH balanced, 100% fragrance- and dye-free liquid soap.
- After washing, rinse the device thoroughly with warm water (**Fig. 8**), and pat dry on both sides with a lint-free cloth.
- The device can also be machine washed (40°C) with a mild detergent and at low spin speed. Fabric softeners, bleaches, and other products/cleaning solutions may result in device damage and should not be used.
- For a device with a textile cover: The device shall be dried overnight before reuse.  
For a device without a textile cover (with Easy Glide): The device can be used directly after washing and being pat dry on both sides with a lint-free cloth.
- Always return the device to its neutral state immediately after cleaning. Do not store inside-out as device damage may result.

**Warning:** Excessive perspiration may compromise stability and/or suspension. Remove the device and dry off both the device and the residual limb when appropriate.

**Caution:** Skin health should be monitored daily. In case of abnormal symptoms, discontinue use of the device and contact the healthcare professional immediately.

**Caution:** Be aware that common household or bath products may cause or contribute to skin irritation, e.g., soaps, deodorants, perfumes, abrasive cleaners, aerosol- or alcohol sprays.

**Caution:** After contact with salt and/or chlorinated water, wash the device with tap water.

**Note:** Initially the inner surface of the Iceross Dermo Liners and Protect Liners should be moist to the touch. Therefore, the device should not be left inside out for an extended period or exposed to excessive heat or sunlight.

### ***Environmental Conditions***

The device is safe to use with caution with waterproof components.

**Caution:** Frequent exposure to salt and/or chlorinated water might affect the durability of the device.

Only applicable for locking liner:

**Warning:** After submerging the device in water, it needs to be removed, and the residual limb and the skin contact layer dried before ambulation.

**Warning:** The suspension is compromised if water gets between the skin and the device.

Only applicable for Relax Transfemoral, Relax Locking, and Relax 3C/6C Cushion.

**Caution:** Prolonged exposure to salt and/or chlorinated water compromises the shielding of the textile cover against electromagnetic influences.

## **SUPPLEMENTARY DEVICES**

The following accessories can be used in conjunction with Össur cushion and locking liners.

Refer to Össur catalog for more information, e.g., sizing.

### **Iceross® Distal Cup**

The device is a prosthetic interface intended to be used as part of a system that replaces a missing lower limb.

The device provides further cushioning for distally conical limbs with poor soft tissue coverage. It conforms to the residual limb and can therefore be used as a filling element to ensure full contact between the silicone liner and the residual limb.

To correctly fit the device:

1. Measure the circumference of the residual limb 4 cm above the distal end (**Fig. 1**).
2. Choose a device one or two sizes smaller than the measurement.
3. Roll the device directly over the limb.
4. When using a liner in combination with a Distal Cup, the liner measurement should be made with the Distal Cup in place (i. e., over the Distal Cup).
5. Choose the correct liner size based on the measurement obtained with the Distal Cup in place.

**Note:** For daily care, refer to Cleaning and Care instructions.

### **Iceross® Pads**

The device is a prosthetic interface intended to be used as part of a system that replaces a missing lower limb.

The device can be used for prosthetic volume management to temporarily accommodate lost distal end volume as well as to assess distal contact while fitting a check socket.

Place the pad either in the distal end of the liner or in the distal end of the socket.

**Note:** For daily care, refer to Cleaning and Care instructions.

### **Iceross® Interface Sheath**

The device is intended to be used for prosthetic socket donning.

The device can be worn over the liner to allow easier donning/doffing of the socket by reducing the friction between the socket and liner.

Pull the device over the liner. When used in conjunction with a locking liner ensure the metal pin and distal connector fully penetrates the knitted hole to ensure the proper function of the lock.

Machine wash at max 40°C with light-colored fabric. Do not tumble dry.

### **Iceross® Socks**

The device is intended to be used for accommodating residual limb volume reduction.

The device can be worn over the liner for volume management.

If the patient experiences compromised suspension caused by residual limb volume fluctuations, one or more sock can be pulled over the liner to restore socket fit.

Pull the device over the liner. When used in conjunction with a locking liner ensure the metal pin and distal connector fully penetrates the knitted hole to ensure the proper function of the lock.

Machine wash at max 40°C with light-colored fabric. Do not tumble dry.

## **REPORT OF SERIOUS INCIDENT**

Any serious incident in relation to the device must be reported to the manufacturer and relevant authorities.

## **DISPOSAL**

The device and packaging must be disposed of in accordance with respective local or national environmental regulations.

## **LIABILITY**

Össur prosthetic devices are designed and verified to be safe and compatible in combination with each other and custom-made prosthetic sockets with Össur adapters, and when used in accordance with their intended use.

Össur does not assume liability for the following:

- Device not maintained as instructed by the instructions for use.
- Device assembled with components from other manufacturers.
- Device used outside of recommended use condition, application, or environment.

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