ENGLISH

INTENDED USE

The device is intended for plantar surface offloading and immobilization. The device must be fitted and adjusted by a healthcare professional.

Indications for use

- Plantar foot ulcers and wound prevention.

Contraindications

This device is not intended for use on patients with chronic Charcot feet with a rocker bottom deformity.

Warnings and Cautions:

- The device must be prescribed by a healthcare professional.
- Weight limit: The patient's maximum weight is 136 kg (300 lbs).
 - Hexplugs should only be removed or altered by a healthcare professional.
 - Wounds under treatment should be monitored closely for maintenance of pressure relief and healing progress.
 - The insole should not be used in conjunction with any other footwear device.
 - In case of emergency, the compliance strap can be cut with scissors for the removal of the device.
 - Before reapplying the device inspect:
 - The foot for any rash, sores, new ulcers, or other indications of skin irritation.
 - The boot, bladders, and liners for any damage.
 - If the device has become loose or uncomfortable to wear, readjust the straps as necessary and contact your healthcare professional if discomfort is not relieved.
 - Patient should use caution when walking on slippery or wet surfaces to avoid injury.
 - Do not remove the device unless instructed by a healthcare professional.
 - This device is only intended for treating or preventing plantar ulcers of the foot.

GENERAL SAFETY INSTRUCTIONS

The healthcare professional should inform the patient about everything in this document that is required for safe use of this device.

Any serious incident in relation to the device must be reported to the manufacturer and relevant authorities.

The patient should immediately contact a healthcare professional:

- If there is a change or loss in device functionality, or if the device shows signs of damage or wear hindering its normal functions.
 - If any pain, skin irritation, or unusual reaction occurs with the use of the device.
- The device is for single patient – multiple use.

FITTING INSTRUCTIONS

While carrying out the following instructions, please refer to the overview figure for locating components mentioned in the text (Fig. 1).

Device Application

Note: Do not remove Pull Tabs (A) until further instructed.

- Remove Liner (B) by opening Contact Closure Straps (C) and lifting it out of the device (Fig. 2).
- Place lower leg and foot into the Liner and inner Bootie (D). Close the inner Bootie and then Liner snugly over the foot and leg (Fig. 3, 4).
- Insert foot into the device by gently spreading the Struts (E) outwards (Fig. 5). Ensure Struts are centered over the ankle bone (malleoli). The lower leg should be kept at a 90° angle with the foot and ensure that no gap is in between the Liner heel and the Strut (Fig. 6).
- Remove Pull Tabs located in between the Struts and the Liner. This allows the Liner to be secured into position. Pulling the tabs in a forward/downward direction makes them easier to remove (Fig. 7).

- Remove foot from the device. Cover wound with material such as a transparent film dressing, tape, gauze, etc. Mark material with a transferable mark using iodine, lipstick, indelible ink pencil, etc. (Fig. 8).
- Place foot inside the device. Ensure a proper location of the patient's foot before transferring the mark to the insole's surface (Fig. 9). Close the inner Bootie and then Liner snugly over the foot and leg. Bear weight to ensure transfer of relief zone mark. After the patient returns to a seated position, remove the foot from the device.
- Remove marked Hexplugs (F) from the insole (Fig. 10).
Caution: This device should not be used if the relieved area is too large and inhibits pressure relief by allowing the ulcer or wound to bottom out.
- Scythe the edges of Hexplugs surrounding the relief zone with a cutting tool such as a scalpel or scissors (Fig. 11).
Note: To provide additional support, a transparent film dressing can be placed over the relief zone, or a low-density foam can be used as filler inside the relief zone, no more than 0.5" (1 cm) thick.
- Fasten the Contact Closure Straps starting at the ankle, foot and then working upwards (Fig. 12). Strap adjustments may be necessary with changes in lower limb volume.

Device Removal

Undo the straps in any order and attach the straps back to themselves. This will make it easier when fitting the device next time and will prolong the life of the straps. Open the Liner to lift the foot and lower leg out of the device.

Accessories and Replacement Parts

Please refer to the Össur catalog for a list of available replacement parts or accessories.

USAGE

Cleaning and care

Washing Instructions

Device plastic can be washed with a damp cloth and mild soap.

Remove liner from the device before washing.

- Hand-wash using mild detergent and rinse thoroughly.
- Air dry.

Note: Do not machine-wash, tumble dry, iron, bleach, or wash with fabric softener.

Note: Avoid contact with salt water or chlorinated water. In case of contact, rinse with fresh water and air dry.

DISPOSAL


The device and packaging must be disposed of in accordance with respective local or national environmental regulations.

LIABILITY

Össur does not assume liability for the following:

- Device not maintained as instructed by the instructions for use.
- Device assembled with components from other manufacturers.
- Device used outside of recommended use condition, application, or environment.

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