















Instructions for Use
















---

DIRECT SOCKET TF







		3
EN	Instructions for Use	6
DE	Gebrauchsanweisung	14
FR	Notice d'utilisation	22
ES	Instrucciones para el uso	30
IT	Istruzioni per l'uso	38
NO	Bruksanvisning	46
DA	Brugsanvisning	54
SV	Bruksanvisning	62
EL	Οδηγίες Χρήσης	70
FI	Käyttöohjeet	78
NL	Gebruiksaanwijzing	86
PT	Instruções de Utilização	94
PL	Instrukcja użytkowania	102
CS	Návod k použití	110
TR	Kullanım Talimatları	118
RU	Инструкция по использованию	126
JA	取扱説明書	134
ZH	中文说明书	142
KO	사용 설명서	149

Medical Device

## PRODUCT DESCRIPTION

The Direct Socket contains all necessary materials to fabricate a single Direct Socket for prosthetic use.

The Direct Socket TF consists of a Brim and a Material Kit with the following parts:

- Fiber braids with a distal attachment
- Silicone insulation sheets
- Silicone insulation caps
- Direct Socket resin with static mixer
- Plastic protection sheets
- Trimline strip
- Spacer for Icelock 600\* (clear)
- Spacer for Icelock 562 Hybrid\* (blue)

\*Only supplied with 4" fiber braid version

## INDICATIONS FOR USE

- Lower limb amputation.
- Lower limb deficiency.

## CONTRAINDICATIONS

None known.

## INTENDED USE

The device is a modular component for a lower limb prosthesis that replaces the function and aesthetics of a missing limb. Suitability of this component for the prosthesis must be evaluated by a healthcare professional.

The Direct Socket TF is intended for creating a transfemoral socket (a device) using one of these suspension methods:

- Seal-In®
- Pin Locking

The device is for single patient use.

The device is designed for everyday use of low to high active users.

The device is not to be used by patients with the following conditions:

- Residual limb with extreme volume fluctuations when Seal-In liners are used.
- Össur does not recommend prescribing the device for patients that do not have sufficient hand strength to don the device properly, and/ or do not have a caregiver to assist with donning.

## GENERAL SAFETY INSTRUCTIONS

Inform the patient about everything in this document that is required for safe use of this device.

**WARNING:** If the device shows signs of damage or wear hindering its normal functions, the patient should stop using the device and contact a healthcare professional.

**WARNING:** In cases of functional change or functional loss, the patient should stop using the device and contact a healthcare professional.

## REQUIRED COMPONENTS

For creating a socket the following is required:

Direct Socket Tool Kit

Casting Liner

Suspension components based on suspension type selected.

Bottom Seal (vacuum suspension only).

If using a 4" fiber braid and

Icelock 562: Put the blue spacer between the lock and distal attachment.

Icelock 600: Put the clear spacer between the lock and the modular spacer.

## DEVICE SELECTION

### Material Kit

Please refer to the table below to verify device selection.

**WARNING:** Do not exceed weight limit. Risk of device failure.

The load level specification is shown in the tables below.

Braid size	Weight Limit	Label Text
4"	100 kg / 220 lbs	ISO 10328-P5-100kg
5"	166 kg / 366 lbs	ISO 10328-P7-166kg
7"	166 kg / 366 lbs	ISO 10328-P7-166kg
9"	166 kg / 366 lbs	ISO 10328-P7-166kg

### Brim

1. Measure the circumference 1–2 cm (0.4" to 0.8") below the level of the perineum (**Figure 1**).

2. Select the correct brim size from the table below.

**NOTE:** If the circumference is between two sizes, select the larger size. If you will put the brim high on the lateral side, think of selecting one size larger.

Proximal measurement	Recommended Brim size:
39 – 42 cm	40
42 – 44 cm	42,5
44 – 47 cm	45
47 – 49 cm	47,5
49 – 52 cm	50
52 – 55 cm	52,5
55 – 57 cm	55
57 – 59 cm	57,5
59 – 62 cm	60
62 – 64 cm	62,5
64 – 67 cm	65
67 – 69 cm	67,5
69 – 72 cm	70

## ASSEMBLY INSTRUCTIONS

Direct Socket production is divided into three phases: preparation, casting, and finishing.

For some steps of the fabrication, the assistance of a technician is recommended.

**CAUTION:** Use protective clothing, gloves, mask and glasses, during fiber braid handling and socket grinding process.

**CAUTION:** Read the Material Safety Data Sheet (MSDS) before using resin. Refer to the MSDS in case of emergency.

**WARNING:** Avoid direct contact with resin. Make sure everyone involved wears safety glasses and gloves. Direct contact with the resin can cause serious eye irritation and skin irritation.

**WARNING:** Do not use a different type of resin. Socket may break.

**WARNING:** Do not cut fiber braids without air extraction. Fiber particles can cause skin, eye, or respiratory irritation.

## **Preparation**

### **Relief Pads**

Localize bony prominences that require pressure relief. If required, select the appropriate relief pad(s) and place on the residual limb.

**NOTE:** Too many relief pads or relief pads put on top of each other should be avoided, to prevent cracks in the silicone liner and a bad socket fit. The relief pads are for single use.

### **Patient Position During Socket Casting**

- Patients with bad muscle response must lie on their back during socket casting (**Figure 2**).
- Patients with good muscle response must lie on their side during socket casting (**Figure 2**).

### **Attaching the Injection Valve**

Turn the injection valve into the injection channel. An arrow on the distal connector shows you the injection channel (**Figure 3**).

**WARNING:** Make sure you attach the injection valve correctly. If the injection valve is not in the injection channel, the resin will erupt when you inject it. Resin eruption can cause safety risks.

### **Taping Ring**

Use the screws to attach the taping ring to the distal connector (**Figure 4**).

### **Liner for Casting**

1. Determine the correct size. Measure the circumference of the residual limb at 4 cm from the distal end.
2. Choose the device size which corresponds to the measurement, or the one which is immediately below. Example: If the residual limb measures 24,5 cm at the indicated location, choose liner size 23,5 instead of size 25.
3. Roll on the liner used for casting (**Figure 5**).
4. Use the lanyard cord to examine soft tissue stabilization (**Figure 6**). If you are not sure, try a casting liner one size above or below.

**WARNING:** The casting liner must be cleaned and disinfected before each use. See Casting Liner section below.

**NOTE:** A wrong liner size can lead to an uncomfortable socket. If you make an uncomfortable socket, you will have to make a new socket.

### **Distal Attachment Pin**

1. Turn the distal attachment pin into the casting liner umbrella.
2. Make sure the pin and umbrella are correctly aligned with the patient's limb.

### **Using the Plastic Protection Sheet**

You need to put some of the plastic protection sheet under the liner to protect the patient. Do as follows:

1. Make a hole in the middle of the plastic protection sheet.
2. Roll 1/3 of the liner down.
3. Pull the plastic protection sheet up over the residual limb (**Figure 7**).
4. Roll the liner up over the plastic protection sheet.

### **Silicone Insulation Sheet, First Layer**

1. Pull the first length of silicone insulation sheet over the casting liner (**Figure 8**).

**NOTE:** The textile layer must be on the inner side. If it is not on the inner side, you will laminate the socket to the silicone insulation sheet.

2. Make sure to pull the silicone insulation sheet over the plastic protection sheet (**Figure 9**).
3. Use the sheet straps to keep the silicone insulation sheet in position during socket production (**Figure 10**).

### Attaching the brim

1. Put the brim on the residual limb. The thick part of the brim must be on the medial side (**Figure 11**).
2. Position the brim 1 to 2cm (0.4" to 0.8") below the level of the perineum.
3. Use the brim tape to attach the proximal end of the brim to the silicone insulation sheet. This will prevent resin flow under the brim (**Figure 12**).

### If using Icelock Expulsion Valve 552, do as follows:

1. Put a piece of double-sided adhesive tape onto the back of the dummy.
2. Stick the dummy onto the silicone insulation sheet on the residual limb (**Figure 13**).

**NOTE:** If you want to use the Icelock 544 Plate Unity as an alternative, you must use the Icelock 544 Direct Socket Adapter Plate.

### Putting on the Silicone Insulation Cap

1. Use scissors to cut away silicone insulation sheet that is not necessary (**Figure 14**). Make sure the insulation sheet ends a small distance above the distal part of the casting liner.
2. Put the silicone insulation cap on the distal part of the casting liner (**Figure 15**).

### Attach the Fiber Braids with Distal Attachment

1. Fiber braids with small or medium distal attachment (**Figure 16**):
  - Position the injection channel anteriorly (1 o'clock). This will improve resin flow.Fiber braids with large distal attachment:
  - The orientation of the attachment offset will determine injection channel position. If the injection channel not positioned anteriorly, pay more attention to distributing the resin evenly.
2. Turn the distal attachment nut onto the distal attachment pin and tighten nut firmly by hand.  
**NOTE:** Make sure no fiber braids are caught between the residual limb and the distal connector.
3. Pull the first layer of fiber braids up over the silicone insulation sheet. Pull the second layer of fiber braids up over the first layer. Repeat for all layers.
4. Check the fiber braids for folds and creases. Make the fiber braids smooth where necessary.
5. Mark a cutting line 25mm above the expected trimline (**Figure 17**).
6. Remove the fiber braids with distal attachment.
7. Cut the fiber braids along the cutting line with scissors. Do not cut along expected trimline!
8. Assemble again. Repeat steps 1 to 4.

### Silicone Insulation Sheet, Second Layer

Use the sheet applicator to roll the second length of silicone insulation sheet over the fiber braids (**Figure 18**).

**NOTE:** The textile layer must be on the outer side. If it is not on the outer side, you will laminate the direct socket to the silicone insulation sheet.

### O-rings

1. Take the large O-ring and put it on the top edge of the distal connector (**Figure 19**).  
**NOTE:** When using small distal attachment, do not apply the large O-ring.
2. Take the small O-ring and put it into the groove on the taping ring (**Figure 20**).  
**CAUTION:** Make sure there are no wrinkles and air between the O-rings. Wrinkles and air can cause the resin to leak.
3. Tape over the O-rings to fix the second layer of silicone insulation sheet to the distal connector and the taping ring (**Figure 21**).  
**NOTE:** When using small distal attachment, do not tape above the distal connector.  
**NOTE:** Do not tape above the proximal O-ring. Tape above the proximal O-ring can stop the resin flow. You will not be able to cast a usable socket if the resin does not flow and saturate the fiber braids.
4. Use scissors to cut away silicone insulation sheet that is not necessary below the taping ring (**Figure 22**).
5. Place the Air Escape Tube at the proximal end laterally between the first and second layer of silicone insulation sheet. (**Figure 23**). This will allow air to escape during resin injection. Both layers of silicone insulation sheet must fully cover the residuum and fiber braids.

## Casting

### Preparing the Resin Injection Tool

1. Cover the working area and the floor with protective sheeting.
2. Remove the parts from the packaging.
3. Shake the resin cartridge before use.
4. Lightly turn off the cap from the resin cartridge. The cap will lift off the green safety pin.
5. Remove the metal safety ring from the green safety pin.
6. Remove the green safety pin and discard it.
7. Put the cap onto the static mixer (**Figure 24**).
8. Push the static mixer 1 cm / 0.4 inches deep into the injection tube (**Figure 25**).
9. Turn the cap tight to safely attach the static mixer to the resin cartridge (**Figure 25**).
10. Keep the resin cartridge pointed up and insert it into the resin injection tool.
11. Lightly press the handle to secure the resin cartridge until it softly connects to bottom of the cartridge.
12. Keep the resin injection tool with the cartridge pointed up and sit in front of the patient.
13. The technician fully attaches the injection tube to the injection valve (**Figure 26**).

### Resin Injection

**WARNING:** Make sure the injection tube is fully secured to the static mixer and injection valve. If the tube is not fully secured, the resin can erupt when it is injected.

**NOTE:** Resin shall be injected gradually and spread evenly towards the proximal end. Excessive resin injected at the distal end can cause the resin to get hot and cure faster. This will make it harder to inject more resin.

1. Keep the resin injection tool pointed up while you start to inject the resin slowly (**Figure 27**).
2. When resin has saturated the distal area, injection speed can be increased.
3. The technician slowly pulls the resin into the interior and upwards, pushing out the air during resin injection.
4. Control the resin flow to the rear side of the residual limb while you maintain a slow build-up.

**NOTE:** Make sure to spread and saturate the braids before the resin starts to cure.

5. Use the lanyard cord to push the resin around the socket. Make sure to saturate the braids with resin up to and above the expected trimline (**Figure 28**).
6. Stop the resin flow when the fiber braids are saturated. To stop, push the pressure release lever on the grip of the resin injection tool (**Figure 29**).
7. Use pliers on crimp rings to shut off resin flow in the injection tube (**Figure 30**).
8. Place side cutter between the two crimp rings and cut through the injection tube.

### Resin Curing

1. Tell the patient to slowly extend the residual limb while the resin is still warm.
2. Scars and muscle attachments can change the shape of the residual limb during hip extension. Stretch the outer layer of silicone insulation sheet up when you mold the distal end to hold the fiber braids in place.
3. Do not allow the resin to build up at the distal end. If necessary, massage the resin of the socket when it cures to correctly shape.

### Casting Disassembly

After 10 minutes do as follows:

1. Remove the tape from the taping ring.
2. Remove the two O-rings.
3. Remove the air escape tube
4. Remove the outer layer of silicone insulation sheet.
5. Remove the distal attachment nut and the tape ring.
6. Remove the socket from the residual limb.
7. Pull the inner layer of silicone insulation sheet and the silicone insulation cap out of the socket.
8. For vacuum suspension only: break off the injection valve to seal the injection channel. Screwing out the valve can cause air leak and compromise the vacuum suspension.
9. For other suspension methods: Unscrew the injection valve. Do not break off the valve. The Icelock 214 can not be assembled if valve is not removed.

## **Finishing**

1. Use scissors to cut away fiber braids that are not saturated (**Figure 31**).
  2. Machine grind and sand the edge (**Figure 32**).  
**NOTE:** Do not grind the textile layer of the brim.
  3. Water sand the edge (**Figure 33**). Make sure the edge is rounded and not sharp.
  4. Remove unwanted resin from the inside of the socket.
  5. If used, grind the Icelock Expulsion Valve 552 dummy.
- Refer to the Instructions for Use of the Icelock Expulsion Valve 552 for the assembly process.
6. For vacuum suspension, insert Bottom Seal (**Figure 34**).
  7. Check socket fit for the following:
    - The user can securely don the Direct Socket, with or without assistance.
    - Proximal fit is good.
    - No pistoning in the socket.
    - Visible movement is minimal when pulling on prosthesis.
    - Rotational stability is good.
    - Vacuum suspension: Full movement is achieved without compromising the distal vacuum.
  8. Adjust socket if needed or cast a new one if required.

## **USAGE**

### **Donning**

**NOTE:** Do not pull on the brim when socket is donned. Pulling on brim can damage the socket.

1. Fold the brim over the socket (**Figure 35**).
2. Fully step into the socket (**Figure 35**).
3. Fold the brim back up (**Figure 35**).

### **Doffing**

1. Fold the brim over the socket (**Figure 36**).
2. Release the suspension and pull the socket off (**Figure 36**).
3. Fold the brim back up (**Figure 36**).

**NOTE:** Keep the brim folded up when socket is not in use.

### **Cleaning and care**

Clean with a damp cloth and a mild soap.

## **ENVIRONMENTAL CONDITIONS**

Device:

- Operating temperature: -15°C to 50°C (5°F to 122°F)
- Storage temperature: 0°C to 50°C (32°F to 122°F)

Resin:

- Optimal usage temperature: 25°C (77°F)
- Storage temperature: 15°C to 25°C (59°F to 77°F)

## **SUPPLEMENTARY PRODUCTS**

### **Casting Liner**

The liner can be reused if reprocessed according to method in IFU.

Because the casting liners are to be applied on a number of amputees, it is extremely important that each liner is cleaned and disinfected after each and every use.

#### **Cleaning and disinfecting the Iceross Casting Liner:**

1. Turn the liner inside out and wipe the surface with alcohol impregnated wipes until completely wet. Allow to dry.
2. Turn the liner back to normal shape and repeat the procedure on the outer surface.

#### **CAUTION:**

- Only use ethanol or isopropanol solutions for manual cleaning of the liner.
- Only use a washing disinfectant if the liner will be sterilized afterwards. If the liner will not be sterilized, only clean and disinfect manually.

### **Liner Inspection**

1. Check the liner for visible contamination and tears in the silicone layer.
2. Dispose liner if any tears are visible in the silicone layer.
3. If any contamination is visible, repeat manual cleaning and disinfection (see above).

### **WARNING:**

Casting liner shall be disposed after use on amputees that have diagnostic open wound(s) with MRSA bacteria or similar superbugs, e.g.:

- Streptococcus pneumoniae (*S. pneumoniae*)
- Clostridium difficile (*C. diff.*)
- Campylobacter
- Neisseria gonorrhoeae (gonorrhea)
- Salmonella
- Methicillin-resistant Staphylococcus aureus (MRSA)

### **WARNING:**

For patients with open wounds on their amputated limb, all open wounds shall be covered before donning the casting liner to prevent contamination. Thick bandages can affect the socket fit.

**NOTE:** The liner should not be exposed to glass, basalt, or carbon fibers or other foreign particles. Such substances can become embedded in the silicone causing aggravation of the skin. Washing alone may not be sufficient to eliminate the problem.

**NOTE:** Casting liners are for casting purpose only and should not be used as an everyday liner.

### **Direct Socket Tool kit**

The Direct Socket tool kit includes the following parts. Some are also available separately.

- Icecast Pump
- Icecast Bladder
- Resin Injection Tools (200ml and 400ml)
- Distal Attachment Tool
- Tape Ring Kits (Standard and Large)
- Lanyard Cord
- Relief Pads EVA
- Sheet straps
- Brim Tape

**NOTE:** Icecast Bladder is sensitive to puncturing, especially when inflated and inverted. It can be easily damaged by improper handling of sharp objects such as scissors, knives and cured braids. Fingernails may also cause damage during donning and doffing. Ensure Icecast Bladder is not inverted when not in use to minimize risk of damage.

**NOTE:** Clean the Icecast Bladder with a damp cloth and a mild soap.

### **REPORT A SERIOUS INCIDENT**

#### **Important notice to users and/or patients established in Europe:**

The user and/or patient must report any serious incident that has occurred in relation to the device to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

### **DISPOSAL**

The device and packaging must be disposed of in accordance with respective local or national environmental regulations.

### **LIABILITY**

Össur does not assume liability for the following:

- Device not maintained as instructed by the instructions for use.
- Device assembled with components from other manufacturers.
- Device used outside of recommended use condition, application, or environment.

### **Compliance**

This device has been tested according to ISO 10328 standard to three million load cycles. Depending on patient activity, this may correspond to 3-5 years of use.

ISO 10328 - "P" - "m"kg \*)

\*) Body mass limit not to be exceeded!

For specific conditions and limitations of use see manufacturer's written instructions on intended use!

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