























Instructions for Use

DIRECT SOCKET TT







		3
EN	Instructions for Use	5
DE	Gebrauchsanweisung	12
FR	Notice d'utilisation	19
ES	Instrucciones para el uso	26
IT	Istruzioni per l'uso	33
NO	Bruksanvisning	40
DA	Brugsanvisning	47
SV	Bruksanvisning	54
EL	Οδηγίες Χρήσης	61
FI	Käyttöohjeet	68
NL	Gebruiksaanwijzing	75
PT	Instruções de Utilização	82
PL	Instrukcja użytkowania	89
CS	Návod k použití	96
TR	Kullanım Talimatları	103
RU	Инструкция по использованию	110
JA	取扱説明書	118
ZH	中文说明书	125
KO	사용 설명서	131

Direct Socket pressure casting	New amputee	Fleshy	Muscular	Bony
mmHg	40-50	60	70	80

Medical Device

PRODUCT DESCRIPTION

The Direct Socket contains all necessary materials to fabricate a single Direct Socket for prosthetic use.

The Direct Socket TT consists of a Material Kit with the following parts:

- Fiber braids with a distal attachment
- Silicone insulation sheets
- Silicone insulation caps
- Direct Socket resin with static mixer
- Plastic protection sheets
- Trimline strip
- Spacer for Icelock 600* (clear)
- Spacer for Icelock 562 Hybrid* (blue)

*Only supplied with 4" fiber braid version

INDICATIONS FOR USE

- Lower limb amputation.
- Lower limb deficiency.

CONTRAINDICATIONS

None known.

INTENDED USE

The device is a modular component for a lower limb prosthesis that replaces the function and aesthetics of a missing limb. Compatibility of components and the overall prosthesis must be evaluated by a medical professional.

The device is intended for creating a transtibial socket using one of these suspension methods:

- Seal-In® (distal vacuum)
- Seal-In® X
- Pin Locking
- Cushion
- Seal-In® X Locking

The device is for single patient use.

The device is designed for everyday use of low to high active users.

The device is not to be used by patients with the following conditions:

- Residual limb with extreme volume fluctuations when Seal-In liners are used.
- Ultra-short transtibial amputation performed directly under the tuberosity of the tibia.
- Long residual limb exceeding bladder length.

GENERAL SAFETY INSTRUCTIONS

Inform the patient about everything in this document that is required for safe use of this device.

WARNING: If the device shows signs of damage or wear hindering its normal functions, the patient should stop using the device and contact a healthcare professional.

WARNING: In cases of functional change or functional loss, the patient should stop using the device and contact a healthcare professional.

REQUIRED COMPONENTS

For creating a socket the following is required:

Direct Socket Tool Kit

Icross Direct Casting Liner

Suspension components based on suspension type selected.

Bottom Seal (vacuum suspension only)

If using a 4" fiber braid and

Icelock 562: Put the blue spacer between the lock and distal attachment.

Icelock 600: Put the clear spacer between the lock and the modular spacer.

DEVICE SELECTION

Please refer to the table below to verify device selection.

WARNING: Do not exceed weight limit. Risk of device failure.

Braid size	Weight Limit	Label Text
4"	100 kg / 220 lbs	ISO 10328-P5-100 kg
5"	166 kg / 366 lbs	ISO 10328-P7-166 kg
7"	166 kg / 366 lbs	ISO 10328-P7-166 kg
9"	166 kg / 366 lbs	ISO 10328-P7-166 kg

ASSEMBLY INSTRUCTIONS

Direct Socket production is divided into three phases: preparation, casting, and finishing.

For some steps of the fabrication, the assistance of a technician is recommended.

CAUTION: Use protective clothing, gloves, mask and glasses, during fiber braid handling and socket grinding process.

CAUTION: Read the Material Safety Data Sheet (MSDS) before using resin. Refer to the MSDS in case of emergency.

WARNING: Avoid direct contact with resin. Make sure everyone involved wears safety glasses and gloves. Direct contact with the resin can cause serious eye irritation and skin irritation.

WARNING: Do not use a different type of resin. Socket may break.

WARNING: Do not cut fiber braids without air extraction. Fiber particles can cause skin, eye, or respiratory irritation.

Preparation

Relief Pads

Localize bony prominences that require pressure relief. If required, select the appropriate relief pad(s) and place on the residual limb.

NOTE: Too many relief pads or relief pads put on top of each other should be avoided, to prevent cracks in the silicone liner and a bad socket fit. The relief pads are for single use.

Patient Position During Socket Casting

Patients should be seated during socket casting.

Attaching the Injection Valve

Turn the injection valve into the injection channel. An arrow on the distal connector shows you the injection channel (**Figure 1**).

WARNING: Make sure you attach the injection valve correctly. If the injection valve is not in the injection channel, the resin will erupt when you inject it. Resin eruption can cause safety risks.

Taping Ring

Use the screws to attach the taping ring to the distal connector (**Figure 2**).

Liner for Casting

1. Determine the correct size. Measure the circumference of the residual limb at 4 cm from the distal end.
2. Choose the device size which corresponds to the measurement, or the one which is immediately below. Example: If the residual limb size measures 24,5 cm at the indicated location, choose liner size 23,5 instead of size 25.
3. Roll on the liner used for casting (**Figure 3**).
4. Use the lanyard cord to examine soft tissue stabilization (**Figure 4**). If you are not sure, try a casting liner one size above or below.

WARNING: The casting liner must be cleaned and disinfected before each use. See Casting Liner section below.

NOTE: A wrong liner size can lead to an uncomfortable socket. If you make an uncomfortable socket, you will have to make a new socket.

Distal Attachment Pin

1. Turn the distal attachment pin into the casting liner umbrella.
2. Make sure the pin and umbrella are correctly aligned with the patient's limb.

Silicone Insulation Sheet, First Layer

Pull the first length of silicone insulation sheet over the casting liner (**Figure 5**).

NOTE: The textile layer must be on the inner side. If it is not on the inner side, you will laminate the socket to the silicone insulation sheet.

Putting on the Silicone Insulation Cap

1. Use scissors to cut away silicone insulation sheet that is not necessary (**Figure 6**). Make sure the insulation sheet ends a small distance above the distal part of the casting liner.
2. Put the silicone insulation cap on the distal part of the casting liner (**Figure 7**).

Attach the Fiber Braids with Distal Attachment

1. Fiber braids with small or medium distal attachment (**Figure 8**):
Position the injection channel anteriorly (1 o'clock). This will improve resin flow.

Fiber braids with large distal attachment:

The orientation of the attachment offset will determine injection channel position.

If the injection channel not positioned anteriorly, pay more attention to distributing the resin evenly.

2. Turn the distal attachment nut onto the distal attachment pin and tighten nut firmly by hand.
NOTE: Make sure no fiber braids are caught between the residual limb and the distal connector.
3. Pull the first layer of fiber braids up over the silicone insulation sheet. Pull the second layer of fiber braids up over the first layer. Repeat for all layers.
4. Check the fiber braids for folds and creases. Make the fiber braids smooth where necessary.
5. Mark a cutting line 25mm above the expected trimline (**Figure 9**).
6. Remove the fiber braids with distal attachment.
7. Cut the fiber braids along the cutting line with scissors. Do not cut along expected trimline! Assemble again.
Repeat steps 1 to 4.

Silicone Insulation Sheet, Second Layer

Use the sheet applicator to roll the second length of silicone insulation sheet over the fiber braids (**Figure 10**).

NOTE: The textile layer must be on the outer side. If it is not on the outer side, you will laminate the direct socket to the silicone insulation sheet.

O-rings

1. Take the large O-ring and put it on the top edge of the distal connector (**Figure 11**).
NOTE: When using small distal attachment, do not apply the large O-ring.
2. Take the small O-ring and put it into the groove on the taping ring (**Figure 12**).
CAUTION: Make sure there are no wrinkles and air between the O-rings. Wrinkles and air can cause the resin to leak.
3. Tape over the O-rings to fix the second layer of silicone insulation sheet to the distal connector and the taping ring (**Figure 13**).
NOTE: When using small distal attachment, do not tape above the distal connector.
NOTE: Do not tape above the proximal O-ring. Tape above the proximal O-ring can stop the resin flow. You will not be able to cast a usable socket if the resin does not flow and saturate the fiber braids.
4. Use scissors to cut away silicone insulation sheet that is not necessary below the taping ring (**Figure 14**).
5. Place the Air Escape Tube at the proximal end laterally between the first and second layer of silicone insulation sheet. (**Figure 15**). This will allow air to escape during resin injection. Both layers of silicone insulation sheet must fully cover the residual limb and fiber braids.

Casting

Preparing the Resin Injection Tool

1. Cover the working area and the floor with protective sheeting.
2. Remove the parts from the packaging.
3. Shake the resin cartridge before use.
4. Lightly turn off the cap from the resin cartridge. The cap will lift off the green safety pin.
5. Remove the metal safety ring from the green safety pin.
6. Remove the green safety pin and discard it.
7. Put the cap onto the static mixer (**Figure 16**).
8. Push the static mixer 1 centimeter/ 0.4 inches deep into the injection tube (**Figure 17**).
9. Turn the cap tight to safely attach the static mixer to the resin cartridge.
10. Keep the resin cartridge pointed up and insert it into the resin injection tool (**Figure 17**).
11. Lightly press the handle to secure the resin cartridge until it softly connects to bottom of the cartridge.
12. Keep the resin injection tool with the cartridge pointed up and sit in front of the patient.
13. The technician fully attaches the injection tube to the injection valve (**Figure 18**).

Resin Injection

WARNING: Make sure the injection tube is fully secured to the static mixer and injection valve. If the tube is not fully secured, the resin can erupt when it is injected.

NOTE: Resin shall be injected gradually and spread evenly towards the proximal end. Excessive resin injected at the distal end can cause the resin to get hot and cure faster. This will make it harder to inject more resin .

1. Keep the resin injection tool pointed up while you start to inject the resin slowly (**Figure 19**).
2. When resin has saturated the distal area, injection speed can be increased.
3. The technician slowly pulls the resin into the interior and upwards, pushing out the air during resin injection (**Figure 20**). The lanyard cord can be used to push the resin around the socket
4. Control the resin flow to the posterior side of the residual limb while you maintain a slow build-up.
NOTE: Make sure to spread and saturate the braids up to and above the expected trimline, before the resin starts to cure.
5. Stop the resin flow when the fiber braids are saturated. To stop, push the pressure release lever on the grip of the resin injection tool (**Figure 21**).
6. Use pliers on crimp rings to shut off resin flow in the injection tube (**Figure 22**).
7. Place side cutter between the two crimp rings and cut through the injection tube.
8. Make sure the patient keeps the residual limb straight and in position at this stage.
9. Attach Icecast Pump to Icecast Bladder. Inflate Icecast Bladder to 20mmHg. Roll back Icecast Bladder and attach to Distal Attachment (**Figure 23**).
10. Keep knee and pump aligned straight
11. Push on Icecast Bladder, until covering the entire cast. Make sure resistance is as low as possible. If necessary, lower resistance by releasing air from the Icecast Bladder.
12. Use pump and inflate to predetermined pressure level (**Figure 24**).

Bladder pressure	New amputee	Fleshy	Muscular	Bony
mmHg	40-50	60	70	80

NOTE: Do not inflate to more than 100 mmHg.

13. Disconnect Icecast Pump. Ask the patient to relax the quadriceps muscle to create normal knee flexion. Add support. Allow resin to cure for 10 minutes.

Casting Disassembly

1. Remove the tape from the taping ring.
2. Remove the two O-rings.
3. Remove the air escape tube
4. Remove the outer layer of silicone insulation sheet.
5. Mark the trimline on the finished socket (**Figure 25**).
6. Remove the distal attachment nut and the tape ring.
7. Remove the socket from the residual limb.
8. Pull the inner layer of silicone insulation sheet and the silicone insulation cap out of the socket.

9. For vacuum suspension only: break off the injection valve to seal the injection channel. Screwing out the valve can cause air leak and compromise the vacuum suspension.
10. For other suspension methods: Unscrew the injection valve. Do not break off the valve. The Icelock 214 can not be assembled if valve is not removed.

Finishing

1. Use a jigsaw to cut along the trimline and shape proximal end of socket (**Figure 26**).
2. Insert lanyard again into casting liner, for easy donning of socket. Don the socket and test the trimline on the residual limb with the casting liner. Find the tendons on the back of the knee. Make a relief to obtain 90° flexion (**Figure 27**).
3. Machine grind and sand the edge (**Figure 28**).
4. Water sand the edge (**Figure 29**). Make sure the edge is rounded and not sharp.
5. Cover the socket edges with trimline strip.
6. For vacuum suspension, insert Bottom Seal (**Figure 30**).
7. Check socket fit for the following:
 - The user can securely don the Direct Socket, with or without assistance.
 - Proximal fit is good.
 - No pistoning in the socket.
 - Visible movement is minimal when pulling on prosthesis.
 - Rotational stability is good.
 - Vacuum suspension: Full movement is achieved without compromising the distal vacuum.
8. Adjust socket if needed or cast a new one if required.

Cleaning and care

Clean with a damp cloth and a mild soap.

ENVIRONMENTAL CONDITIONS

Device:

- Operating temperature: -15°C to 50°C (5°F to 122°F)
- Storage temperature: 0°C to 50°C (32°F to 122°F)

Resin:

- Optimal usage temperature: 25°C (77°F)
- Storage temperature: 15°C to 25°C (59°F to 77°F)

SUPPLEMENTARY PRODUCTS

Casting Liner

Icross Direct Casting Liner is intended to be used with Direct Socket during casting. The liner is not intended to be used by the patient for any other purpose.

The casting liner is a multiple patient reusable device. The liner can be reused if reprocessed according to method in IFU.

Because the casting liners are to be applied on a number of amputees, it is extremely important that each liner is cleaned and disinfected after each and every use.

Cleaning and disinfecting the Icross Casting Liner:

1. Turn the liner inside out and wipe the surface with alcohol impregnated wipes until completely wet. Allow to dry.
2. Turn the liner back to normal shape and repeat the procedure on the outer surface.

CAUTION:

- Only use ethanol or isopropanol solutions for manual cleaning of the liner.
- Only use a washing disinfectant if the liner will be sterilized afterwards. If the liner will not be sterilized, only clean and disinfect manually.

Liner Inspection

1. Check the liner for visible contamination and tears in the silicone layer.
2. Dispose liner if any tears are visible in the silicone layer.

3. If any contamination is visible, repeat manual cleaning and disinfection (see above).

WARNING:

Casting liner shall be disposed after use on amputees that have diagnostic open wound(s) with MRSA bacteria or similar superbugs, e.g.:

- Streptococcus pneumoniae (S. pneumoniae)
- Clostridium difficile (C. diff.)
- Campylobacter
- Neisseria gonorrhoeae (gonorrhea)
- Salmonella
- Methicillin-resistant Staphylococcus aureus (MRSA)

WARNING:

For patients with open wounds on their amputated limb, all open wounds shall be covered before donning the casting liner to prevent contamination. Thick bandages can affect the socket fit.

NOTE: The liner should not be exposed to glass, basalt, or carbon fibers or other foreign particles. Such substances can become embedded in the silicone causing aggravation of the skin. Washing alone may not be sufficient to eliminate the problem.

NOTE: Casting liners are for casting purpose only and should not be used as an everyday liner.

Direct Socket Tool kit

The Direct Socket tool kit includes the following parts. Some are also available separately.

- Icecast Pump
- Icecast Bladder
- Resin Injection Tools (200ml and 400ml)
- Distal Attachment Tool
- Tape Ring Kits (Standard and Large)
- Lanyard Cord
- Relief Pads EVA
- Sheet straps
- Brim Tape

NOTE: Icecast Bladder is sensitive to puncturing, especially when inflated and inverted. It can be easily damaged by improper handling of sharp objects such as scissors, knives and cured braids. Fingernails may also cause damage during donning and doffing. Ensure Icecast Bladder is not inverted when not in use to minimize risk of damage.

NOTE: Clean the Icecast Bladder with a damp cloth and a mild soap.

REPORT A SERIOUS INCIDENT

Important notice to users and/or patients established in Europe:

The user and/or patient must report any serious incident that has occurred in relation to the device to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

DISPOSAL

The device and packaging must be disposed of in accordance with respective local or national environmental regulations.

LIABILITY

Össur does not assume liability for the following:

- Device not maintained as instructed by the instructions for use.
- Device assembled with components from other manufacturers.
- Device used outside of recommended use condition, application, or environment.

Compliance

This device has been tested according to ISO 10328 standard to three million load cycles. Depending on patient activity, this may correspond to 3-5 years of use. The load level specification is shown in the tables below.

ISO 10328 - "P" - "m"kg *) 
*) Body mass limit not to be exceeded!  For specific conditions and limitations of use see manufacturer's written instructions on intended use!

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