






Instructions for Use

FLEX[®]

		3
EN	Instructions for Use	5
DE	Gebrauchsanweisung	7
FR	Notice d'utilisation	10
ES	Instrucciones para el uso	12
IT	Istruzioni per l'uso	15
NO	Bruksanvisning	17
DA	Brugsanvisning	19
SV	Bruksanvisning	22
EL	Οδηγίες Χρήσης	24
FI	Käyttöohjeet	27
NL	Gebruiksaanwijzing	29
PT	Instruções de Utilização	32
CS	Návod k použití	34
SK	Indikácie použitia	36

INTENDED USE

The device is intended for external support, stabilization and protection of the knee.

The device must be fitted and adjusted by a healthcare professional.

Indications for use

- ACL, MCL, LCL, PCL, rotary, and combined instabilities.
- Long-term use for post-op, rehab, and functional support.

No known contraindications.

Warnings and Cautions:

- The use of the device may increase the risk of deep vein thrombosis and pulmonary embolism.

GENERAL SAFETY INSTRUCTIONS

The healthcare professional should inform the patient about everything in this document that is required for safe use of this device.

Any serious incident in relation to the device must be reported to the manufacturer and relevant authorities.

The patient should stop using the device and contact a healthcare professional:

- If there is a change or loss in device functionality, or if the device shows signs of damage or wear hindering its normal functions.
- If any pain, skin irritation, or unusual reaction occurs with the use of the device.

The device is for single patient – multiple use.

This device is a supportive brace only and is not intended or guaranteed to prevent injury.

FITTING INSTRUCTIONS

Device Application

1. Unfasten all straps from the posterior side of the device and loosen Anterior Tibial Strap (A).
2. Bend the knee at a 90° angle, with the foot flat on the floor (**Fig. 2**).
3. Place the device on the leg such that the patella is centered between the Condyle Pads (B).

Ensure proper alignment of the device on the leg:

4. Height positioning: Align the center of the Hinge (C) slightly above the middle of the patella.
5. The device can be adjusted for individual knee width: To widen or narrow the width at the knee, the gel Condyle Pad can be removed and replaced with thinner neoprene pads [3 mm (1/8") or 1.5 mm (1/16")] or you can install a shim behind the gel Condyle Pad to make the fit more snug. Any combination of pads can be used to achieve a more comfortable fit (**Fig 6, 7**). Secure neoprene pads in place with covers, which are also included (**Fig. 8, 9**).

Fasten the straps in the following order. All straps can be adjusted to length.

6. Fit the Lower Cruciate Strap (D) first, located directly below the knee (**Fig. 3**). This strap is located above the calf muscle and aids in the correct positioning of the device as it helps create the device-to-bone contact required to help keep the device in place.
 7. Next, snugly fasten the Bottom Strap (E) by inserting it through the D-Ring.
 8. Continue with the Upper Cruciate Strap (F), located just above the knee (**Fig. 4**).
 9. Fasten the Top Strap (G). Make sure that the strap hook tab does not overlap the elastic insert on the strap.
- Note:** Overtightening the top two straps can cause the device to migrate.
10. Last, snugly fasten the Anterior Tibial Strap (A) (**Fig. 5**).
 - If the device is too far anterior, loosen the anterior tibial strap and tighten the lower and bottom straps.
 - If the device is too far posterior, loosen the lower and bottom straps and tighten the tibial control strap.

Device Adjustments

Range of Motion (ROM) Adjustment

The amount of flexion and extension can be adjusted by changing the size of the flexion and extension stops. Changing stops should be done under the advice of your healthcare professional.

Note: If you change the stop on one hinge, you must change the stop to the same size on the opposite hinge.

Extension Stop Sizes: 10°, 20°, 30°, 40°, lockout (full extension).

Flexion Stop Sizes: 45°, 60°, 75°, 90°.

To adjust flexion and extension stops:

1. Remove the screw securing stop cover. Note the EXT and FLX markings on side of the stop cover indicating the placement of extension and flexion stops (**Fig. 10**).
2. Lift the stop cover and insert desired extension stop into the top half of the groove (**Fig. 11**).
3. Make sure the notch in the extension stop faces towards the metal bumper located in the middle of the groove (**Fig. 12**).
4. Insert desired flexion stop into the lower half of the groove with the notch facing towards the metal bumper (**Fig. 13**).
5. Close stop cover and replace the screw.

Device Removal

1. Sit down and bend the knee at an 80° angle with the foot flat on the floor.
2. Unfasten all posterior straps of the device and attach the straps back to themselves. This will make it easier to put the device on next time and will also prolong the life of the straps.

Accessories and Replacement Parts

Please refer to the Össur catalog for a list of available replacement parts or accessories.

USAGE

Cleaning and care

Washing the device with the soft goods detached allows for more thorough cleaning.

Washing Instructions

- Hand-wash using mild detergent and rinse thoroughly.
- Air dry.

Note: Do not machine-wash, tumble dry, iron, bleach, or wash with fabric softener.

Note: If the device is used in salt water or chlorinated water, rinse with fresh water and air dry.

Hinge

- Remove foreign materials (e.g., dirt or grass) and clean using fresh water.

DISPOSAL

The device and packaging must be disposed of in accordance with respective local or national environmental regulations.

LIABILITY

Össur does not assume liability for the following:

- Device not maintained as instructed by the instructions for use.
- Device assembled with components from other manufacturers.
- Device used outside of recommended use condition, application, or environment.

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