ENGLISH

DESCRIPTION OF SYMBOLS

MD Medical Device

USES

Neck support

CAUTION

- Patients should be instructed to inspect their neck and this product before and after each use.
- Modifications to this product are the sole responsibility of the prescribing Physician.
- Adjustments may be necessary for function and comfort.

PHOTO CAPTIONS

- V-notch
- Contact Closure
- Foam Brace

FITTING INSTRUCTIONS

- Detach contact closure strap. Center the V-notch of the brace underneath the chin
- Wrap brace around neck.
- Fasten contact closure and secure to a snug fit.

WASHING INSTRUCTIONS

- The brace may be cleaned using a mild facial soap and cold water. Do not wash in washing machine.
- Lay brace out flat on a towel and air dry. Do not squeeze or wring out the brace. Do not place in dryer. Do not place near heater. Do not use bleach or strong detergents.

GENERAL SAFETY INFORMATION

WARNING: In case of functional change or functional loss, the patient should stop using the device and contact a healthcare professional.

WARNING: If the device shows signs of damage or wear hindering its normal functions, the patient should stop using the device and contact a healthcare professional.

REPORT A SERIOUS INCIDENT

Important notice to users and/or patients established in Europe:

The user and/or patient must report any serious incident that has occurred in relation to the device to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

DISPOSAL

The device and packaging must be disposed of in accordance with respective local or national environmental regulations.

LIABILITY TEXT

Össur does not assume liability for the following:

- Device not maintained as instructed by the instructions for use.
- Device assembled with components from other manufacturers.
- Device used outside of recommended use condition, application, or environment.

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