Note: Do not machine-wash, tumble dry, iron, bleach, or wash with fabric softener.
Note: Avoid contact with salt water or chlorinated water. In case of contact, rinse with fresh water and air dry.
Note: Bleach can cause the device to discolor and decrease elasticity. If elasticity is decreased, the device will not work as intended.

DISPOSAL

The device and packaging must be disposed of in accordance with respective local or national environmental regulations.

LIABILITY

Össur does not assume liability for the following:

- Device not maintained as instructed by the instructions for use.
- Device assembled with components from other manufacturers.
- Device used outside of recommended use condition, application, or environment.

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SYMBOLS

Medical Device

Caution: This device contains natural rubber latex which may cause allergic reactions.

INTENDED USE

The device is intended to provide gross immobilization of the lumbar spine and/or compression

Indications for use

Lower back conditions that may benefit from gross motion restriction, such as:

- Muscular low back pain
- Lumbosacral Strains
- Chronic low back pain
- Sciatica

Contraindications

- Conditions requiring more rigid or higher immobilization (> D12) or a surgical procedure, unless medically indicated.
- Pelvic floor dysfunction.

Warnings and Cautions:

- Do not use the device if your skin is prone to skin lesions from pressure or rubbing. Compression from the device can cause a skin lesion.
- Do not use the device if you have a skin lesion in the lower back region where the device is worn. Compression from the device can make a skin lesion worse.
- The device should not be used for extended periods of time and removed during long periods of rest, e.g., sleep.

GENERAL SAFETY INSTRUCTIONS

Read these instructions carefully before use. Keep them for future reference.

Any serious incident in relation to the device must be reported to the manufacturer and relevant authorities.

The patient should stop using the device and contact a healthcare professional:

- If there is a change or loss in device functionality, or if the device shows signs of damage or wear hindering its normal functions.
- If any pain, skin irritation, or unusual reaction occurs with the use of the device.

The device is for single patient – multiple use.

USAGE

Device Application

1. Unfasten the elastic outer straps and secure them to the device so that there is no tension on the outer straps (Fig. 1).
2. Align the device behind the patient with the lumbar pad centered on the lower back (Fig. 2).
3. Wrap the arms around the waist in the order shown (Fig. 3, 4)
4. Pull the ends of the outer straps forward equally and fasten the two ends onto the ends of the inner belt (Fig. 5, 6)
5. Air version only: Turn the valve clockwise until the inflate icon is on top and press the air pump until the desired inflation is achieved (Fig. 7, 8). Do not over-turn the valve.

Device Removal

Unfasten the outer straps and then unfasten the belt arms to remove the device.

Air version only: Turn the valve counterclockwise until the deflate icon is on top and press the air pump to deflate the bladder.

Accessories and Replacement Parts

Please refer to the Össur catalog for a list of available replacement parts or accessories.

Cleaning and care

- Hand-wash using mild detergent and rinse thoroughly.
- Air dry, do not expose to direct heat.

Össur Americas
27051 Towne Centre Drive
Foothill Ranch, CA 92610, USA
Tel: +1 (949) 382 3883
Tel: +1 800 233 6263
ossurusa@ossur.com

UK REP
Össur UK Ltd
Unit No 1, S:Park
Hamilton Road
Stockport SK1 2AE, UK
Tel: +44 (0) 8450 065 065
ossuruk@ossur.com

Össur Europe BV
De Schakel 70
5651 GH Eindhoven
The Netherlands
Tel: +800 3539 3668
Tel: +31 499 462840
info-europe@ossur.com

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©Copyright Össur 2023-05-25