## ENGLISH

**MD** Medical Device

### INTENDED USE

The device is intended to support the ankle.

### Indications for use

Situations requiring inversion/eversion restriction. These may include.

- Mild to severe ankle sprains
- Chronic ankle instability
- Post-cast and post walker support
- Ankle injury rehabilitation and/or post-surgical use

No known contraindications.

### Warnings and Cautions:

- Always wear an absorbent sock under the device.
- Care should be taken not to overtighten the device.
- Patients with already compromised blood circulation in their extremities, like those with diabetes or peripheral vascular disease, need to be extra careful and are advised to consult with a healthcare professional before using the device. If you experience pain, swelling, sensation changes, or if your extremity shows signs of insufficient blood flow (turns blue, white, or cold) while using this device, discontinue the use immediately and consult with a healthcare professional.

### GENERAL SAFETY INSTRUCTIONS

Read these instructions carefully before use. Keep them for future reference.

Any serious incident in relation to the device must be reported to the manufacturer and relevant authorities. Never wear the device directly over an open wound.

The patient should stop using the device and contact a healthcare professional:

- If there is a change or loss in device functionality, or if the device shows signs of damage or wear hindering its normal functions.
- If any pain, skin irritation, or unusual reaction occurs with the use of the device.

The device is for single patient – multiple use.

### USAGE

#### Device Application

This device is designed to be worn over an athletic sock. For safe use and optimal support wear the device inside a laced athletic shoe.

1. Undo the contact closure straps.
2. Spread the side stirrups and step into the device. Center the sides along the lower leg and ankle (**Fig. 1**).
3. Hold the stirrups up against the leg while wrapping the lower contact closure strap snugly around the device (**Fig. 2**). Secure contact closure.
4. If required, adjust the heel width by peeling back the bottom layer and reapply by engaging the hook and loop where the brace is snug (**Fig. 3**). Position the bottom of the stirrups to be at or slightly above the bottom of the foot.
5. Wrap the contact closure top strap snugly around device (**Fig. 4**). Secure contact closure.
6. Apply shoe and tighten laces (**Fig. 5**). Readjust the device as necessary for maximum comfort and support.

#### Device removal

Unfasten the straps to remove the device.

#### Cleaning and care

- Wipe with a damp cloth and mild soap.
- Air dry

**Note:** Do not machine-wash, tumble dry, iron, bleach, or wash with fabric softener.

**Note:** Avoid contact with salt water or chlorinated water. In case of contact, rinse with fresh water and air dry.

### DISPOSAL

The device and packaging must be disposed of in accordance with respective local or national environmental regulations.

### LIABILITY

Össur does not assume liability for the following:

- Device not maintained as instructed by the instructions for use.
- Device assembled with components from other manufacturers.
- Device used outside of recommended use condition, application, or environment.

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