







Instructions for Use

---

FORMFIT<sup>®</sup> POST-OP KNEE  
16"

		3
EN	Instructions for Use	4
DE	Gebrauchsanweisung	6
FR	Notice d'utilisation	9
ES	Instrucciones para el uso	11
IT	Istruzioni per l'uso	14
NO	Bruksanvisning	16
DA	Brugsanvisning	19
SV	Bruksanvisning	21
EL	Οδηγίες Χρήσης	23
FI	Käyttöohjeet	26
NL	Gebruiksaanwijzing	29
PT	Instruções de Utilização	31
PL	Instrukcja użytkowania	34
CS	Návod k použití	36
TR	Kullanım Talimatları	39
RU	Инструкция по использованию	41
JA	取扱説明書	44
ZH	中文说明书	46
KO	사용 설명서	48
SK	Návod na použitie	51
ET	Kasutusjuhend	53
RO	Indicații pentru Utilizare	56
AR	تعليمات الاستخدام	58

# ENGLISH

---

## SYMBOLS

Medical Device

## INTENDED USE

The device is intended to immobilize the knee

The device must be fitted and adjusted by a healthcare professional.

### *Indications for use*

- Post-op/post-injury ROM control or immobilization.
- Ligament and fracture-related surgeries and when a controlled range of motion is necessary.

No known contraindications.

### **Warnings and Cautions:**

- The use of the device may increase the risk of deep vein thrombosis and pulmonary embolism.
- Care should be taken not to overtighten the device.
- Regular healthcare professional supervision is recommended for patients with peripheral vascular disease, neuropathy, and sensitive skin.
- The device is not designed nor intended to prevent injury or withstand substantial loads during impacts such as falling or jumping.
- Make sure to apply the stop locks over the flexion and extension buttons to secure the desired ROM setting. Not using the stop locks may result in an unintentional change to both flexion and extension angles.
- Any adjustments to the device should be performed by a healthcare professional.
- Proper rehabilitation and activity modification are also essential parts of a safe treatment program.
- Prolonged immobilization may result in decreased strength and decreased range of motion around immobilized joint.
- Open wounds or damaged skin should be covered with a bandage or other appropriate covering to prevent direct contact between the wound and device.

## GENERAL SAFETY INSTRUCTIONS

The healthcare professional should inform the patient about everything in this document that is required for safe use of this device.

Any serious incident in relation to the device must be reported to the manufacturer and relevant authorities.

The patient should immediately contact a healthcare professional:

- If there is a change or loss in device functionality, or if the device shows signs of damage or wear hindering its normal functions.
- If any pain, skin irritation, excessive pressure, or unusual reaction occurs with the use of the device.

The device is for single patient – multiple use.

## FITTING INSTRUCTIONS

### *Device Application*

**Warning:** Make sure to keep the leg straight when you slide the foam pads and straps under it.

1. Unfasten the leg straps.

2. Fully pull the open ends of the straps through the double D-rings to remove them from the left strut.
3. Remove the foam pads from the two struts.
4. Fully open the foam pads.
5. Wrap the foam calf and thigh pads tightly around the leg and fasten at the front. **(Fig. 2)**.
6. Put the struts on each side of the leg. If necessary, use bending tools to adapt the brace to the patient's leg shape **(Fig. 3)**.  
**Note:** Lower strut should extend no lower than calf/Achilles tendon junction. Upper strut should extend no higher than 2" or 5 cm from groin level.
7. Attach the struts to the foam pads on each side of the leg.  
**Note:** Align hinge center to the top of the patella and ensure both hinge center and struts are aligned with the mid-line of the leg **(Fig. 4)**.
8. Carefully slide the straps under the foam wraps.
9. Put the straps back through the double D-rings.
10. Pull and fasten the strap below the knee first; followed by the strap above the knee. Complete the remaining straps alternating below and above the knee until all straps are fastened **(Fig. 5)**.  
**Caution:** Do not overtighten the straps. If the straps are too tight, the blood circulation can be cut off.

### **Device Adjustments**

Use the drop lock to fully immobilize the knee. It can be set at -10°, 0°, 10°, 20°, 30° of flexion. An arrow shows the angle. The angle is shown by the scale on the strut, below the hinge.

Adjust immobilization as follows **(Fig. 7)**:

1. Push the blue drop-lock button up to unlock the hinge.
2. Move the brace until the arrow and the correct angle agree.
3. Push the blue drop-lock button down to lock the hinge angle.

Adjust the ROM as follows **(Fig. 6)**:

Extension: Pull the black button out and slide it to the desired angle.

Extension can be limited at: -10°, 0°, 10°, 20°, 30°, 40°, 50°, 60°, 70°, 80°, 90°.

Flexion: Pull the white button out and slide it to the desired angle. Flexion can be limited at: 0°, 10°, 20°, 30°, 40°, 50°, 60°, 70°, 80°, 90°, 100°, 110°, 120°.

Release a button to let it click into place. Put the stop locks onto the buttons to secure them **(Fig. 6)**.

### **Device Removal**

1. Unfasten the straps and the foam pads. Attach the straps back to themselves to make it easier to put the device on next time and to prolong the life of the straps.
2. Remove the device **(Fig. 8)**. You can put it back on again fully assembled.

### **Accessories and Replacement Parts**

Please refer to the Össur catalog for a list of available replacement parts or accessories.

## **USAGE**

### **Cleaning and care**

- Hand-wash using mild detergent and rinse thoroughly.
- Air dry.

**Note:** Do not machine-wash, tumble dry, iron, bleach, or wash with fabric softener.

**Note:** Avoid contact with salt water or chlorinated water. In case of contact, rinse with fresh water and air dry.

## DISPOSAL

The device and packaging must be disposed of in accordance with respective local or national environmental regulations.

## LIABILITY

Össur does not assume liability for the following:

- Device not maintained as instructed by the instructions for use.
- Device assembled with components from other manufacturers.
- Device used outside of recommended use condition, application, or environment.

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