





Instructions for Use

GENU SLEEVE
PROTECTOR SLEEVE

		3
EN	Instructions for Use	4
CS	Návod k použití	7
SK	Indikácie použitia	10
HR	Upute za uporabu	16
HU	Használati javallat	19
BG	Указания за употреба	22
SL	Indikácie pre použitie	25

Medical Device

DESCRIPTION

The device is a sleeve intended to be used as a suspension mechanism for Össur cushion liners.

INTENDED USE

The device is a prosthetic interface with suspension properties intended to be used as part of a system that replaces a missing lower limb. Suitability of the device for the prosthesis and the patient must be evaluated by a healthcare professional.

The device must be fitted and adjusted by a healthcare professional.

Indications For Use and Target Patient Population

- Lower limb loss, amputation, or deficiency
- No known contraindications

The device is for low to moderate impact use, e.g., walking.

GENERAL SAFETY INSTRUCTIONS

Warning: Use of a lower limb prosthetic device carries an inherent risk of falling which may lead to injury.

The healthcare professional should inform the patient about everything in this document that is required for safe use of this device.

Warning: If there is a change or loss in device functionality, or if the device shows signs of damage or wear hindering its normal functions, the patient should stop using the device and contact a healthcare professional.

Warning: Misuse of the device may result in loss of suspension.

Caution: Potential side-effects are skin irritation/reaction problems.

Caution: The device must not come in contact with glass, carbon fibers, or other foreign particles. Washing alone may not be sufficient to eliminate the problem. If the device is exposed to foreign substances or chemicals, it should be returned to the healthcare professional.

The device is for single patient use.

DEVICE SELECTION

To determine the correct size:

Measure the circumference of the residual limb 15 cm above mid-patella (**Fig. 1**) and choose the corresponding device size. See the sizing chart below.

Size	Measurement (cm)
S	30-40 cm
M	38-50 cm
L	48-64 cm

FITTING INSTRUCTIONS

Trimming Instructions

The device may be trimmed at either end to improve the fit to the residual limb and/or the prosthetic socket.

Caution: When cutting the device, ensure a clean cut is made. Any sharp edges or notches can lead to tears in the material and/or skin irritation.

Socket Fitting

Caution: The socket must be airtight in order to maintain suction suspension when the sleeve is applied.

Caution: The proximal brim of the socket must be smooth and without sharp edges that could damage the device.

USAGE

Donning

1. Turn the device inside out, halfway, with the inner surface of the distal part half exposed.
2. Stretch the device over the brim of the socket and roll the distal end completely down (**Fig. 2**).
3. Roll the upper part down until the socket brim is completely exposed.
4. Align the prosthesis and step into the socket.
5. Roll the device up with both hands. Do not pull on the edges as this may place strain on the skin (**Fig. 3**).

The device should be put on with the knee bent slightly to avoid high pressure on the patellar tendon.

At its proximal end, the device shall contact the intact skin directly to obtain correct suspension performance.

Caution: Open wounds or damaged skin should be covered with a bandage or other appropriate covering to prevent direct contact between the wound and device.

Caution: Do not apply lotion to the residual limb immediately before donning the device. The residual limb should be clean and dry to prevent device damage.

Cleaning and care

Daily cleaning of the residual limb is essential. Use of a mild, pH balance, 100% fragrance- and dye-free liquid soap is recommended. If the skin is dry, apply a pH balance, 100% fragrance- and dye-free lotion to nourish and soften the skin.

Clean the device daily after use and before first use.

- Turn the device inside out.
- Wash with pH balanced, 100% fragrance- and dye-free liquid soap.
- After washing, rinse the device thoroughly with warm water and pat dry on both sides with a lint-free cloth (**Fig. 4**).
- Do not machine wash.
- Always return the device to its neutral state immediately after cleaning. Do not store inside-out as device damage may result.
- The device should not be exposed to excessive heat or sunlight.
- Take care not to scrub the surface of the device. That can roughen the surface, which can cause skin irritation.

Note: Please insert the foam supplied when storing the device to prevent it from sticking together.

Warning: Excessive perspiration may compromise stability and/or suspension. Remove the device and dry off both the device and the residual limb when appropriate.

Caution: Skin health should be monitored daily. In case of abnormal symptoms, discontinue use of the device and contact the healthcare professional immediately.

Caution: Be aware that common household or bath products may cause or contribute to skin irritation, e.g., soaps, deodorants, perfumes, abrasive cleaners, aerosol- or alcohol sprays.

Caution: After contact with salt and/or chlorinated water, wash the device with tap water.

Environmental Conditions

The device is safe to use with caution with waterproof components.

Warning: After submerging the device in water, it needs to be removed, and the residual limb and the skin contact layer dried before ambulation.

Warning: The suspension is compromised if water gets between the skin and the device.

Caution: Frequent exposure to salt and/or chlorinated water might affect the durability of the device.

SUPPLEMENTARY DEVICES

Genu Sleeve Protector

The device is intended to protect a prosthetic sleeve a part of a prosthetic system.

The sleeve can be used with the Protector to prolong the use of the sleeve. The Protector shall be put on the socket before the sleeve, with the adhesive strap downwards (**Fig. 5**).

REPORT OF SERIOUS INCIDENT

Any serious incident in relation to the device must be reported to the manufacturer and relevant authorities.

DISPOSAL

The device and packaging must be disposed of in accordance with respective local or national environmental regulations.

LIABILITY

Össur prosthetic devices are designed and verified to be safe and compatible in combination with each other and custom-made prosthetic sockets with Össur adapters, and when used in accordance with their intended use.

Össur does not assume liability for the following:

- Device not maintained as instructed by the instructions for use.
- Device assembled with components from other manufacturers.
- Device used outside of recommended use condition, application, or environment.

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