









Instructions for Use

ICEFORM® CUSHION
ICEFORM® LOCKING

		3
EN	Instructions for Use	5
DE	Gebrauchsanweisung	10
FR	Notice d'utilisation	16
ES	Instrucciones para el uso	22
IT	Istruzioni per l'uso	27
NO	Bruksanvisning	33
DA	Brugsanvisning	38
SV	Bruksanvisning	43
EL	Οδηγίες Χρήσης	48
FI	Käyttöohjeet	54
NL	Gebruiksaanwijzing	59
PT	Instruções de Utilização	65
PL	Instrukcja użytkowania	71
TR	Kullanım Talimatları	76
RU	Инструкция по использованию	81
JA	取扱説明書	87
ZH	中文说明书	92
KO	사용 설명서	96

DESCRIPTION

The device is a transtibial TPE liner that covers the residual limb and provides an interface between the residual limb and the socket of the prosthesis. The device provides a secure suspension to the prosthesis, protection of the residual limb, and comfort.

A locking liner has a distal attachment, to be used with a mechanical lock. A cushion liner has no distal attachment and needs to be used together with a sleeve that holds vacuum.

INTENDED USE

The device is a prosthetic interface with suspension properties intended to be used as part of a system that replaces a missing lower limb. Suitability of the device for the prosthesis and the patient must be evaluated by a healthcare professional.

The device must be fitted and adjusted by a healthcare professional.

Indications For Use and Target Patient Population

- Lower limb loss, amputation, or deficiency
- No known contraindications

The device is for low to moderate impact use, e.g., walking.

GENERAL SAFETY INSTRUCTIONS

Warning: Use of a lower limb prosthetic device carries an inherent risk of falling which may lead to injury.

The healthcare professional should inform the patient about everything in this document that is required for safe use of this device.

Warning: If there is a change or loss in device functionality, or if the device shows signs of damage or wear hindering its normal functions, the patient should stop using the device and contact a healthcare professional.

Warning: Misuse of the device may result in loss of suspension.

Caution: Potential side-effects are skin irritation/reaction problems.

Caution: The device must not come in contact with glass, carbon fibers, or other foreign particles. Washing alone may not be sufficient to eliminate the problem. If the device is exposed to foreign substances or chemicals, it should be returned to the healthcare professional.

The device is for single patient use.

DEVICE SELECTION

To determine the correct size:

- Measure the circumference of the residual limb 4 and 30 cm from the distal end (**Fig. 1**)
- Choose the corresponding device size. If the measurement taken does not have a corresponding device size, round the measurement down (**Fig. 2**).

Note: Do not choose a device size which exceeds the measurement, i.e., do not round up.

Caution: A too small/tight device may result in pistoning, numbness, or blistering.

Caution: A too large/loose device may result in increased perspiration and movement of the limb inside the device, which may lead to blisters and rashes.

If any of the above symptoms, or any other indications of improper fit are experienced, the patient should contact the healthcare professional immediately.

FITTING INSTRUCTIONS

Caution: The device must be thoroughly inspected for defects before fitting.

Trimming Instructions

The device may be trimmed to increase range of motion, or for patient preference (**Fig. 3**). The Iceross Curvemaster is a recommended tool, specifically designed to produce a smooth, rounded edge when trimming.

Caution: Do not cut the device below the level of the matrix or below the proximal trim-lines of the socket. Excessive trimming may reduce the suction between the device and limb, and compromise suspension.

Socket Fitting

Achieving the intended function of the device is dependent on the fit of the socket. Appropriate socket design is essential to achieve suspension, control, and comfort.

The following checklist should be used when checking the socket fit:

- Patient can don the socket easily without the use of lubricant aids.
- Full contact between the liner and the socket inner wall.
- No air pockets develop when the patient loads the prosthesis.
- Minimized pistoning in the socket thus, no visible movement when pulling on the prosthesis.
- Full movement possible without compromising the suspension.

Caution: If feedback indicates any of the above conditions are not met, the socket will need to be adjusted or remade accordingly.

Caution: For a cushion liner, the socket must be airtight in order to maintain vacuum suspension when a sleeve is applied.

Custom Shape Capture

Due to the material properties of the device, the unique shape of the limb will be automatically captured during normal use. To accelerate shape capture, “Cold Capture” or “Hot Capture” techniques can be utilized. Both techniques require a cast/scan of the residual limb with no liner and a subsequent positive model.

Cold Capture

1. Add a protective layer to the clean positive model.
2. Apply device to positive model.
3. Leave for a minimum time of 60 minutes.

Hot Capture

Caution: Direct physical contact with heating elements will permanently damage the device.

Caution: Use heat-protective gloves at all times during this process.

1. Add a protective layer to the clean, dry positive model.
2. Place the positive model into a preheated oven at 65°C (150°F) for 60 minutes.
3. Apply the device to the positive model and place it into the oven.
4. Remove from the oven after a maximum of 10 minutes.
5. Allow to cool.

Note: Make sure the material of the positive model and protective layer can withstand 65°C (150°F).

Note: The device must not be put into the oven unless it is on a positive model.

Locking Liner

Warning: Failure to adhere to the following instructions may cause the pin to detach from the liner resulting in loss of suspension, which may lead to injury.

It is recommended that the device is used in conjunction with a compatible Össur lock.

Caution: Do not use a flangeless pin with the device (**Fig. 4**). A flangeless pin can push through the distal attachment.

Fixed distal attachment

If a lanyard suspension system is used: Turn the connecting screw into the distal attachment of the liner. Refer to the instructions for use for the lanyard suspension system for correct instructions.

If a locking pin is used: Apply Loctite 410/411 or an equivalent compound to the threads of the pin and tighten the attachment pin to the fixed distal attachment, torque to 4 Nm. Refer to the instructions for use for the lock used for correct instructions.

USAGE

Donning

Caution: Open wounds or damaged skin should be covered with a bandage or other appropriate covering to prevent direct contact between the wound and device.

Caution: Ensure the inside of the device is clean, dry, and free from any foreign objects that may cause skin irritation.

Caution: Do not apply lotion to the residual limb immediately before donning the device. The residual limb should be clean and dry to prevent device damage.

Liner Donning

1. Grip the top of the device from the inside and slide it over the hand until the inner surface is fully exposed (**Fig. 5**). Take care not to damage the device with fingernails.
2. After exposing as much of the distal end of the device as possible, position it against the distal end of the residual limb (**Fig. 6**).
Note: Be careful that the pin on the locking liner is positioned centrally on the distal end of the residual limb (**Fig. 8**).
3. Bend the knee slightly.
4. With light compression, roll the device all the way upwards onto the limb (**Fig. 7**).
Caution: Do not tug or pull the upper end of the liner when rolling the device upwards onto the limb. This can result in tension on the skin, causing blisters or rashes.
5. Check that no air pockets are present. If air pockets are found, reapply the device.
Note: The device must be orientated with the external labeling positioned at the front of the residual limb (**Fig. 9**).

Socket Donning

Cushion Liner: Use together with a sleeve. Refer to the instructions for use for the sleeve for correct use.

Step into the socket using body weight to achieve full engagement before using the sleeve to maintain vacuum suspension.

Locking Liner: Step into the socket using bodyweight to fully engage the distal attachment into the lock until a secure connection is achieved.

Warning: Ensure that the connection between the distal attachment and lock is secure and without play.

Cleaning and care

Daily cleaning of the residual limb is essential. Use of a mild, pH balance, 100% fragrance- and dye-free liquid soap is recommended. If the skin is dry, apply a pH balance, 100% fragrance- and dye-free lotion to nourish and soften the skin.

Clean the device daily after use and before first use.

- Turn the device inside out (**Fig. 5**).
- Wash with pH balanced, 100% fragrance- and dye-free liquid soap.
- After washing, rinse the device thoroughly with warm water (**Fig. 10**), and pat dry on both sides with a lint-free cloth.
- The device shall be dried overnight before reuse.
- Always return the device to its neutral state immediately after cleaning. Do not store inside-out as device damage may result.

Warning: Excessive perspiration may compromise stability and/or suspension. Remove the device and dry off both the device and the residual limb when appropriate.

Caution: Skin health should be monitored daily. In case of abnormal symptoms, discontinue use of the device and contact the healthcare professional immediately.

Caution: Be aware that common household or bath products may cause or contribute to skin irritation, e.g., soaps, deodorants, perfumes, abrasive cleaners, aerosol- or alcohol sprays.

Caution: After contact with salt and/or chlorinated water, wash the device with tap water.

Environmental Conditions

The device is safe to use with caution with waterproof components.

Caution: Frequent exposure to salt and/or chlorinated water might affect the durability of the device.

Only applicable for locking liner:

Warning: After submerging the device in water, it needs to be removed, and the residual limb and the skin contact layer dried before ambulation.

Warning: The suspension is compromised if water gets between the skin and the device.

SUPPLEMENTARY DEVICES

The following accessories can be used in conjunction with Össur cushion and locking liners.

Refer to Össur catalog for more information, e.g., sizing.

Iceross® Interface Sheath

The device is intended to be used for prosthetic socket donning.

The device can be worn over the liner to allow easier donning/doffing of the socket by reducing the friction between the socket and liner.

Pull the device over the liner. When used in conjunction with a locking liner ensure the metal pin and distal connector fully penetrates the knitted hole to ensure the proper function of the lock.

Machine wash at max 40°C with light-colored fabric. Do not tumble dry.

Iceross® Socks

The device is intended to be used for accommodating residual limb volume reduction.

If the patient experiences insecurity or decreased contact with the socket due to volume fluctuations, one or more sock can be pulled over the liner to restore socket fit.

Pull the device over the liner. When used in conjunction with a locking liner ensure the metal pin and distal connector fully penetrates the knitted hole to ensure the proper function of the lock.

Machine wash at max 40°C with light-colored fabric. Do not tumble dry.

REPORT OF SERIOUS INCIDENT

Any serious incident in relation to the device must be reported to the manufacturer and relevant authorities.

DISPOSAL

The device and packaging must be disposed of in accordance with respective local or national environmental regulations.

LIABILITY

Össur prosthetic devices are designed and verified to be safe and compatible in combination with each other and custom-made prosthetic sockets with Össur adapters, and when used in accordance with their intended use.

Össur does not assume liability for the following:

- Device not maintained as instructed by the instructions for use.
- Device assembled with components from other manufacturers.
- Device used outside of recommended use condition, application, or environment.

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