

Instructions for Use

ICELOCK® 544
















		3
EN	Instructions for Use	5
CS	Návod k použití	9
SK	Indikácie použitia	14
RO	Indicații de utilizare	19
HR	Upute za uporabu	24
HU	Használati javallat	29
BG	Указания за употреба	34
SL	Indikácie pre použitie	39

1

	A	B	C	D	E	F	G	H	I	J	K	L	M
L-544611	●	●			●		●				●		
L-544621	●	●			●		●		●				
L-544612	●		●		●			●		●			
L-544613	●	●			●			●			●		
L-544641	●	●			●		●					●	●
L-544620	●		●		●								
L-544640	●			●		●							

ENGLISH

MD Medical Device

PRODUCT DESCRIPTION

The Icelock® 544, referred to as the device in the following document, consists of a variety of component configurations to support the use of elevated and passive vacuum that help to provide and/or facilitate secure coupling of the residual limb to the prosthesis specifically in non-weight-bearing situations, i.e. swing phase during walking.

The different components can be configured as shown in **Figure 1** and the corresponding table. The names of the components are shown in the following table.

A	Socket Adapter
B	Seal
C	Duckbill valve
D	Mesh Filter
E	Expulsion Housing
F	E2 Connector
G	Large Seal
H	O-Ring
I	Icelock® 544 Expulsion Plate
J	Icelock® 544 EP Plate
K	Icelock® 544 Reservoir Plate / Icelock® 544 Unity Plate
L	Icelock® 544 LL Plate
M	LimbLogic Pump Seal
X ₁	Housing dummy
X ₂	Lamination dummy

INDICATIONS FOR USE

- Lower limb amputation and/or congenital deficiency
- No known contraindications

INTENDED USE

The device is intended to connect and release a prosthetic system that replaces a missing lower limb.

The device contains a series of locks systems that can be configured in different ways, making it compatible with many different suspension systems.

The device is designed to facilitate suspension of a prosthesis, specifically by:

- Creating differential air pressure in prosthetic socket;
- Permitting healthcare professional to locate and secure a locking mechanism into socket during fabrication.

Note: The valves are positioned in a housing that seals securely against the inner surface of the hard socket near the distal end. The valves expel air under positive pressure when entering a socket. For configurations L-544611, L-544621, L-544613 and L-544641 the push button needs to be pressed during removal of the prosthesis to allow air back into the hard socket.

Warning: Never use configurations with an expulsion valve (L-544612, L-544620, L-544640) or other valves without a functional release button when combined with Iceross Seal-In® or other solutions where liner is sealed against the inside wall of a socket. Some means of allowing air back in the socket must be in place.

Note: The device must be fitted and adjusted by a healthcare professional.

Note: Intended users of the device are healthcare professionals and patients.

GENERAL SAFETY INFORMATION

Caution: The device is for single patient use.

Warning: In case of functional change or functional loss, the patient should stop using the device and contact a healthcare professional.

Warning: If the device shows signs of damage or wear hindering its normal functions, the patient should stop using the device and contact a healthcare professional.

Note: This device has been tested according to ISO 10328 standard to three million load cycles. Depending on patient activity, this may correspond to 3–5 years of use.


Warning: If combined with products with different maximum user weight, the maximum loading corresponds to that of the weakest product.

Caution: Össur products and components are designed and tested according to the applicable official standards or an in-house defined standard when no official standard applies. Compatibility and compliance with these standards are achieved only when Össur products are used with other Össur components.

PRODUCT SELECTION

The device is for moderate impact use, e.g. walking.

The weight limit for the device is 136 kg.

ASSEMBLY

Icelock L-544 Series can either be set up on a thermoplastic socket for temporary test fit purposes or laminated directly into a socket.

Thermoplastic Socket Procedure

1. Drape a thermoplastic socket with conventional methods. Roughen distal area.
Note: Distal area must be clean of dirt, dust, and grease, as otherwise the bonding connection may not work.
2. Using a designated 2-component PU adhesive putty, bond the adapter (A) assembly (**Fig. 2**) to desired location on the socket (**Fig. 3**).
3. Wrap thermoplastic socket and adapter firmly with glass reinforced casting tape, enabling test fitting of the prosthesis (**Fig. 4**).

4. Drill through housing dummy into cast using a sharp 3 mm drill bit. Take care to align drill with centerline (**Fig. 5**).
5. Remove housing dummy (X_1) and replace with seal (B) / duckbill valve (C) / mesh filter (D) and then screw in the expulsion housing (E) / E2 connector (F) (**Fig. 6**).
Note: When inserting the housing/connector, make sure that the seal/valve/mesh is not crooked, otherwise the housing may not hold vacuum properly.
6. Place large seal (G) / O-ring (H) below housing. Complete assembly with the plate (I/J/K/L), applying medium strength thread locker to 4-hole screws and torque twice to 10 Nm (**Fig. 7**).
Note: Make sure threads are clean and provide full engagement of screws without resistance. Only use screws supplied with kit.
7. [For L-544611, L-544612 and L-544613] Secure Unity tube to the connector on the plate (J/K).

Laminated Socket Procedure

1. Drape a layer of ~3 mm PETG or equivalent rigid thermoplastic sheet over cast positive (**Fig. 8**). Roughen distal area.
Note: Distal area must be clean of dirt, dust, and grease, as otherwise the bonding connection may not work.
2. Screw housing dummy (X_1) into socket adapter (A) (**Fig. 9**). Using a designated 2-component PU adhesive putty, bond the adapter assembly to desired location directly on thermoplastic layer (**Fig. 10**). Take care to fill the center area for a bubble free bond between housing dummy and thermoplastic layer.
3. Remove housing dummy (X_1). Place lamination dummy (X_2) in adapter (**Fig. 11**).
4. Tie reinforcement material in groove on socket adapter and do a second tie between socket and adapter. Apply the reinforcement material layup to agree with the expected stress (**Fig. 12**).
5. Apply PVA bag and impregnate laminate with conventional methods (**Fig. 13**).
6. Tie around lamination dummy with electrical tape to push out the resin between the PVA bag and dummy. This will produce clean edge and minimize finishing work (**Fig. 13**).
7. Once resin has set, drill through lamination dummy (X_2) into cast using a sharp 3 mm drill bit. Take care to align drill with centerline (**Fig. 14**).
8. Remove excess resin, PVA bag and lamination dummy (X_2). Remove plaster from inside of socket.
9. Place seal (B) / duckbill valve (C) / mesh filter (D) and then screw in the expulsion housing (E) / E2 connector (F) (**Fig. 15**).
Note: When inserting the housing/connector, make sure that the seal/valve/mesh is not crooked, otherwise the housing may not hold vacuum properly.
10. Place large seal (G) / O-ring (H) below housing. Complete assembly with the plate (I/J/K/L), applying medium strength thread locker to 4-hole screws and torque twice to 10 Nm (**Fig. 16**).
Note: Make sure threads are clean and provide full engagement of screws without resistance. Only use screws supplied with kit.
11. [For L-544611, L-544612 and L-544613] Secure Unity tube to the connector on the plate (J/K).

USAGE

Icelock 544 series components that are equipped with a release button (I/J/K/L) require the patient to press the release button when doffing the prosthesis in case of an internal seal, like Iceross Seal-In products.

Note: Pressing the release button upon donning will also facilitate faster and more complete engagement and is strongly recommended when using Iceross Seal-In liners. Capability to follow this process is essential.

CLEANING

Disassembly and cleaning

1. Unscrew the plate (I/J/K/L) and expulsion housing (E) / E2 connector (F) and remove seal (B) / duckbill valve (C) / mesh filter (D).
2. Clean under running water with mild detergent and clear airways with an air hose. Hold down the release button while clearing airway.

Reassembly

1. Place seal (B) / duckbill valve (C) / mesh filter (D) and then screw in the expulsion housing (E) / E2 connector (F) (**Fig. 15**).

Note: When inserting the housing, make sure that the seal/valve/mesh is not crooked, otherwise the housing may not hold vacuum properly.

2. Place large seal (G) / O-ring (H) below housing. Complete assembly with the plate (I/J/K/L), applying medium strength thread locker to 4-hole screws and torque twice to 10 Nm (**Fig. 16**).

Note: Make sure threads are clean and provide full engagement of screws without resistance. Use screws supplied with kit only.

3. [For L-544611, L-544612 and L-544613] Secure Unity tube to the connector on the plate (J/K).

ENVIRONMENTAL CONDITIONS

The device can be used in temperatures between -15 °C (5 °F) to 50 °C (122 °F).

The device is Weatherproof.

Caution: A weatherproof device allows use in a wet and/or humid environment but does not allow submersion. Fresh water splashing against the enclosure from any direction shall have no harmful effect. The device should be dried thoroughly after contact with fresh water. Fresh water: Includes tap water. Excludes salt and chlorinated water.

Note: For a weatherproof system, weatherproof components must be selected for use with the device.

Warning: Do not submerge a weatherproof product in water. This can cause corrosion. A corroded product can break.

Caution: Rinse the product with fresh water and then dry thoroughly after use in wet and/or humid environment.

REPORT A SERIOUS INCIDENT

Any serious incident in relation to the device must be reported to the manufacturer and relevant authorities.

DISPOSAL

All components of the product and packaging should be disposed of in accordance with respective national environmental regulations.

LIABILITY

Össur does not assume liability for the following:

- Device not maintained as instructed by the instructions for use.
- Device assembled with components from other manufacturers.
- Device used outside of recommended use condition, application, or environment.

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