







Instructions for Use







ICELOCK 562 HYBRID











L-562000, L-562400, L-562500, L-562410

		3
EN	Instructions for Use	6
CS	Návod k použití	16
SK	Indikácie použitia	26
RO	Indicații de utilizare	36
HR	Upute za uporabu	46
HU	Használati javallat	56
BG	Указания за употреба	66
SL	Indikácie pre použitie	76

ENGLISH

Medical Device

DESCRIPTION

The device uses two different locking mechanisms to provide secure suspension, mechanical and vacuum locking. It is offered in two versions, with either a passive or active vacuum suspension.

Both locking mechanisms are released by pressing a single button, making it easy to doff a prosthesis.

Components

- A. Central Screw
- B. Fixation Dummy
- C. Distal Connector
- D. Nails
- E. Manufacturing Dummy 1
- F. Manufacturing Dummy 2
- G. Aluminum Dummy
- H. Spacers
- I. Flat-head Screws
- J. Attachment Pin for Icelock 562
- K. Lock Body
- L. Large Square Ring
- M. Large O-ring
- N. Small Square Ring

INTENDED USE

The device is intended to connect and release a prosthetic system that replaces a missing lower limb.

Suitability of the device for the prosthesis and the patient must be evaluated by a healthcare professional.

The device must be fitted and adjusted by a healthcare professional.

Indications for Use and Target Patient Population

- Lower limb amputation and/or congenital deficiency
- None known contraindications

The device is for high impact use, e.g., walking and occasional running.

The weight limit for the device is 166 kg.

GENERAL SAFETY INFORMATION

The healthcare professional should inform the patient about everything in this document that is required for safe use of this device.

Warning: If there is a change or loss in device functionality, or if the device shows signs of damage or wear hindering its normal functions, the patient should stop using the device and contact a healthcare professional.

The device is for single patient use.

REQUIRED COMPONENTS

Use in combination with an Iceross® Seal-In X Locking liner, Iceross® Seal-In X Seal and airtight socket. Please refer to the Instructions for Use supplied with Iceross® Seal-In X Locking liner and Iceross® Seal-In X Seals.

FITTING INSTRUCTIONS

Applicable to all socket fabrications:

Caution: Make sure no debris contaminates the lock body during grinding (**Fig. II**). The lock body will not operate if it is contaminated with debris.

Socket Construction

Prepare plaster negative as usual.

Plaster Positive

1. Model the plaster positive to agree with the residual limb. Mold the distal end of the liner as precisely as possible.
 2. Flatten the distal end with a plaster rasp until the diameter is the same size as the Manufacturing Dummy 1+2 (**Fig. 1**). For liner size 22 and above, use the two Manufacturing Dummies. For liner sizes smaller than size 22, only use the Manufacturing Dummy 1.
 3. Use one of the two following methods to ensure air to flow to the distal connector at the distal end:
 - Drill a hole through the whole model from the distal center to the proximal end before fixing the Manufacturing Dummy.
- OR
- Apply a thin stocking to the plaster.
4. Use the Nails supplied to attach the correct Manufacturing Dummy/Dummies to the plaster model (**Fig. 2**).

Note: If using a moist plaster model, isolate properly.

THERMOPLASTIC SOCKETS

1. Put the Small square ring onto the cut-out on the Distal Connector. Use the Central Screw and the Fixation Dummy to secure the Distal Connector on the Manufacturing Dummy/Dummies and the plaster model (**Fig. 3**).

Note:

- Orientate the distal connector to the line of progression.
 - For liners below size 22, replace the Large Square Ring with the Large O-ring before you drape the material (**Fig. 10a**).
2. Select the applicable material depending on the expected stresses and drape as usual (**Fig. 4**).
 3. Tightly wind a strap into the groove of the Distal Connector during draping. This will help the socket become airtight.
 4. After fully hardening, grind and remove the material at the distal end until completely exposing the head of the Central Screw (**Fig. 5**).
 5. Remove the Central Screw.
 6. Continue grinding, until fully exposing the Fixation Dummy.
 7. Remove the Fixation Dummy and the Small Square Ring.
 8. Grind carefully until the material is flush with the surface of the Distal Connector (**Fig. 6**).
 9. Put the Aluminum Dummy on the socket (**Fig. 7a**).
 10. Align the four holes of the Aluminum Dummy with the four holes on the Distal connector.
 11. Secure the Aluminum Dummy with the Central screw.
 12. Drill \varnothing 6,5 mm holes through the four holes of the Aluminum Dummy, until you reach the heads of the Set Screws.
 13. Remove the Set Screws (**Fig. 8a**).
 14. Remove the socket from the plaster.
 15. Carefully push the Distal Connector out of the socket (**Fig. 9a**).
 16. Replace the Large Square Ring with the Large O-ring (**Fig. 10a**).
 17. Use a \varnothing 10 mm drill to drill again into the \varnothing 6,5 mm holes in the PETG (**Fig. 11a**).

Note: Do not damage the cylindrical walls on the inside of the socket, then the sealing won't work.

18. Remove all dust, dirt or similar contamination from the socket and the Distal Connector.
19. Put the four Spacers into the \varnothing 10 mm holes (**Fig. 12a**).
20. Assemble the Distal Connector and Lock Body.
21. Apply medium strength threadlocker to the threads of the Flat-head Screws to fasten them.
22. Torque the Flat-head Screws to 10 Nm.

Note: Make sure the four holes stay aligned.

Caution:

- Make sure not to manipulate the connector area when adjusting the socket with heat during the fitting process. Fix this area with rigid tape before heating.
- Do not let threadlocker or resin enter the lock body. The lock body will not operate if contaminated.

LAMINATED SOCKETS

Note: To make sure the definitive socket is airtight, use PETG for a thin inner layer.

1. Do the steps given in the “*Preparation*” section.
2. To drape, do steps 1 through 8 given in the “*Thermoplastic Sockets*” section
3. Put the Small Square Ring back on again.
Note: It prevents resin leaks towards the distal connector.
4. Align the four holes of the Aluminum Dummy with the four holes on the Distal Connector.
5. Secure the Aluminum Dummy with the Central Screw (**Fig. 7b**).

6. The thin PETG replaces the lower PVA-bag. Seal the transition between the PETG and plaster at the proximal end.
7. Make the outer surface rough.
8. Fill the head of the Central Screw with plasticine or similar.
9. Attach reinforcement material between the Distal Connector and Aluminum Dummy (**Fig. 8b**).
10. Drape unidirectional carbon tape between the Distal Connector and Aluminum Dummy.
11. Apply the reinforcement material layup to agree with the expected stress.
12. Put the upper PVA bag above the material and Aluminum Dummy (**Fig. 9b**).
13. Before the resin gels, fasten a cord around the PVA bag, distally to the Aluminum Dummy.
14. Drill four \varnothing 6,5 mm holes through the Aluminum Dummy until you reach the heads of the Set Screws (**Fig. 10b**).
15. Remove the Central Screw, Aluminum Dummy and the four Set Screws (**Fig. 11b**).
16. Remove the socket from the plaster.
17. Carefully push the Distal Connector out of the socket (**Fig. 12b**).
18. Exchange the Large Square Ring with the Large O-ring (**Fig. 13b**).
19. Remove all dust, dirt or similar contamination from the socket and the distal connector.
20. Assemble the Distal Connector and the Lock Body (**Fig. 14b**).
Note: Make sure the four holes stay aligned.
21. Apply medium strength threadlocker to the threads of the Flat-head Screws to fasten them.
22. Torque the Flat-head Screws to 10 Nm.

USE WITH THE DIRECT SOCKET SYSTEM

Icelock® 562 Hybrid (L-562400) and (L-562410) are compatible with the Direct Socket system.

Refer to Instructions for use on Direct Socket to make the socket.

1. Mount the Lock Body with the M6X14 screws supplied in the Lock kit for Direct Socket. It must be fully engaged with the lamination ring.
2. Fasten the Flat-head Screws with medium strength threadlocker and torque to 10 Nm.

Caution: Do not let Loctite or resin enter the lock body. The lock body will not operate if contaminated with threadlocker or resin.

ASSEMBLY INSTRUCTIONS

Retrofitting the Lock Body between vacuum suspension, the replacement parts need to be bought separately.

- 1a. If changing to a passive vacuum system, replace the 90° Unity® tube connector an outlet port (**Fig. I**).
- 1b. If changing to the Unity® vacuum system, remove one of the outlet ports and replace with the 90° Unity® tube connector.
2. Position the release button on the lateral/medial side, as required by the amputee.
3. Shorten the plastic release button as necessary.

Attachment Pin

To attach and secure the Attachment Pin for Icelock® 562 to the Iceross® Seal-In X Locking liner.

1. Apply Loctite 410/411 or an equivalent to the threads of the Attachment Pin.
2. Turn the Attachment Pin into the liner umbrella.
3. Torque the Attachment Pin to **4 Nm**.

USAGE

Cleaning and care

Clean with damp cloth and a mild soap. Dry with a cloth after cleaning

Environmental Conditions

The device is Waterproof.

A Waterproof device can be used in a wet or humid environment and submerged in up to 1-meter-deep fresh water for a maximum of 30 minutes. No contact with salt water or chlorinated water is allowed.

Dry with a cloth after contact with fresh water or humidity. Clean with fresh water in case of accidental exposure to other liquids, chemicals, sand, dust, or dirt and dry with a cloth.

MAINTENANCE

The device and the overall prosthesis should be examined by a healthcare professional. Interval should be determined based on patient activity.

Note: Do not disassemble the 90° Unity® tube connector in its components, as this may result in damage of internal sealings and thus loss of elevated vacuum.

Caution: The locking pin should be exchanged at least once a year. Check the release button for ease of use. If the release button is hard to use, the mechanism is damaged. Exchange the device.

Check if the device is airtight. If the device is not airtight, clean it as follows:

1. Examine the pin hole. If necessary, remove any dust, dirt, or similar contamination.
2. Trap the pin hole with e.g., tape or plasticine.
3. Use low air pressure to clean the socket.

If noise or leakage occurs, do as follows:

1. If using an active vacuum suspension, use a 3 mm hex key to remove the 90° Unity® tube connector.
2. Remove the outlet port with an 8 mm wrench.
3. Use low air pressure to clean the air channels.
4. Clean or exchange the duckbill valves and assemble again correctly (Fig. I).

REPORT A SERIOUS INCIDENT

Any serious incident in relation to the device must be reported to the manufacturer and relevant authorities.

DISPOSAL

All components of the product and packaging should be disposed of in accordance with respective national environmental regulations.

LIABILITY

Össur does not assume liability for the following:

- Device not maintained as instructed by the instructions for use.
- Device assembled with components from other manufacturers.
- Device used outside of recommended use condition, application, or environment.

Compliance

This device has been tested according to ISO 10328 standard to three million load cycles

Depending on patient activity, this may correspond to 3-5 years of use.

ISO 10328 - P7 - 166 kg *)

*) Body mass limit not to be exceeded!

For specific conditions and limitations of use
see manufacturer's written instructions on
intended use!

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