


Instructions for Use






---

ICEROSS SEAL-IN<sup>®</sup> X SEALS

		3
EN	Instructions for Use	5
DE	Gebrauchsanweisung	9
FR	Notice d'utilisation	13
ES	Instrucciones para el uso	17
IT	Istruzioni per l'uso	21
NO	Bruksanvisning	25
DA	Brugsanvisning	29
SV	Bruksanvisning	33
FI	Käyttöohjeet	37
EL	Οδηγίες Χρήσης	41
NL	Gebruiksaanwijzing	45
PT	Instruções de Utilização	49
PL	Instrukcja użytkowania	53
TR	Kullanım Talimatları	57
RU	Инструкция по использованию	61
JA	取扱説明書	65
ZH	中文说明书	69
KO	사용 설명서	72

### DESCRIPTION

The device consists of a silicone seal ring and a textile donning aid. The anti-friction coating makes it easy to don and doff a prosthetic socket without any lubricant aids such as alcohol spray. When the device is correctly positioned onto a donned liner from the **Iceross Seal-In X** family, it will be held in tension against the silicone bands on the exterior surfaces of the liner and is then compressed against the interior surface of the prosthetic socket. This creates a distal vacuum chamber that reliably suspends the residual limb in the prosthetic leg.

The device is offered in 3 variants:

- Seal-In X-Grip
- Seal-In X-Volume
- Seal-In X-Classic

### INTENDED USE

The device is a prosthetic interface with suspension properties intended to be used as part of a system that replaces a missing lower limb.

Suitability of the device for the prosthesis and the patient must be evaluated by a healthcare professional.

The device must be fitted and adjusted by a healthcare professional.

### *Indications For Use and Target Patient Population*

- Lower limb loss, amputation, or deficiency
- No known contraindications

For a transtibial device: The device is for low to high impact use, e.g., walking and occasional running.

For a transfemoral device: The device is for low to extreme impact use, e.g., walking, running, and sports.

### GENERAL SAFETY INSTRUCTIONS

**Warning:** Use of a lower limb prosthetic device carries an inherent risk of falling which may lead to injury.

The healthcare professional should inform the patient about everything in this document that is required for safe use of this device.

**Warning:** If there is a change or loss in device functionality, or if the device shows signs of damage or wear hindering its normal functions, the patient should stop using the device and contact a healthcare professional.

**Warning:** Misuse of the device may result in loss of suspension.

**Warning:** The device shall not be used by patients with drastic volume changes within a short period of time.

- Grip: Volume changes exceeding 3 ply.
- Classic: Volume changes exceeding 6 ply.
- Volume: Volume changes exceeding 10 ply.

**Caution:** The device shall only be used in combination with a suitable liner from the Iceross Seal-In® X family.

The device is for single patient use.

## DEVICE SELECTION

**Caution:** A too small/tight device can lead to excessive pressure on the residual limb which can affect residual limb health.

**Note:** A too small/tight device can lead to overstretching and permanent deformation of the device, affecting device performance and/or cause premature wear and tear.

**Note:** A too big/loose device can lead to compromised suspension, as the connection between the seal wall and silicone bands may break, and/or a crease may form when donning the socket.

One of the following methods can be used to determine the correct device size.

### Method 1: Seal Ring Kits

The preferred method to determine the correct device size is by using an **Iceross Seal-In X Kit** (please refer to catalog for further details)

1. Place a device from the kit at the desired location on top of the liner, taking into account the recommended distance between the lowest socket edge and seal lip (**Fig. 1**).
2. Assess deformation of the device and residual limb. The device should sit securely in place without excessive compression (**Fig. 2**).
3. If the device does not fit at the desired location, select another size, or relocate to obtain proper tightness.

### Method 2: Measurements

Legend for **Fig. 3**:

- X: Seal size
- Y: Circumference (cm)
- A: Tight fit
- B: Loose fit

Measure the circumference (in cm) at the desired location directly on the skin. Refer to the graph in **Fig. 3** to estimate the correct device size.

Examples:

- Circumference of 40 cm results in seal ring sizes 44 or 47.
- Circumference of 26 cm results in seal ring sizes 28 or 30.

**Warning:** The device shall not be used by patients with an extremely short residual limb. Refer to the instructions for use for the suitable liner of the Iceross Seal-In X family for the minimum residual limb length.

## FITTING INSTRUCTIONS

**Caution:** The device must be thoroughly inspected for defects before fitting.

Achieving the intended function of the device is dependent on the fit of the socket. Appropriate socket design is essential to achieve suspension, control, and comfort. Refer to the instructions for use of the suitable liner of the Iceross Seal-In X family for socket construction guidelines and fitting.

### Iceross Seal-In® Socks

If the patient experiences compromised suspension caused by residual limb volume fluctuations, one or more sock can be pulled over the liner to restore socket fit.

The sock shall be placed below the seal lip (**Fig. 1A and 8**).

**Note:** This only applies to the variants Classic and Volume.

**Caution:** The sock shall never be placed between the liner and the seal (**Fig. 8**).

**Caution:** If volume fluctuation cannot be adequately compensated for by this means, the socket will need to be remade to suitably reduced dimensions.

## USAGE

### *Donning*

**Caution:** The device must not be worn directly on an exposed residual limb.

1. The device must be clean and free from dirt before use. If needed, rinse with water and wipe dry with a lint-free cloth.
2. Use the fabric end of the device to pull it on and adjust to the desired location (**Fig. 4 and 5**).

**Note:** The uppermost sealing lip of the seal ring should not be closer to the socket brim than shown in **Fig. 1**.

**Caution:** The silicone portion of the device must not extend proximal to the socket brim (**Fig. 6**)

**Caution:** The device must be aligned with the silicone bands on the liner to achieve proper seal (**Fig. 6**).

### **Socket Donning**

With the liner and seal ring correctly configured, step into the socket, expelling all air until fully engaged into the distal end (**Fig. 7**).

**Caution:** Do not use creams or lotions to aid socket donning, as this will cause the device to soften and expand over time, affecting its ability to be tensioned on the liner.

**Note:** Do not use alcohol spray as a donning aid for the device as this can destroy the low-friction coating, making socket donning difficult.

### *Cleaning and care*

**Caution:** Skin health should be monitored daily. In case of abnormal symptoms, discontinue use of the device and contact the healthcare professional immediately.

**Caution:** After contact with salt and/or chlorinated water, wash the device with tap water.

Clean the silicone portion of the device daily using a damp cloth or rinse with water.

**Note:** Do not use soap or any other cleaning aids with the device.

Inspect the device daily for holes and before first use. The presence of any holes will affect the seal performance.

Store the device with the low-friction coating facing out. Do not store inside out.

### *Environmental Conditions*

The device is safe to use with caution with waterproof components.

**Warning:** After submerging the device in water, it needs to be removed, and the residual limb and the liner skin contact layer and socket are dried before ambulation.

**Warning:** The suspension is compromised if water gets between the skin and the device or the seal and socket.

**Caution:** Frequent exposure to salt and/or chlorinated water might affect the durability of the device.

## REPORT OF SERIOUS INCIDENT

Any serious incident in relation to the device must be reported to the manufacturer and relevant authorities.

## **DISPOSAL**

The device and packaging must be disposed of in accordance with respective local or national environmental regulations.

## **LIABILITY**

Össur prosthetic devices are designed and verified to be safe and compatible in combination with each other and custom-made prosthetic sockets with Össur adapters, and when used in accordance with their intended use.

Össur does not assume liability for the following:

- Device not maintained as instructed by the instructions for use.
- Device assembled with components from other manufacturers.
- Device used outside of recommended use condition, application, or environment.

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