








Instructions for Use





---

ICERROSS DERM O SEAL-IN<sup>®</sup>  
ICERROSS SEAL-IN<sup>®</sup>  
TRANSFEMORAL

		3
EN	Instructions for Use	5
DE	Gebrauchsanweisung	11
FR	Notice d'utilisation	18
ES	Instrucciones para el uso	24
IT	Istruzioni per l'uso	30
NO	Bruksanvisning	36
DA	Brugsanvisning	41
SV	Bruksanvisning	47
EL	Οδηγίες Χρήσης	53
FI	Käyttöohjeet	59
NL	Gebruiksaanwijzing	65
PT	Instruções de Utilização	71
PL	Instrukcja użytkownika	77
TR	Kullanım Talimatları	83
RU	Инструкция по использованию	89
JA	取扱説明書	95
ZH	中文说明书	100
KO	사용 설명서	105

### DESCRIPTION

The device is a silicone liner that covers the residual limb and provides an interface between the residual limb and the socket of the prosthesis. The device forms a suspension system where a distal vacuum chamber is created below the level of the seal that reliably suspends the residual limb in the prosthesis.

### INTENDED USE

The device is a prosthetic interface with suspension properties intended to be used as part of a system that replaces a missing lower limb. Suitability of the device for the prosthesis and the patient must be evaluated by a healthcare professional.

The device must be fitted and adjusted by a healthcare professional.

### *Indications For Use and Target Patient Population*

- Lower limb loss, amputation, or deficiency
- No known contraindications

Iceross Dermo Seal-In is a transtibial liner.

The device is for low to high impact use, e.g., walking and occasional running.

Iceross Transfemoral Seal-In is a transfemoral liner.

The device is for low to extreme impact use, e.g., walking, running, and sports.

### GENERAL SAFETY INSTRUCTIONS

**Warning:** Use of a lower limb prosthetic device carries an inherent risk of falling which may lead to injury.

The healthcare professional should inform the patient about everything in this document that is required for safe use of this device.

**Warning:** If there is a change or loss in device functionality, or if the device shows signs of damage or wear hindering its normal functions, the patient should stop using the device and contact a healthcare professional.

**Warning:** Misuse of the device may result in loss of suspension.

**Warning:** The device shall not be used by patients with drastic volume changes within a short period of time, exceeding 6 ply.

**Caution:** Potential side-effects are skin irritation/reaction problems.

**Caution:** The device must not come in contact with glass, carbon fibers, or other foreign particles. Washing alone may not be sufficient to eliminate the problem. If the device is exposed to foreign substances or chemicals, it should be returned to the healthcare professional.

The device is for single patient use.

### DEVICE SELECTION

For a transtibial liner, the correct size is determined by bending the knee slightly, allowing the soft tissue to hang down.

To determine the correct size:

- Measure the circumference of the residual limb 4 cm from the distal end (**Fig. 1**)
- Choose the corresponding device size. If the measurement taken does not have a corresponding device size, round the measurement down.

**Note:** Do not choose a device size which exceeds the measurement, i.e., do not round up.

**Caution:** A too small/tight device may result in pistoning, numbness or blistering, or produce excessive ring tension to the residual limb from the seal.

**Caution:** A too large/loose device may result in increased perspiration and movement of the limb inside the device, which may lead to blisters and rashes.

If any of the above symptoms, or any other indications of improper fit are experienced, the patient should contact the healthcare professional immediately.

**Note:** For a transtibial liner, make sure that the liner matrix covers at least 4 cm of bony structure. Do not allow the matrix to extend higher than the fibula head.

**Note:** A slight tension must be around the proximal edge of the seal to ensure optimal suspension. If the seal is loose, a crease may form when donning the socket. This may be prevented by selecting the next size below.

**Warning:** The device shall not be used by patients with an extremely short residual limb (refer to the Table below).

Transtibial Liner		Transfemoral Liner	
Device Size	Minimum Residual Limb Length (cm)	Device Size	Minimum Residual Limb Length (cm)
16-23.5	11	23.5-28	16
25-30	12	30-34	17
32-36	13	36-40	18
40	14	45	19
		50	20
		55	21

Residual limb length is measured as follows:

- Transtibial: From mid-patella tendon to distal end.
- Transfemoral: From ischial tuberosity to distal end, with the tissue hanging down.

## FITTING INSTRUCTIONS

**Caution:** The device must be thoroughly inspected for defects before fitting.

### *Trimming Instructions*

The device may be trimmed to increase range of motion, or for patient preference (**Fig. 2**). The Icross Curvemaster is a recommended tool, specifically designed to produce a smooth, rounded edge when trimming.

**Caution:** Do not cut the device below the level of the matrix, through the seal(s), or below the proximal trim-lines of the socket. Excessive trimming may reduce the suction between the device and limb, and compromise suspension.

### *Socket Construction Guidelines*

Recommended socket materials to ensure long-term seal function:

- All PETG materials
- Laminated sockets
- Thermoplastic inner socket materials with a smooth surface

The inner surface of the socket must be airtight and as smooth as possible, irrespective of the socket material.

To achieve the smoothest surface possible on the plaster positive model, it is recommended to use fine abrasive for final smoothing, followed by silicone spray.

Use of soft/smooth socket lining material will prolong the device sealing function.

A one-way valve with a release function such as Icelock 500 or 544 series is essential for the proper function of seal suspension. The socket valve must be positioned as distally as possible, at least below the seal.

The device is Unity® compatible.

### **Socket Fitting**

Achieving the intended function of the device is dependent on the fit of the socket. Appropriate socket design is essential to achieve suspension, control, and comfort.

The following checklist should be used when checking the socket fit:

- Full contact between the seal lips and the socket inner wall.
- No air pockets develop when the patient loads the prosthesis.
- No pistoning in the socket
- No visible movement when pulling on the prosthesis.
- Full movement possible without compromising the suspension.

**Caution:** If feedback indicates any of the above conditions are not met, the socket will need to be adjusted or remade accordingly.

## **USAGE**

### **Donning**

**Caution:** Do not apply lotion to the residual limb immediately before donning the device. The residual limb should be clean and dry to prevent device damage.

**Caution:** Open wounds or damaged skin should be covered with a bandage or other appropriate covering to prevent direct contact between the wound and device.

**Caution:** Ensure the inside of the device is clean, dry, and free from any foreign objects that may cause skin irritation.

### **Liner Donning**

1. Grip the top of the device from the inside and slide it over the hand until the inner surface is fully exposed (**Fig. 3**). Take care not to damage the device with fingernails.
2. After exposing as much of the distal end of the device as possible, position it against the distal end of the residual limb (**Fig. 4**).
3. With light compression, roll the device all the way upwards onto the limb (**Fig. 5**).

**Caution:** Do not tug or pull the upper end of the liner when rolling the device upwards onto the limb. This can result in tension on the skin, causing blisters or rashes.

4. Check that no air pockets are present. If air pockets are found, reapply the device.
5. Ensure the seal sits correctly and horizontally around its full circumference (**Fig. 6**).
6. Ensure the seal is correctly positioned by running your fingers around the seal and assess the tension.

## Socket Donning

1. Lightly spray the seal around its full circumference and inside the socket with Iceross® Clean & Simple lubricant spray or similar for smooth entry (**Fig. 7**).

**Warning:** Do not apply the lubricant spray to the skin. Avoid contact with eyes.

2. Align the prosthesis and push into the socket expelling air (**Fig. 8**).

**Note:** If difficulties are experienced while donning the prosthesis, more lubricant may be needed, or the valve is blocked.

**Note:** If stretching of the skin is experienced to be more to one side, the lubrication might not be evenly distributed around the full circumference of the seal.

## Cleaning and care

Daily cleaning of the residual limb is essential. Use of a mild, pH balance, 100% fragrance- and dye-free liquid soap is recommended. If the skin is dry, apply a pH balance, 100% fragrance- and dye-free lotion to nourish and soften the skin.

Clean the device daily after use and before first use.

- Turn the device inside out (**Fig. 3**).
- Wash with pH balanced, 100% fragrance- and dye-free liquid soap.
- After washing, rinse the device thoroughly with warm water (**Fig. 9**), and pat dry on both sides with a lint-free cloth.
- The device can also be machine washed (40°C) with a mild detergent and at low spin speed. Fabric softeners, bleaches, and other products/cleaning solutions may result in device damage and should not be used.
- The device shall be dried overnight before reuse.
- Always return the device to its neutral state immediately after cleaning. Do not store inside-out as device damage may result.
- The device should not be exposed to excessive heat or sunlight.
- The seal must be clean and free from dirt before use.

**Warning:** Excessive perspiration may compromise stability and/or suspension. Remove the device and dry off both the device and the residual limb when appropriate.

**Caution:** Skin health should be monitored daily. In case of abnormal symptoms, discontinue use of the device and contact the healthcare professional immediately.

**Caution:** Be aware that common household or bath products may cause or contribute to skin irritation, e.g., soaps, deodorants, perfumes, abrasive cleaners, aerosol- or alcohol sprays.

**Caution:** After contact with salt and/or chlorinated water, wash the device with tap water.

## Environmental Conditions

The device is safe to use with caution with waterproof components.

**Warning:** After submerging the device in water, it needs to be removed, and the residual limb and the liner skin contact layer and socket are dried before ambulation.

**Warning:** The suspension is compromised if water gets between the skin and the device or the seal and socket.

**Caution:** Frequent exposure to salt and/or chlorinated water might affect the durability of the device.

## SUPPLEMENTARY DEVICES

The following accessories can be used in conjunction with Össur seal liners.

Refer to Össur catalog for more information, e.g., sizing.

### Iceross® Distal Cup

The device is a prosthetic interface intended to be used as part of a system that replaces a missing lower limb.

The device provides further cushioning for distally conical limbs with poor soft tissue coverage. It conforms to the residual limb and can therefore be used as a filling element to ensure full contact between the silicone liner and the residual limb.

To correctly fit the device:

1. Measure the circumference of the residual limb 4 cm above the distal end (**Fig. 1**).
2. Choose a device one or two sizes smaller than the measurement.
3. Roll the device directly over the limb.
4. When using a liner in combination with a Distal Cup, the liner measurement should be made with the Distal Cup in place (i.e., over the Distal Cup).
5. Choose the correct liner size based on the measurement obtained with the Distal Cup in place.

**Note:** For daily care, refer to Cleaning and Care instructions.

### Iceross® Pads

The device is a prosthetic interface intended to be used as part of a system that replaces a missing lower limb.

The device can be used for prosthetic volume management to temporarily accommodate lost distal end volume as well as to assess distal contact while fitting a check socket.

Place the pad either in the distal end of the liner or in the distal end of the socket.

**Note:** For daily care, refer to Cleaning and Care instructions.

### Iceross Seal-In® Socks

The device is intended to be used for accommodating residual limb volume reduction.

If the patient experiences compromised suspension caused by residual limb volume fluctuations, one or more sock can be pulled over the liner to restore socket fit.

**Caution:** If volume fluctuation cannot be adequately compensated for by this means, the socket will need to be remade to suitably reduced dimensions.

- Invert the seal completely in the distal direction (**Fig. 10**).
- Pull the device over the inverted seal until the distal end of the sock is level with the proximal border of the seal (**Fig. 11**).
- Reflect the seal over the end of the sock (**Fig. 12**).
- Don the prosthesis as normal.

**Caution:** Loss of suspension may be encountered if a sock is placed over the seal.

Machine wash at max 40°C with light-colored fabric. Do not tumble dry.

## REPORT OF SERIOUS INCIDENT

Any serious incident in relation to the device must be reported to the manufacturer and relevant authorities.

## **DISPOSAL**

The device and packaging must be disposed of in accordance with respective local or national environmental regulations.

## **LIABILITY**

Össur prosthetic devices are designed and verified to be safe and compatible in combination with each other and custom-made prosthetic sockets with Össur adapters, and when used in accordance with their intended use.

Össur does not assume liability for the following:

- Device not maintained as instructed by the instructions for use.
- Device assembled with components from other manufacturers.
- Device used outside of recommended use condition, application, or environment.

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