






Instructions for Use

ICEROSS POST-OP TT & TF

| | | |
|----|-----------------------------|-----|
| | | 3 |
| EN | Instructions for Use | 4 |
| DE | Gebrauchsanweisung | 12 |
| FR | Notice d'utilisation | 21 |
| ES | Instrucciones para el uso | 29 |
| IT | Istruzioni per l'uso | 37 |
| NO | Bruksanvisning | 45 |
| DA | Brugsanvisning | 53 |
| SV | Bruksanvisning | 61 |
| EL | Οδηγίες Χρήσης | 69 |
| FI | Käyttöohjeet | 78 |
| NL | Gebruiksaanwijzing | 86 |
| PT | Instruções de Utilização | 94 |
| PL | Instrukcja użytkownika | 102 |
| CS | Návod k použití | 110 |
| TR | Kullanım Talimatları | 118 |
| RU | Инструкция по использованию | 126 |
| JA | 取扱説明書 | 134 |
| ZH | 中文说明书 | 141 |
| KO | 사용 설명서 | 148 |

Please continue with [Treatment Chart](#)

IFU – INSTRUCTIONS FOR USE

***NOTE:** Nomenclature and the definition of use in this IFU refers to “ISO 21065:2017. *Prosthetics and orthotics — Terms relating to the treatment and rehabilitation of persons having a lower limb amputation.*”

Immediate post-operative care (normally up to three days): Treatment provided after surgery where the focus is recovery from the procedure, achieving medical stability, preventing complications and initiating mobility

Early rehabilitation: Treatment, provided once the person is medically stable, to achieve optimum conditions for rehabilitation with or without a prosthesis

Temporary (preparatory) prosthesis: Individually manufactured prosthesis, to permit early mobilization and gait re-education prior to provision of a definitive prosthesis

PRODUCT DESCRIPTION

The Iceross Post-Op liner is a reusable liner used in post-operative therapy after trans-tibial or trans-femoral amputation.

The liner is available in two different versions. One version is for trans-tibial patients. It has an umbrella and a matrix. The other version is for trans-femoral patients. This version of the liner is a sleeve only. It does not have an umbrella or a matrix.

The two versions are made of silicone and have an EasyGlide coating. They are used for compression treatment of the residual limb after immediate postoperative care.

Compression treatment can be performed under medical supervision by care providers (i. e. doctors, physiotherapists, nurses, certified prosthetists). Therapy should begin 5–7 days post-operatively, ideally after an Össur Rigid Dressing (ORD) has been applied to the residual limb immediately following amputation surgery and used for 5–7 days.

Use of the liner should commence at such an early stage, as the largest change in residual limb volume occurs during the first weeks after amputation. Compression treatment is the most effective during this period. The liner is introduced gradually, ensuring the amputee and residual limb become accustomed to compression and the new environment. Compression therapy is done until the optimal goal of wound healing is obtained, and the patient is prepared for the next step, i.e. rehabilitation with or without a prosthesis.

The trans-tibial version of the liner can be included as a liner in a temporary (preparatory) prosthesis* interface in the later part of early rehabilitation, when the residual limb is healed and stable.

INDICATIONS FOR USE

Lower extremity amputation.

CONTRAINDICATIONS

Do not use Iceross Post-Op liners for compression therapy on amputees with the following conditions:

- PAOD (Doppler index below 0.6 in the non-amputated area)
- Decompensated heart failure
- Septic phlebitis
- Phlegmasia cerulea dolens
- Oozing dermatosis
- Severe sensory disturbances
- Intolerance and allergic reactions to used materials

INTENDED USE

Iceross Post-Op liners are used as part of early rehabilitation treatment of a residual limb after trans-tibial or trans-femoral amputation.

First, for compression therapy in the early stage of treatment.

Second, after trans-tibial amputation as a part of a temporary prosthetic interface in a later stage of treatment.

GENERAL SAFETY INSTRUCTIONS

Inform the patient concerning proper handling and care of the liner as stated in this instruction manual.

WARNING: Do not use the trans-femoral version as part of a temporary (preparatory) prosthesis interface.

PRODUCT MATERIAL LIST

Silicone, polyamide 6.6, elastane, polyester, steel

LINER SIZE SELECTION

Follow the steps listed below to provide the patient with a first experience of using a liner and avoid the risk of an uncomfortable degree of compression:

1. For trans-tibial amputees, determine the correct liner size by bending the knee slightly, allowing the soft tissue to hang down.
2. Measure 4 cm above the distal end of the residual limb. Make a small mark on the dressing/ residual limb.
3. Measure the circumference of the residual limb at the mark to determine the correct liner size (**Fig. 1**).
4. Select a liner matching the actual circumference measured, or the next available smaller size.

Example:

If the residual limb measures 24,5 cm at the indicated location, choose liner size 23,5 instead of size 25.

COMPRESSION THERAPY

Wound Dressing

Place a bandage over the wound area. An occlusive dressing is favorable for effective wound healing. Use an absorbent dressing (**Fig. 2**).

Trimming

Trim the liner according to user preference. Use an Iceross Curvemaster or a similar tool to produce a smooth and rounded edge, minimizing the risk of skin irritation and/or the liner tearing.

Donning and Doffing

Always communicate with the patient to ensure no discomfort is being experienced. If the patient is not able to actively participate in applying the liner, always make sure an assistant holds the limb in the correct position during application.

PRECAUTION: Do not use alcohol spray as a donning aid, as this can destroy the Easy Glide coating.

Donning

1. Remove the Iceross Post-Op liner from its packaging.
2. Place one hand inside the liner and turn it inside out, taking care not to damage it with fingernails (**Fig. 3**).
3. Ensure the liner is fully inverted exposing as much of the bottom as possible. Position the bottom of the liner directly against the distal residual limb (**Fig. 4**).
4. For transtibial amputees: Bend the knee slightly.
5. Roll the liner all the way up the limb (**Fig. 5**).
6. Feel along the length of the liner to ensure no air pockets are present; if air pockets are found reapply the liner.

Doffing

1. Before removing the Iceross liner, record the time it has been worn on the limb in the treatment chart.
2. Roll off the liner.
3. Following removal, invert the liner back to its normal position.

Daily Use

Don Iceross Post-Op liner twice daily, as following:

Day 1: don the liner for one hour in the morning and the afternoon.

Day 2: don the liner for two hours in the morning and the afternoon.

Day 3: don the liner for three hours in the morning and the afternoon.

Day 4: don for four hours in the morning and the afternoon.

Continue to don for four hours in the morning and the afternoon on each following day, until end of therapy.

| Time of day | Day 1 | Day 2 | Day 3 | Day 4 | Day 5 | Day 6 | Day 7 |
|-------------|--------|---------|---------|---------|---------|---------|---------|
| Morning | 1 hour | 2 hours | 3 hours | 4 hours | 4 hours | 4 hours | 4 hours |
| Afternoon | 1 hour | 2 hours | 3 hours | 4 hours | 4 hours | 4 hours | 4 hours |

NOTE:

- choose a larger size liner if the patient experiences discomfort or pain caused by compression of the residual limb (see section “*Liner Size Selection*”).
- don liner when the patient is being transferred e.g. from bed to wheelchair, and when going to the toilet. This will provide required compression during movement. When the liner is not in use, apply a soft dressing to maintain the effect of the compression therapy. Use a rigid dressing (ORD) overnight to protect the residual limb and to maintain contracture prophylaxis.

- Make sure of keeping a one-hour break between morning and afternoon sessions.

Size Changes During Compression Therapy

Record daily measurements for the residual limb in the treatment chart template to follow treatment progress. Fill in other information in the treatment chart template as applies.

As oedema decreases, change the liner size to provide continuous compression treatment. Measure the residual limb as stated above and choose proper liner size.

The patient may require 3–4 different liner sizes during the first 3 weeks of compression therapy.

Depending on success of compression therapy provided by the liner, the residual limb may be ready for prosthetic fitting.

In case of trans-tibial amputation, the last liner used for compression therapy can be used for prosthetic fitting and the following interim supply.

NOTE:

- If the size change leads to notable discomfort for the patient, use the previous larger size instead until the next day.
- Treatment is optimal if compression time is identical when changing from liner to liner. The patient should always decide on the length of the treatment, if able to do so.

CLEANING AND DISINFECTION WITHOUT STERILIZATION

Clean, disinfect and inspect the liner before first use and twice a day, after removing the liner from the amputee.

Manual Cleaning and Disinfection

1. Turn the liner inside out and wipe the surface with alcohol-impregnated wipes until completely wet. Allow to dry.
2. Turn the liner back to normal shape and repeat the procedure on the outer surface.

PRECAUTIONS:

- Only use ethanol or isopropanol solutions for manual cleaning of the liner.
- Only use a washing disinfectant if the liner will be sterilized afterwards. If the liner will not be sterilized, only clean and disinfect manually.

Liner Inspection

1. Check the liner for visible contamination and tears in the silicone layer.
2. Dispose liner if any tears are visible in the silicone layer.
3. If any contamination is visible, repeat manual cleaning and disinfection (see above).

STERILIZATION

WARNING: Always sterilize the liner before use on another patient.

PRECAUTIONS:

- Sterilization is limited to max. 10 cycles.
- Always wash and disinfect the liner before sterilization.

Preparation for Sterilization

Track the number of sterilization cycles on a chart. Mark the number directly on the liner. Identify the liner using the serial number at the distal end.

Prepare the liner for sterilization by performing only one of the two following methods:

- a. Manual cleaning and disinfection (see “*Manual Cleaning and Disinfection*” below)

OR

- a. By use of a washing disinfectant (see “*Washing Disinfectant*” below).

Inspect the liner (see “*Liner Inspection*” below) after performing a) or b).

Manual Cleaning and Disinfection

1. Turn the liner inside out and wipe the surface with alcohol-impregnated wipes until completely wet. Allow to dry.
2. Turn the liner back to normal shape and repeat the procedure on the outer surface.

Washing Disinfectant

1. Clean and disinfect using a washing disinfectant at 90 °C and an exposure time of 5 min.
2. Place the liner in the washing disinfectant, with the distal end at the top (**Fig. 6**). Make sure the stand is long enough to keep the proximal end open.
3. Add cleaning agent and start the washing disinfectant.

PRECAUTION: After using a washing disinfectant the liner is required to be sterilized before reuse (see “*Sterilization Process*” below).

Liner Inspection

1. Check the liner for visible contamination and tears in the silicone layer.
2. Dispose liner if any tears are visible in the silicone layer.
3. If any contamination is visible, repeat cleaning and disinfection (see “*Preparation for Sterilization*” above).
4. If neither tears nor contamination are visible, sterilize the liner (see “*Sterilization Process*” below).

Sterilization Process

1. Place the liner in a sterilization film and seal both ends.
2. Apply the necessary labeling.
3. Sterilize in an autoclave with fractionated pre-vacuum and moist heat at 134 °C and an exposure time of at least 5 min.

PRECAUTION: Do not use other media for sterilization, as silicone can absorb the substances it is exposed to.

USE WITH A TEMPORARY (PREPARATORY) PROSTHESIS

WARNING: Only use the liner as part of a temporary (preparatory) prosthesis for trans-tibial patients.

PRECAUTION: Monitor the performance of the liner. If there are changes in performance, stop use of the liner and contact the practitioner.

Liner Size Selection for a Temporary (Preparatory) Prosthesis

Determine the correct liner size (see the section “*Liner Size Selection*”).

Make sure the liner matrix covers 4 cm of bony structure. Do not allow the matrix to extend higher than the fibula head.

PRECAUTION: If the liner is too tight, it can cause pistoning and/ or numbness. If the liner is too loose, it can increase perspiration or let the limb move. The two can cause blisters and rashes.

Trimming

Trim the liner according to user preference and to increase range of motion. Use an Iceross Curvemaster or a similar tool to produce a smooth and rounded edge, minimizing the risk of skin irritation and/or the liner tearing.

PRECAUTION: do not cut the liner below the proximal trim-lines of the socket. Excessive trimming may lead to discomfort and may reduce the suction between the liner and limb and compromise suspension.

Locking suspension

If a lanyard suspension system is used, turn the connecting screw into the distal attachment of the liner. Refer to the instructions for use of the lanyard suspension system for the correct instructions.

If a locking pin is used, apply Loctite 410/411 or an equivalent compound to the threads. Then tighten the attachment pin to the fixed distal attachment. Össur recommends to torque to 4 Nm (3 ft lb). Refer to the instructions for use of the attachment pin for the correct instructions.

WARNINGS:

- Do not use flangeless pins with locking liners. A flangeless pin can cause damage to the distal attachment.
- If excessive movement or product wear is detected between the attachment pin and the distal attachment at any time, the user should be instructed to immediately discontinue use of the liner and consult with the practitioner. Please follow these recommendations to prevent health consequences.

Socket Fitting

Fit the socket to the user.

Follow the checklist to control the fit of the socket:

- The user can securely don socket without using lubricants
- Liner and inner wall of socket maintain full contact.
- Proximal fit is good.
- Minimized pistoning in the socket.
- If the user experiences insecurity or decreased contact with the socket due to volume fluctuations, one or more Iceross Socks can

be added to restore the socket fit.

PRECAUTION: If feedback indicates any of the above conditions are not met, the socket will need to be adjusted or remade accordingly.

Liner Care

Wash and clean liner daily on both sides following use. Remove the liner, turn the liner inside out and wash with a pH balanced, 100% fragrance- and dye-free soap.

PRECAUTION: Do not use fabric softeners, bleaches and other products/ cleaning solutions. Using these may result in liner damage.

Rinse the liner thoroughly with water after cleaning and wipe both sides dry. The liner can be used directly after washing. It does not need to dry overnight.

Always return the liner to its normal state with the distal attachment facing outwards directly after cleaning. Always check for damage or wear. Any damage may weaken the effectiveness of the liner and should be reported to the user's prosthetist immediately. Misuse may result in potential loss of suspension.

Skin Care

Wash and clean residual limb daily. Wash with a mild liquid, pH balanced, 100% fragrance and dye-free soap. If dry skin is noted, apply a pH balanced, 100% fragrance- and dye-free lotion to nourish and soften the skin.

PRECAUTION:

- Do not apply lotion to the residual limb immediately prior to donning the liner. To prevent damage to the liner, the residual limb should be clean and dry.
- Inspect the residual limb and skin each day. Make sure the two are healthy and continued use of the liner is safe.
- Do not use common household or bath products, including soaps, deodorants, perfumes, aerosol or abrasive cleaners. These can cause or contribute to skin irritation.

COMPLIANCE

Össur products and components are designed and tested according to the applicable official standards or an in-house defined standard when no official standard applies. Compatibility and compliance with these standards are only achieved if Össur products are used with other recommended Össur components.

We recommend carrying out regular yearly safety checks.

DISPOSAL

Return to practitioner for disposal.

All components of the product and packaging should be disposed of in accordance with respective national environmental regulations. Users should contact their local governmental office for information on how these items can be disposed of in an environmentally sound manner.

LIABILITY

The manufacturer recommends using the device only under the conditions specified and for the intended purposes. The device must be maintained according to the instructions for use.

The manufacturer is not liable for damage caused by component combinations that were not authorized by the manufacturer.

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EN – Caution: Össur products and components are designed and tested according to the applicable official standards or an in-house defined standard when no official standard applies. Compatibility and compliance with these standard is achieved only when Össur products are used with other recommended Össur components. If un-usual movement or product wear is detected in a structural part of a device at any time, the patient should be instructed to immediately discontinue use of the device and consult his/her clinical specialist. This product has been designed and tested based on single patient usage. This device should NOT be used by multiple patients. If any problems occur with the use of this product, immediately contact your medical professional.

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