
Instructions for Use

---

i-DIGITS<sup>®</sup> QUANTUM  
i-DIGITS<sup>®</sup> ACCESS  
i-DIGITS<sup>®</sup>

		3
EN	Instructions for Use	4
DE	Gebrauchsanweisung	13
FR	Notice d'utilisation	22
ES	Instrucciones para el uso	31
IT	Istruzioni per l'uso	40
NO	Bruksanvisning	49
DA	Brugsanvisning	58
SV	Bruksanvisning	67
EL	Οδηγίες Χρήσης	76
FI	Käyttöohjeet	85
NL	Gebruiksaanwijzing	94
PT	Instruções de Utilização	103
PL	Przewodnik techniczny	112
TR	Kullanım Talimatları	121
JA	取扱説明書	139
ZH	中文说明书	148
KO	사용 설명서	156
AR	تعليمات الاستخدام	165

## DESCRIPTION

The device is a microprocessor-controlled partial hand prosthesis with powered digits.

The wristband houses the microprocessor and removable batteries (**Fig. 1**).

Device labels (**Fig. 2**) can be found at the rear of the device, on the power supply and on the internal face of the battery pack.

This document is for patients (intended device operator) and healthcare professionals.

## INTENDED USE

The device is intended as a part of a prosthetic system that replaces the function of a missing upper limb.

Suitability of the device for the prosthesis and the patient must be evaluated by a healthcare professional.

The device must be fitted and adjusted by a healthcare professional.

### *Indications For Use and Target Patient Population*

- Upper limb loss, amputation, or deficiency

No known contraindications.

For weight limit and impact level, see **Table 3**.

## GENERAL SAFETY INSTRUCTIONS

The healthcare professional should inform the patient about everything in this document that is required for safe use of this device.

**Warning:** If there is a change or loss in device functionality, or if the device shows signs of damage or wear hindering its normal functions, the patient should stop using the device and contact a healthcare professional.

**Warning:** Improper handling and/or adjustment of device can cause malfunction.

**Caution:** The device does not provide sensation, such as temperature or moisture.

**Warning:** Do not disassemble or modify device in any way except for changing digit covers.

**Caution:** Do not expose to high forces or vibration.

**Caution:** Do not carry objects using only the tips of the digits. Carry objects by evenly distributing weight across the digits, close to the knuckles and palm of the hand (**Fig. 3**).

**Warning:** Do not use to operate machinery with moving parts that may cause personal injury or damage.

**Warning:** Do not use device for extreme activities or in hazardous environments that may cause injuries to a natural hand.

**Caution:** Electrodes may contain small amounts of nickel. If skin irritation or discomfort arises, stop using the device immediately and contact your healthcare professional.

**Caution:** Patient must comply with local regulations on the operation of automobiles, aircraft, sailing vessels of any kind, and any other motorized vehicle or device. It is entirely the user's responsibility to seek confirmation that they are physically and legally able to drive using the device and to the fullest extent permitted by law.

**Caution:** Do not use in an environment where exposure to highly electrical and/or magnetic fields can occur (e. g. near electrical transformers, high-power radio/TV transmitters, RF surgical equipment, CT and MRI scanners) as interference with EMG signal is likely to occur.

**Caution:** Do not use your device while undergoing specific investigations or treatments near other medical electrical equipment.

**Note:** Do not apply lotions or creams to the residual limb immediately before donning the device. The residual limb should be clean and dry to limit interference with electrodes.

**Caution:** Do not expose electrode to dirt or liquids other than sweat. Limit the impact of sweat by cleaning regularly using the recommended cleaning agents.

**Caution:** The device must be thoroughly inspected for defects before fitting.

**Note:** Electrodes are classified as Type BF applied part.

**Note:** Wristband is classified as Type B applied part.

The device is for single patient use.

## REQUIRED COMPONENTS

The device is delivered with a power supply.

### **Device Control App**

Some device functions can be accessed and configured through the My i-Limb® (patient) and Biosim® (healthcare professional) mobile apps.

The healthcare professional can use the Biosim app for patient assessment, configuration of the device, system checks, calibrations, and updates.

The patient can use the My i-Limb app to check the status of the device, access training programs, direct control of the device (e. g. grip selection) and view instructions for use.

Depending on the device, the apps are either only available in Apple App Store or the Google Play Store.

For more detailed device information, see **Table 1**.

### **DEVICE SELECTION**

See Össur catalog for a list of recommended components.

<b>Table 1 - Feature Comparison</b>			
<b>Control Options</b>	<b>i-Digits® Quantum</b>	<b>i-Digits® (Only Available in France)</b>	<b>i-Digits® Access</b>
App Control	Yes	Yes	Yes
Available in Apple App Store	Yes	Yes	Yes
Available in Google Play Store	Yes	Yes	No
Muscle Control	Yes	Yes	Yes
Proximity Control	Yes	Yes	No
Gesture Control	Yes	No	No
Grips Available	20	20	12
My Grips	12	12	No

### **Connecting to the App**

Each device has a number for identifying the device, located on the wristband underneath the battery.

When connecting to the My i-Limb app, the number will be displayed on the connection screen. Selecting the number connects the app to the device. Once connected to the app, the device number is displayed in the 'about' section.

### **Controlling the Device**

Control options available vary between different devices or models.

For more detailed device information, see **Table 1**.

### **App Control**

An automated grip can be accessed at the touch of an icon on the Biosim and My i-Limb app. These are called quick grips. The device will exit the grip when the icon is tapped again or by selecting another grip icon to enter.

### **Muscle Control**

Triggers are specific muscle signals that can be used to access an automated grip. There are 4 potential muscle triggers: hold open, double impulse, triple impulse and co-contraction. Use the app to turn on muscle control.

### **Proximity Control (available on i-Digits® Quantum and i-Digits®)**

Grip Chips are small Bluetooth devices that activate grips on the prosthesis when the prosthesis is moved close to them or by tapping the Grip Chip.

Grip Chips activate My Grips through proximity or tapping. Program using the Biosim or My i-Limb app. Before the first use, pull the plastic tab to activate the battery. Ensure the device is not connected to the Biosim or My i-Limb app on your phone.

### **Proximity Activation**

Fully open the device near (15 cm/ 6") to a Grip Chip

Relax the muscles and wait up to 3 seconds for activation.

Give a long open hand position signal to exit a grip accessed via Proximity Control.

## Double-tap Activation

Double-tap a Grip Chip to activate the grip quickly.  
The light on the Grip Chip flashes once.

**Note:** Wait 3 seconds between each double tap to avoid incorrect detection.

**Note:** The device must maintain a fully extended and open hand position with digits stalled open to achieve Proximity Control.

## Grip Chip Battery Replacement

Release the holding catch to remove the battery. Insert the new battery under the retaining clip, ensuring the Össur logo faces outwards.

The battery is expected to last 6 months.

Battery life is reduced if it is stored close to the device.

The battery is not rechargeable.

**Battery type:** CR1616. Replacement batteries are available from customer support upon request.

## Gesture Control (available on i-Digits® Quantum)

Gesture control enables the patient to quickly access multiple grips, with only the necessity to learn one technique. Grips are assigned to gesture control in the My i-Limb or Biosim apps.

To access Gesture Control:

1. The device must maintain a fully extended and open hand position with digits stalled open.
2. Hold arm parallel to the ground (elbow bent to 90°).
3. Maintain an open signal until the finger twitches.
4. Move the hand within 1 second in a direction set to access a desired grip, the device will adopt that grip.

**Note:** Gesture control can also be accessed using co-contraction, depending on the set-up of the device through the healthcare professional.

**Note:** It is recommended to practice activating Gesture Control in different positions within the activation range.

## My Grips (available on i-Digits® Quantum and i-Digits®)

My Grips are personalized grips that can be customized in My i-Limb App. Each digit can be positioned and set individually to be active or inactive.

To ensure your grip works well for a particular object or task, use the "Test my Grip" button. Follow these steps:

1. Close the digit around the object.
2. The digits will then operate according to your programmed setup.

Remember to manually rotate the thumb to the required position.

## USAGE

### Device operation

#### Turning the Device On

1. Push and hold the power button on the Wristband.
2. The status indicator display will light up and then turn off.

#### Turning the Device OFF

1. Push and hold the power button on the Wristband.
2. The status indicator display will light up and then turn off.

### Battery

Batteries are placed into the battery case on the wristband, where they are secured by the holding catch located towards the wrist. Release the holding catch to remove the battery (**Fig. 4**).

Each battery is equipped with a light that indicates the charge level of the battery. When the battery charge level is low, the light will light up red. The light will stay lit until the battery has been charged to an adequate level. When the device's battery is fully depleted, the light will turn off.

From a single charge of the supplied 800mA batteries, the device can be opened and closed 600 times. This will vary depending on battery age and device use.

After the batteries are drained, swap them out for the backup set.

**Note:** Electronic devices connected to a power outlet might impact the function of the device.

**Warning:** Do not bend or exert excessive pressure on the battery.

**Warning:** Do not expose to high temperatures or to flames.

**Warning:** Do not disassemble, crush, heat above 70 °C or incinerate battery pack

**Warning:** Do not short-circuit the battery.

### Charging

1. Remove the batteries from the wristband and insert into the Charger Base unit.
2. Connect the Charger Base with the power supply.
3. Plug the power supply into the wall socket.

Charging time: 90 minutes to 3 hours.

### Battery Charge Level Indication

The charging state is indicated on the back of the base unit (**Fig. 5**):

- 3rd light on: Charger is powered.
- 2nd and 5th lights flashing green: Batteries are charging.
- 2nd and 5th lights solid green: Batteries are charged.
- 1st and 4th red lights on: Battery fault, unplug and try again. If lights continue to illuminate, contact your healthcare professional.

### Covers

The cover must be worn on digits and/or thumb that have been assembled.

If the cover is damaged, replace it by removing the damaged cover and sliding the new, correctly sized cover onto the digit and/or thumb.

**Caution:** The device must always be worn with a cover(s).

**Note:** Make sure covers are full to the tip of each digit.

**Caution:** Always wear an undamaged cover supplied by Össur to avoid electrostatic build-up and discharge.

### Cleaning and care

It is recommended to maintain the device in good outward condition by cleaning the surface of the device regularly.

Clean with a damp cloth and a mild soap. Dry with a cloth after cleaning.

Clean the device regularly with 70 - 90% isopropyl alcohol to help with disinfection.

### Environmental Conditions

The device should not come in contact with fresh water, salt water or chlorinated water.

Dry with a cloth after contact with fresh water or humidity.

Use of device outside of the specified environmental conditions can result in device malfunction.

See which environmental conditions apply to the device in **Table 2**.

Table 2 - Environmental Conditions			
	Use	Shipping	Extended storage
Temperature	0°C to 40°C (32°F to 104°F)	-40°C to 70°C (-40°F to 158°F)	-25°C to 70°C (-13°F to 158°F)
Relative humidity	0% to 95%	0% to 95%	10% to 95%
Atmospheric pressure	700 hPa to 1060 hPa	700 hPa to 1060 hPa	700 hPa to 1060 hPa

### MAINTENANCE

Examine the device regularly for signs of damage. If you find signs of damage, consult your healthcare professional.

Device should be serviced by a qualified healthcare professional annually.

**Caution:** Batteries are to be replaced annually.

**Warning:** Do not attempt device repairs. Risk of injury or damage to device. Send device to Össur for maintenance or repair.

**Caution:** Do not service or perform maintenance while device is in use.

## SPECIFICATIONS

Table 3 - Device Specifications	
Voltage	7.4 V (nominal)
Max. current	5 A
Battery capacity	Rechargeable lithium polymer 7.4 A (nominal); 800 mAh capacity
Max. device load	20 kg/ 44 lbs
Finger carry load (static limit)	5 kg/ 11 lbs
Time to full open to full close	0.8 seconds
Expected service life	5 years
Impact Level	The device is for low to moderate impact use

Table 4 - Power Supply Specifications	
Manufacturer	FRIWO
Model no.	
Device	FW8030M/24
Charging Base	PL069653
Input	100–240 VAC, 50–60 Hz, 0.6–0.3 A
Input/Output	24 VDC, 1.25A

### ELECTROMAGNETIC COMPATIBILITY

**Warning:** Use of this device adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this device and the other equipment should be observed to verify that they are operating normally.

**Warning:** Use of accessories and cables other than those specified or provided by the manufacturer of this device could result in increased electromagnetic emissions or decreased electromagnetic immunity of this device and result in improper operation.

The device is suitable for use in any environment except in flammable/explosive atmospheres, or where exposure to high electrical and/or magnetic fields can occur (e. g., electrical transformers, high-power radio/TV transmitters, RF surgical equipment, CT and MRI scanners).

Table 5 - Electromagnetic Emissions Compliance Information		
Emissions test	Standard	Compliance level
Conducted and radiated RF emissions	CISPR 11	Group 1 - Class B
Harmonics emissions	IEC 61000-3-2	Group 1 - Class B
Voltage fluctuations / Flicker emissions	IEC 61000-3-3	Group 1 - Class B

**Table 6 - Electromagnetic Immunity Compliance Information**

Immunity Test	Standard	Compliance Level	
		Device	Charging base
Electrostatic discharge	IEC 61000-4-2	±8 kV contact ±15 kV air	
Electrical fast transient/burst	IEC 61000-4-4	Not applicable	± 2 kV power supply lines ± 1 kV input/output lines
Surge	IEC 61000-4-5	Not applicable	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth
Voltage dips	IEC 61000-4-11	Not applicable	0 % UT; 0.5 cycle At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315°
			0 % UT; 1 cycle and 70 % UT; 25/ 30 cycles Single phase: at 0°
Voltage interruptions	IEC 61000-4-11	Not applicable	0 % UT; 250/ 300 cycle
Power frequency	IEC 61000-4-8	Not applicable	30 A/m
Conducted disturbances included by RF fields	IEC 61000-4-6	Not applicable	3V 0.15 MHz – 80 MHz 6V in ISM and amateur bands between 0.15 MHz and 80 MHz 80 % AM at 1 kHz
Radiated RF IEC	IEC 61000-4-3	12 V/m 26 MHz to 1GHz 10 V/m 1000 MHz to 2.7GHz 80% AM at 1kHz	10 V/m 80 MHz – 2.7 GHz 80 % AM at 1 kHz
Proximity magnetic fields	IEC 61000-4-39	Not applicable	30kHz at 8 A/m EN 60601- 1-2
Immunity to proximity fields from RF wireless communications equipment	IEC 61000-4-39	134.2kHz at 65 A/m 13.56MHz at 7.5 A/m 30kHz at 6 A/m	

Table 7 - RF Wireless Communications Compliance Information				
Test frequency (MHz)	Band (MHz)	Service	Modulation	Immunity test level (V/m)
385	380-390	TETRA 400	Pulse modulation 18 Hz	27
450	430-470	GMRS 460, FRS 460	Pulse modulation 18 Hz	28
710	704-787	LTE Band 13, 17	Pulse modulation 217 Hz	9
745				
780				
810	800-960	GSM 800 /900, TETRA 800, iDEN 820, CDMA 850, LTE Band 5	Pulse modulation 18 Hz	28
870				
930				
1720	1700-1990	GSM 1800, CDMA 1900, GSM 1900, DECT, UMTS LTE Band 1, 3, 4, 25	Pulse modulation 217 Hz	28
1845				
1970				
2450	2400-2570	Bluetooth, WLAN, 802.11 b/g/n, RFID 2450, LTE Band 7	Pulse modulation 217 Hz	28
5240	5100-5800	WLAN 802.11 a/n	Pulse modulation 217 Hz	9
5500				
5785				

This device contains the following radio frequency transmitters:

Table 8 - Wireless Modules			
Model	Regulatory certificates	Type and frequency characteristics	Effective radiated power
Bluetooth Low Energy Dual Mode Module Model BR-LE4.0-D2A	<b>FCC</b> Contains FCC ID: XDULE40-D2	(Dual Mode) Version V2.1 +ED (GFSK + $\pi/4$ DQPSK + 8DPSK) 2402 -2480 MHz Version V4.0 (GFSK) 2402 -2480 MHz	Adjustable Power (-23 dBm to 10.5 dBm) short to long range
	<b>Canada</b> Contains IC: 8456A-LE4D2		
	<b>Japan</b> Contains transmitter with certificate number  R 205-160268		

#### USA-Federal Communications Commission (FCC)

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses, and can radiate radio frequency energy. If not installed and used in accordance with the instructions, it may cause harmful interference to radio communications. However, there is no ensured specification that interference will not occur in a particular installation.

If this equipment does cause harmful interference to radio or television reception, which can be determined by tuning the equipment off and on, the user is encouraged to try and correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.

- Increase the distance between the equipment and the receiver.
- Connect the equipment to outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

Any changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

**Caution:** Exposure to Radio Frequency Radiation – This device must not be co-located or operating in conjunction with any other antenna or transmitter.

### Canada – Industry Canada (IC)

This device complies with RSS 210 of Industry Canada.

Operation is subject to the following two conditions:

- this device may not cause interference, and
- this device must accept any interference, including interference that may cause undesired operation of this device.

**Caution:** Exposure to Radio Frequency Radiation – The installer of this radio equipment must ensure that the antenna is located or pointed such that it does not emit RF field in excess of Health Canada limits for the general population; consult Safety Code 6, obtainable from Health Canada's website [www.hc-sc.gc.ca/rpb](http://www.hc-sc.gc.ca/rpb)

Table 9 - Recommended separation distance between portable and mobile RF communications equipment and the Össur Myoelectric Prosthetic Devices			
Rated maximum output power of transmitter in Watt	Separation distance according to frequency of transmitter in meters		
	150 kHz to 80 MHz d = 1.2 √P	80 MHz to 800 MHz d = 1.2 √P	800 MHz to 2.7GHz d = 2.3 √P
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

**Note:** At 80MHz and 800MHz, the separation distance for the higher frequency range applies.

**Note:** Guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people.

### REPORT OF SERIOUS INCIDENT

Any serious incident in relation to the device must be reported to the manufacturer and relevant authorities.

### DISPOSAL

The device and packaging must be disposed of in accordance with respective local or national environmental regulations.

### LIABILITY

Össur does not assume liability for the following:

- Device not maintained as instructed by the instructions for use.
- Device assembled with components from other manufacturers.
- Device used outside of recommended use condition, application, or environment.

## Compliance

The device has been designed and tested according to the applicable international standards or in-house defined standards when no international standard exists or applies.

The device has been tested by an independent accredited laboratory and found to comply to:

- EN60601-1
- EN60601-1-2

## SYMBOLS

Medical Device

BF Applied Part

Consult instructions for use

Device contains electronic components and/or batteries that should not be disposed of in regular waste

**IP22**

Protected against solid foreign objects of 12.5mm diameter and greater and protected against splashing water.

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