







Instructions for Use





---

MIAMI TLSO™ 456

		3
EN	Instructions for Use	5
DE	Gebrauchsanweisung	7
FR	Notice d'utilisation	10
ES	Instrucciones para el uso	13
IT	Istruzioni per l'uso	16
NO	Bruksanvisning	19
DA	Brugsanvisning	22
SV	Bruksanvisning	24
EL	Οδηγίες Χρήσης	27
FI	Käyttöohjeet	30
NL	Gebruiksaanwijzing	33
PT	Instruções de Utilização	36
PL	Instrukcja użytkowania	38
CS	Návod k použití	41
TR	Kullanım Talimatları	44
RU	Инструкция по использованию	47
ZH	中文说明书	50

## SYMBOLS

Medical Device

## INTENDED USE

The device is intended to provide gross immobilization of the thoracic and lumbar spine

The device must be fitted and adjusted by a healthcare professional.

### *Indications for use*

Indications that require gross thoracic and lumbar spine immobilization in all three planes, in which the following conditions may be present:

- Herniated disc
- Degenerative Spinal Pathologies
- Post-Surgical Stabilization
- Stable, non-displaced spinal fractures
- Spinal stenosis
- Spondylolisthesis
- Spondylolysis

### *Contraindications*

- Unstable, displaced fractures

### **Warnings and Cautions:**

- Do not over tighten the compression system to the point where it causes discomfort or difficulty breathing. For sitting, the patient may find it desirable to slightly loosen the compression system.
- This is a supportive device only and is not intended or guaranteed to prevent spinal injury.
- Do not machine wash or tumble dry.
- Össur accepts no responsibility for any injury which might occur while using the orthosis.
- Make sure the hinged adjustment mechanisms are tightened sufficiently to lock them into position.
- In addition to increasing the risk of infection, used product components may contain small unnoticeable defects that may lead to non-optimal performance, and even mechanical failure.
- The braces are not designed to be worn while bathing.
- The device is a non-invasive, adjustable brace for single patient use.

## GENERAL SAFETY INSTRUCTIONS

The healthcare professional should inform the patient about everything in this document that is required for safe use of this device.

Any serious incident in relation to the device must be reported to the manufacturer and relevant authorities.

The patient should stop using the device and contact a healthcare professional:

- If there is a change or loss in device functionality, or if the device shows signs of damage or wear hindering its normal functions.
- If any pain, skin irritation, or unusual reaction occurs with the use of the device.

The device is for single patient – multiple use.

## FITTING INSTRUCTIONS

While carrying out the following instructions, please refer to the overview figure for locating components mentioned in the text (**Fig. 1**).

### Device Application

1. Before application, ensure the Belt Compression System (A) is stretched to its full width and remove the Overlap Closures (B) from both ends of the belt.
2. Log roll the patient onto their side. Omit this step if fitting in standing position. If standing, patient must stand with back straight and chest out.
3. Position the Posterior Panel and PTE (C) centered over the spine, with the bottom of the Posterior Panel at approximately the sacrococcygeal joint (**Fig. 2**).  
**Note:** Panels can be removed from their fabric sleeve and customized by heat-molding, grinding, or trimming, to achieve correct patient fit.
4. Log roll the patient back into the supine position. Omit this step if fitting in standing position.
5. Wrap the lumbar belt around the waist and through the Össur Fit Tool (**Fig. 3**). Adjust the length of the belt arms to position tool centrally. Secure the belt ends to desired sizing (**Fig. 4**).
6. Remove Össur Fit Tool from the belt.  
**Note:** If an Össur Fit Tool is not available, lay the belt flat with Overlap Closures removed and compression system loosened. Fold each belt arm equally to achieve the desired waist circumference and secure the Belt Ends (**Fig. 5**). Follow steps 2-4.  
**Note:** If the patient's waist exceeds 127cm (50"), an extension belt can be added to each belt arm to reach a maximum belt circumference of 178cm (70").
7. Place Belt Ends into the Overlap Closures (**Fig. 6**). To obtain the lowest profile fit, the belt arm ends can be trimmed if later adjustments are not required.
8. Close the belt around the waist.
9. Slide fingers through the holes in the two compression system handles and pull until the belt is at the appropriate tightness (**Fig. 7**). Adjust the compression system cord length if necessary. (**Fig. 8**)
10. Remove Strap Anchor (D) from belt, place the Shoulder Straps (E) over the shoulder, and re-fasten the Strap Anchor onto the belt (**Fig. 9**).  
**Note:** Place the Strap Anchor on a smooth area of the belt as posterior to the patient as possible without introducing axilla discomfort.
11. Pull webbing to tighten Shoulder Straps.
12. Adjust the chest strap position by sliding the D-Ring (F) (**Fig. 10**).
13. Fasten Chest Buckle (G) and tighten chest strap.
14. Upon sitting, standing or ambulation, the brace may require fine tuning for comfort, fit and stabilization. This may include adjusting or tightening straps, compression system, and/or rigid panels.

### Device Adjustments

#### Panel & Belt

- To heat mold the panels, use a heat gun (175 °C/350 °F).
- The belt arms can be trimmed anteriorly where the Belt Ends are folded, or posteriorly by detaching the belt arms from the Hook and Loop Strip adjacent to the Belt Compression System (**Fig. 11**). If trimming posteriorly, ensure the Belt Ends are retained and folded (**Fig. 6**).

## **Device Removal**

1. Detach Chest Strap Buckle.
2. Detach compression system handles and loosen. Affix them onto the belt after loosening.
3. Detach the Overlap Closure and remove device.
4. Before reapplying the device and to ensure a proper fit, be sure the compression system is stretched to its full width by pulling the ends of the belt apart.

## **Accessories and Replacement Parts**

Please refer to the Össur catalog for a list of available replacement parts or accessories.

## **USAGE**

### **Cleaning and care**

Washing the device with the soft goods detached allows for more thorough cleaning.

- Hand-wash using mild detergent and rinse thoroughly.
- Air dry.

**Note:** Do not machine-wash, tumble dry, iron, bleach, or wash with fabric softener.

**Note:** Avoid contact with salt water or chlorinated water. In case of contact, rinse with fresh water and air dry.

## **DISPOSAL**

The device and packaging must be disposed of in accordance with respective local or national environmental regulations.

## **LIABILITY**

Össur does not assume liability for the following:

- Device not maintained as instructed by the instructions for use.
- Device assembled with components from other manufacturers.
- Device used outside of recommended use condition, application, or environment.

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