
















Instructions for Use





MIAMI TLSO™ 464

		3
EN	Instructions for Use	6
DE	Gebrauchsanweisung	10
FR	Notice d'utilisation	14
ES	Instrucciones para el uso	18
IT	Istruzioni per l'uso	22
NO	Bruksanvisning	26
DA	Brugsanvisning	30
SV	Bruksanvisning	34
EL	Οδηγίες Χρήσης	38
FI	Käyttöohjeet	42
NL	Gebruiksaanwijzing	46
PT	Instruções de Utilização	50
PL	Instrukcja użytkownika	54
TR	Kullanım Talimatları	58
RU	Инструкция по использованию	62
JA	取扱説明書	66
ZH	中文说明书	70
KO	사용 설명서	74

ENGLISH

SYMBOLS

Medical Device

INTENDED USE

The device is intended to provide gross immobilization of the thoracic and lumbar spine

The device must be fitted and adjusted by a healthcare professional.

External fixation for fractures or soft tissue injuries

Indications for use

Indications that require gross thoracic and lumbar spine immobilization in all three planes, in which the following conditions may be present:

- Herniated disc
- Degenerative Spinal Pathologies
- Post-Surgical Stabilization
- Stable, non-displaced spinal fractures
- Spinal stenosis
- Spondylolisthesis
- Spondylolysis

Contraindications

- Unstable, displaced fractures

Warnings and Cautions:

- Do not over tighten the compression system to the point where it causes discomfort or difficulty breathing. For sitting, the patient may find it desirable to slightly loosen the compression system.
- This is a supportive device only and is not intended or guaranteed to prevent spinal injury.
- Do not machine wash or tumble dry.
- Össur accepts no responsibility for any injury which might occur while using the orthosis.
- Make sure the hinged adjustment mechanisms are tightened sufficiently to lock them into position.
- In addition to increasing the risk of infection, used product components may contain small unnoticeable defects that may lead to non-optimal performance, and even mechanical failure.
- The braces are not designed to be worn while bathing.
- The device is a non-invasive, adjustable brace for single patient use.

GENERAL SAFETY INSTRUCTIONS

The healthcare professional should inform the patient about everything in this document that is required for safe use of this device.

Any serious incident in relation to the device must be reported to the manufacturer and relevant authorities.

The patient should stop using the device and contact a healthcare professional:

- If there is a change or loss in device functionality, or if the device shows signs of damage or wear hindering its normal functions.
- If any pain, skin irritation, or unusual reaction occurs with the use of the device.

The device is for single patient – multiple use.

FITTING INSTRUCTIONS

While carrying out the following instructions, please refer to the overview figure for locating components mentioned in the text (**Fig. 1**).

Preparation

No preparation is needed for over-the-shoulder strapping configuration.

For axillary strapping configuration:

- Remove strap from the Posterior Thoracic Extension (PTE) (A)
- Reattach the free end of the strap to the strap anchor (**Fig. 2**)

For over-the-shoulder and axillary strapping configuration (requires purchase of additional kit):

- Remove Pectoral Disc assemblies (B) from the ball joints on the Yoke (C) (**Fig. 3**)
- Attach secondary buckles from accessory kit over the ball joints into the groove (**Fig. 4**)
- Reattach the Pectoral Disc assemblies onto the ball joints
- Attach the free end of the strap from the accessory kit to the strap anchors (**Fig. 2**)

Device Application

1. Before application, ensure the Belt Compression System (D) is stretched to its full width and remove the Overlap Closures (E) from both ends of the belt.
2. Log roll the patient onto their side. Omit this step if fitting in standing position. If standing, patient must stand with back straight and chest out.
3. Position the belt by centering the Posterior Panel (PTE) over the spine, with the bottom of the Posterior Panel approximately at the sacrococcygeal joint. (**Fig. 5**).

Note: Panels can be removed from their fabric sleeve and customized by heat-molding, grinding, or trimming, to achieve correct patient fit.

4. Log roll the patient back into the supine position. Omit this step if fitting in standing position.
5. Wrap the lumbar belt around the waist and through the Össur Fit Tool (**Fig. 6**). Adjust the length of the belt arms to position tool centrally. Secure the Belt Ends to desired sizing (**Fig. 7**).
6. Remove Össur Fit Tool from belt.

Note: If an Össur Fit Tool is not available, lay the belt flat with Overlap Closures removed and compression system loosened. Fold each belt arm equally to achieve the desired waist circumference and secure the Belt Ends (**Fig. 8**). Follow steps 2-4.

Note: If patient's waist exceeds 127cm (50"), an extension belt can be added to each belt arm to reach a maximum belt circumference of 178cm (70").

7. Place Belt Ends into the Overlap Closures (**Fig. 9**). To obtain the lowest profile fit, the belt arm ends can be trimmed if later adjustments are not required.
8. Place the Anterior Thoracic Extension (ATE) (F) onto the patient's front, centering the anterior panel on the abdomen with the bottom edge just above the symphysis pubis, while still allowing the patient to sit comfortably (**Fig. 10**).
9. Adjust the height of the ATE by pushing the button on the anterior panel and sliding the Strut (G) to the appropriate length. The superior edge of the pectoral pads should be located approximately one inch below the sternal notch (**Fig. 11**).
10. Using the provided tool, loosen the screw by one full turn and then adjust the angle on the hinged adjustment mechanism to accommodate sternal contours (**Fig. 12**).

Note: Make sure the hinged adjustment mechanism is tightened enough to lock it into position.

11. Fasten the non-labeled belt end onto the ATE panel.
12. Fasten the retaining straps onto the ATE panel soft good, over the non-labeled belt end (**Fig. 13**). Close the belt.
13. Connect both quick-release Buckles (H), tighten the shoulder straps, and position the shoulder pad for optimal comfort (**Fig. 14**).
14. If desired, remove alligator, trim the strap, and replace alligator.
15. Slide fingers through the holes in the two compression system handles and pull until the belt is at the appropriate tightness (**Fig. 17**). If necessary, adjust the compression system cord length (**Fig. 18**).
16. The brace may require fine tuning for comfort, fit and stabilization when sitting, standing, or ambulating. This may include adjusting or tightening straps, compression system, rigid panels, and/or ATE/PTE.

Device Adjustments

Axillary Strap Only Option

- Replace step 13 with the following: fasten the strap anchor onto the belt, attach the axillary strap to the buckles and then tighten strap (**Fig. 16**).
- Resume steps 14 to 16.

Dual Strap Option

- After step 13, attach axillary strap to secondary buckles (**Fig. 15**)
- Fasten the strap anchor onto the belt and then tighten strap (**Fig. 16**)
- Resume steps 14 to 16.

Hinge/Length

If device needs to be taller (**Fig. 19**):

1. Remove hinge screw carefully with provided tool; do not lose the preload spring.
2. Attach extension piece to lower stay.
3. Attach Yoke to extension piece.
4. Using the provided tool, adjust the angle on the hinged adjustment mechanism to accommodate sternal contours.

Note: Make sure the hinged adjustment mechanism is tightened sufficiently to lock it into position.

Panel & Belt

- To heat mold the panels, use a heat gun (175 °C/350 °F).
- The belt arms can be trimmed anteriorly where the Belt Ends are folded, or posteriorly by detaching the belt arms from the Hook and Loop Strip adjacent to the Belt Compression System (**Fig. 20**). If trimming posteriorly, ensure the Belt Ends are retained and folded (**Fig. 9**).

Device Removal

1. Detach strap Buckle(s) on right shoulder.
2. Detach compression system handles and loosen. Fix them onto the belt after loosening.
3. Detach the Overlap Closure and remove device.
4. Before reapplying the device and to ensure a proper fit, be sure the compression system is stretched to its full width by pulling the ends of the belt apart. Reattach right strap Buckle before reapplying the device.

Accessories and Replacement Parts

Please refer to the Össur catalog for a list of available replacement parts or accessories.

USAGE

Cleaning and care

Washing the device with the soft goods detached allows for more thorough cleaning.

- Hand-wash using mild detergent and rinse thoroughly.
- Air dry.

Note: Do not machine-wash, tumble dry, iron, bleach, or wash with fabric softener.

Note: Avoid contact with salt water or chlorinated water. In case of contact, rinse with fresh water and air dry.

DISPOSAL

The device and packaging must be disposed of in accordance with respective local or national environmental regulations.

LIABILITY

Össur does not assume liability for the following:

- Device not maintained as instructed by the instructions for use.
- Device assembled with components from other manufacturers.
- Device used outside of recommended use condition, application, or environment.

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