






Instructions for Use

MAUCH[®] KNEE
MAUCH[®] KNEE PLUS

		3
EN	Instructions for Use	4
DE	Gebrauchsanweisung	9
FR	Notice d'utilisation	14
ES	Instrucciones para el uso	19
IT	Istruzioni per l'uso	24
NO	Bruksanvisning	29
DA	Brugsanvisning	34
SV	Bruksanvisning	39
EL	Οδηγίες χρήσης	44
FI	Käyttöohjeet	49
NL	Gebruiksaanwijzing	54
PT	Instruções de Utilização	59
PL	Instrukcja użytkowania	64
TR	Kullanım Talimatları	69
RU	Инструкция по использованию	74
JA	取扱説明書	79
ZH	中文说明书	84
KO	사용 설명서	89

DESCRIPTION

The Mauch Knees are single axis hydraulic knee systems with swing and stance control (SNS). The devices are designed to provide variable hydraulic resistance to flexion and extension during the swing phase based on changes in the patient's walking speed and yielding stance control. A mode selector switch allows for manual locking and free swing functions.

Mauch Knee

The device has integrated proximal and distal Male Pyramid Adapters.

Mauch Knee Plus

The device has either threaded or through hole proximal and distal 4-Hole Connectors.

INTENDED USE

The device is intended as a part of a prosthetic system that replaces knee function of a missing lower limb.

Suitability of the device for the prosthesis and the patient must be evaluated by a healthcare professional.

The device must be fitted and adjusted by a healthcare professional.

Indications For Use and Target Patient Population

- Lower limb loss, amputation, or deficiency
- No known contraindications

The device is for moderate to high impact use, e.g., walking and occasional running.

Mauch Knee

The weight limit for the device is 136kg.

Mauch Knee Plus

The weight limit for the device is 166 kg.

GENERAL SAFETY INSTRUCTIONS

Warning: Use of a lower limb prosthetic device carries an inherent risk of falling which may lead to injury.

The healthcare professional should inform the patient about everything in this document that is required for safe use of this device.

Warning: If there is a change or loss in device functionality, or if the device shows signs of damage or wear hindering its normal functions, the patient should stop using the device and contact a healthcare professional.

Warning: Avoid placing hands or fingers near moving joints.

Caution: Do not adjust other screws than described in these instructions. The device is for single patient use.

ASSEMBLY INSTRUCTIONS

Adapter Options

Warning: Risk of structural failure. Components from other manufacturers have not been tested and may cause excessive load on the device.

Mauch Knee

Caution: For the device's proximal connection, using a 4-hole adapter for a long residual limb can limit knee flexion, because it may hit the knee frame and cover before the full range of knee flexion is achieved. In this case, it is recommended to use either 3 or 4 Prong Socket Adapter, please refer to Össur catalogue.

Mauch Knee Plus

Caution: Care should be taken with some component arrangements and/or alignment. It can limit knee flexion and cause damage to the knee frame or the cylinder.

ALIGNMENT INSTRUCTIONS

Bench Alignment (Fig. 1)

Alignment Goal

Alignment reference line (B) should:

- pass through midpoint of socket at the ischial tuberosity level (D)
- pass through knee axis (A)
- fall at the 1/3 mark on the inside of the foot cover.

Note: Prioritize knee alignment over foot alignment if there is a mismatch.

Alignment Instructions

1. Position the foot so that the alignment reference line (B) falls at the 1/3 mark on the inside of the foot cover (with the foot cover and shoe on). Consider the external rotation of the foot.
2. Use the applicable adapters to connect the knee to the foot and establish the correct knee center height.
3. Position the knee so that the alignment reference line passes through knee axis (A)
4. On the lateral side of the socket, make a first mark at the midpoint of the socket at the ischial tuberosity level (D). Make a second mark at the midpoint of the socket distally (E). Draw a line through both marks.
5. Position the socket so the alignment reference line (B) goes through first mark at the midpoint of the socket at the ischial tuberosity level (D).
6. Adjust the socket flexion to 5° in addition to the existing position (i. e., hip flexion contracture) and set the height of the full prosthesis.
7. Use the applicable adapters to connect the knee to the socket.

Warning: After adjustments, all screws must be secured with a medium strength threadlocker and tightened with the correct torque.

Note: The proximal connector is shifted 6 mm in an anterior direction from the knee axis.

Note: If the alignment reference line (B) is too far anterior to the knee axis (A), the patient may require more effort to create a hyperextension moment that will disengage the stance control and allow for fluid initiation of the knee flexion.

Static Alignment

- Make sure the patient stands with equal weight on both legs.
- Check for correct prosthesis length.
- Check internal / external rotation.
- Check for correct load on toe and heel.

Dynamic Alignment

- Ensure the patient is familiar with the functioning of the device.

- Ensure that the patient can carry out the procedure for engaging the Swing Only mode and the Manual Lock mode, please refer to the following section: **Dynamic Alignment - U-Lever/Mode Selector Switch.**

Swing Control

The variable hydraulic resistance to flexion and extension during swing phase can be adjusted by rotating the cap atop the control cylinder, which is accessible on the posterior side of the device.

Flexion resistance (Fig. 2)

To adjust the flexion resistance, rotate the 11-point "•HYDRAULIK•" scale on top of the cap relative to the Indicator Mark (2A). Aligning the Indicator Mark with the following letters or dots on the scale sets the flexion resistance to:

- Minimum (2B), at the letter "H".
- Medium (2C), at the letter "A".
- Maximum (2D), at the letter "K".

Note down the selected flexion resistance setting.

Extension resistance (Fig. 3)

To engage the extension resistance adjustment, rotate the cap all the way in one direction until the Indicator Mark moves relative to the numeric scale on the side of the cylinder. Aligning the Indicator Mark with the following values on the numbered scale sets the extension resistance to:

- Minimum (3A), at 1.
- Medium (3B), at 5 - 6.
- Maximum (3C), at 10.

Caution: After setting the extension resistance, remember to re-set the flexion resistance to the previously selected setting.

Stance phase stability (Fig. 4)

The device's stance control system is designed to provide hydraulic yielding resistance to knee flexion when the prosthesis is loaded. This assists the user in descending stairs foot over foot and in stumble recovery. The stance resistance or yielding rate can be adjusted from the posterior side of the device using the Stance Control Adjustment Screw (4A) near the U-Lever (4B) on the piston rod. The Stance Control Adjustment Screw only adjusts 1/3 of a rotation:

- Counterclockwise rotation from the 12 o'clock position to the 8 o'clock position creates maximum yield rate (minimum stance resistance) (4C).
- Clockwise rotation from the 8 o'clock position to the 12 o'clock position creates minimum yield rate (maximum stance resistance) (4D).

Caution: Do not exceed the 1/3 rotation.

U-Lever/Mode Selector Switch (Fig. 5: A)

The device is designed to provide the user with a Manual Lock and a Swing Only mode of operation. These functions are engaged and disengaged through the U-Lever or mode selection switch located near the top of the piston rod on the device's posterior side.

U-Lever in the upward position (5B):

- Manual Lock or Swing Only mode is engaged.

U-Lever in the downward position (5C):

- Manual Lock or Swing Only mode is disengaged.

Procedure for engaging Swing Only mode:

1. Hyperextend the knee by putting weight on the tip of the prosthetic foot (**Fig. 6: A**).
2. Raise the U-Lever (5B).

The device's stance resistance is now eliminated.

Procedure for engaging Manual Lock mode:

1. Do not hyperextend the knee. Keep the knee in the desired locked position, either bent or straight. (**Fig. 6: B**).
2. Raise the U-Lever (5B).

The device is now locked against flexion.

USAGE

Cleaning and care

Wipe the device with a soft cloth. Do not use solvents.

Environmental Conditions

The device should not come in contact with fresh water, salt water or chlorinated water.

Caution: The device must be stored upright. Storing the device on its side will cause the cylinder to draw in air.

Caution: Device should not be used in dusty environment. Exposure to sand, talcum, or similar should be avoided.

The device can be used in temperatures between -15°C to 50°C.

MAINTENANCE

The device and the overall prosthesis should be examined by a healthcare professional. Interval should be determined based on patient activity. Recommended interval is every 6 months.

Check for damage, excessive wear and dirt.

Caution: Do not disassemble the device or attempt your own repairs. The cylinder cannot be replaced in the field. For cylinder replacement and/or repairs or if any problem should arise, send the device in for service.

Air can be bled from the device by engaging Swing Only mode and cycling the cylinder until all of the air has been eliminated.

Wipe the device with a soft cloth moistened with small amount of general purpose oil or sewing machine oil.

Caution: Do not use compressed air to clean device. Air forces pollutants into bearings and may cause malfunctions and wear.

Lubricate piston rod with silicone spray or similar.

Caution: Do not use talcum powder for lubrication.

REPORT OF SERIOUS INCIDENT

Any serious incident in relation to the device must be reported to the manufacturer and relevant authorities.

DISPOSAL

The device and packaging must be disposed of in accordance with respective local or national environmental regulations.

LIABILITY

Össur prosthetic devices are designed and verified to be safe and compatible in combination with each other and custom-made prosthetic sockets with Össur adapters, and when used in accordance with their intended use.

Össur does not assume liability for the following:

- Device not maintained as instructed by the instructions for use.
- Device assembled with components from other manufacturers.
- Device used outside of recommended use condition, application, or environment.

Compliance

This device has been tested according to ISO 10328 standard to three million load cycles.

Depending on patient activity, this may correspond to 3-5 years of use.

Mauch Knee

The load level specification is ISO10328 - P6 - 136 kg

Mauch Knee Plus

The load level specification is ISO10328 - P6 - 166 kg

ISO 10328 - "P" - "m"kg *)

*) Body mass limit not to be exceeded!

For specific conditions and limitations of use
see manufacturer's written instructions on
intended use!

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