












Instructions for Use

MIAMI LSO™

		3
EN	Instructions for Use	5
DE	Gebrauchsanweisung	7
FR	Notice d'utilisation	9
ES	Instrucciones para el uso	12
IT	Istruzioni per l'uso	14
NO	Bruksanvisning	17
DA	Brugsanvisning	19
SV	Bruksanvisning	21
EL	Οδηγίες Χρήσης	24
FI	Käyttöohjeet	26
NL	Gebruiksaanwijzing	28
PT	Instruções de Utilização	31
PL	Instrukcja użytkowania	33
CS	Návod k použití	35
TR	Kullanım Talimatları	38
JA	取扱説明書	40
ZH	中文说明书	42
KO	사용 설명서	44

ENGLISH

Medical Device

INTENDED USE

The device is intended to provide gross immobilization of the lumbar spine and/or compression

The device must be fitted and adjusted by a healthcare professional.

Indications for use

Indications that require gross immobilization of the lumbar region of the spine, this may include stable, non-displaced spinal fractures (L1-L5), spinal stenosis, herniated discs, degenerative spinal pathologies, spondylolisthesis, spondylolysis.

Contraindications

- Unstable, displaced fractures.

Warnings and Cautions:

- This is a supportive device only and is not intended or guaranteed to prevent spinal injury.
- The patient should be instructed how to tighten and loosen the compression system while at home. Do not overtighten the compression system to the point where it causes discomfort or difficulty breathing.
- For sitting, the patient may find it desirable to slightly loosen the compression system.

GENERAL SAFETY INSTRUCTIONS

The healthcare professional should inform the patient about everything in this document that is required for safe use of this device.

Any serious incident in relation to the device must be reported to the manufacturer and relevant authorities.

The patient should stop using the device and contact a healthcare professional:

- If there is a change or loss in device functionality, or if the device shows signs of damage or wear hindering its normal functions.
- If any pain, skin irritation, or unusual reaction occurs with the use of the device.

The device is for single patient – multiple use.

FITTING INSTRUCTIONS

Device Application

Össur's recommended application technique is with the patient in a supine position.

1. Before application, ensure the Belt Compression System (A) is stretched to its full width and remove the Overlap Closures (B) from both ends of the belt.
2. Secure Posterior Panel (C) to the inside of the Adjustable Belt (D) via Hook-and-Loop tabs and close the overlap flap (**Fig. 2**).
3. Log roll the patient onto their side. Position the Posterior Panel centered over the spine, with the bottom of the Posterior Panel at approximately the sacrococcygeal joint.

Note: Panels can be removed from their fabric sleeve and customized by heat-molding, grinding, or trimming, to achieve correct patient fit (**Fig. 3**).

4. Log roll the patient back into the supine position.
5. Wrap the Adjustable Belt around the waist and insert the Össur Fit Tool (E) (sold separately) (**Fig. 4**). Adjust the length of the belt arms to position the tool centrally. Secure the Belt Ends (F) to desired sizing (**Fig. 5**).
Note: If the patient's waist exceeds 127 cm (50"), an optional Belt Extension (G) can be added to the belt arms to reach a maximum belt circumference of 178 cm (70").
6. Remove Össur Fit Tool from belt (**Fig. 6**).
Note: If an Össur Fit Tool is not available, lay the belt flat with Overlap Closures removed and compression system loosened. Fold each belt arm equally to achieve the desired waist circumference and secure the Belt Ends (**Fig. 7**).
Follow steps 2–4.
7. Place Belt Ends into the Overlap Closures (**Fig. 8**).
8. Once Belt Ends are secured, the Overlap Closures will serve as belt ends (**Fig. 9**).
Note: To obtain the lowest profile fit, the Belt Ends can be trimmed if later adjustments are not required.
9. If required, the Anterior Panel (H) should be centered on the abdomen with the bottom edge just above the symphysis pubis, while still allowing the patient to sit comfortably (**Fig. 10**). The Anterior Panel will attach via Hook-and-Loop to the inside of the Overlap Closure.
10. Close the belt securing the Overlap Closures. Slide thumbs through the holes in the two Compression System Handles (I) and pull until the device is at the appropriate tightness (**Fig. 11**). Adjust the Compression System Cord (J) length, if necessary (**Fig. 12**).

Device Adjustments

The Belt Compression System also allows adjustments by detaching the Belt Ends from the Hook-and-Loop Strip posteriorly (**Fig. 13**).

- Trimming the belts at an angle can give the opportunity to accommodate different anatomies.
- If trimming posteriorly, ensure the Belt Ends are retained and folded (**Fig. 8**).

To heat-mold the panels, use a heat gun (175°C/350°F).

Device Removal

1. Detach Compression System Handles slowly from the belt to loosen and re-attach them to the sides of the Adjustable Belt.
2. Detach the Overlap Closures and remove the device.
3. To ensure a proper fit, ensure the Belt Compression System is stretched to its full width before re-applying the device.

Accessories and Replacement Parts

Please refer to the Össur catalog for a list of available replacement parts or accessories.

USAGE

Cleaning and care

Washing the device with the soft goods detached allows for more thorough cleaning.

- Hand-wash using mild detergent and rinse thoroughly.
- Air dry.

Note: Do not machine-wash, tumble dry, iron, bleach, or wash with fabric softener.

Note: Avoid contact with salt water or chlorinated water. In case of contact, rinse with fresh water and air dry.

DISPOSAL

The device and packaging must be disposed of in accordance with respective local or national environmental regulations.

LIABILITY

Össur does not assume liability for the following:

- Device not maintained as instructed by the instructions for use.
- Device assembled with components from other manufacturers.
- Device used outside of recommended use condition, application, or environment.

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