












Instructions for Use


MIAMI JTO[®] THORACIC
EXTENSION

		3
EN	Instructions for Use	6
DE	Gebrauchsanweisung	9
FR	Notice d'utilisation	12
ES	Instrucciones para el uso	15
IT	Istruzioni per l'uso	18
NO	Bruksanvisning	21
DA	Brugsanvisning	24
SV	Bruksanvisning	27
EL	Οδηγίες Χρήσης	30
FI	Käyttöohjeet	33
NL	Gebruiksaanwijzing	36
PT	Instruções de Utilização	39
PL	Instrukcja użytkowania	42
CS	Návod k použití	45
ZH	中文说明书	48
SK	Indikácie použitia	51
SL	Indikácie pre použitie	54

ENGLISH

MD

Medical Device

MR

Magnetic Resonance (MR) safe

INTENDED USE

The device is intended to restrict movement in the low cervical and high thoracic spine

X-ray and CT lucent.

The device must be fitted and adjusted by a healthcare professional.

Indications for use

For conditions requiring gross immobilization of the low cervical and/or high-thoracic spine:

- Low C-Spine and high T-Spine precaution for trauma patients.
- Immobilization for pre and post low C-spine and high T-spine surgery.
- Other conditions requiring gross immobilization of the low-cervical and/or high-thoracic.

No known contraindications.

Warnings and Cautions:

Regular cleaning of the device and the skin beneath along with inspection of the skin for any signs of irritation is required to reduce the risk of skin ulceration. Bedridden patients are at an increased risk for skin ulceration.

The Chin Strut may be removed (with a healthcare professional's approval) for sleep or when the patient is supine.

Patients under C-Spine precautions require supine application of the device.

GENERAL SAFETY INSTRUCTIONS

The healthcare professional should inform the patient about everything in this document that is required for safe use of this device.

Any serious incident in relation to the device must be reported to the manufacturer and relevant authorities.

The patient should immediately contact a healthcare professional:

- If there is a change or loss in device functionality, or if the device shows signs of damage or wear hindering its normal functions.
- If any pain, skin irritation, or unusual reaction occurs with the use of the device.

The device is for single patient – multiple use.

FITTING INSTRUCTIONS

Size Selection

The device is designed exclusively for use with the Miami J Collar, Miami Jr collar (Size P3), and Miami J Select (Setting 1-5). Fits adult patients with a waist circumference greater than 24 inches / 61 cm.

Device Application

The device should be worn over light clothing.

Ensure that correct spinal protocol is followed.

1. Remove sternal pad from collar (**Fig. 1**) and Chin Strut from the Front Panel of the device (**Fig. 2**).

2. Squeeze collar to expose chin rivet. Attach Chin Strut. The connection on the top of the Chin Strut will slide around the chin rivet and snap into place (**Fig. 3**).
3. Attach feet of Chin Strut to the front panel in the lowest position.
 - The collar's sternal plate attaches to the hook and loop at the top of the front panel. Make sure the padding of the collar is behind the front panel padding (**Fig. 4**).
4. Apply the front panel to the patient (**Fig. 5**).
 - Adjust the height of the chin strut on the front panel to suit the patient's anatomy. Align the Rivet holes in the Chin Strut with corresponding holes on the front panel to ensure symmetrical positioning (**Fig. 6**).
 - Adjust the height of the front panel by pushing the button on the front plate and sliding the strut to the highest possible position, it should be above the costal margin (**Fig. 7**).
5. Slide the collar back piece behind the patient's neck and center it, noting the arrow on the panel pointing up. Ensure the padding extends beyond the edge of the plastic. **Note:** Long hair should be placed outside of the plastic.
6. Flare the sides of the collar's front piece out, slide it up the chest wall and scoop it up under the chin. The sides of the front piece should be oriented up, off the trapezius, and toward the ears (**Fig. 8**).
7. While holding the front piece securely, curl the ends snugly against the patient's neck. With the front piece of the device inside the back piece, apply the Hook and Loop Strap and secure the opposite side of the device in the same fashion (**Fig. 9A**). Tighten straps alternately to an equal length on both sides. Hook and Loop Straps must be aligned to the front adhesive sections. Excess may be trimmed.
8. Log roll the patient onto their side (**Fig. 9B**).
9. Position the back panel by centering it over the spine, with the top edge of the panel approximately at T1/T2 (**Fig. 10A**). Proper positioning of the back panel will allow the scapula to move freely.
10. Log roll the patient back into the supine position (**Fig. 10B**).
11. Attach the front panel to the back panel using the waist buckles. Tighten the waist straps equally to center the front panel (**Fig. 11, 12A, 12B**). Excess may be trimmed.
12. Attach the shoulder buckles and tighten the shoulder straps equally (**Fig. 13**). Excess may be trimmed.
13. Re-adjust all straps as necessary. When the final fit has been confirmed. Attach the rivets on both sides to secure the chin strut to the front panel (**Fig. 14**).

Final Fitting Checklist

A properly applied device will look like (**Fig. 15**).

- Chin is centered comfortably in Chin Support. Chin should not extend over the edge of the pad, nor fall inside the collar.
- Hook and Loop Straps are equal length and aligned to front adhesive sections
- Front of collar angled up toward ears. The lower plastic edge should not be resting on the patient's clavicles nor digging into the trapezius.
- Front and back panels are centered on the body.
- No plastic touching skin. Pads extend beyond all plastic edges.
- No slack or gaps in any of the straps.
- Chin Strut set at lowest permissible treatment position.

Device Removal

Ensure that correct spinal protocol is followed.

Undo the buckles of the front panel and straps of collars front piece, remove the front panel and slide the back piece out from behind the head. Removing the back panel may require a log roll.

Accessories and Replacement Parts

Please refer to the Össur catalog for a list of available replacement parts or accessories.

USAGE

Cleaning and care

- Wipe with a damp cloth and mild soap.
- Air dry.

Note: Do not machine-wash, tumble dry, iron, bleach, or wash with fabric softener.

Note: If the device is used in salt water or chlorinated water, rinse with fresh water and air dry.

DISPOSAL

The device and packaging must be disposed of in accordance with respective local or national environmental regulations.

LIABILITY

Össur does not assume liability for the following:

- Device not maintained as instructed by the instructions for use.
- Device assembled with components from other manufacturers.
- Device used outside of recommended use condition, application, or environment.

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