











Instructions for Use

MIAMI J® SELECT


		3
EN	Instructions for Use	6
DE	Gebrauchsanweisung	10
FR	Notice d'utilisation	14
ES	Instrucciones para el uso	18
IT	Istruzioni per l'uso	23
NO	Bruksanvisning	27
DA	Brugsanvisning	31
SV	Bruksanvisning	35
EL	Οδηγίες Χρήσης	38
FI	Käyttöohjeet	43
NL	Gebruiksaanwijzing	47
PT	Instruções de Utilização	51
PL	Instrukcja użytkowania	56
TR	Kullanım Talimatları	60
RU	Инструкция по использованию	64
JA	取扱説明書	68
ZH	中文说明书	72
KO	사용 설명서	75
SL	Indikácie pre použitie	79
RO	Indicații de utilizare	83

12

a


b


c


d


e


f

g

h

ENGLISH

SYMBOLS

Medical Device

Magnetic Resonance (MR) safe

Product overview (Fig. 1):

- a. Front
- b. Height Indicator Marks
- c. Chin Support
- d. Hook Landing Area
- e. Tracheal Opening
- f. Patient Compliance Lockout
- g. Height Adjustment Button
- h. Sternum Relief Knob
- i. Sternum Contact
- j. Back Panel
- k. Repeatable Fit Tabs (sold separately)
- l. Strap

INTENDED USE

The device is intended to provide gross immobilization to the cervical spine

X-ray and CT lucent.

The device must be fitted and adjusted by a healthcare professional.

Indications for use

Conditions requiring gross immobilization of the cervical spine. This may include:

- C-Spine precaution for trauma patients
- Immobilization for pre and post c-spine surgery
- Other conditions requiring gross immobilization of the mid-cervical spine

Contraindications

- Patients with a compromised airway or known spinal deformities such as ankylosing spondylitis.
- Patients with penetrating trauma injuries.

Warnings and Cautions:

Warning: If an unstable fracture is suspected or unknown, with or without a sustained trauma, ensure additional spinal precautions are implemented to immobilize the spine.

Warning: Use of a cervical collar may increase intracranial pressure (ICP) through jugular venous compression.

Warning: Use of a cervical collar may increase complexity of airway management.

Warning: Cervical spine immobilization, including use of a cervical collar, has been associated with:

- Impaired respiratory effort and forced expiratory volume
- Pneumonia
- Aspiration
- Worsening of existing cervical spine injury
- Severe neurological deterioration in patients with ankylosing spondylitis

- Triggering of non-compliance or agitation
- Delay to definitive treatment
- Impaired physical examination/secondary survey

Caution: Cervical spine immobilization, including use of a cervical collar, has been associated with:

- Increased pain and discomfort, which may lead to increased spinal movement
- Impaired swallowing

Caution: Perfumes and harsh cleaning agents should not be used on or under the collar, as they may compromise the integrity of materials.

Caution: Use of a cervical collar for cervical spine immobilization is not recommended in situations where patients are awake, alert, not intoxicated, without neck tenderness or pain and without abnormal sensory or motor findings on examination.

Caution: Regular cleaning of the collar, pads and the skin beneath along with inspection of the skin for any signs of irritation is required to reduce the risk of skin ulceration. Bedridden patients are at an increased risk for skin ulceration.

Precautions:

- At least two persons are needed for initial collar placement: one to maintain the patient's head and neck in proper alignment, the second to fit the collar.
- Collar removal or any adjustments to the collar should be done with physician permission only and according to the physician's instructions.
- Unless otherwise specified by the physician, the patient should not remove the collar except to wash under it and change the pads.
- The patient will need a second person to help when cleaning the collar.

These instructions do not supersede hospital protocol and/or direct orders of the patient's healthcare professional.

GENERAL SAFETY INSTRUCTIONS

The healthcare professional should inform the patient about everything in this document that is required for safe use of this device.

Any serious incident in relation to the device must be reported to the manufacturer and relevant authorities.

The patient should immediately contact a healthcare professional:

- If there is a change or loss in device functionality, or if the device shows signs of damage or wear hindering its normal functions.
- If any pain, skin irritation, or unusual reaction occurs with the use of the device.

The device is for single patient – multiple use.

FITTING INSTRUCTIONS

Size Selection

- Miami J Select is an adjustable collar suitable for most adult and adolescent patients; some patient phenotypes (**Fig. 2**) and anatomical variations may not be accommodated.
 - The Miami J Select is suitable for Stout (**Fig. 2b**), Tall (**Fig. 2d**), Regular (**Fig. 2e**), and Short (**Fig. 2f**) patient phenotypes.
 - The Miami J Select has not been tested in Super Short (**Fig. 2a**) and XSmall (**Fig. 2c**) patient phenotypes and should therefore not be used on these patients.

- Adjust to the height that fits most comfortably and maintains desired treatment position.
 - Height settings 2-5 will be used most frequently. Height settings 1, 6, and 7 will be required less frequently.
- Once patient can sit/stand following initial application, check spine is in desired treatment position and adjust height setting if necessary.
- To ensure proper sizing, it is a good idea to fit operative patients prior to surgery.
 - Collar size may be incorrectly identified in the presence of postoperative dressings and swelling.

Device Application

Ensure that correct spinal protocol is followed.

1. Ensure the padding extends beyond the edge of the plastic and the Sternum Relief Knob is oriented vertically as shown in **Fig. 4**.
 2. Slide the Back Panel behind the patient's neck and center it, noting arrow on panel pointing up (**Fig. 3**).
- Note:** Long hair should be placed outside of the plastic.
3. Slide the Front of collar up the chest wall until the bottom of the tracheal opening is at the level of sternal notch. Chin Support may not contact chin at this point (**Fig. 4**). Sides of the collar Front should be oriented up, off the trapezius, and toward the ears.
 4. Holding the collar against the chest with the bottom of the tracheal opening at the level of the sternal notch, depress the Height Adjustment Button and manually adjust the Chin Support to achieve the desired height (**Fig. 5**).
 5. When desired height is attained, release Height Adjustment Button to lock in place.
 6. Verify that desired height is achieved on both sides of the device. Depress Height Adjustment Button and make micro-adjustments as necessary.
 7. While holding the Front securely, place the sides of the collar Front within the sides of the Back Panel. Apply the Straps of the back panel to the Hook Landing Areas on the front. Tighten Straps alternately to an equal length on both sides. (**Fig. 6**)
 8. If further height adjustment is required, depress Height Adjustment Button, and manually move to position. Use the Height Indicator Marks to record patient's collar height.
 9. When desired position is attained, engage the Patient Compliance Lockout located behind the Height Adjustment Button by moving the lever to the left. (**Fig. 7**)
 10. Straps must be aligned with the Hook Landing Areas. When patient is properly fit, there should be equal amounts of excess straps overhanging the Hook Landing Area. These may be trimmed (**Fig. 8**).

Device Adjustments

1. Optional Repeatable Fit Tabs can be applied to trimmed straps to provide a contact point for the patient to don and doff the collar to a repeatable position. To use, thread the strap end through the tab and attach to the hook molded into the tab (**Fig. 9**).
2. The blue Sternum Relief Knob adds a degree of freedom for the patient so that activities such as eating or talking can be performed with greater comfort. The knob turns 90 degrees to a horizontal position. (**Fig. 10**)
3. Upon healthcare professionals recommendation, collar corners and edges may be modified or trimmed to relieve discomfort or pressure.

Final Fitting Checklist

A properly applied device will look like **Fig. 11**.

- Collar extends from mandible to just below sternal notch.
- Chin is centered comfortably in Chin Support. Chin should not extend over edge of Sorbatex™ pad, nor fall inside collar.
- Straps are aligned to the Hook landing area and are an equal length.
- Sides of Back should overlap sides of Front.
- Front of collar angled up toward ears. Lower plastic edge should not be resting on patient's clavicles nor digging into trapezius.
- No plastic touching skin. Sorbatex™ pads extend beyond all plastic edges.
- No slack or gaps in any of the straps.
- Tracheal opening and posterior vent are midline.
- Collar sits away from the neck. There should be a "finger" gap between the tracheal opening and the neck.

If the collar is fit too closely/tightly to the neck, size up to the next taller size.

Accessories and Replacement Parts

Please refer to the Össur catalog for a list of available replacement parts or accessories.

USAGE

Cleaning and care

Remove the pads from the plastic shells (**Fig. 12**)

- Hand-wash using mild detergent and rinse thoroughly.
- Wring out the excess water and squeeze in a towel. Lay flat to air dry.

Note: Do not machine-wash, tumble dry, iron, bleach, or wash with fabric softener.

Note: Avoid contact with salt water or chlorinated water. In case of contact, rinse with fresh water and air dry.

Device plastic can be washed with a damp cloth and mild soap.

- Replace the pads by attaching the grey/dull side to the hook on the inside of the plastic shells.

DISPOSAL

The device and packaging must be disposed of in accordance with respective local or national environmental regulations.

LIABILITY

Össur does not assume liability for the following:

- Device not maintained as instructed by the instructions for use.
- Device assembled with components from other manufacturers.
- Device used outside of recommended use condition, application, or environment.

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