




Instructions for Use


---


MIAMI LUMBAR<sup>®</sup> POSTEO

		3
EN	Instructions for Use	6
DE	Gebrauchsanweisung	9
FR	Notice d'utilisation	12
ES	Instrucciones para el uso	15
IT	Istruzioni per l'uso	18
NO	Bruksanvisning	21
DA	Brugsanvisning	24
SV	Bruksanvisning	27
FI	Käyttöohjeet	30
NL	Gebruiksaanwijzing	33
PT	Instruções de Utilização	36
PL	Instrukcje Użytkowania	39
JA	取扱説明書	42
ZH	中文说明书	45
KO	사용 설명서	48

## INTENDED USE

The device is intended to provide gross immobilization of the thoracic and lumbar spine. The device must be fitted and adjusted by a healthcare professional.

### Indications for use

- Back pain caused by stable vertebral compression fracture (VCF).
- Post-operative immobilization

### Contraindications

- Unstable, displaced fractures.

### Warnings and Cautions:

- Do not over tighten the compression system to the point where it causes discomfort or difficulty breathing. For sitting, the patient may find it desirable to slightly loosen the compression system.
- These instructions do not supersede the specific orders of your healthcare professional.
- This is a supportive device only and is not intended or guaranteed to prevent spinal injury.
- The device should be worn over light clothing.
- The device is not meant to be worn during sleep, nor while bathing.
- Do not alter the personalized configurations and settings created by your healthcare professional, unless instructed to do so.

## GENERAL SAFETY INSTRUCTIONS

The healthcare professional should inform the patient about everything in this document that is required for safe use of this device.

Any serious incident in relation to the device must be reported to the manufacturer and relevant authorities.

The patient should immediately contact a healthcare professional:

- If there is a change or loss in device functionality, or if the device shows signs of damage or wear hindering its normal functions.
- If any pain, skin irritation, or unusual reaction occurs with the use of the device.

The device is for single patient – multiple use.

## FITTING INSTRUCTIONS

### Size Selection

Measure the back length from the top of the sacrum up to a finger width above the bony protrusion at the base of the neck (C7 vertebrae) - follow the spinal curvature. Choose the appropriate size by following the guidance in the table below

Size	Back Length
X-Small	14.5" - 17" (37 - 44 cm)
Small	17" - 20" (44 - 51 cm)
Medium	19" - 22" (49 - 56 cm)
Large	21" - 24" (55 - 61 cm)

### Device Adjustments

**Adjusting Belt length:** Measure the patient's waist circumference. Following the guidelines on the chart below and the belt arms (B), trim both belt arms to the proper size (**Fig. 3**).

**Note:** Be careful not to cut through the cord on the opposite side of the brace while trimming.

Remove and discard the cut portion of the belt and reattach both belts using the belt attachments (C) (**Fig. 4**).

Please note that the attachment of the pulley panel hook and loop and the belt arms is not designed to be readjusted more than 50 times.

Size	Waist circumference
XS	26" - 30" (66 - 76 cm)
S	30" - 34" (76 - 86 cm)
M	34" - 38" (86 - 96 cm)
L	38" - 42" (96 - 106 cm)
XL	42" - 46" (106 - 116 cm)
2XL	46" - 50" (116 - 127 cm)
3XL	50" - 54" (127 - 137 cm)
4XL	54" - 58" (137 - 147 cm)
5XL	58" - 62" (147 - 157 cm)
6XL	62" - 66" (157 - 167 cm)

**Adjusting the back panel:** Hold the back panel (A) against the patient's back so that the lower edge of the device is level with the coccyx (Fig. 5). The upper edge should reach to approximately 1.18 – 1.97 inches / 3 – 5 cm below C7. The aluminium strut (D) contained inside the device is cold-formable and can be shaped without a tool. Start by shaping at the lordosis and then work upwards (Fig. 6).

#### **Device Application**

Before application, ensure the Belt Compression System is stretched to its full width.

Device should be applied while standing with back in corrected position.

1. Position the posterior panel centered over the spine, with the bottom of the posterior panel at approximately the sacrococcygeal joint.
2. Wrap the belt arms around the waist and secure with the overlapping closures (I) (Fig. 7).
3. Slide fingers through the holes in the two compression system handles (H) and pull outwards until the belt is at the appropriate tightness. Attach handles to the oval area of the overlapping closures (Fig. 8).
4. Remove the axillary strap handles from the external button hooks (J), place the axillary straps (E) over the shoulders (Fig. 9) and refasten to the internal button hooks (K) (Fig. 10).
5. Slide Axillary Strap pads downwards to access the alligator clips (M). Remove alligator clips, adjust the length of the straps pull the straps downwards until a snug fit is achieved (Fig. 11).

**Note:** The straps should not be too tight, digging into the skin or axilla. A finger width gap should be left. Trim the strap to the appropriate length, if needed, and replace alligator clip on strap. (Fig. 12).

6. Reposition the axillary strap pads (F) for comfort (Fig. 13).

#### **Final Fitting Checklist**

A properly applied device will look like (Fig. 14).

- The back portion should fit the patient's spinal curvature .
- The belt arms run snugly around top of the pelvis / waist .
- The axillary straps are flat against the body.
- Compression system handles for lower back support are placed in the front oval area.
- Axillary strap handles (L) are clipped on the internal button hooks.
- To avoid the possibility of the brace slipping, the belt arms need to sufficiently grip the patient's waist, and have a higher level of tightness than the shoulder straps.

#### **Device removal**

Clip the axillary strap handles back to the external button hooks (Fig. 15). Detach compression system handles slowly from the belt to loosen, and attach them at the original position – on the side of the belt arms (Fig. 16).

Detach the overlapping front closure (**Fig. 17**). Remove axillary straps from shoulders to remove the device from the body.

### ***Accessories and Replacement Parts***

Please refer to the Össur catalog for a list of available replacement parts or accessories.

### **USAGE**

#### ***Cleaning and care***

- Hand-wash using mild detergent and rinse thoroughly.
- Air dry.

**Note:** Do not machine-wash, tumble dry, iron, bleach, or wash with fabric softener.

**Note:** Avoid contact with salt water or chlorinated water. In case of contact, rinse with fresh water and air dry.

### **DISPOSAL**

The device and packaging must be disposed of in accordance with respective local or national environmental regulations.

### **LIABILITY**

Össur does not assume liability for the following:

- Device not maintained as instructed by the instructions for use.
- Device assembled with components from other manufacturers.
- Device used outside of recommended use condition, application, or environment.

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