


















Instructions for Use

MIAMI TLSO™ 464


		3
EN	Instructions for Use	6
DE	Gebrauchsanweisung	10
FR	Notice d'utilisation	14
ES	Instrucciones para el uso	18
IT	Istruzioni per l'uso	22
NO	Bruksanvisning	26
DA	Brugsanvisning	30
SV	Bruksanvisning	34
EL	Οδηγίες Χρήσης	38
FI	Käyttöohjeet	42
NL	Gebruiksaanwijzing	46
PT	Instruções de Utilização	50
PL	Instrukcja użytkowania	54
CS	Návod k použití	58
TR	Kullanım Talimatları	62
RU	Инструкция по использованию	66
JA	取扱説明書	70
ZH	中文说明书	74
KO	사용 설명서	78


19

20

Medical Device

PRODUCT DESCRIPTION

Miami TLSO 464 is a modular rigid thoraco-lumbo-sacral spine orthosis designed to support and provide motion restriction to a patient's spine during routine daily use.

The device is a non-invasive, adjustable brace for single patient use. Device length can be adjusted to patient morphology requirements. A pulley system is used to tighten the belt and belt panels to the patient.

Product Overview (Fig. 1)

- a. Yoke
- b. Strut
- c. Overlap Closure/Belt End
- d. Buckle
- e. PTE
- f. Belt
- g. Pectoral Disc
- h. ATE

INDICATIONS FOR USE

Indications that require gross thoracic and lumbar spine immobilization in all three planes, in which the following conditions may be present:

- Herniated disc
- Degenerative spinal pathologies
- Post-surgical stabilization
- Stable, non-displaced spinal fractures
- Spinal stenosis
- Spondylolisthesis
- Spondylolysis

CONTRAINDICATIONS

Unstable, displaced fractures.

GENERAL SAFETY INSTRUCTIONS / WARNINGS AND PRECAUTIONS

Warning: In cases of functional change or functional loss, the patient should stop using the device and contact a healthcare professional.

Warning: If the device shows signs of damage or wear hindering its normal functions, the patient should stop using the device and contact a healthcare professional.

- Do not over tighten the compression system to the point where it causes discomfort or difficulty breathing. For sitting, the patient may find it desirable to slightly loosen the compression system. If brace causes discomfort or irritation, contact your healthcare provider for further instructions and discontinue.
- These are supportive braces only and are not intended or guaranteed to prevent spinal injury. They should be applied by a licensed healthcare practitioner and as directed by a physician or other qualified medical.
- Do not machine wash or place in a tumble dryer.

- Össur accepts no responsibility for any injury which might occur while using the orthosis.
- The braces are not designed to be worn while bathing.
- The device is a non-invasive, adjustable brace for single patient use.

FITTING INSTRUCTIONS

BRACE PREPARATION

No preparation is needed for over-the-shoulder strapping configuration.

For axillary strapping configuration:

- Remove strap from the posterior thoracic extension (PTE)
- Reattach the free end of the strap to the strap anchor (**Fig. 2**)

For over-the-shoulder and axillary strapping configuration (requires purchase of additional kit):

- Remove pectoral disc assemblies from the ball joints on the yoke (**Fig. 3**)
- Attach secondary buckles from accessory kit over the ball joints into the groove (**Fig. 4**)
- Reattach the pectoral disc assemblies onto the ball joints
- Attach the free end of the strap from the accessory kit to the strap anchors (**Fig. 2**)

OVER-THE-SHOULDER BRACE APPLICATION

1. Before application, ensure the Belt Compression System is stretched to its full width and remove the Overlap Closures from both ends of the belt.
2. Log roll the patient onto their side. Omit this step if fitting in standing position. If standing, patient must stand with back straight and chest out.
3. Position the belt by centering the posterior panel (PTE) over the spine, with bottom of posterior panel approximately at the sacrococcygeal joint. (**Fig. 5**).
NOTE: Panels can be removed from their fabric sleeve and customized by heat-molding, grinding, or trimming, to achieve correct patient fit.
4. Log roll the patient back into the supine position. Omit this step if fitting in standing position.
5. Wrap the lumbar belt around the waist and through Össur Fit Tool (**Fig. 6**). Adjust length of belt arms to position tool centrally. Secure the belt ends to desired sizing (**Fig. 7**).
6. Remove Össur Fit Tool from belt.
NOTE: If an Össur Fit Tool is not available, lay the belt flat with Overlap Closures removed and compression system loosened. Fold each belt arm equally to achieve the desired waist circumference and secure the Belt Ends (**Fig. 8**). Follow steps 2-4.
NOTE: If patient's waist exceeds 50" (127 cm), an extension belt can be added to each belt arm to reach a maximum belt circumference of 70" (178 cm).
7. Place belt ends into the overlap closures (**Fig. 9**). To obtain the lowest profile fit, the belt arm ends can be trimmed if later adjustments are not required.
8. Place the anterior thoracic extension (ATE) onto the patient's front, centering the anterior panel on the abdomen with the bottom edge just above the symphysis pubis, while still allowing the patient to sit comfortably (**Fig. 10**).

- Adjust the height of the ATE by pushing the button on the anterior panel and sliding the strut to the appropriate length. The superior edge of the pectoral pads should be located approximately one inch below the sternal notch (**Fig. 11**).
- Using the provided tool, loosen the screw by one full turn and then adjust the angle on the hinged adjustment mechanism to accommodate sternal contours (**Fig. 12**).
NOTE: Make sure the hinged adjustment mechanism is tightened enough to lock it into position.
- Fasten the non-labeled belt end onto the ATE panel.
- Fasten the retaining straps onto the ATE panel soft good, over the non-labeled belt end (**Fig. 13**). Close the belt.
- Connect both quick-release buckles, tighten the shoulder straps, and position the shoulder pad for optimal comfort (**Fig. 14**).
- If desired, remove alligator, trim the strap, and replace alligator.
- Slide fingers through the holes in the two compression system handles and pull until the belt is at the appropriate tightness (**Fig. 17**). If necessary, adjust the compression system cord length (**Fig. 18**).
- The brace may require fine tuning for comfort, fit and stabilization when sitting, standing, or ambulating. This may include adjusting or tightening straps, compression system, rigid panels, and/or ATE/PTE.

AXILLARY STRAP ONLY OPTION / ADJUSTMENTS

- Replace step 13 with the following: fasten the strap anchor onto the belt, attach the axillary strap to the buckles and then tighten strap (**Fig. 16**).
- Resume steps 14 to 16.

DUAL STRAP OPTION / ADJUSTMENTS

- After step 13, attach axillary strap to secondary buckles (**Fig. 15**)
- Fasten the strap anchor onto the belt and then tighten strap (**Fig. 16**)
- Resume steps 14 to 16.

HINGE/LENGTH ADJUSTMENTS

If brace needs to be taller (**Fig. 19**):

- Remove hinge screw carefully with provided tool; do not lose the pre-load spring.
- Attach extension piece to lower stay.
- Attach yoke to extension piece.
- Using the provided tool, adjust the angle on the hinged adjustment mechanism to accommodate sternal contours.
NOTE: Make sure the hinged adjustment mechanism is tightened sufficiently to lock it into position.

BRACE REMOVAL

- Detach strap buckle(s) on right shoulder.
- Detach compression system handles and loosen. Fix them onto the belt after loosening.
- Detach the overlap closure and remove brace.
- Before reapplying the brace and to ensure a proper fit, be sure the compression system is stretched to its full width by pulling the ends of the belt apart. Reattach right strap buckle before reapplying the brace.

PANEL & BELT ADJUSTMENTS

- To heat mold the panels, use a heat gun (350°F/175°C).
- The belt arms can be trimmed anteriorly where the Belt Ends are folded, or posteriorly by detaching the belt arms from the Hook and Loop Strip adjacent to the Belt Compression System (**Fig. 20**). If trimming posteriorly, ensure the Belt Ends are retained and folded per **Fig 9**.

MAINTENANCE / PRODUCT CARE / CLEANING AND CARE

- Miami TLSO is designed for minimum maintenance and care
- The padded covers may be removed from the rigid panels for washing
- Hand-wash using mild soap and rinse thoroughly
- Air dry flat
- Do not machine wash or place in tumble dryer

REPORT A SERIOUS INCIDENT

Important notice to users and/or patients established in Europe:

The user and/or patient must report any serious incident that has occurred in relation to the device to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

DISPOSAL

The device and packaging must be disposed of in accordance with respective local or national environmental regulations.

LIABILITY

Össur does not assume liability for the following:

- Device not maintained as instructed by the instructions for use.
- Device assembled with components from other manufacturers.
- Device used outside of recommended use condition, application, or environment.

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