

## Instructions for Use




 **ÖSSUR**<sup>®</sup>  
LIFE WITHOUT LIMITATIONS

**Navii<sup>®</sup>**

3

|    |                             |     |
|----|-----------------------------|-----|
| EN | Instructions for Use        | 4   |
| DE | Gebrauchsanweisung          | 16  |
| FR | Notice d'utilisation        | 29  |
| ES | Instrucciones para el uso   | 42  |
| IT | Istruzioni per l'uso        | 55  |
| NO | Bruksanvisning              | 68  |
| DA | Brugsanvisning              | 80  |
| SV | Bruksanvisning              | 93  |
| EL | Οδηγίες Χρήσης              | 105 |
| FI | Käytöohjeet                 | 119 |
| NL | Gebruiksaanwijzing          | 132 |
| PT | Instruções de Utilização    | 145 |
| PL | Instrukcja użytkowania      | 158 |
| TR | Kullanım Talimatları        | 171 |
| RU | Инструкция по использованию | 183 |
| JA | 取扱説明書                       | 196 |
| ZH | 中文说明书                       | 209 |
| KO | 사용 설명서                      | 222 |
| CS | Návod k použití             | 234 |
| SK | Indikácie použitia          | 247 |
| RO | Indicații de utilizare      | 260 |
| HR | Upute za uporabu            | 273 |
| HU | Használati javallat         | 286 |
| BG | Указания за употреба        | 299 |
| SL | Indikácie pre použitie      | 312 |
| UA | Інструкція з використання   | 325 |

## DESCRIPTION

This document is in regard to Navii® (RKA01).

The device is a microprocessor-controlled prosthetic knee. The device provides structural support to the patient in the stance phase of gait activities and uses a magneto-rheological actuator to resist movement.

The device consists of the following main components (**Fig. 1**):

1. Proximal pyramid adapter
2. Magnetorheologic Actuator
3. Mechanical lock
4. Charging port
5. Power button and Charge status indicator
6. Device label
7. Distal pyramid adapter
8. Protective cover
9. Magnetic seam
10. Cover ridges
11. Cover nob
12. Recess

Device labels (**Fig. 2**) can be found at the rear of the device, on the power supply and on the internal face of the battery pack.

This document is for patients (intended device operator) and healthcare professionals.

## INTENDED USE

The device is intended as a part of a prosthetic system that replaces knee function of a missing lower limb.

Suitability of the device for the prosthesis and the patient must be evaluated by a healthcare professional.

The device must be fitted and adjusted by a healthcare professional.

### ***Indications For Use and Target Patient Population***

- Lower limb loss, amputation, or deficiency
- No known contraindications

The device is for low to high impact use, e.g., walking and occasional running.

The weight limit for the device is from 45 kg (99 lbs) to 136 kg (300 lbs) for low to moderate impact level and 110 kg (243 lbs) for high impact level.

## GENERAL SAFETY INSTRUCTIONS

**Warning:** Use of a lower limb prosthetic device carries an inherent risk of falling which may lead to injury.

The healthcare professional should inform the patient about everything in this document that is required for safe use of this device.

**Warning:** If there is a change or loss in device functionality, or if the device shows signs of damage or wear hindering its normal functions, the patient should stop using the device and contact a healthcare professional.

**Warning:** Avoid placing hands or fingers near moving joints.

**Warning:** A case of unexpected power failure or battery depletion during use will cause the device to stop functioning normally and may subject the patient to the risk of falling.

**Warning:** Unexpected power loss increases the risk of falling.

**Warning:** The device's misinterpretation of patient activity increases the risk of falling.

**Warning:** If device emits a warning signal, stop using device and refer to **Table 1** for more information. Note that certain errors can cause the device to shut down.

**Note:** In case of power failure, malfunction, or warnings it is recommended to engage the mechanical lock. When locked, the device can support full weight-bearing.

**Note:** High intensity uses or operation in high ambient temperature will cause the device actuator to heat up and become hot to the touch.

**Note:** Improper handling and/or adjustment of the device can cause malfunction. Patient should avoid:

- Direct impact to the device;
- Excessive shocks or vibrations;
- High-impact activities, sports, excessive loading, and heavy-duty use.

**Warning:** Not suitable for use in the presence of flammable anesthetic mixture with air, oxygen, or nitrous oxide.

**Note:** Do not attempt to modify device in any way.

**Warning:** Risk of structural failure. Components from other manufacturers have not been tested and may cause excessive load on the device.

The device is for single patient use.

## REQUIRED COMPONENTS

### Össur Logic App

For the initial fitting, the healthcare professional must use the Össur Logic app. It is available on Apple App Store. See Össur Logic on the App Store for mobile device compatibility.

The healthcare professional can use the app to adjust the knee to the user's physiology (weight, strength, etc.), gait style, and personal preferences.

The patient can use the Össur Logic app to check the status of the device, monitor the battery charge level, count steps, access training programs, and view instructions for use.

See instructions on how to connect to the app in Usage section.

## DEVICE SELECTION

See Össur catalog for a list of recommended components.

Use Össur titanium adapters to achieve waterproof performance.

**Note:** Adapter installation and set-screw tightening shall be performed according to instructions provided in their respective accompanying documents.

## ALIGNMENT INSTRUCTIONS

### Bench Alignment (Fig. 4)

#### Alignment Goal

Alignment reference line (B) should:

- pass through midpoint of socket at the ischial tuberosity level (D)
- pass 0-5 mm anterior or posterior to knee axis (A)
- fall at the 1/3 mark on the inside of the foot cover.

**Note:** Prioritize knee alignment over foot alignment if there is a mismatch.

#### Alignment Instructions

1. Position the foot so that the alignment reference line (B) falls at the 1/3 mark on the inside of the foot cover (with the foot cover and shoe on). Consider the external rotation of the foot.
2. Use the applicable adapters to connect the knee to the foot and establish the correct knee center height.
3. Position the knee so that the alignment reference line passes 0-5 mm anterior or posterior to knee axis (A).
4. On the lateral side of the socket, make a first mark at the midpoint of the socket at the ischial tuberosity level (D). Make a second mark at the midpoint of the socket distally (E). Draw a line through both marks.
5. Position the socket so the alignment reference line (B) goes through first mark at the midpoint of the socket at the ischial tuberosity level (D).
6. Adjust the socket flexion to 5° in addition to the existing position (i.e., hip flexion contracture) and set the height of the full prosthesis.
7. Use the applicable adapters to connect the knee to the socket.

**Caution:** At maximum flexion, ensure that a minimum distance of 3 mm (1/8") remains between the device and the socket (Fig. 5). If contact between device and socket cannot be avoided due to socket volume, make sure that the point of contact is on the device frame directly below the user interface panel. Keep the contact area flat and cushioned for pressure distribution.

#### Static Alignment

Turn the device on.

- Make sure the patient stands with equal weight on both legs.
- Check for correct prosthesis length.
- Check internal / external rotation.
- Check for correct load on toe and heel.

#### Dynamic Alignment

1. Keep device powered ON.
2. Ensure the patient is familiar with how the device functions.
3. Tell the patient to walk between parallel bars with usual gait. This lets the patient get used to the behavior of the device.
4. Adjust dynamic alignment as required.
  - Make sure that the motion in swing and stance phases stays in the line of progression.
  - Make sure to achieve step-length symmetry.

5. The Össur Logic App must be used for setting device parameters during dynamic alignment.

**Note:** Use the “Initial Setup” function in the Össur Logic app to adjust the device for basic use. Use the “Advanced Settings” function to fine tune the device.

If stability modifications are desired after patient assessment, the reference line can be shifted up to 5 mm (3/16") anterior or 5 mm (3/16") posterior to the knee axis.

**Caution:** As the alignment reference line is positioned more posterior to the knee center, the patient will experience more stance flexion during loading response. Additional voluntary control will be required to maintain knee stability when the device is powered OFF.

## After Alignment

After alignment, the healthcare professional should:

- Have the patient load the heel and allow the knee to flex while it is loaded. Ensure the patient experiences the resistance provided by the knee.
- Have the patient experience the roll-over of the foot. Encourage the patient to apply hip extension.
- Ask the patient to sit down using the stance resistance of the device.
- Show the function of the manual lock as this could be helpful when the knee is turned OFF.
- Ask the patient to walk both short and slightly longer distances under supervision to evaluate the setup of the device fully.

The healthcare professional should instruct the patient on:

- How to operate the device.
- How to understand the warning signals that the device may issue (see. Device Warning Signals section).
- How to recognize and activate the different states of the device (see. Supported Activities section).

## USAGE

### Device operation

#### Turning the Device On

1. Push and hold the power button for longer than 1 seconds.
2. The device will play three ascending tones and vibrate.
3. When the status indicator flashes green, the device is ready to use.

#### Turning the Device Off

1. Push and hold the power button for longer than 1 seconds.
2. The device will play three descending tones and turn off.

### Connecting to the mobile app

1. Turn on the device.
2. Open the app and follow the instructions.

**Note:** The required PIN code can be found on the device label at the back of the device.

### Battery

**Warning:** A damaged power supply can lead to electric shock.

**Warning:** Always disconnect the power supply from the device before donning.

**Caution:** The Battery Pack charge connector is Nickel-plated which may cause allergic reactions.

### Charging

1. Plug the power supply into the wall socket.
2. Plug the power supply into the charging port on the back of the device.
3. After a short while, the battery indicator will show the battery charge level with solid green lights and indicate charging with a flashing green light. If the power supply is disconnected during charging, the battery indicator may continue flashing for a short while.
4. After charging, disconnect the power supply from the charging port.
5. Disconnect the power supply from the wall socket to save energy.

Battery performance is affected by temperature below -10°C (14°F) or over 45°C (113°F), which may result in battery warnings or prevent device to perform as expected.

**Note:** It takes approximately 3 hours to charge the device battery pack when fully depleted.

**Note:** A fully charged battery pack can last for 48 to 72 hours, depending on activity level.

**Note:** The battery pack cannot be overcharged.

**Note:** The battery pack can only be charged in temperatures between 10°C and 35°C. If charge is attempted outside these limits, it will be delayed until the battery pack temperature is within the specified range.

## Battery Charge Level Indication

Push the battery indicator button. The green lights on the battery indicator show approximate battery level.

| Charge indicator | State of charge |
|------------------|-----------------|
|                  | 0% - 20%        |
|                  | 21% - 40%       |
|                  | 41% - 60%       |
|                  | 61% - 80%       |
|                  | 81% - 100%      |

## Device mechanical functions

### Mechanical Lock

To engage mechanical lock, press mechanical lock button from the left (**Fig. 6**). To unlock, press the lock button from the right.

### Protective Cover

The knee joint is fully functional without the protective cover. However, it is recommended to use the protective cover during general use to protect the knee joint from superficial scratches and clothing from wear and tear.

### Use in open water

It is recommended to remove the protective cover when using the device in open water (i. e. lake, river, ocean) to reduce risk of loss.

### Removing and attaching covers

The cover reaches around the device and is attached with a magnetic seam in the back.

To remove covers, pull the cover apart at the magnetic seam and remove from the knee (**Fig. 7a**).

To attach the cover, wrap it around the knee and close the magnetic seam. Ensure that the ridges inside the top part fit correctly into the grooves on the actuator (**Fig. 7b**) and cover knob fits into the recess.

### Supported activities

The device automatically adapts to activities such as walking, standing up and more. It has several different states that the patient must learn to recognize and activate.

### Standing State

The device is in standing state by default and it will revert to standing state if it does not recognize a movement pattern.

The device will provide support when the patient applies weight to the prosthesis. The knee will release/swing freely when the prosthesis is unloaded/no weight is applied to it. The patient can take small steps and turns in a natural and physiological way.

### Standing and Locking

The device can be locked in several positions. When locked, the device can support full loading.

To lock the device in full extension:

1. Press the mechanical lock button.
2. Fully extend the device.

To unlock the device:

1. Push the knee in hyperextension.
2. Press the mechanical lock button.

To lock the device in flexion at 10° or 20°:

1. Flex the device beyond the desired lock position and press the mechanical lock button.
2. Extend the device to the desired lock position.

## **Walking on Level Ground**

The device provides adaptive support in stance phase to ensure consistent support and controlled swing motion. The device automatically adapts to the patient's walking speed and style.

To walk on level ground:

1. Before the first step, put weight onto the device to trigger the standing state.
2. Take the first step with either the sound limb or prosthetic limb and continue walking at a preferred pace.

If using a cane or a crutch, the patient must take care to put weight on the device during stance to maintain support.

## **Sitting Down**

The device automatically detects sitting down motion and provides yielding resistance, allowing the patient to distribute weight equally to both legs and control speed while sitting down.

To sit down:

1. Stand in front of the chair.
2. Put weight on the device.
3. Lean back slightly and start flexing the knee until it yields.
4. Load the device and use the bodyweight to maintain a continuous knee flexion until seated.

Once in a sitting position, the device is unlocked and moves freely.

## **Standing Up**

To stand up from a seated position:

1. Ensure that the prosthetic foot is placed directly under the knee.
2. Apply weight onto prosthetic foot. Placing hands on knees helps to maintain equal weight distribution.
3. Begin to stand up. The device will allow smooth extension all the way to the standing state.

## **Walking Up Stairs**

**Caution:** Always use the banister or handrail when walking up stairs.

Leading with the prosthetic side:

1. Stop in front of the first step of the stairs.
2. Using residual limb muscles, raise the prosthesis up and position the foot flat on the step.
3. Put weight on the prosthesis. Squeeze the gluteal and residual limb muscles to extend the hip and move over the foot and knee.
4. Place the sound-side foot onto the next step and start extending the hip.
5. Use hip muscles on the prosthetic side to bring the prosthesis to the next step. The knee flexes to provide toe clearance.
6. Following swing extension, the device maintains the shank segment in vertical alignment to allow correct prosthetic foot placement on the next step.

Leading with the sound side:

1. Stop in front of the first step of the stairs.
2. Place the sound-side foot onto the first step, extending the hip on the prosthetic side.
3. Flex the hip on the prosthetic side to bring the prosthesis to the next step. The knee flexes to provide toe clearance.
4. Following swing extension, the device maintains the shank segment in vertical alignment to allow correct prosthetic foot placement on the next step.
5. Put weight on the prosthesis. Squeeze the gluteal and residual limb muscles to extend the hip and move over the foot and knee.

At the top of the stairs:

- a. If the prosthetic foot lands first, continue walking.
- b. If the sound foot lands first, wait for a moment and don't drive the residual limb upward like on the other steps. Gently flick the prosthesis backward to clear the step lip and then continue to step forward onto the prosthesis.

## **Walking Up Ramps**

Walking up ramps requires no change from level ground walking.

## **Walking Down Stairs and Ramps**

**Caution:** Always use the handrail when you walk down ramps and stairs.

The knee will give support while it flexes when walking down ramps and stairs.

To walk down stairs or ramps:

1. Take the first step down with the prosthesis.
2. Put weight on the prosthesis. Lean slightly backwards to flex the knee. The knee will flex and give support.
3. Lower the other foot onto the next step or ramp.
4. Continue to walk down the stairs or ramp. The knee will adapt the support to the walking speed.

**Note:** Leaning backwards keeps from falling forward if balance is lost.

**Note:** When walking down stairs, ensure to place half of the foot on the step to ensure adequate support.

## **Kneeling**

The device automatically detects kneeling motion. The device will provide support until the knee reaches the ground.

To kneel:

1. Put weight onto the prosthetic limb.
2. Take a step forward with the sound leg.
3. Lightly flex the hip on the prosthetic side and hold the pressure until the knee yields.
4. Lower to a kneeling position.

## **Cycling**

**Note:** Always rely on sound limb when getting on and off the bicycle.

1. Stand on sound limb, swing the prosthetic limb over the bicycle, and place the prosthetic foot on the pedal.
2. Push off with sound limb and start to pedal. Keep a minimal load on the device during the first two rotations to allow cycling mode activation. The device indicates that it is in cycling mode with one long beep.
3. Cycle with a steady and smooth rotation.
4. When stepping off the bike, stand on sound limb and swing the prosthetic limb over the bicycle. Fully extend the device to return to normal knee function, which is confirmed by one short beep.

## **Running**

Device function automatically adjusts to running whenever running gait is recognized.

## **Extension hold**

This function keeps the knee extended while lifted.

To enable Extension hold:

1. Flick the knee into extension by swinging it forward while snapping back with the residual limb.
2. The knee will remain extended for up to 3 seconds.

## **Device warning signals**

The device may give warnings in the form of audible beeps and vibration. In case of warnings, refer to **Table 1** for a list of possible causes and take action.

In case of warnings, stop using device and troubleshoot as follows:

1. Check battery status and recharge if necessary.
2. Check if the device is hot and allow to cool down.

If this does not help, contact your healthcare professional.

**Table 1 - Warning Signals**

| Warning Type       | Device State                     | Device Feedback   | Action  |
|--------------------|----------------------------------|---|---|
| Low battery        | Low battery (5% charge)          | Beeps 10 times every 7 seconds.<br>Vibration 3 times every 10 seconds.<br>Device status indicator turns red.<br>One battery indicator flashing green. | Charge battery.                               |
|                    | Critical low battery (2% charge) | Automatic shutdown.   | Charge battery.                               |
| Device temperature | Critical high device temperature | Beeps 10 times every 7 seconds.<br>Vibration 3 times every 10 seconds.<br>Device status indicator turns red.<br>Automatic shutdown after 2 minutes.   | Stop using device. Send to Össur for service. |
| Critical error     | Critical hardware error          | Beeps 10 times every 7 seconds.<br>Vibration 3 times every 10 seconds.<br>Device status indicator turns red.  | Stop using device. Send to Össur for service. |
|                    | Forced shutdown error            | Automatic shutdown.   | Stop using device. Send to Össur for service. |

**Automatic shutdown**

When the battery is depleted, the device encounters a critical error or high temperature, the device will emit a long steady beep with pulsing vibration for 30 seconds while the status indicator turns red.

The device will then emit three descending tones, vibrate three times, and shut down.

**Cleaning and care**

It is recommended to maintain the device in good outward condition by cleaning the surface of the device regularly.

Remove protective cover before cleaning (**Fig. 7**).

- Wipe with a damp cloth and mild soap.

**Caution:** Do not use solvents stronger than isopropyl alcohol as it might degrade the material.

**Caution:** Do not use compressed air to clean device.

**Note:** Device is not delivered in a sterilized state or intended to be sterilized.

**Environmental Conditions**

The device is Waterproof.

The device can be used in a wet or humid environment and submerged in up to 2-meter-deep water for a maximum of 30 minutes.

It can tolerate contact with: Salt water, chlorinated water, perspiration, urine, and mild soaps.

It can also tolerate occasional exposure to sand, dust, and dirt. Continuous exposure is not allowed.

**Note:** If sand or dirt gets into the device, range of motion may be affected.

**Cleaning after exposure to sand, salt water, or chlorinated water**

1. Turn the device OFF.
2. Remove the protective cover.
3. Rinse the device and the protective cover thoroughly with fresh water. Flex and extend the knee joint while rinsing to clean all surfaces. Engage and disengage the lock to remove all debris.
4. Dry the device and protective cover thoroughly with a soft cloth.

See which environmental conditions apply to the device in **Table 2**.

If the device is used outside the allowed temperature range, device performance may be unexpected and battery warnings may occur (See. **Table 1**). If the device is stored below 0°C (32°F) for an extended time, device will not perform as expected. When the device warms up, performance will be restored.

**Table 2 - Environmental Conditions**

|                      | Use                                  | Charging                           | Shipping                              | Extended storage                   |
|----------------------|--------------------------------------|------------------------------------|---------------------------------------|------------------------------------|
| Temperature          | -10 °C to 45 °C<br>(14 °F to 113 °F) | 10 °C to 35 °C<br>(50 °F to 95 °F) | -25 °C to 70 °C<br>(-13 °F to 158 °F) | 0 °C to 45 °C<br>(32 °F to 113 °F) |
| Relative humidity    | 0% to 100%                           | 0% to 100% (non-condensing)        | 10% to 100%                           | 10% to 90%                         |
| Atmospheric pressure | 700 hPa to 1060 hPa                  | 700 hPa to 1060 hPa                | 700 hPa to 1060 hPa                   | 700 hPa to 1060 hPa                |

## MAINTENANCE

The device and the overall prosthesis should be examined by a healthcare professional. Interval should be determined based on patient activity. Recommended interval is every 40 months.

Examine the device regularly for signs of damage. If you find signs of damage, consult your healthcare professional.

**Warning:** Do not attempt device repairs. Risk of injury or damage to device. Send device to Össur for maintenance or repair.

Expected life duration under normal usage is 6 million steps. Depending on patient activity, this corresponds to a duration of use of 4 to 6 years.

### Re-Packaging for shipping

The device shall always be shipped in the original packaging.

**Note:** If the device is shipped without the cover, the extra foam piece included in the packaging shall be placed at the bottom of the device to prevent damage during shipping.

## SPECIFICATIONS

**Table 3 - Device Specifications**

|   |  |
|---|--|
| Device weight   | 1.8 kg (4.0 lbs)   |
| Device build height (Fig. 3)                                  | 236 mm (9.29")   |
| Distance from actuator center to proximal connection (Fig. 3) | 35 mm (1.38)   |
| Device width (Fig. 3)   | 80 mm (3.15)   |
| Battery pack specification                                    | Lithium-Ion / 3500 mAh / 50.89 Wh                        |
| Energy to charge the battery pack                             | 141 VAh at 230 V AC, 50 Hz<br>119 VAh at 110 V AC, 60 Hz |
| Primary packaging weight                                      | 2.4 kg (5.3 lbs)   |
| Primary packaging materials                                   | Polypropylene (PP 5) plastic and Polyurethane Foam       |

**Caution:** Only use the provided power supply. Do not use a different power supply.

**Table 4 - Power Supply Specifications**

|              |                                  |
|--------------|----------------------------------|
| Manufacturer | FRIWO                            |
| Model no.    | FW8030M/24                       |
| Input        | 100–240 VAC, 50–60 Hz, 0.6–0.3 A |
| Output       | 24 VDC, 1.25A                    |

## ELECTROMAGNETIC COMPATIBILITY

**Warning:** Use of this device adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this device and the other equipment should be observed to verify that they are operating normally.

**Warning:** Use of accessories and cables other than those specified or provided by the manufacturer of this device could result in increased electromagnetic emissions or decreased electromagnetic immunity of this device and result in improper operation.

**Warning:** Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the device, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result. The device is suitable for use in the home healthcare environment except in flammable/explosive atmospheres, or where exposure to high electrical and/or magnetic fields can occur (e. g., electrical transformers, high-power radio/TV transmitters, RF surgical equipment, CT and MRI scanners).

The device may be susceptible to electromagnetic interference from portable and mobile RF communications equipment such as mobile (cellular) telephones or other equipment, even if that other equipment complies with CISPR EMISSION requirements.

EM-disturbances may lead to the device shutting down, resetting, freezing or adopting an erratic behavior.

**Table 5 - Electromagnetic Emissions Compliance Information**

| Emissions test                           | Standard      | Compliance level |
|--|---------------|------------------|
| Conducted and radiated RF emissions      | CISPR 11      | Group 1Class B   |
| Harmonics emissions                      | IEC 61000-3-2 | Class A          |
| Voltage fluctuations / Flicker emissions | IEC61000-3-3  | Complies         |

**Table 6 - Electromagnetic Immunity Compliance Information**

| <b>Immunity Test</b>                        | <b>Standard</b> | <b>Compliance Level</b>   |
|---|-----------------|---|
| Electrostatic discharge                     | IEC 61000-4-2   | ± 8 kV contact<br>± 15 kV air   |
| Electrical fast transient/burst             | IEC 61000-4-4   | ± 2 kV for power supply lines<br>± 1 kV for input/output lines  |
| Surge                                       | IEC 61000-4-5   | ± 1 kV line(s) to line(s)<br>± 2 kV line(s) to earth  |
| Voltage dips                                | IEC 61000-4-11  | 0 % UT; 0.5 cycle<br>At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315°                                     |
|   |                 | 0 % UT; 1 cycle and<br>70 % UT; 25/30 cycles<br>Single phase: at 0°                                       |
| Voltage interruptions                       | IEC 61000-4-11  | 0 % UT; 250/300 cycle   |
| Power frequency (50/60 Hz) magnetic field   | IEC 61000-4-8   | 30 A/m  |
| Conducted disturbances induced by RF fields | IEC 61000-4-6   | 3 V<br>0,15 MHz 80 MHz<br>6 V in ISM and amateur bands<br>between 0,15 MHz and 80 MHz<br>80 % AM at 1 kHz |
| Radiated RF EM fields                       | IEC 61000-4-3   | 10 V/m<br>80 MHz 2,7 GHz<br>80 % AM at 1 kHz  |
| Immunity to proximity magnetic fields       | IEC 61000-4-39  | 134.2 kHz @ 65A/m<br>3.56 MHz @ 7.5A/m<br>30 kHz @ 8A/m   |

**Table 7 - RF wireless communications Compliance Information**

| <b>Test frequency (MHz)</b> | <b>Band (MHz)</b> | <b>Service</b>   | <b>Modulation</b>       | <b>Immunity test level (V/m)</b> |
|-----------------------------|-------------------|--|-------------------------|----------------------------------|
| 385                         | 380-390           | TETRA 400  | Pulse modulation 18 Hz  | 27                               |
| 450                         | 430-470           | GMRS 460, FRS 460  | Pulse modulation 18 Hz  | 28                               |
| 710                         | 704-787           | LTE Band 13, 17  | Pulse modulation 217 Hz | 9                                |
| 745                         |                   |  |                         |                                  |
| 780                         |                   |  |                         |                                  |
| 810                         |                   |  |                         |                                  |
| 870                         | 800-960           | GSM 800/900,<br>TETRA 800, iDEN 820,<br>CDMA 850, LTE Band 5         | Pulse modulation 18 Hz  | 28                               |
| 930                         |                   |  |                         |                                  |
| 1720                        |                   |  |                         |                                  |
| 1845                        | 1700-1990         | GSM 1800, CDMA 1900,<br>GSM 1900, DECT, UMTS<br>LTE Band 1, 3, 4, 25 | Pulse modulation 217 Hz | 28                               |
| 1970                        |                   |  |                         |                                  |
| 2450                        | 2400-2570         | Bluetooth, WLAN,<br>802.11 b/g/n,<br>RFID 2450,<br>LTE Band 7        | Pulse modulation 217 Hz | 28                               |
| 5240                        | 5100-5800         | WLAN 802.11 a/n  | Pulse modulation 217 Hz | 9                                |
| 5500                        |                   |  |                         |                                  |
| 5785                        |                   |  |                         |                                  |

## Wireless modules regulated information

This device contains the following radio frequency transmitters:

**Table 8 - Wireless Modules**

| Model                                     | Regulatory certificates  | Type and frequency characteristics            | Effective radiated power |
|---|--|---|--------------------------|
| WiFi Module Model:<br>ATWINC1500- MR210PB | <p><b>FCC</b><br/>Contains FCC ID:<br/>2ADHKATWINC1500</p> <p><b>Canada</b><br/>Contains transmitter module IC:<br/>20266-WINC1500PB</p> <p><b>Japan</b><br/>Contains transmitter with<br/>certificate number<br/> 005-101762</p> <p><b>Korea</b><br/>Contains transmitter with<br/>certificate number<br/> R-CRM-mcp-WINC1500MR210P</p> <p><b>Taiwan</b><br/>Contains transmitter with<br/>certificate number<br/> CCAN18LP0321T2</p> <p><b>China</b><br/>CMIIT ID: 2018DJ1305</p> <p><b>Brazil</b><br/>Este product contém a placa Modelo ATWINC1500-MR210PB<br/>código de homologação ANATEL08497-18-08759</p> | 24122462 MHz (5 MHz<br>Interval, 11 Channels) | 4,72 mW/MHz              |
| Bluetooth 4.2 Module Model:<br>BMD-300    | <p><b>FCC</b><br/>Contains FCC ID: 2AA9B04</p> <p><b>Canada</b><br/>Contains IC: 12208A-04</p> <p><b>Japan</b><br/>Contains transmitter with<br/>certificate number<br/> R 210-106799</p> <p><b>Korea</b><br/>Contains transmitter with<br/>certificate number<br/> R-CRM-Rgd-BMD-300</p> <p><b>Brazil</b><br/>Contains ANATEL approved module # 00820-21-05903</p> <p><b>Mexico</b><br/>Este equipo contiene el módulo con IFT #: NYCE/CT/0146/17/TS</p>  | F1D 2402-2480 MHz                             | 0.35 mW/MHz              |

## USA-Federal Communications Commission (FCC)

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses, and can radiate radio frequency

energy. If not installed and used in accordance with the instructions, it may cause harmful interference to radio communications. However, there is no ensured specification that interference will not occur in a particular installation.

If this equipment does cause harmful interference to radio or television reception, which can be determined by tuning the equipment off and on, the user is encouraged to try and correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the distance between the equipment and the receiver.
- Connect the equipment to outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

Any changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

**Caution:** Exposure to Radio Frequency Radiation This device must not be co-located or operating in conjunction with any other antenna or transmitter.

### **Canada/Industry Canada (IC)**

This device complies with RSS 210 of Industry Canada.

Operation is subject to the following two conditions:

- this device may not cause interference, and
- this device must accept any interference, including interference that may cause undesired operation of this device.

**Caution:** Exposure to Radio Frequency Radiation The installer of this radio equipment must ensure that the antenna is located or pointed such that it does not emit RF field in excess of Health Canada limits for the general population; consult Safety Code 6, obtainable from Health Canada's website [www.hc-sc.gc.ca/rpb](http://www.hc-sc.gc.ca/rpb)

### **REPORT OF SERIOUS INCIDENT**

Any serious incident in relation to the device must be reported to the manufacturer and relevant authorities.

### **DISPOSAL**

The device and packaging must be disposed of in accordance with respective local or national environmental regulations.

### **LIABILITY**

Össur prosthetic devices are designed and verified to be safe and compatible in combination with each other and custom-made prosthetic sockets with Össur adapters, and when used in accordance with their intended use.

Össur does not assume liability for the following:

- Device not maintained as instructed by the instructions for use.
- Device assembled with components from other manufacturers.
- Device used outside of recommended use condition, application, or environment.

### **Compliance**


The device has been designed and tested according to the applicable international standards or in-house defined standards when no international standard exists or applies.

The device complies with Directive 2011/65/EU of the European Parliament and of the Council of 8 June 2011 on the restriction of the use of certain hazardous substances in electrical and electronic equipment (RoHS3).

The device complies with Directive 2014/53/EU of the European Parliament and of the Council of 16 April 2014 on the making available on the market of radio equipment.

This device has been tested according to ISO 10328 standard to three million load cycles.

The load level specification is ISO10328 – P6 – 136 kg

## SYMBOLS

| Symbol  | Description  |
|---|--|
|    | Medical Device   |
|    | Type B Applied Part  |
|     | Consult instructions for use   |
|    | Device contains electronic components and/or batteries that should not be disposed of in regular waste |
|     | Serial number  |
|     | See operating instructions   |
|    | Fragile, handle with care  |
|     | Keep Dry   |
|     | For indoor use only  |
|    | ON / OFF   |
|   | Alternating current  |
|   | Direct current   |
|   | QR Code containing device specific wireless communication PIN Code                                     |
|  | This way up  |
|   | Recyclable symbol  |

| Symbol   | Description   |
|--|---|
|   | Class II Equipment  |
|  | Energy Efficiency Level VI  |
|  | Manufacturer address  |
|  | Year of manufacture   |
| <b>IP68</b>  | Protected against harmful effects of dust and continuous immersion in water |

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