


ENGLISH

Medical Device

INTENDED USE

The device is intended for external support or stabilization of the knee.

Indications for use

For patients with unstable or injured knees that may benefit from additional knee support, such as: post-operation rehabilitation, tendonitis, bursitis, arthritis, sprains and strains.
No known contraindications.

Warnings and Cautions:

- Inspect the skin under the device's contact area before and after each use.
- The use of the device may increase the risk of deep vein thrombosis and pulmonary embolism.

GENERAL SAFETY INSTRUCTIONS

Read these instructions carefully before use. Keep them for future reference.

Any serious incident in relation to the device must be reported to the manufacturer and relevant authorities.

The patient should stop using the device and contact a healthcare professional:

- If there is a change or loss in device functionality, or if the device shows signs of damage or wear hindering its normal functions.
- If any pain, skin irritation, or unusual reaction occurs with the use of the device.

The device is for single patient – multiple use.

USAGE

Device Application

1. Open the device and place it behind the knee such that the wider section of the device faces upward.
2. With the leg in full extension wrap the device around the leg and fasten the contact closures over the front. Position the device, such that the kneecap opening and the middle of the hinge are centered on the kneecap (**Fig. 3**).
3. Insert the straps through the corresponding D-ring (**Fig. 4**) and fasten the straps at the calf and thigh (**Fig. 5**).
4. Stand up and verify the correct position and fit of the device. Adjust the device positioning and strap pressure, if needed.

Device removal

Undo the straps and contact closures and unwrap the device.

Cleaning and care

- Hand-wash using mild detergent and rinse thoroughly.
- Air dry.

Note: Do not machine-wash, tumble dry, iron, bleach, or wash with fabric softener.

Note: Avoid contact with salt water or chlorinated water. In case of contact, rinse with fresh water and air dry.

DISPOSAL

The device and packaging must be disposed of in accordance with respective local or national environmental regulations.

LIABILITY

Össur does not assume liability for the following:

- Device not maintained as instructed by the instructions for use.
- Device assembled with components from other manufacturers.
- Device used outside of recommended use condition, application, or environment.

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