Any serious incident in relation to the device must be reported to the manufacturer and relevant authorities.

The patient should stop using the device and contact a healthcare professional:

- If there is a change or loss in device functionality, or if the device shows signs of damage or wear hindering its normal functions.
- If any pain, skin irritation, or unusual reaction occurs with the use of the device.

The device is for single patient – multiple use.

FITTING INSTRUCTIONS

Device Application

1. Before application, ensure the Belt Compression System is stretched to its full width.
 2. If required, secure the Posterior Panel to the inside of the Adjustable Belt using the posterior panel attachments. The attachment straps on the belt should Hook-and-Loop directly into the panels.
 3. Log roll the patient onto their side. Position the Posterior Panel centered over the spine, with the bottom of the Posterior Panel at approximately the sacrococcygeal joint (**Fig. 1**).
- Note:** Panels can be removed from their fabric sleeve and customized by heat-molding, grinding, or trimming, to achieve correct patient fit.
4. Log roll the patient back into the supine position.
 5. If required, the Anterior Panel should be centered on the symphysis pubis, while still allowing the patient to sit comfortably. Attach the Anterior Panel to the inside of the adjustable belt using the anterior panel attachments.
 6. Wrap both belt arms around and to the front of the patient and secure with the overlapping closures.
 7. Slide thumbs through the holes in the two compression system handles and pull until the device is at the appropriate tightness (**Fig. 2**). Please note that the best placement for the handles on the belt front is within the oval area designated in the middle (**Fig. 3**).

Device Adjustments

To heat mold the panels, use a heat gun (350°F/175°C). The belt arms can be trimmed by detaching the belt arms from the Hook-and-Loop Strip posteriorly (**Fig. 4**). Trimming the belts at an angle can give the opportunity to accommodate different anatomies.

Device Removal

1. Detach compression system handles slowly from the belt to loosen and re-attach them to the sides of the Adjustable Belt.
2. Detach the Overlap Closures and remove the device.
3. To ensure a proper fit, ensure the Belt Compression System is stretched to its full width before re-applying the device.

Accessories and Replacement Parts

Please refer to the Össur catalog for a list of available replacement parts or accessories.

USAGE

Cleaning and care

Washing the device with the soft goods detached allows for more thorough cleaning.

- Hand-wash using mild detergent and rinse thoroughly.
- Air dry.

Note: Do not machine-wash, tumble dry, iron, bleach, or wash with fabric softener.

Note: Avoid contact with salt water or chlorinated water. In case of contact, rinse with fresh water and air dry.

DISPOSAL

The device and packaging must be disposed of in accordance with respective local or national environmental regulations.

LIABILITY

Össur does not assume liability for the following:

- Device not maintained as instructed by the instructions for use.
- Device assembled with components from other manufacturers.
- Device used outside of recommended use condition, application, or environment.

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MD Medical Device

INTENDED USE

The device is intended to provide gross immobilization of the lumbar spine and/or compression

The device must be fitted and adjusted by a healthcare professional.

Indications for use

Indications that require gross immobilization of the lumbar region of the spine, this may include stable, non-displaced spinal fractures (L1-L5), spinal stenosis, herniated discs, degenerative spinal pathologies, spondylolisthesis, spondylolysis.

Contraindications

- Unstable, displaced fractures.

Warnings and Cautions:

- This is a supportive device only and is not intended or guaranteed to prevent spinal injury.
- The patient should be instructed how to tighten and loosen the compression system while at home. Do not overtighten the compression system to the point where it causes discomfort or difficulty breathing.
- For sitting, the patient may find it desirable to slightly loosen the compression system.

GENERAL SAFETY INSTRUCTIONS

The healthcare professional should inform the patient about everything in this document that is required for safe use of this device.

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© Copyright Össur 2023-05-24 IFU 0600 1397_001 Rev. 7