









Instructions for Use

OCCIAN[®] BACK

		3
EN	Instructions for Use	5
DE	Gebrauchsanweisung	8
FR	Notice d'utilisation	12
ES	Instrucciones para el uso	16
IT	Istruzioni per l'uso	20
NO	Bruksanvisning	23
DA	Brugsanvisning	26
SV	Bruksanvisning	29
EL	Οδηγίες Χρήσης	32
FI	Käyttöohjeet	36
NL	Gebruiksaanwijzing	39

ENGLISH

SYMBOLS

Medical Device

Magnetic Resonance (MR) safe

DESCRIPTION

These Instructions for Use (IFU) provide additional information that enables the usage of an accessory with a particular medical device. Refer to the applicable IFU for the usage of the medical device.

INTENDED USE

The device is intended to provide gross immobilization to the cervical spine

The device must be fitted and adjusted by a healthcare professional. X-ray and CT lucent.

Indications for use

- Conditions requiring gross immobilization of the mid cervical spine
- Patients not awake and alert within 24-48 hours after admission
- Patients with pre-existing skin injuries and/or lacerations to the head/neck
- Burn patients who must wear a collar
- Elective surgery patients desiring nighttime post-op comfort

No known contraindications.

Warnings and Cautions:

Regular cleaning of the device and the skin beneath along with inspection of the skin for any signs of irritation is required to reduce the risk of skin ulceration. Bedridden patients are at an increased risk for skin ulceration.

GENERAL SAFETY INSTRUCTIONS

The healthcare professional should inform the patient about everything in this document that is required for safe use of this device.

Any serious incident in relation to the device must be reported to the manufacturer and relevant authorities.

The patient should immediately contact a healthcare professional:

- If there is a change or loss in device functionality, or if the device shows signs of damage or wear hindering its normal functions.
- If any pain, skin irritation, or unusual reaction occurs with the use of the device.

The device is for single patient – multiple use.

FITTING INSTRUCTIONS

Size Selection

The Occian® Back is a universally sized device to fit all adult Miami J® front collars as well as the Miami J® Select front collar and the Miami Jr® P3 front collar. The Occian® Back Jr. is designed to fit the Miami Jr® P2 front collar. Refer to the Miami J® Size Selection Instructions to select the appropriate size Miami J® collars.

Supine Application

Ensure that correct spinal protocol is followed throughout the application.

1. Remove any pillows from behind the patient's head. Position patient with arms to the side, shoulders down and head centrally aligned (**Fig. 1**).
2. Remove the original collar back piece and slide the Occian Back piece behind the patient's neck and center it (**Fig. 2**). Ensure the padding extends beyond all plastic edges of the collar and ensure long hair be placed outside of the Occian Back.

Note: Retain the original collar back for use when the patient becomes ambulatory.

3. Flare the sides of the collar front out, slide it up the chest wall and scoop it up under the chin (**Fig. 3**). The sides of the collar front should be oriented up, off the trapezius, and toward the ears. Place the plastic of the collar front in between the back piece padding and plastic; this will ensure that the padding of the back piece is against the skin.
4. While holding the front collar securely, curl the ends snugly against the patient's neck. Apply the hook-and-loop strap and secure the opposite side of the collar in the same fashion. Tighten straps alternately, one at a time, to an equal length on both sides (**Fig. 4**).
5. Hook-and-loop Straps must be aligned to the front collar adhesive sections. Excess may be trimmed (**Fig. 5**).

Note: Collar and Occian Back application must be secure to ensure proper fit, maintain alignment, and prevent the chin from slipping inside. If the patient can slip his/her chin inside the collar, it is a clear indication that it is not snug enough.

Device Removal

Ensure that correct spinal protocol is followed throughout the device removal.

Before removing the collar, note where the ends of the hook-and-loop straps are. When the collar is replaced, the straps should be in the same position. You may use a pen to mark the positions on the white plastic shell.

1. Position patient flat in bed without a pillow. Keep head in a neutral position. Do not allow the neck to bend forward, backward or from side to side. Position patient with arms to the side, shoulders down and head aligned centrally (**Fig. 5**).
2. While maintaining c-spine precautions, gently slide the Occian Back piece from behind the patient's neck and then remove the collar front (**Fig. 8, 9**).

USAGE

Cleaning and care

Keeping the device and the skin beneath clean is an important part of treatment. Daily cleaning will help to prevent skin irritation. Inspect and assess the skin beneath the device for breakdown once per shift. Clean the device once per day or more often as needed for patients with excessive oral drainage or diaphoresis resulting in dampness of the padding.

With healthcare professional approval only:

1. Gently roll the patient onto his/her side. (Place a small pillow or folded towel under his/her cheek to prevent the head from tilting sideways) (**Fig. 6**).

2. Undo the hook-and-loop strap on the top side of the collar and fold the Occian Back down to clean and check the back of the patient's head for redness or irritation (**Fig. 7**). Visualize and feel skin, scalp and under hair for any irritation or breakdown. Make sure the Occian Back is replaced and securely re-fastened after skin check.
3. Refasten the collar and return the patient to a supine position.

Note: Having a beard may cause skin irritation. Redness on the chin may indicate improper use of a pillow. If permitted, pillows must always be placed behind the head and shoulders.

If unable to log-roll:

1. Undo the hook-and-loop straps. Remove the Occian Back, sliding out from under the neck if necessary (**Fig. 8**).
2. While maintaining c-spine precautions, remove the Miami J[®] front. Clean the collar front following the instructions in the Miami J[®] Instructions for Use.
3. Palpate the back of the head and wash the front and back of the neck with mild facial soap and water and completely dry the skin. Observe for any redness or irritation under the collar, especially over bony areas like the chin, or clavicles.

The Occian Collar Back / Occian Jr. may be easily cleaned now.

Peel the soiled Occian pad off and perform the following on both the pad and the outer shell (**Fig. 10,11**):

- Wipe with a damp cloth and mild soap.
- Air dry.

Note: Do not machine-wash, tumble dry, iron, bleach, or wash with fabric softener.

Note: Avoid contact with salt water or chlorinated water. In case of contact, rinse with fresh water and air dry.

Dry the collar and the surface of the pad thoroughly with a towel and re-attach the pad. Fold the pads in half with the logo printed side next to the hook-and-loop; center the pad in the white shell (**Fig. 12**).

DISPOSAL

The device and packaging must be disposed of in accordance with respective local or national environmental regulations.

LIABILITY

Össur does not assume liability for the following:

- Device not maintained as instructed by the instructions for use.
- Device assembled with components from other manufacturers.
- Device used outside of recommended use condition, application, or environment.

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