




Instructions for Use

Philadelphia[®] Tracheotomy Collar

		3
EN	Instructions for Use	4
DE	Gebrauchsanweisung	6
FR	Notice d'utilisation	9
ES	Instrucciones para el uso	12
IT	Istruzioni per l'uso	14
NO	Bruksanvisning	17
DA	Brugsanvisning	19
SV	Bruksanvisning	22
EL	Οδηγίες Χρήσης	25
FI	Käyttöohjeet	27
NL	Gebruiksaanwijzing	30
PT	Instruções de Utilização	33
PL	Instrukcja użytkowania	35
JA	取扱説明書	38
ZH	中文说明书	41
KO	사용 설명서	43
RO	Indicații de utilizare	45
HR	Upute za uporabu	48

SYMBOLS

Medical Device

Magnetic Resonance (MR) safe

INTENDED USE

The device is intended to provide gross immobilization to the cervical spine

The device must be fitted and adjusted by a healthcare professional.

X-ray and CT lucent.

Indications for use

Conditions requiring gross immobilization of the cervical spine. This may include:

- C-Spine precaution for trauma patients
- Immobilization for pre and post c-spine surgery
- Other conditions requiring gross immobilization of the mid-cervical spine

Contraindications

- Patients with a compromised airway or known spinal deformities such as ankylosing spondylitis.
- Patients with penetrating trauma injuries.

Warnings and Cautions:

Warning: If an unstable fracture is suspected or unknown, with or without a sustained trauma, ensure additional spinal precautions are implemented to immobilize the spine.

Warning: Use of a cervical collar may increase intracranial pressure (ICP) through jugular venous compression.

Warning: Use of a cervical collar may increase complexity of airway management.

Warning: Cervical spine immobilization, including use of a cervical collar, has been associated with:

- Impaired respiratory effort and forced expiratory volume
- Pneumonia
- Aspiration
- Worsening of existing cervical spine injury
- Severe neurological deterioration in patients with ankylosing spondylitis
- Triggering of non-compliance or agitation
- Delay to definitive treatment
- Impaired physical examination/secondary survey

Caution: Cervical spine immobilization, including use of a cervical collar, has been associated with:

- Increased pain and discomfort, which may lead to increased spinal movement
- Impaired swallowing

Caution: Perfumes and harsh cleaning agents should not be used on or under the collar, as they may compromise the integrity of materials.

Caution: Use of a cervical collar for cervical spine immobilization is not recommended in situations where patients are awake, alert, not intoxicated, without neck tenderness or pain and without abnormal sensory or motor findings on examination.

Caution: Regular cleaning of the collar, pads and the skin beneath along with inspection of the skin for any signs of irritation is required to reduce the risk of skin ulceration. Bedridden patients are at an increased risk for skin ulceration.

Precautions:

- At least two persons are needed for initial collar placement: one to maintain the patient's head and neck in proper alignment, the second to fit the collar.

- Collar removal or any adjustments to the collar should be done with physician permission only and according to the physician's instructions.
- Unless otherwise specified by the physician, the patient should not remove the collar except to wash under it and change the pads.
- The patient will need a second person to help when cleaning the collar.

These instructions do not supersede hospital protocol and/or direct orders of the patient's healthcare professional.

GENERAL SAFETY INSTRUCTIONS

The healthcare professional should inform the patient about everything in this document that is required for safe use of this device.

Any serious incident in relation to the device must be reported to the manufacturer and relevant authorities.

The patient should immediately contact a healthcare professional:

- If there is a change or loss in device functionality, or if the device shows signs of damage or wear hindering its normal functions.
- If any pain, skin irritation, or unusual reaction occurs with the use of the device.

The device is for single patient – multiple use.

FITTING INSTRUCTIONS

Size Selection

Proper sizing is critical to achieve comfort and proper immobilization.

Ensure that correct spinal protocol is followed.

Height: Measure the vertical distance from the top of the patient's shoulder to the tip of the chin (**Fig. 1-a**). The measurement corresponds to the height of the device. If the patient's measurement falls between two consecutive sizes, apply the smaller size first.

Circumference: Measure around the patient's neck (**Fig. 2-b**). The measurement corresponds to the size circumference of the device.

Measurement A - Neck Height				
Infant	1.25" (3.2 cm)			
Pediatric	1.75" (4.4 cm)			
Adult	2.25" (5.7 cm)	3.25" (8.3 cm)	4.25" (10.8 cm)	5.25" (13.3 cm)

Measurement B - Neck Circumference	
Size	Circumference
Infant	6" - 8" (15.2 - 20.3 cm)
Pediatric	8" - 11" (20.3 - 27.9 cm)
Adult Small	10" - 13" (25.4 - 33.0 cm)
Adult Medium	13" - 16" (33.0 - 40.6 cm)
Adult Large	16" - 19" (40.6 - 48.3 cm)
Adult Extra Large	19"+ (48.3 cm+)

Supine Application

- Slide the back piece behind the patient's neck and center it, noting arrow on panel pointing up (**Fig. 3**).

Note: Long hair should be placed outside of the back piece.

- Apply the front piece with the chin secured in the recess and center it to secure neutral alignment. The front piece overlaps the back piece to ensure effective immobilization and comfort (**Fig. 4**). The front arrow should point upward.
- Apply the Hook and Loop strap on both sides. Tighten straps alternately to an equal length on both sides (**Fig. 5**). Excess may be trimmed.
- The Plastazote material can be trimmed or heat moulded as necessary depending on the patient's anatomy. Ensure that all edges are rounded off smoothly.

Note: Device application must be secure to ensure proper fit, maintain alignment, and prevent the chin from slipping inside.

Device Removal

Ensure that correct spinal protocol is followed.

Undo the straps, remove the front piece and slide the back piece out from behind the head.

Accessories and Replacement Parts

Please refer to the Össur catalog for a list of available replacement parts or accessories.

USAGE

Cleaning and care

- Wipe with a damp cloth and mild soap.
- Air dry.

Note: Do not machine-wash, tumble dry, iron, bleach, or wash with fabric softener.

Note: Avoid contact with salt water or chlorinated water. In case of contact, rinse with fresh water and air dry.

DISPOSAL

The device and packaging must be disposed of in accordance with respective local or national environmental regulations.

LIABILITY

Össur does not assume liability for the following:

- Device not maintained as instructed by the instructions for use.
- Device assembled with components from other manufacturers.
- Device used outside of recommended use condition, application, or environment.

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