ENGLISH

SYMBOLS

MD Medical Device

LATEX **Caution:** This device contains natural rubber latex which may cause allergic reactions.

The device is for single patient – multiple use.

USAGE

Size Selection

Measure the ankle circumference. The sizing chart is shown on the package. If between sizes always select the smaller size for optimal fit.

• If any pain, skin irritation, or unusual reaction occurs with the use of the device.

Remove the device if you experience pain, swelling, sensation changes, or if your extremity shows signs of insufficient blood flow (turns blue, white, or cold) while using this device, discontinue the use immediately and consult with a healthcare professional before using the device. If you experience pain, swelling, sensation changes, or if your extremity shows signs of insufficient blood flow (turns blue, white, or cold) while using this device, discontinue the use immediately and consult with a healthcare professional.

• The device should be used directly on intact skin and never over an open wound.

Warnings and Cautions:

- The device should not be used for extended periods of time and removed during long periods of rest, e.g., sleep.
- Patients with already compromised blood circulation in their extremities, like those with diabetes or peripheral vascular disease, need to be extra careful and are advised to consult with a healthcare professional before using the device. If you experience pain, swelling, sensation changes, or if your extremity shows signs of insufficient blood flow (turns blue, white, or cold) while using this device, discontinue the use immediately and consult with a healthcare professional.

• The device should be used directly on intact skin and never over an open wound.

GENERAL SAFETY INSTRUCTIONS

Read these instructions carefully before use. Keep them for your future reference. Any serious incident in relation to the device must be reported to the manufacturer and relevant authorities. The patient should stop using the device and contact a healthcare professional:

- If there is a change or loss in device functionality, or if the device shows signs of damage or wear hindering its normal functions.

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Ossur Deutschland GmbH