2. With use of this type of device, it is possible to receive various types of skin irritation, including rashes and pressure points. If you experience any unusual reactions while using this product, consult your health care provider immediately.
3. If you have any questions as to how this knee brace applies to your own specific knee condition, consult your health care provider.
4. If brace causes discomfort, irritation, or is damaged, discontinue use immediately and contact your health care provider.
5. This brace is single patient use only. If you do not agree to the above conditions, contact your brace provider before using this product.



WARNING: If there is a change or loss in device functionality or if the device shows signs of damage or wear hindering its normal functions, the patient should stop using the device and contact a healthcare professional.

FITTING INSTRUCTIONS

Brace Application

1. Unfasten straps #1, #2, #3 and #4. Loosen strap #5 (Figure 2).
2. With the knee bent in approximately 90 degrees, place the brace over your knee such that the patella is centered between the condyle pads. (Figure 2) **NOTE:** This places the Ossur hinge alignment guide (center of the hinge) slightly above the joint line so the brace can properly seat itself when you stand. (Figure 3)
3. For All straps that are trimmed or adjusted to length.
3. Fasten the lower cruciate strap just below the knee. Fasten the strap snugly above the calf muscle as this will aid in the brace suspension. (Figure 4)
4. Snugly fasten the bottom strap by inserting it through the D-ring.
5. Snugly fasten the upper cruciate strap above the knee.
6. Snugly fasten the top strap. **NOTE:** Over-tightening the top two straps can cause the brace to migrate.
7. Tightly fasten tibial strap #5, while making sure it fits comfortably. After fitting, the tibial strap will only need to be readjusted to improve brace comfort/suspension (Figure 5).

For more helpful information including fitting videos, scan the unique QR code located on the outside of the brace frame.

Brace Adjustments

Adjusting Pad and Strap Length: For comfort and comfort, ensure the strap pads are centered on the leg when tightening straps. If a strap pad is too long, remove the pad from the strap and use scissors to shorten the pad. To adjust the length of a strap, detach the alligator clip at the end of the strap. Adjust the length as necessary and reattach the alligator tab to the strap.

Adjusting the Brace Contouring:

The Rebound DUAL brace can be adjusted to individual anatomies. Use bending irons to cold form the aircraft grade aluminum frame to fit different thigh and calf shapes.

Adjusting the Varus / Valgus angle of the Brace:

The brace can be adjusted for individual varus or valgus angles and to create an offloading force to either the medial or lateral compartment.

- To adjust the angle, loosen the screws on the desired side of the brace and slide the upper frame to the desired position with a steady force.
- Adjust the medial or lateral upper frame to the desired position or adjust the lateral upper frame to offload the lateral compartment. The slide bar has a lettered scale and may be moved to increase or decrease the amount of varus or valgus force on the leg. (Figure 6). The lettered scale does not relate to degrees of offloading and should be used as a guide.
- Tighten the screws until snug (recommended torque value of 3 in-lb).
- All adjustments should be determined by a healthcare provider, in accordance with the patient's comfort level. Each time the brace is adjusted, the healthcare provider should note the letter for a record of progression.

Adjusting the Height of the Brace:

The brace has 0" extension stops attached to the hinges. If needed, range of motion can be adjusted using the extension and flexion stop kits.

- The brace height can be adjusted by loosening the screws on both sides of the upper frame and pulling the frame up with a steady force (Figure 6). Ensure the brace is properly aligned by checking that both hinge arms are located on the same position on the medial and lateral hinge arm.
- Retighten the screws until snug.

Range of Motion adjustment (ROM) (optional accessory, not included): The brace has 0" extension stops attached to the hinges. If needed, range of motion can be adjusted using the extension and flexion stop kits.

1. Remove the condyle pads from the medial and lateral hinges. Remove retaining screws.
2. Using a short Phillips head screwdriver, slide the 0" stop out of the hinge on both sides.
3. Re-apply the desired extension stop and gently secure it with the retaining screw.
4. Repeat steps to change the flexion stop.

NOTE: Do not overtighten retaining screw.

PCL Strap Attachment

- (optional accessory, not included): The PCL kit is designed to counteract posterior subluxation of the tibia. Trim and attach PCL strap for proper fit.
- Attaching PCL Strap – Acculign version
1. Note height position on medial and lateral hinge arms.
 2. Remove hook-and-loop fastener to access connection area.
 3. Remove Acculign hinge screws.
 4. Remove Acculign hinge d-ring.
 5. Position PCL Acculign hinge d-ring.
 6. Position PCL Acculign hinge inserts with raised portion into recessed pocket.
 7. Return frame to initial positions on hinge arms.
 8. Tighten Acculign screws.
 9. Reapply hook-and-loop fastener.
 10. Attach PCL strap.

ENGLISH

Medical Description

PRODUCT DESCRIPTION

Rebound DUAL Basic knee brace is a functional support for knee ligament instabilities.

List of Parts (Figure 1):

1. Top Strap
2. Upper Cruciate Strap
3. Lower Cruciate Strap
4. Bottom Strap
5. Anterior Tibial Strap
6. Upper Frame
7. Lower Frame
8. Acculign™ Hinge Arm
9. Acculign™ Hinge Arm Screw
10. D-Ring
11. Condyle pad

INDICATIONS FOR USE

The brace is indicated for combined and non-combined instabilities of the: Anterior cruciate ligament (ACL), Medial collateral ligament (MCL), Lateral collateral ligament (LCL), Posterior cruciate ligament (PCL); & conditions requiring unicompartmental load control.

CONTRAINDICATIONS

None known.

INTENDED USE

The device is intended for external support or stabilization of the knee.

GENERAL SAFETY INSTRUCTIONS

This knee brace is a supportive knee brace only. It is not intended or guaranteed to prevent knee injury.

The purchaser agrees to the following:

1. Ossur accepts no responsibility for any injury which might occur while using Ossur knee bracing products.

REPLACEMENT PARTS

- Condyle Pads
- Liner Kits
- Strap Kits

OPTIONAL ACCESSORIES

- Range-of-motion (ROM) Kit
- PCL Kit

TO PERFORM A SERIOUS INCIDENT

Important notice to users and/or patients established in Europe: The user and/or patient must report any serious incident that has occurred in relation to the device to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

DISPOSAL

The device and packaging must be disposed of in accordance with respective local national environmental regulations.

LIABILITY

Ossur does not assume liability if the device is not maintained as instructed by the instructions for use, if the device is assembled with components from other manufacturers, or if the device is used outside of recommended use condition, application, or environment.

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