
Instructions for Use




---

PROPRIO FOOT<sup>®</sup>

PSX02

		3
EN	Instructions for Use	4
DE	Gebrauchsanweisung	17
FR	Notice d'utilisation	32
ES	Instrucciones para el uso	46
IT	Istruzioni per l'uso	60
NO	Bruksanvisning	74
DA	Brugsanvisning	87
SV	Bruksanvisning	101
EL	Οδηγίες Χρήσης	114
FI	Käyttöohjeet	128
NL	Gebruiksaanwijzing	142
PT	Instruções de Utilização	156
PL	Instrukcja użytkowania	170
TR	Kullanım Talimatları	184
RU	Инструкция по использованию	197
JA	取扱説明書	211
ZH	中文说明书	224
KO	사용 설명서	237

## DESCRIPTION

This document is in regard to PROPRIO FOOT® (PSX02).

The device is a microprocessor controlled prosthetic ankle and foot system.

The device consists of the following components (**Fig. 1**):

1. Charging port;
2. Actuator;
3. Pyramid connector;
4. Power button;
5. Status indicator;
6. Unit Serial Number;
7. Carbon blades.

This device must be used with a Foot Cover and a Spectra Sock.

This document is for patients (intended device operator) and healthcare professionals.

## INTENDED USE

The device is intended as a part of a prosthetic system that replaces the foot and ankle function of a missing lower limb.

Suitability of the device for the prosthesis and the patient must be evaluated by a healthcare professional.

The device must be fitted and adjusted by a healthcare professional.

Healthcare professional must be authorized to fit the device by Össur after completing corresponding training.

### *Indications For Use and Target Patient Population*

- Unilateral transtibial, transfemoral, and knee disarticulation amputation;
- Bilateral transtibial, transfemoral, and knee disarticulation amputation.

No known contraindications.

The device is for low to moderate impact use, e.g., walking.

The weight limit for the device is 125 kg.

## GENERAL SAFETY INSTRUCTIONS

**Warning:** Use of a lower limb prosthetic device carries an inherent risk of falling which may lead to injury.

The healthcare professional should inform the patient about everything in this document that is required for safe use of this device.

**Warning:** If there is a change or loss in device functionality, or if the device shows signs of damage or wear hindering its normal functions, the patient should stop using the device and contact a healthcare professional.

**Warning:** Avoid placing hands or fingers near moving joints.

**Warning:** Not suitable for use in the presence of flammable anesthetic mixture with air, oxygen, or nitrous oxide.

**Warning:** In case of an unexpected power loss the device will behave like a mechanical foot. Note that the ankle alignment may not be optimal.

**Warning:** Unexpected power loss increases the risk of falling.

**Note:** Improper handling and/or adjustment of the device can cause malfunction. Patient should avoid:

- Direct impact to the device;
- Excessive shocks or vibrations;
- High-impact activities, sports, excessive loading, and heavy-duty use.

**Note:** Do not attempt to remove covers or modify device in any way.

**Caution:** Footwear that restricts movement of the ankle area e.g. due to height or stiffness of material will impair device function.

The device is for single patient use.

## REQUIRED COMPONENTS

The device is delivered with a power supply.

## Össur Logic App

For the initial fitting, the healthcare professional must use the Össur Logic app. It is available on Apple App Store. See Össur Logic on the App Store for mobile device compatibility.

The patient can use the Össur Logic app to connect to the device and access some functions.

## DEVICE SELECTION

See Össur catalog for a list of recommended components.

Use Össur titanium adapters to achieve waterproof performance.

**Note:** Adapter installation and set-screw tightening shall be performed according to instructions provided in their respective accompanying documents.

Verify that selected variant of the device is suitable for the impact level and weight limit according to the following table.

**Warning:** Do not exceed weight limit. Risk of device failure.

Incorrect category selection may also result in poor device function.

Table 1

Weight kg	45-52	53-59	60-68	69-77	78-88	89-100	101-116	117-125
Weight lbs	99-115	116-130	131-150	151-170	171-194	195-220	221-256	257-275
Low Impact Level	1	1	2	3	4	5	6	7
Moderate Impact Level	1	2	3	4	5	6	7	8

## ASSEMBLY INSTRUCTIONS

### Foot Cover and Spectra Sock

**Caution:** To avoid pinching fingers, always use a shoehorn.

1. Put the foot into the Spectra Sock.
2. Remove the foot cover attachment.
3. Use the applicable (straight) end of a shoehorn to put the foot with the Spectra Sock into the Foot Cover.
4. Move the shoehorn up to fully push the foot into the Foot Cover.
5. Put the foot cover attachment back onto the Foot Cover.
6. After alignment is complete, fix the Spectra Sock to the prosthesis to seal against dust and dirt.

**Note:** The Spectra Sock must be pulled up to prevent it from interfering with moving parts of the foot.

If required doff the Foot Cover as follows:

1. Remove the foot cover attachment.
2. Insert the applicable (edged) end of a shoehorn behind the foot.
3. Push the shoehorn down and pull the foot out of the Foot Cover.
4. Fully remove the Spectra Sock.

### Prosthesis

Assemble prosthesis with applicable devices.

**Warning:** Risk of structural failure. Components from other manufacturers have not been tested and may cause excessive load on the device.

**Warning:** Ensure proper attachment by following the applicable device assembly instructions.

**Caution:** Do not use heel wedges. Heel wedges will cause premature wear and tear.

## ALIGNMENT INSTRUCTIONS

Before alignment and setup the device should be powered ON and connected to the Össur Logic App.

### Initial Setup

Turn the device on.

After turning the device ON, two orange lights on the LED Indicator Ring will show initial set-up is required.

**The setup consists of five steps.**

1. Ankle Alignment
2. Bench Alignment

3. Dynamic Alignment
4. Auto-Adjustment
5. Advanced Set-Up

Device indicator shows the setup progress (**Table 2**).

**Ankle Alignment (without wearing prosthesis)**

Ankle Alignment adjusts the angle of the ankle to allow optimal posture with shoes of various heel heights.

**Note:** The device must be fitted with suitable footwear before performing Ankle Alignment.

**Note:** Heel height adjustment range is 0-50 mm (2"). If heel height exceeds 50 mm, the device will emit an audio signal and not adjust to height.

1. Fit the device with suitable footwear and ensure free ankle motion.
2. Place the device on a level surface; make sure device is fully-charged device turned ON.
3. Initiate Ankle Alignment by pressing the power button quickly twice successively or in the app.
4. The device will beep once confirming initiation, and the indicator will light up, showing the alignment progress (**Table 2**).
5. The actuator will move the ankle into appropriate position.
6. If successful, the device will emit an ascending audio tone, and the indicator will flash white. If unsuccessful, the indicator will flash yellow, showing process needs to be repeated.

**Bench Alignment (Fig. 3)**

**Alignment Goal**

Alignment reference line (B) should:

- pass through midpoint of socket at either the patellar tendon level or at the ischial tuberosity level (D).
- fall at the 1/3 mark on the inside of the foot cover.

**Note:** Prioritize knee alignment over foot alignment if there is a mismatch.

**Alignment Instructions**

1. Position the foot so that the alignment reference line (B) falls at the 1/3 mark on the inside of the foot cover (with the foot cover and shoe on). Consider the external rotation of the foot.
2. Use the applicable adapters to connect either the socket or the knee to the foot and establish the correct knee center height.
3. If using a prosthetic knee: Position knee according to knee alignment instructions
4. On the lateral side of the socket, make a first mark (alignment reference point) at the midpoint of the socket at either patellar tendon level or the ischial tuberosity level (D). Make a second mark at the midpoint of the socket distally (E). Draw a line through both marks.
5. Position the socket so the alignment reference line (B) goes through the alignment reference point (D).
6. Set the correct socket angles for flexion/extension and abduction/adduction.
7. If using a prosthetic knee: Use the applicable adapters to connect the socket to the prosthesis.

**Static Alignment**

- Make sure the patient stands with equal weight on both legs.
- Check for correct prosthesis length.
- Check internal / external rotation.
- Check for correct load on toe and heel.

**Dynamic Alignment**

**Note:** Alignment should aim to minimize medial-lateral motion during swing phase, as excessive side motion may affect the initiation of dorsiflexion during swing phase.

The heel to toe action can be influenced by:

- Heel Stiffness
- Anterior-Posterior positioning of device
- Dorsi-Plantarflexion
- Shoe Characteristics

Consider the following actions if needed:

#### Symptoms

- Device comes to flat position too early (patient feels he/she is sinking into a hole)
- Climbing over the toe requires extra energy
- Toe feels too stiff
- Knee hyperextends

#### Action

- Check Wedge instructions (if applicable)
- Shift socket anterior (or device posterior)
- Consider dorsiflexion
- Check heel of the shoe and shoe performance

#### Symptoms

- Rapid heel to toe movement
- Poor control over prosthesis at initial contact
- Minimal energy return feeling
- Too little push off from the toe
- Knee becomes unstable

#### Action

- Check Wedge instructions (if applicable)
- Shift socket posterior (or device anterior)
- Consider plantarflexion
- Check heel of the shoe and shoe performance

### Auto Adjustment

Auto adjustment calibrates the device to fit the patient's gait.

Before performing the Auto Adjustment

- Ensure the patient has access to a level ground walkway long enough to complete at least 15 adequate strides without interruption at a consistent self-selected speed.
- Ensure the walkway surface is completely flat. It is recommended to perform Auto Adjustment indoors.

**Note:** The device will not raise the toe in level walking during the Auto-Adjustment procedure.

### To perform Auto Adjustment

1. Ask the patient to stand still.
2. Enable Auto-Adjustment using the Össur Logic app.
3. A single white light on the status indicator will light up to confirm Auto-Adjustment has been initiated (**Table 2**)
4. Ask the patient to walk in a straight line on level ground at a self-selected walking speed, avoiding turns or walking in a semi-circle. The indicator will light up white in a clockwise direction, showing successful step registration.  
**Note:** If the walkway is not long enough for completing 15 adequate strides, the patient should be instructed to stop, turn around and then continue to walk in a straight line.
5. If successful, the device will emit an ascending audio tone, and the indicator will flash white. If unsuccessful, the indicator will flash yellow, showing process needs to be repeated.

### Troubleshooting

If the Auto Adjustment does not complete successfully:

- Ensure the alignment is correct.
- Ensure there are no obvious gait deviations, such as significant circumduction or rotation, as this may prevent successful completion of the Auto Adjustment procedure.
- Ensure the patient is walking at consistent speed.
- Ensure the patient is not walking in a circle or turning.
- If Auto Adjustment is not confirmed after 15 adequate strides, it must be repeated until successful. The device will not operate normally until the Auto Adjustment is successfully completed.
- Following changes to the alignment, the Auto Adjustment must be repeated.
- If LEDs are flashing yellow, the procedure must be repeated.

## Advanced Setup

Advanced functions like stair and ramp adaptation can be finetuned in the Össur Logic app.

Table 2 - Initial setup

Procedure	Indicator appearance	Device Status
Initial setup		Initial Setup required
		Auto-Adjustment required
		Heel Height Adjustment/Ankle Alignment required
Ankle Alignment Auto Adjustment		Procedure Initiation Confirmed
		Procedure in progress
		Procedure Completed
		Procedure Failed

## USAGE

### Device operation

#### Turning the Device On

1. Push and hold the power button for longer than 1 seconds.
2. The device will play three ascending tones and vibrate.
3. When the status indicator flashes green, the device is ready to use.

#### Turning the Device Off

1. Push and hold the power button for longer than 4 seconds.
2. The device will play three descending tones and turn off.

#### Connecting to the App

1. Turn on the device.
2. Open the app and follow the instructions in the app.

**Note:** The required PIN code can be found on the device label at the back of the device (see. Descriptions of Symbols section).

### Battery

**Warning:** A damaged power supply can lead to electric shock.

**Warning:** Always disconnect the power supply from the device before donning.

**Caution:** Only use the provided power supply. Do not use a different power supply.

**Caution:** The Battery Pack charge connector is Nickel-plated which may cause allergic reactions. The device can be charged whether it is turned ON or OFF.

## Charging

1. Plug the power supply into the wall socket.
2. Plug the power supply into the charging port on the back of the device.
3. After a short while, the battery indicator will show the battery charge level with solid green lights and indicate charging with a flashing green light. If the power supply is disconnected during charging, the battery indicator may continue flashing for a short while.
4. After charging, disconnect the power supply from the charging port.
5. Disconnect the power supply from the wall socket to save energy.

Device will turn ON automatically when the power supply is disconnected.

**Note:** It takes approximately 4 hours to charge the device battery pack when fully depleted.

Charge the battery daily for optimal performance.

**Note:** The battery pack cannot be overcharged.

**Note:** A fully charged battery pack can last for 24 to 48 hours, depending on activity level.

**Note:** The battery pack can only be charged in temperatures between 0°C and 45°C. If charge is attempted outside these limits, it will be delayed until the battery pack temperature is within the specified range.

## Battery Charge Level Indication

Push the battery indicator button. The green lights on the battery indicator show approximate battery level.

Charge Level	0% - 15%	16% - 29%	30% - 39%	40% - 49%	50% - 59%	60% - 69%	70% - 79%	80% - 89%	90% - 94%	95% - 100%
Indicator appearance										

## Ankle Alignment (while wearing prosthesis)

Ankle Alignment adjusts the angle of the ankle to allow optimal posture with shoes of various heel heights.

1. Patient should be in a seated position.
2. Ensure that the foot is positioned on a flat surface, with the sole from heel to toe on the ground.
3. Initiate Ankle Alignment by pressing the power button quickly twice successively or in the app.
4. The device will beep once confirming initiation, and the indicator will light up, showing the alignment progress (Table 2).
5. When device beeps again, lift the foot from the ground to allow for ankle adaptation.
6. The actuator will move the ankle into appropriate position.
7. If successful, the device will emit an ascending audio tone, and the indicator will flash white. If unsuccessful, the indicator will flash yellow, showing process needs to be repeated.

## Supported activities

If device is unable to move to the desired position, it will stop moving.

**Note:** The device can only move while unloaded in swing phase. Short swing phase duration (under 0.4 sec), heavy shoes, interrupted swing cycle, or foot contact with the ground can result in incomplete dorsiflexion movement.

## Walking on Level Ground

- The device automatically adapts to the patients walking speed and style.

**Note:** Circumduction or any abnormal side movement may hinder the initiation of dorsiflexion during swing phase.

## Inclined Ascent

- On inclines, the device gradually raises the toe to the degree of the slope.
- During swing phase, the device will raise the toe to provide additional ground clearance.

## Inclined Descent

- On declines, the device lowers the toe.

**Note:** It is not necessary to lead with a particular limb when ascending or descending slopes.

**Warning:** Ascending or descending inclines of more than 15° may trigger stair response.

### Walking Up Stairs

- During stair ascent, the device will raise the toe.
- Leading with the prosthetic side upstairs will provide stair adaptation after the second prosthetic step on the stairs.
- Leading with the sound side upstairs will provide stair adaptation after the first prosthetic step.
- When finishing the stair ascent and stepping onto level ground, it will take one stride for the device to adapt to the level ground position.

### Walking Down Stairs

- During stair descent, the device will raise the toe.
- Place the entire foot on the step.
- Leading with the prosthetic side when going down will provide stair adaptation after the second prosthetic step on the stairs.
- When reaching the bottom of the steps, it will take one stride for the device to adapt to the neutral level ground.

**Note:** If socket pressure is felt during the stair response, the patient should contact their healthcare professional, who may alter the stair response to decrease this socket pressure.

### Relax Mode

- In RELAX MODE the forefoot is fully lowered to the ground to enable a more relaxed position.
- To return device to a neutral state, tap the foot on the ground.
- RELAX MODE is disabled by default. It can be configured in the Össur Logic Mobile App.
- When sitting it is recommended to activate RELAX MODE.

### Chair Exit Mode

- When CHAIR EXIT MODE is triggered the device will move to dorsiflexion, enabling the patient to flex the knee and place the foot beneath the chair, providing more controlled transition from sitting to standing. The ankle will move back to neutral after the first step.
- CHAIR EXIT MODE function can be configured using Össur Logic Mobile App.
- If CHAIR EXIT MODE is enabled, exiting RELAX MODE will cause the device to move into CHAIR EXIT MODE.
- If RELAX MODE is disabled, but CHAIR EXIT MODE is enabled, the patient must first extend the leg (and virtually enter RELAX MODE) before being able to enter CHAIR EXIT MODE.
- If RELAX MODE is disabled, but CHAIR EXIT MODE is enabled, the patient must first extend the leg (and virtually enter RELAX MODE) before being able to enter CHAIR EXIT MODE.

### Standing State

- While standing, the device operates like a normal energy storing and returning foot.
- While standing, flexing the hip and placing the prosthetic foot in front of the sound foot will cause the device to enter RELAX MODE and move to plantar flexion.
- Extending the hip from the flexed position causes the ankle to move back to neutral position.

### Kneeling

- When kneeling, the device will enter RELAX MODE, allowing for a more relaxed position with less excessive knee flexion.

### Driving

When driving a vehicle, the device must be turned off.

**Caution:** Inadvertent movement of the ankle may interfere with the brake and/or accelerator pedals during driving.

### Device Feedback

See **Table 4** for a list of possible device feedback and what it means.

**Table 4 - Additional feedback**

Indicator appearance	Audio Feedback	Device Status
	N/A	<b>Normal</b> Device is ON and operational

Indicator appearance	Audio Feedback	Device Status
	3 beeps	<b>Low Battery</b> Battery charge level below 15%
	N/A	<b>Error</b> Restart required, return to Össur for servicing if error persists
	10 beeps	<b>Critical Error</b> Device non-operational, return to Össur for servicing
	2 ascending tones	<b>Connecting</b> Connecting to Össur Logic app
N/A	2 descending tones	<b>Disconnecting</b> Disconnecting from Össur Logic app
	N/A	<b>Connected</b> Device is ON and operational and connected to Össur Logic app

### Cleaning and care

It is recommended to maintain the device in good outward condition by cleaning the surface of the device regularly. Clean with a damp cloth and a mild soap. Dry with a cloth after cleaning.

**Caution:** Do not use solvents stronger than isopropyl alcohol as it might degrade the material.

**Caution:** Do not use compressed air or pressurized water jets to clean device.

**Note:** Device is not delivered in a sterilized state or intended to be sterilized.

### Environmental Conditions

The device can be used in a wet or humid environment and submerged in up to 1 meter-deep water for a maximum of 30 minutes.

It can tolerate contact with: Freshwater, perspiration, urine, and mild soaps.

It can also tolerate occasional exposure to salt water, sand, dust, and dirt. Continuous exposure is not allowed.

**Caution:** Risk of damage to device if exposed to chlorinated water.

Dry with a cloth after contact with fresh water or humidity.

Clean with fresh water after exposure to other liquids, chemicals, sand, dust, or dirt and dry with a cloth.

See which environmental conditions apply to the device in **Table 5**.

Battery performance is affected by temperature below -10°C (14°F) or over 45°C (113°F), which may result in battery warnings or prevent device to perform as expected.

**Note:** Prolonged exposure to temperatures below 0°C can result in device malfunction.

**Table 5 Environmental Conditions**

	Use	Charging	Shipping	Extended storage
<b>Temperature</b>	-10°C to 45°C (14°F to 113°F)	10°C to 40°C (50°F to 105°F)	-25°C to 70°C (-13°F to 158°F)	0°C to 45°C (32°F to 113°F)
<b>Relative humidity</b>	0% to 100%, condensing	0% to 100%, relative humidity (non-condensing)	10% to 100%, condensing	10% to 90%, non- condensing
<b>Atmospheric pressure</b>	700 hPa to 1060 hPa	700 hPa to 1060 hPa	700 hPa to 1060 hPa	700 hPa to 1060 hPa

### MAINTENANCE

The device and the overall prosthesis should be examined by a healthcare professional. Interval should be determined based on patient activity.

The patient should inspect the device regularly. If unusual movement or product wear is detected in the device at any time, the patient should be instructed to immediately discontinue use of the device and consult healthcare professional.

**Warning:** Do not attempt device repairs or remove covers. Risk of injury or damage to device. Send device to Össur for maintenance or repair.

**Noise from Foot**

Noise may occur if sand or debris is present in device. In that case, the healthcare professional should doff the foot, clean it with the help of compressed air and replace the Spectra sock if it is damaged.

**Note:** Do not use compressed air to clean the actuator housing.

**SPECIFICATIONS**

**Table 6 - Power Supply Specifications**

<b>Manufacturer</b>	FRIWO
<b>Model no.</b>	FW8030M/24
<b>Input</b>	100–240 VAC, 50–60 Hz, 0.6–0.3 A
<b>Output</b>	24 VDC, 1.25A

**Power supply specifications**

**Table 7 - Device Specifications**

<b>Device weight</b>	1.3 kg (2.9 lbs) (Size 24 Category 2) 1.5 kg (3.3 lbs) (Size 27 Category 5) 1.6 kg (3.5 lbs) (Size 30 Category 8)
<b>Device build height (Fig. 2)</b>	174 mm (6 7/8") (Size 24 Category 2) 180 mm (7 1/16") (Size 27 Category 5) 188mm (7 3/8") (Size 30 Category 8)
<b>Ankle range of motion</b>	33°
<b>Heel height adjustment range</b>	0-50 mm (Size 27 Category 5)
<b>Battery pack specification</b>	Lithium-Ion / 3500 mAh / 25.2Wh
<b>Battery pack weight</b>	521 g (1.15 lbs)
<b>Energy to charge the battery pack</b>	69.9 VAh at 230 VAC / 50 Hz 59.9 VAh at 110 VAC / 60 Hz
<b>Primary packaging weight</b>	3.1 kg (6.8 lbs)
<b>Primary packaging materials</b>	Polypropylene (PP 5) plastic and Polyurethane Foam

**Note:** The primary packaging is made from recyclable materials.

**ELECTROMAGNETIC COMPATIBILITY**

**Warning:** Use of this device adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this device and the other equipment should be observed to verify that they are operating normally.

**Warning:** Use of accessories and cables other than those specified or provided by the manufacturer of this device could result in increased electromagnetic emissions or decreased electromagnetic immunity of this device and result in improper operation.

**Warning:** Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the device, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.

The device is suitable for use in any environment except in flammable/explosive atmospheres, or where exposure to high electrical and/or magnetic fields can occur (e. g., electrical transformers, high-power radio/TV transmitters, RF surgical equipment, CT and MRI scanners).

The device may be susceptible to electromagnetic interference from portable and mobile RF communications equipment such as mobile (cellular) telephones or other equipment, even if that other equipment complies with CISPR EMISSION requirements.

**Table 8 - Electromagnetic Emissions Compliance Information**

Emissions test	Standard	Compliance level
Conducted and radiated RF emissions	CISPR 11	Group 1 – Class B
Harmonics emissions	IEC 61000-3-2	Class A
Voltage fluctuations / Flicker emissions	IEC61000-3-3	Complies

**Table 9 - Electromagnetic Immunity Compliance Information**

Immunity Test	Standard	Compliance Level
Electrostatic discharge	IEC 61000-4-2	± 8 kV contact ± 15 kV air
Electrical fast transient/burst	IEC 61000-4-4	± 2 kV for power supply lines ± 1 kV for input/output lines
Surge	IEC 61000-4-5	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth
Voltage dips	IEC 61000-4-11	0 % UT; 0.5 cycle At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315°
		0 % UT; 1 cycle and 70 % UT; 25/30 cycles Single phase: at 0°
Voltage interruptions	IEC 61000-4-11	0 % UT; 250/300 cycle
Power frequency (50/60 Hz) magnetic field	IEC 61000-4-8	30 A/m
Conducted disturbances induced by RF fields	IEC 61000-4-6	3 V 0,15 MHz – 80 MHz 6 V in ISM and amateur bands between 0,15 MHz and 80 MHz 80 % AM at 1 kHz
Radiated RF EM fields	IEC 61000-4-3	10 V/m 80 MHz – 2,7 GHz 80 % AM at 1 kHz

**Table 10 - RF wireless communications Compliance Information**

Test frequency (MHz)	Band (MHz)	Service	Modulation	Immunity test level (V/m)
385	380-390	TETRA 400	Pulse modulation 18 Hz	27
450	430-470	GMRS 460, FRS 460	Pulse modulation 18 Hz	28
710	704-787	LTE Band 13, 17	Pulse modulation 217 Hz	9
745				
780				
810	800-960	GSM 800/900, TETRA 800, iDEN 820, CDMA 850, LTE Band 5	Pulse modulation 18 Hz	28
870				
930				
1720	1700-1990	GSM 1800, CDMA 1900, GSM 1900, DECT, UMTS LTE Band 1, 3, 4, 25	Pulse modulation 217 Hz	28
1845				
1970				

Test frequency (MHz)	Band (MHz)	Service	Modulation	Immunity test level (V/m)
2450	2400-2570	Bluetooth, WLAN, 802.11 b/g/n, RFID 2450, LTE Band 7	Pulse modulation 217 Hz	28
5240	5100-5800	WLAN 802.11 a/n	Pulse modulation 217 Hz	9
5500				
5785				

### Wireless modules regulated information

This device contains the following radio frequency transmitters:

Table 11 - Wireless Modules

Model	Regulatory certificates	Type and frequency characteristics	Effective radiated power
WiFi Module Model: ATWINC1500- MR210PB	<b>FCC</b> Contains FCC ID: 2ADHKATWINC1500	2412 - 2462 MHz (5 MHz Interval, 11 Channels)	
	<b>Canada</b> Contains transmitter module IC: 20266-WINC1500PB		
Bluetooth 4.2 Module Model: BMD-300	<b>FCC</b> Contains FCC ID: 2AA9B04	FID 2402-2480 MHz	0.35 mW/MHz
	<b>Canada</b> Contains IC: 12208A-04		
	<b>Japan</b> Contains transmitter with certificate number  R210-106799		

### USA-Federal Communications Commission (FCC)

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses, and can radiate radio frequency energy. If not installed and used in accordance with the instructions, it may cause harmful interference to radio communications. However, there is no ensured specification that interference will not occur in a particular installation.

If this equipment does cause harmful interference to radio or television reception, which can be determined by tuning the equipment off and on, the user is encouraged to try and correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the distance between the equipment and the receiver.
- Connect the equipment to outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

Any changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

**Caution:** Exposure to Radio Frequency Radiation – This device must not be co-located or operating in conjunction with any other antenna or transmitter.

### Canada – Industry Canada (IC)

This device complies with RSS 210 of Industry Canada.

Operation is subject to the following two conditions:

- this device may not cause interference, and

- this device must accept any interference, including interference that may cause undesired operation of this device.

**Caution:** Exposure to Radio Frequency Radiation – The installer of this radio equipment must ensure that the antenna is located or pointed such that it does not emit RF field in excess of Health Canada limits for the general population; consult Safety Code 6, obtainable from Health Canada’s website [www.hc-sc.gc.ca/rpb](http://www.hc-sc.gc.ca/rpb)

**REPORT OF SERIOUS INCIDENT**

Any serious incident in relation to the device must be reported to the manufacturer and relevant authorities.

**DISPOSAL**

The device and packaging must be disposed of in accordance with respective local or national environmental regulations.

**LIABILITY**

Össur prosthetic devices are designed and verified to be safe and compatible in combination with each other and custom-made prosthetic sockets with Össur adapters, and when used in accordance with their intended use.

Össur does not assume liability for the following:

- Device not maintained as instructed by the instructions for use.
- Device assembled with components from other manufacturers.
- Device used outside of recommended use condition, application, or environment.

**Compliance**

The device has been designed and tested according to the applicable international standards or in-house defined standards when no international standard exists or applies.

Compatibility and compliance with these standards are achieved only when Össur products are used with other recommended Össur components.

The device has been tested by an independent accredited laboratory and found to comply to:

- EN60601-1:2006/A1:2013
- EN60601-1-2:2015

This device has been tested according to ISO 10328 standard to two million load cycles.

Depending on patient activity this may correspond to 2-3 years of use.

ISO 10328 - Label		
Category	Weight (Kg) Moderate Impact	Label Text
1	52	ISO10328 – P3* – 52 kg
2	59	ISO10328 – P3* – 59 kg
3	68	ISO10328 – P3 – 68 kg
4	77	ISO10328 – P4 – 77 kg
5	88	ISO10328 – P4 – 88 kg
6	100	ISO10328 – P5 – 100 kg
7	116	ISO10328 – P5 – 116 kg
8	125	ISO10328 – P5 – 125 kg
*Actual test loads reflect maximum body mass		

**ISO 10328 - “P” - “m”kg \*)** 

---

\*) Body mass limit not to be exceeded!

 For specific conditions and limitations of use see manufacturer’s written instructions on intended use!

## SYMBOLS

Medical Device

Caution Symbol

Consult instructions for use

See operating instructions

BF Applied Part

ON / OFF

Direct current

Device contains electronic components and/or batteries that should not be disposed of in regular waste

**IP67**

Protected against harmful effects of dust and temporary immersion in water

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