













Instructions for Use

PAPOOSE[®]

		3
EN	Instructions for Use	5
DE	Gebrauchsanweisung	8
FR	Notice d'utilisation	12
ES	Instrucciones para el uso	15
IT	Istruzioni per l'uso	19
NO	Bruksanvisning	23
DA	Brugsanvisning	26
SV	Bruksanvisning	29
FI	Käyttöohjeet	33
NL	Gebruiksaanwijzing	36
PL	Instrukcja użytkowania	40
KO	사용 설명서	44

SYMBOLS

Medical Device

Magnetic Resonance (MR) safe

INTENDED USE

The device is intended to provide gross immobilization to the cervical spine.

The device must be fitted and adjusted by a healthcare professional. X-ray and CT lucent.

Indications for use

Conditions requiring gross immobilization of the infant cervical spine.

This may include:

- C-Spine precaution
- Immobilization for pre and post c-spine surgery

Contraindications

- Patients with a compromised airway or known spinal deformities such as ankylosing spondylitis.
- Patients with penetrating trauma injuries.

Warnings and Cautions:

Warning: If an unstable fracture is suspected or unknown, with or without a sustained trauma, ensure additional spinal precautions are implemented to immobilize the spine.

Warning: Use of a cervical collar may increase intracranial pressure (ICP) through jugular venous compression.

Warning: Use of a cervical collar may increase complexity of airway management.

Warning: Cervical spine immobilization, including use of a cervical collar, has been associated with:

- Impaired respiratory effort and forced expiratory volume
- Pneumonia
- Aspiration
- Worsening of existing cervical spine injury
- Severe neurological deterioration in patients with ankylosing spondylitis
- Triggering of non-compliance or agitation
- Delay to definitive treatment
- Impaired physical examination/secondary survey

Caution: Cervical spine immobilization, including use of a cervical collar, has been associated with:

- Increased pain and discomfort, which may lead to increased spinal movement
- Impaired swallowing

Caution: Perfumes and harsh cleaning agents should not be used on or under the collar, as they may compromise the integrity of materials.

Caution: Regular cleaning of the collar, pads and the skin beneath along with inspection of the skin for any signs of irritation is required to reduce

the risk of skin ulceration. Bedridden patients are at an increased risk for skin ulceration.

Precautions:

- At least two persons are needed for initial collar placement: one to maintain the patient's head and neck in proper alignment, the second to fit the collar.
- Collar removal or any adjustments to the collar should be done with physician permission only and according to the physician's instructions.
- Unless otherwise specified by the physician, the patient should not remove the collar except to wash under it and change the pads.
- The patient will need a second person to help when cleaning the collar.

These instructions do not supersede hospital protocol and/or direct orders of the patient's healthcare professional.

GENERAL SAFETY INSTRUCTIONS

The healthcare professional should inform the patient about everything in this document that is required for safe use of this device.

Any serious incident in relation to the device must be reported to the manufacturer and relevant authorities.

The patient should immediately contact a healthcare professional:

- If there is a change or loss in device functionality, or if the device shows signs of damage or wear hindering its normal functions.
- If any pain, skin irritation, or unusual reaction occurs with the use of the device.

The device is for single patient – multiple use.

FITTING INSTRUCTIONS

Size Selection

The device's size is based upon Broselow-Luten color-code for the age 0-3 months.

At least two people are needed for initial device placement: to maintain the patient's head, neck, and airway in proper alignment; and one to apply and fit the brace.

Ensure that correct spinal protocol is followed.

Supine Application

1. Using proper c-spine precautions, position infant with arms to the side, shoulders down, and head aligned centrally (nose over belly button) (**Fig. 1**).
2. Carefully log roll the infant and center the device behind the head and torso (**Fig. 2**). Ensure the padding extends beyond all edges of the plastic.
3. Apply the head strap. Feed the strap through the slot on the opposite side of the patient's forehead and then back upon itself to attach the end to the hook and loop fastener (**Fig. 3**).
4. The length of the leg strap is adjustable to accommodate varying body lengths. To adjust, lift up the blue pad at bottom of the device and unbutton the bottom strap. Refasten using the appropriate hole for the patient's body length (**Fig. 4**).
5. Secure the leg and waist straps in three steps:
 - (1) Hook and loop Waist Strap;
 - (2) Leg Strap: Attach hook and loop to hook on Waist Strap;
 - (3) Hook and loop Waist Strap. The scratchy hook should always face away from the patient's skin (**Fig. 5**).

For additional stability, the device may be used in conjunction with the Miami Jr.® P0 size Front.

6. To apply the Miami Jr front, Flare the sides of the front out, slide it up the chest wall and scoop it up under the patient's chin. Sides of the collar front should be oriented up, off the trapezius, and toward the ears. **(Fig. 6)**.
7. While holding the Front securely, curl the ends snugly against the patient's neck **(Fig. 7)**. With the front of the collar inside the device, apply the hook and loop strap and secure the opposite side of the collar in the same fashion. Tighten straps alternately to an equal length on both sides.
8. Hook and Loop Straps must be aligned to the Front adhesive sections **(Fig. 8)**. Excess may be trimmed.

Note: Device application must be secure to ensure proper fit, maintain alignment, and prevent the chin from slipping inside.

Final Fitting Checklist

- No plastic touching the skin. Pads extend beyond all plastic edges **(Fig. 12)**.
- No slack or gaps in any of the straps.

If the device is used in conjunction with the Miami Jr.® P0 size Front **(Fig. 13)**

- The collar extends from the mandible to just below the sternal notch.
- Chin is centered comfortably in Chin Support. Chin should not extend over the edge of the pad, nor fall inside the collar.
- Hook and Loop Straps are equal in length and aligned to Front adhesive sections.
- Sides of Back should overlap sides of Front.
- Front of collar angled up toward ears. The lower plastic edge should not be resting on the patient's clavicles nor digging into the trapezius.
- Tracheal opening and posterior vent are midline.
- Front of collar angled up toward ears. The lower plastic edge should not be resting on the patient's clavicles nor digging into the shoulders.

Device removal

Before removing the device, note or mark the ends of the hook and loop straps with a pen. Then when the device is re-applied, the straps should be returned to the same position.

Undo the collar straps, and remove the collar **(Fig. 9)**.

With the patient on the back, undo the straps in the sequence shown in **(Fig. 10)**.

Hold the device Front in place and gently roll the patient onto the side, and slide the device carefully away behind the patient's neck and back **(Fig. 11)**.

Accessories and Replacement Parts

Please refer to the Össur catalog for a list of available replacement parts or accessories.

USAGE

Cleaning and care

Remove the pads from the plastic shells **(Fig. 14)**

- Hand-wash using mild detergent and rinse thoroughly.
- Air dry.
- Replace the pads by attaching the grey/dull side to the hook on the inside of the plastic shells.

Plastic shells:

- Wipe with a damp cloth and mild soap.

Note: Do not machine-wash, tumble dry, iron, bleach, or wash with fabric softener.

Note: Avoid contact with salt water or chlorinated water. In case of contact, rinse with fresh water and air dry.

DISPOSAL

The device and packaging must be disposed of in accordance with respective local or national environmental regulations.

LIABILITY

Össur does not assume liability for the following:

- Device not maintained as instructed by the instructions for use.
- Device assembled with components from other manufacturers.
- Device used outside of recommended use condition, application, or environment.

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