






Instructions for Use

PASO KNEE

		3
EN	Instructions for Use	5
DE	Gebrauchsanweisung	10
FR	Notice d'utilisation	15
ES	Instrucciones para el uso	20
IT	Istruzioni per l'uso	25
NO	Bruksanvisning	30
DA	Brugsanvisning	34
SV	Bruksanvisning	39
EL	Οδηγίες Χρήσης	44
FI	Käyttöohjeet	49
NL	Gebruiksaanwijzing	54
PT	Instruções de Utilização	59
PL	Instrukcja użytkowania	64
TR	Kullanım Talimatları	69
RU	Инструкция по использованию	74
JA	取扱説明書	79
ZH	中文说明书	83
KO	사용 설명서	87

ENGLISH

MD

Medical Device

DESCRIPTION

The device is a pneumatic polycentric knee joint with a geometric locking system and swing phase control that adapts automatically to walking speed.

The device has a distal tube connector and is supplied with one of these four proximal adapter options (**Fig. 1**). Only replace the adapter with one of these adapter options:

- Pyramid Adapter (A)
- 36 mm Thread (B) for use with 3-Prong Adapter
- Loop Adapter (C)
- IKF Adapter (D)

INTENDED USE

The device is intended as a part of a prosthetic system that replaces knee function of a missing lower limb.

Suitability of the device for the prosthesis and the patient must be evaluated by a healthcare professional.

The device must be fitted and adjusted by a healthcare professional.

Indications For Use and Target Patient Population

- Lower limb loss, amputation, or deficiency
- No known contraindications

The device is for moderate to high impact use, e.g., walking and occasional running.

The weight limit for the device is 136kg.

GENERAL SAFETY INSTRUCTIONS

Warning: Use of a lower limb prosthetic device carries an inherent risk of falling which may lead to injury.

The healthcare professional should inform the patient about everything in this document that is required for safe use of this device.

Warning: If there is a change or loss in device functionality, or if the device shows signs of damage or wear hindering its normal functions, the patient should stop using the device and contact a healthcare professional.

Warning: Avoid placing hands or fingers near moving joints.

Caution: Do not adjust other screws than described in these instructions.

The device is for single patient use.

ASSEMBLY INSTRUCTIONS

Warning: Risk of structural failure. Components from other manufacturers have not been tested and may cause excessive load on the device.

36 mm Thread for use with 3-Prong Socket Adapter

Please refer to 3-Prong Socket Adapter Instructions for use.

After lamination, screw the threaded adapter (**Fig. 1: B**) completely into the 3-prong adapter and tighten the lock screw to 10 Nm. After adjustments, the lock screw must be secured with medium strength threadlocker.

Caution: When prong adapters are used, resin may get under the adapter and hinder full thread engagement. Carefully remove resin without damaging the adapter so that adapter threads can be fully engaged when mounted on the device.

Caution: When carrying out any later adjustments, do not unscrew the threaded adapter more than one full turn. Threaded adapter might break.

Loop Adapter (Fig. 1)

Lamination of adapter to socket:

- Create a dummy disc to replace the adapter disc during lamination. Fix it to the loop adapter (E).
- Laminate the loop adapter (E) to the socket.
- Protect adapter threads during lamination.
- Remove the dummy disc and attach adapter disc (F) before attaching the loop adapter to the knee.

Caution: Fill the central thread of the loop adapter completely with Plastilin, so that no resin can run into the thread during lamination and the central screw can later be completely screwed in.

Caution: Check that the fine serrations on the inner surfaces of the adjustment discs face towards each other.

All Adapters (Fig. 1)

1. Fit the Washers (L and M) onto the Central Screw (I, J or K),
 - match the conical surfaces together.
 - match the fine serrations surfaces together.
2. Mount the adapter on the knee with the Central Screw, using medium strength threadlocker.
3. Tighten the screw to 35 Nm.

Warning: Screws for other adapters may have different lengths. Do not mix the central screws. Check screw for full thread engagement.

Caution: The IKF adapter has reference marks A and P for Anterior and Posterior (G). Make sure to position correctly.

ALIGNMENT INSTRUCTIONS

Bench Alignment (Fig. 2)

Alignment Goal

Alignment reference line (B) should:

- pass through midpoint of socket at the ischial tuberosity level (D)
- pass through midtube (A)
- fall at the 1/3 mark on the inside of the foot cover.

Note: Prioritize knee alignment over foot alignment if there is a mismatch.

Alignment Instructions

1. Position the foot so that the alignment reference line (B) falls at the 1/3 mark on the inside of the foot cover (with the foot cover and shoe on). Consider the external rotation of the foot.
2. Use the applicable adapters to connect the knee to the foot and establish the correct knee center height.
3. Position the knee so that the alignment reference line passes through midtube (A)
4. On the lateral side of the socket, make a first mark at the midpoint of the socket at the ischial tuberosity level (D). Make a second mark at the midpoint of the socket distally (E). Draw a line through both marks.

5. Position the socket so the alignment reference line (B) goes through first mark at the midpoint of the socket at the ischial tuberosity level (D).
6. Adjust the socket flexion to 5° in addition to the existing position (i. e., hip flexion contracture) and set the height of the full prosthesis.
7. Use the applicable adapters to connect the knee to the socket.

Warning: After adjustments, all screws must be secured with a medium strength threadlocker and tightened with the correct torque.

Caution: Adapters used on the distal connection must be cut straight and inserted down to the end stop of the device's tube receiver. No spacer should be used (**Fig. 3**).

Tighten the screws with the following torque:

- Central Screw (**Fig. 1: I, J, or K**): 35 Nm
- Tube Clamp Screw (**Fig. 4: A**): 16 Nm

Static Alignment

- Make sure the patient stands with equal weight on both legs.
- Check for correct prosthesis length.
- Check internal / external rotation.
- Check for correct load on toe and heel.

Dynamic Alignment

Ensure the patient is familiar with the functioning of the device.

Ensure the patient can activate and deactivate the geometric lock.

Adjustment of stance phase stability (**Fig. 5 and 6**)

Stance phase stability can be adjusted by shifting the knee or by using the supplied wedges to change axle geometry according to the following:

To increase stance phase stability:

- shift knee backwards compared to socket (5A) or
- remove wedge (6B)

To allow for a more dynamic gait:

- shift knee forward compared to socket (5B) or
- Add wedge A (6B) or
- Add wedge B (most dynamic)

To change or remove the wedge (6B), fully flex the knee and loosen the screws (6A). After changing the wedge, apply medium strength threadlocker to the screws and tighten by hand.

For IKF adapter only: Stance flexion adjustment

The back bumper (**Fig. 1: H**) can be changed to increase or decrease stance flexion:

Back Bumper (with holes)		
Yellow	Soft	More stance phase damping
Red	Medium	Standard
Brown	Hard	Less stance phase damping

To replace a bumper, unscrew the central screw. Replace the bumper (holes of back bumper facing up) and reassemble according to assembly instructions above (torque to 35 Nm).

USAGE

Knee cap removal and mounting

To remove the knee cap, extend the knee and loosen the two screws inside the bottom of the knee cap and pull straight down (**Fig. 7**). When mounting the knee cap, use medium strength threadlocker to secure screws.

Cleaning and care

Wipe the device with a soft cloth. Do not use solvents.

Environmental Conditions

The device should not come in contact with fresh water, salt water or chlorinated water.

Caution: Device should not be used in dusty environment. Exposure to sand, talcum, or similar should be avoided.

The device can be used in temperatures between -10 °C to 40 °C.

MAINTENANCE

The device and the overall prosthesis should be examined by a healthcare professional. Interval should be determined based on patient activity. Recommended interval is every 6 months.

Check for damage, excessive wear and dirt.

Wipe the device with a soft cloth moistened with small amount of general purpose oil or sewing machine oil.

Lubricate piston rod with silicone spray or similar.

Please change the 3 filter foams of the pneumatic cylinder every 12 months

Caution: Do not use compressed air to clean device. Air forces pollutants into bearings and may cause malfunctions and wear.

Caution: Do not use talcum powder for lubrication.

REPORT OF SERIOUS INCIDENT

Any serious incident in relation to the device must be reported to the manufacturer and relevant authorities.

DISPOSAL

The device and packaging must be disposed of in accordance with respective local or national environmental regulations.

LIABILITY

Össur prosthetic devices are designed and verified to be safe and compatible in combination with each other and custom-made prosthetic sockets with Össur adapters, and when used in accordance with their intended use.

Össur does not assume liability for the following:

- Device not maintained as instructed by the instructions for use.
- Device assembled with components from other manufacturers.
- Device used outside of recommended use condition, application, or environment.

Compliance

This device has been tested according to ISO 10328 standard to three million load cycles.

Depending on patient activity, this may correspond to 3-5 years of use.

ISO 10328 - P6 - 136 kg *) 
*) Body mass limit not to be exceeded!  For specific conditions and limitations of use see manufacturer's written instructions on intended use!

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