• Post-operative, trauma or rehab use  
No known contraindications.

#### Warnings and Cautions:

Do not remove the device unless instructed by a healthcare professional.

- Caution should be used if prescribed for patients with diabetes, vascular deficiency, and neuropathy;
- The skin should be checked regularly
- This device should not be applied on patients who are not able to communicate discomfort.
- Patient should use caution when walking on slippery or wet surfaces to avoid injury.
- Ensure the patient performs circulation checks. If a loss of circulation is felt, or if the walker feels too tight, release contact closure straps and adjust to a comfortable level. If discomfort continues, the patient should contact their healthcare professional immediately.
- Do not remove the inner foam pads which are adhered to the plastic walker boot.

#### GENERAL SAFETY INSTRUCTIONS

The healthcare professional should inform the patient about everything in this document that is required for safe use of this device.

Any serious incident in relation to the device must be reported to the manufacturer and relevant authorities.

The patient should immediately contact a healthcare professional:

- If there is a change or loss in device functionality, or if the device shows signs of damage or wear hindering its normal functions.
- If any pain, skin irritation, or unusual reaction occurs with the use of the device.

The device is for single patient – multiple use.

#### FITTING INSTRUCTIONS

While carrying out the following instructions, please refer to the overview figure for locating components mentioned in the text (Fig. 1).

#### Device Application

**Note:** Do not remove Pull Tabs (A) until further instructed.

1. Remove Liner (B) by opening Contact Closure Straps (C) and lifting it out of the device.
2. Place lower leg and foot into the Liner. Close the Liner snugly over the leg (Fig. 2).
3. Insert foot into the device by gently spreading the Struts (D) outwards (Fig. 3). Ensure Struts are centered over the ankle bone (malleoli). The lower leg should be kept at a 90° angle with the foot and ensure that no gap is in between the Liner heel and the Strut (Fig. 4).
4. Remove Pull Tabs located in between the Struts and the Liner. This allows the Liner to be secured into position. Pulling the tabs in a forward/downward direction makes them easier to remove.

5. To accommodate for rubbing or spacing, additional Foam Pads (E) may be placed inside the Liner, between the Liner and the plastic, or under the Contact Closure Straps for additional comfort (Fig. 5).
6. Fasten the Contact Closure Straps starting at the ankle, foot and working upward. Strap adjustments may be necessary with changes in lower limb volume.

#### Device Removal

Undo the straps in any order and attach the straps back to themselves. This will make it easier when fitting the device next time and will prolong the life of the straps. Open the Liner to lift the foot and lower leg out of the device.

#### USAGE

##### Cleaning and care

Washing the device with the soft goods detached allows for more thorough cleaning.

##### Washing Instructions

- Hand-wash using mild detergent and rinse thoroughly.
- Air dry

**Note:** Do not machine-wash, tumble dry, iron, bleach, or wash with fabric softener.

**Note:** Avoid contact with salt water or chlorinated water. In case of contact, rinse with fresh water and air dry.

#### DISPOSAL

The device and packaging must be disposed of in accordance with respective local or national environmental regulations.

#### LIABILITY

Össur does not assume liability for the following:

- Device not maintained as instructed by the instructions for use.
- Device assembled with components from other manufacturers.
- Device used outside of recommended use condition, application, or environment.

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**MD** Medical Device

#### INTENDED USE

The device is intended for immobilization of the foot and ankle

The device must be fitted and adjusted by a healthcare professional.

#### Indications for use

Foot and ankle conditions that may benefit from motion restriction, such as:

- Grade 2 and 3 ankle sprains
- Stable fractures
- Soft tissue injuries

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